



República de Colombia  
Ministerio de Salud y Protección Social  
Instituto Nacional de Vigilancia de Medicamentos y Alimentos – INVIMA

**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**  
**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

**ANTECEDENTES**

Que mediante escrito número 20211100250 de fecha 24 de mayo de 2021, la Doctora Viviana Carolina Lemus Blanco, actuando en calidad de Representante Legal de la empresa NUVASIVE COLOMBIA S.A.S, solicitó Registro Sanitario para el producto PEDICLE SCREW FIXATION SYSTEM - SISTEMA DE FIJACIÓN DE TORNILLOS PEDICULARES a favor de NUVASIVE COLOMBIA S.A.S, con domicilio en BOGOTÁ D.C., en la modalidad de IMPORTAR y VENDER.

Que mediante escrito No. 20211132611 de fecha 08 de julio de 2021, la Doctora Viviana Carolina Lemus Blanco, actuando en calidad de Representante Legal de la empresa NUVASIVE COLOMBIA S.A.S, allegó anexo al expediente en el sentido de aportar el formulario de solicitud ajustado, incluyendo las indicaciones y usos y complementando las referencias.

Que mediante Auto No. 2021009403 de fecha 19 de julio de 2021, se le informo al interesado que una vez revisada el expediente y para continuar con el trámite solicitado, debe cumplir con los siguientes requerimientos:

- 1. Una vez estudiada la documentación allegada con su solicitud de Registro Sanitario se identifica que pretende agrupar varios dispositivos médicos en un mismo registro sanitario, lo cual puede ser posible siempre y cuando se dé cumplimiento con el artículo 28 del Decreto 4725 de 2005. Lo anterior dado que se identifican varios sistemas independientes (RELINE, ARMADA, VUEPOINT, PRECEPT, X-CORE, entre otros), cada uno con su instrumental quirúrgico asociado, como lo refiere el Certificado de Venta Libre, las Declaraciones de Conformidad y la información técnica aportada. Si bien en el Acta 7 de 2018 de la sala especializada de Dispositivos Médicos, se menciona que es posible agrupar el material de osteosíntesis por áreas anatómicas siendo uno de estas Columna vertebral, se debe mantener las premisas establecidas en el Artículo 28, por sí solo un sistema agrupa dispositivos médicos diferentes para que actúen de forma interdependiente y en este caso los diferentes sistemas no disponen de implantes comunes para que puedan ser agrupados, adicionalmente cada uno dispone de indicaciones de uso específicas. En este sentido, se debe allegar formulario de solicitud de registro sanitario corregido en la sección del nombre del producto, nombre genérico y referencias, en donde se incluya solamente uno de los sistemas junto con sus respectivas referencias del cual se compone y se excluyan las demás referencias que no hacen parte del mismo.*
- 2. Teniendo en cuenta el primer requerimiento, se debe allegar carta emitida por el fabricante en la que se aclare el nombre del producto, por cuanto en el formulario describe "Orthopedic Manual Surgical Implants and instruments - Implantes e instrumentos quirúrgicos manuales ortopédicos" y como nombre genérico "Pedicule Screw fixation- Fijación de tornillos pediculares", lo cual no tiene relación. Se aclara que el nombre genérico debe ser descriptivo del nombre del producto.*
- 3. Teniendo en cuenta el artículo 28 del Decreto 4725 de 2005 y el requerimiento No. 1 del presente auto, se debe allegar declaración emitida por el fabricante en la que aclare y justifique que el producto (sistema a amparar con todas las referencias del sistema) corresponde a un SISTEMA, indicando el nombre del producto y nombre genérico del mismo, especificando los modelos y/o referencias que hacen parte del sistema (con su respectiva descripción). Cabe señalar, que un sistema es el "conjunto de dispositivos INTERDEPENDIENTES que interactúan entre sí para realizar una actividad, que solo pueden adaptarse al producto al que pertenecen y que ninguno de ellos puede ser reemplazado por otro. Cabe anotar que si alguno de los componentes se puede emplear solo o con otro conjunto de productos ya no se considera un sistema".*
- 4. Allegar formulario de solicitud de registro sanitario corregido en el ítem vida útil, toda vez que no se evidencia valor alguno y en la información allegada de los folios 945 a 996 respecto a la vida útil del producto se declaran cinco (5) años para el producto estéril.*
- 5. Allegar el formulario de solicitud de registro sanitario corregido en la sección de indicaciones y usos, donde concuerde con lo descrito en la descripción del dispositivo médico del folio 565 al 592, toda vez que*



República de Colombia  
Ministerio de Salud y Protección Social  
Instituto Nacional de Vigilancia de Medicamentos y Alimentos – INVIMA

**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**  
**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

*no fue diligenciada esta información del formulario. Lo anterior, de conformidad del literal d) del artículo 18 del decreto 4725 del 2005. Tener en cuenta que las indicaciones y usos deben corresponder con el sistema que se desee amparar en el registro sanitario acorde al requerimiento No. 1 del presente auto.*

*6. Allegar el formulario de solicitud de registro sanitario corregido en la sección de componentes y composición, toda vez que fue mal diligenciada la información en el campo "partes que componen el Dispositivo Médico", y en el campo "composición cualitativa", teniendo en cuenta que las partes que componen el Dispositivo Médico se refiere al componente del producto (p. ej. Tornillo) y la composición cualitativa se refiere al material (p. ej. Titanio) con que está elaborado. Tener en cuenta que la información debe corresponder con la descripción del dispositivo médico del folio 565 al 592.*

*7. Allegar documento de descripción del dispositivo médico para el producto que se desea amparar en el registro sanitario en el cual se evidencie, los componentes con su respectiva composición adicionalmente al nombre del producto, indicaciones, contraindicaciones, advertencias, componentes principales, accesorios, relación con pacientes; con traducción al español, para el sistema que se desee amparar acorde al requerimiento No. 1 del presente auto.*

*8. Teniendo en cuenta el requerimiento No. 1, se debe allegar etiquetas del fabricante para el sistema que se desee amparar y Sticker de importador con el nombre del producto que se va a amparar en el Registro sanitario.*

*9. Teniendo en cuenta el requerimiento No. 6 del presente auto, en caso que se mantengan los componentes descritos en el formulario, se debe allegar el desarrollo de las pruebas de evaluación biológica del producto (Citotoxicidad, Sensibilidad, Irritación/Reactividad intercutánea, estudio de toxicidad subcrónica (subaguda) -hemocompatibilidad, genotoxicidad, implantación y demás) que le apliquen, con su respectiva traducción al castellano acorde al artículo 49 del Decreto 4725 de 2005. Cuando sean tecnologías suficientemente demostradas se podrá cumplir con este requisito aportando evidencia científica que se encuentre indexada o publicada en revistas reconocidas a nivel mundial en las que se relacione la evaluación biológica del material específico para el dispositivo médico objeto de la solicitud del registro sanitario. Lo anterior se solicita, por cuanto en la información aportada como evaluación biológica de los folios 1004 a 1049, solo corresponde al material titanio y aleación cromo cobalto, haciendo falta para el material PEEK del cual se pueden componer varias referencias dependiendo además del sistema que sea elegido acorde al requerimiento No. 1, y del cual no se presenta información.*

*10. Teniendo en cuenta el requerimiento No. 1 del presente auto, en caso de no dejar el sistema RELINE se requiere allegar el desarrollo del análisis de riesgos del producto, en donde se describan cuáles son los riesgos puntuales del uso del producto encontrados por el fabricante durante el proceso de diseño y manufactura, sus causas, severidad, ocurrencia y detectabilidad de dichos riesgos, así como las soluciones planteadas por el fabricante para mitigar cada uno de ellos (tabla de evaluación y control del riesgo), acorde al literal j del artículo 18 del Decreto 4725 de 2005 con su respectiva traducción al castellano acorde al artículo 49 del Decreto 4725 de 2005. Lo anterior se solicita, por cuanto en el análisis de riesgos allegado en los folios 1050 a 1173, corresponde únicamente al sistema RELINE.*

*11. Teniendo en cuenta el requerimiento No. 1 del presente auto, se debe allegar documento emitido por el fabricante en el sentido de indicar si el producto dispone de Alertas Sanitarias involucradas con el dispositivo médico objeto de la solicitud del Registro Sanitario y en caso tal que disponga de ellas, los motivos y las gestiones realizadas al respecto, por cuanto en la información aportada no se evidencia si el producto ha presentado o no alertas sanitarias.*

*12. Teniendo en cuenta el requerimiento No. 1 del presente auto, se debe allegar en idioma original con su respectiva traducción al castellano acorde al artículo 49 del Decreto 4725 de 2005, los estudios clínicos (ensayos clínicos/estudios multicéntricos) realizados en pacientes, publicados en revistas indexadas con el mayor nivel de evidencia disponible para demostrar la seguridad y efectividad del producto para dar cumplimiento al literal k) del artículo 18 del decreto 4725 del 2005 para el sistema que se desee amparar.*

Página 2 de 164



República de Colombia  
Ministerio de Salud y Protección Social  
Instituto Nacional de Vigilancia de Medicamentos y Alimentos – INVIMA

**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**  
**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

*Se pueden entregar estudios clínicos publicados de tecnologías similares o equivalentes a la presentada en la solicitud del registro sanitario. Lo anterior se solicita por cuanto la información allegada en los folios 1177 a 1539, corresponden con el sistema PEDICLE SCREW SYSTEM - SISTEMA DE TORNILLOS PEDICULARES.*

*13. Allegar Certificado de Venta Libre - CVL debidamente apostillado o consularizado y legalizado, toda vez que el aportado, no cumple los requisitos que debe reunir el CVL consagrados en el Decreto 4725 de 2005. Lo anterior dado que visualizamos que el Apostille aportado certifica la Firma de Ariani Carrillo en calidad de Notario Público del Estado de Florida, quien da fe de que el CVL que nos ocupa no ha sido adulterado y es una copia exacta, pero no se visualiza que se haya aportado un apostille del CVL radicado al expediente, donde se certifique la firma y las calidades de quien suscribe el CVL, CDR. Cesar A. Pérez, en calidad de Director de Asuntos Regulatorios de la FDA.*

Que mediante escrito No. 20211254942 radicado el 30 de noviembre de 2021, la Doctora Viviana Carolina Lemus Blanco, actuando en calidad de Representante Legal de la empresa NUVASIVE COLOMBIA S.A.S., allego respuesta al Auto No. 2021009403 de fecha 19 de julio de 2021.

#### CONSIDERACIONES

Que una vez realizado el estudio de la respuesta al auto de requerimiento allegada bajo el radicado No. 20211254942 radicado el 30 de noviembre de 2021, se considera SATISFACTORIA por cuanto el interesado allegó los soportes y justificaciones respectivas aclarando que el producto es un solo sistema compuesto por las líneas *RELINE, ARMADA, PRECEPT, VUEPOINT* y *SPHERx*, dando respuesta satisfactoria al punto 1 del auto de requerimientos.

Para dar respuesta al 2, 3 y 4 punto del auto de requerimientos, el interesado allegó certificación y documentos aclaratorios emitida por el fabricante, donde se aclara lo solicitado, dando respuesta SATISFACTORIA a estos puntos.

Adicionalmente, se allega formulario de solicitud de registro sanitario corregido en las indicaciones y usos y componentes y composición, para dar respuesta a los puntos 5 y 6 del auto de requerimientos, considerándose que la respuesta a estos puntos es SATISFACTORIA.

En cuanto al punto 7, 8, 9 y 10, el interesado allegó descripción del dispositivo médico acorde a lo solicitado, etiquetas del fabricante, sticker de importador, estudios de biocompatibilidad y los análisis de riesgo faltantes considerándose que la respuesta a estos puntos es SATISFACTORIA.

Para dar respuesta SATISFACTORIA a los puntos 11 y 12 del auto de requerimientos, se allega por parte del interesado documento expedido por el fabricante donde indica que el producto no dispone de alertas sanitarias y los estudios clínicos realizados, considerándose SATISFACTORIA la respuesta a estos puntos del auto.

En cuanto al punto 13, el interesado allegó Certificado de Venta Libre donde se reúnen los requisitos del Decreto 4725 de 2005.

Que ante este Instituto se ha solicitado la Concesión de un Registro Sanitario con base en la documentación allegada inicialmente y posterior respuesta a Auto, previo estudio técnico y legal de la Dirección de Dispositivos Médicos y Otras Tecnologías del Invima, se emitió concepto favorable en cumplimiento de los requisitos preceptuados en el Decreto 4725 de 2005 modificado por el Decreto 582 de 2017, para la autorización de este Registro Sanitario.

En consecuencia, la Directora de Dispositivos Médicos y Otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos – INVIMA,

#### RESUELVE



República de Colombia  
**Ministerio de Salud y Protección Social**  
 Instituto Nacional de Vigilancia de Medicamentos y Alimentos – INVIMA

**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**  
**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

**ARTICULO PRIMERO.-** Conceder REGISTRO SANITARIO por el término de DIEZ (10) años a  
**PRODUCTO:** PEDICLE SCREW FIXATION SYSTEM - SISTEMA DE FIJACIÓN DE TORNILLOS PEDICULARES,  
**MARCA:** NUVASIVE ®  
**REGISTRO SANITARIO No.:** **INVIMA 2022DM-0025052**  
**TIPO DE REGISTRO:** IMPORTAR Y VENDER  
**TITULAR(ES):** NUVASIVE COLOMBIA S.A.S. con domicilio en BOGOTA - D.C.  
**FABRICANTE(S):** NUVASIVE, INC. con domicilio en ESTADOS UNIDOS DE AMÉRICA  
**IMPORTADOR(ES):** NUVASIVE COLOMBIA S.A.S. con domicilio en BOGOTA - D.C.  
**ACONDICIONADOR(ES):** NUVASIVE COLOMBIA S.A.S. con domicilio en BOGOTA - D.C.  
**TIPO DE DISPOSITIVO** INVASIVO  
**RIESGO:** IIb  
**COMPOSICIÓN:**

PARTES QUE COMPONEN EL DISPOSITIVO MEDICO	COMPOSICION CUALITATIVA
TORNILLOS	TITANIO (Ti-Al-64V) , TITANIO ANODIZADO, ALEACIÓN DE CROMO COBALTO (StarLoyC), ACERO INOXIDABLE
VARILLAS	TITANIO (Ti-Al-64V) , TITANIO ANODIZADO, ALEACIÓN DE CROMO COBALTO (StarLoyC ) , ACERO INOXIDABLE
GANCHOS	TITANIO (Ti-Al-64V) , TITANIO ANODIZADO, ALEACIÓN DE CROMO COBALTO (StarLoyC ) ACERO INOXIDABLE
CONECTORES	TITANIO (Ti-Al-64V) , TITANIO ANODIZADO, ALEACIÓN DE CROMO COBALTO (StarLoyC ) , ACERO INOXIDABLE
ANILLOS	TITANIO (Ti-Al-64V) , TITANIO ANODIZADO, ALEACIÓN DE CROMO COBALTO (StarLoyC ) , ACERO INOXIDABLE

**USOS:** SISTEMA DE IMPLANTES PARA COLUMNA VERTEBRAL: DESTINADO A SER UTILIZADO COMO UN SISTEMA DE FIJACIÓN DE TORNILLO PEDICULAR EN PACIENTES ESQUELÉTICAMENTE MADUROS, ESTÁ DESTINADO A PROPORCIONAR INMOVILIZACIÓN Y ESTABILIZACIÓN DE SEGMENTOS ESPINALES COMO COMPLEMENTO DE LA FUSIÓN EN EL TRATAMIENTO DE LAS SIGUIENTES INESTABILIDADES AGUDAS Y CRÓNICAS O DEFORMIDADES DELA UNIÓN CRANEOCERVICAL, LA REGIÓN POSTERIOR TORÁCICA, LUMBAR Y COLUMNA SACRA:  
 1. ENFERMEDAD DEGENERATIVA DE DISCO (DEFINIDA COMO DOLOR DE ESPALDA DE ORIGEN DISCOGÉNICO CON DEGENERACIÓN DEL DISCO CONFIRMADA POR LA HISTORIA CLÍNICA DEL PACIENTE Y ESTUDIOS RADIOGRÁFICOS)  
 2. ESPONDILOLISTESIS DEGENERATIVA CON LA PRUEBA OBJETIVA DE TRASTORNOS NEUROLÓGICOS  
 3. FRACTURA  
 4. DISLOCACIÓN



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

- 5. ESCOLIOSIS
- 6. CIFOSIS
- 7. TUMOR ESPINAL Y/O
- 8. FUSIÓN PREVIA FALLIDA (SEUDOARTROSIS)

EN EL CASO DE LA FIJACIÓN NO CERVICAL POSTERIOR DE TORNILLOS EN PACIENTES PEDIÁTRICOS, EL SISTEMA ESPINAL ESTÁ INDICADO COMO UN COMPLEMENTO DE FUSIÓN PARA TRATAR LA ESCOLIOSIS IDIOPÁTICA ADOLESCENTE. ADEMÁS, ESTÁ DESTINADO A TRATAR PACIENTES PEDIÁTRICOS DIAGNOSTICADOS CON LAS SIGUIENTES DOLENCIAS: ESPONDILOLISTESIS /ESPONDILOLISIS, Y FRACTURA CAUSADA POR TUMOR Y/O TRAUMA. LA FIJACIÓN PEDIÁTRICA PEDICULAR DE TORNILLOS SE LIMITA A UN MÉTODO POSTERIOR Y ESTÁ DESTINADA A USARSE CON AUTOINJERTO Y/O ALOINJERTO.

EL SISTEMA ESPINAL TAMBIÉN ESTÁ INDICADO PARA EL TRATAMIENTO DE LA ESPONDILOLISTESIS (GRADOS 3 Y 4) DE LA ARTICULACIÓN VERTEBRAL EN PACIENTES CON MADUREZ ESQUELÉTICA QUE RECIBAN TRATAMIENTO DE FUSIÓN MEDIANTE INJERTO ÓSEO AUTÓGENO, EN QUIENES EL DISPOSITIVO SE FIJE O UNA A LA COLUMNA LUMBAR Y LA SACRA, CON EXTRACCIÓN DE LOS IMPLANTES UNA VEZ LOGRADA UNA FUSIÓN SÓLIDA.

CUANDO SE LO UTILIZA COMO SISTEMA ANTEROLATERAL DE TORNILLO NO PEDICULAR EN LA ESPINA TORÁCICA Y LUMBAR, EL SISTEMA ESPINAL TAMBIÉN ESTÁ INDICADO EN LOS CASOS SIGUIENTES:

- 1. ENFERMEDAD DEGENERATIVA DE DISCO (DEFINIDA COMO DOLOR DE ESPALDA DE ORIGEN DISCOGÉNICO CON DEGENERACIÓN DEL DISCO CONFIRMADA POR LA HISTORIA CLÍNICA DEL PACIENTE Y ESTUDIOS RADIOGRÁFICOS)
- 2. ESTENOSIS DE LA COLUMNA
- 3. ESPONDILOLISTESIS
- 4. DEFORMIDADES DE LA COLUMNA
- 5. FRACTURA
- 6. FUSIÓN PREVIA FALLIDA (SEUDOARTROSIS)
- 7. RESECCIÓN DE UN TUMOR

PRESENTACIONES COMERCIALES:  
OBSERVACIONES:

EMPAQUE POR UNIDAD, PAQUETES POR 2, 4, 6, 8, 10 UNIDADES  
BAJO ESTE REGISTRO SANITARIO SE AMPARAN LAS SIGUIENTES REFERENCIAS:

**CÓDIGO, MODELO O REFERENCIA**

**DESCRIPCIÓN**

1594520	MAS PLIF FIXATION SHANK, 4.5X20MM 3Z
1594525	MAS PLIF FIXATION SHANK, 4.5X25MM 3Z
1594530	MAS PLIF FIXATION SHANK, 4.5X30MM 3Z
1594535	MAS PLIF FIXATION SHANK, 4.5X35MM 3Z
1594540	MAS PLIF FIXATION SHANK, 4.5X40MM 3Z
1595020	MAS PLIF FIXATION SHANK, 5.0X20MM 3Z
1595025	MAS PLIF FIXATION SHANK, 5.0X25MM 3Z
1595030	MAS PLIF FIXATION SHANK. 5.0X30MM 3Z
1595035	MAS PLIF FIXATION SHANK, 5 0X35MM 3Z
1595040	MAS PLIF FIXATION SHANK. 5.0X40MM 3Z



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
1595520	MAS PLIF FIXATION SHANK, 5.5X20MM 3Z
1595525	MAS PLIF FIXATION SHANK, 5.5X25MM 3Z
1595530	MAS PLIF FIXATION SHANK, 5.5X30MM 3Z
1595535	MAS PLIF FIXATION SHANK, 5.5X35MM 3Z
1595540	MAS PLIF FIXATION SHANK, 5.5X40MM 3Z
1615030	MAS PLIF FIXATION SHANK, 5.0X30MM FC
1615035	MAS PLIF FIXATION SHANK, 5.0X35MM FC
1615530	MAS PLIF FIXATION SHANK, 5.5X30MM FC
1615535	MAS PLIF FIXATION SHANK, 5.5X35MM FC
1675540	PRECEPT SCREW, 5.5X40MM POLY C.R.
1675545	PRECEPT SCREW, 5.5X45MM POLY C.R.
1675550	PRECEPT SCREW, 5.5X50MM POLY C.R.
1676540	PRECEPT SCREW, 6.5X40MM POLY C.R.
1676545	PRECEPT SCREW, 6.5X45MM POLY C.R.
1676550	PRECEPT SCREW, 6.5X50MM POLY C.R.
1677500	PRECEPT SCREW, 7.5X100MM POLY C.R.
1677540	PRECEPT SCREW, 7.5X40MM POLY C.R.
1677545	PRECEPT SCREW, 7.5X45MM POLY C.R.
1677550	PRECEPT SCREW, 7.5X50MM POLY C.R.
1677590	PRECEPT SCREW, 7.5X90MM POLY C.R.
1678500	PRECEPT SCREW, 8.5X100MM POLY C.R.
1678580	PRECEPT SCREW. 8.5X80MM POLY C.R.
1678590	PRECEPT SCREW. 8.5X90MM POLY C.R.
7070000	MAS PLIF TULIP, MODULAR
7070002	MAS PLIF TI LOCK SCREW
7070003	MAS PLIF FIXATION LOCK SCREW
7074525	MAS PLIF SHANK, 4.5X25MM 3Z BLAST
7074530	MAS PLIF SHANK, 4.5X30MM 3Z BLAST
7074535	MAS PLIF SHANK, 4.5X35MM 3Z BLAST
7074540	MAS PLIF SHANK, 4.5X40MM 3Z BLAST
7075025	MAS PLIF SHANK, 5.0X25MM 3Z BLAST
7075030	MAS PLIF SHANK, 5.0X30MM 3Z BLAST
7075035	MAS PLIF SHANK, 5.0X35MM 3Z BLAST
7075040	MAS PLIF SHANK, 5.0X40MM 3Z BLAST
7075525	MAS PLIF SHANK, 5.5X25MM 3Z BLAST
7075530	MAS PLIF SHANK, 5.5X30MM 3Z BLAST
7075535	MAS PLIF SHANK, 5.5X35MM 3Z BLAST
7075540	MAS PLIF SHANK, 5.5X40MM 3Z BLAST
7076530	MAS PLIF SHANK, 6.5X30MM 3Z BLAST



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
7076535	MAS PLIF SHANK, 6.5X35MM 3Z BLAST
7076540	MAS PLIF SHANK, 6.5X40MM 3Z BLAST
7076545	MAS PLIF SHANK, 6.5X45MM 3Z BLAST
7094525	MAS PLIF SHANK, 4.5X25MM 3Z
7094530	MAS PLIF SHANK, 4.5X30MM 3Z
7094535	MAS PLIF SHANK, 4.5X35MM 3Z
7094540	MAS PLIF SHANK, 4.5X40MM 3Z
7094545	MAS PLIF SHANK, 4.5X45MM 3Z
7095025	MAS PLIF SHANK, 5.0X25MM 3Z
7095030	MAS PLIF SHANK, 5.0X30MM 3Z
7095035	MAS PLIF SHANK, 5.0X35MM 3Z
7095040	MAS PLIF SHANK. 5.0X40MM 3Z
7095045	MAS PLIF SHANK, 5.0X45MM 3Z
7095050	MAS PLIF SHANK, 5.0X50MM 3Z
7095525	MAS PLIF SHANK, 5.5X25MM 3Z
7095530	MAS PLIF SHANK, 5.5X30MM 3Z
7095535	MAS PLIF SHANK, 5.5X35MM 3Z
7095540	MAS PLIF SHANK, 5.5X40MM 3Z
7095545	MAS PLIF SHANK, 5.5X45MM 3Z
7095550	MAS PLIF SHANK, 5.5X50MM 3Z
7095555	MAS PLIF SHANK, 5.5X55MM 3Z
7096525	MAS PLIF SHANK, 6.5X25MM 3Z
7096530	MAS PLIF SHANK, 6.5X30MM 3Z
7096535	MAS PLIF SHANK, 6.5X35MM 3Z
7096540	MAS PLIF SHANK, 6.5X40MM 3Z
7096545	MAS PLIF SHANK, 6.5X45MM 3Z
7096550	MAS PLIF SHANK, 6.5X50MM 3Z
7096555	MAS PLIF SHANK, 6.5X55MM 3Z
7124525	MAS PLIF SCREW, 4.5X25MM POLYAXIAL 3Z
7124530	MAS PLIF SCREW, 4.5X30MM POIYAXIAI 3Z
7124535	MAS PLIF SCREW, 4.5X35MM POLYAXIAL 3Z
7124540	MAS PLIF SCREW, 4.5X40MM POLYAXIAL 3Z
7125025	MAS PLIF SCREW, 5.0X25MM POLYAXIAL 3Z
7125030	MAS PLIF SCREW, 5.0X30MM POLYAXIAL 3Z
7125035	MAS PLIF SCREW, 5.0X35MM POLYAXIAL 3Z
7125040	MAS PLIF SCREW, 5.0X40MM POLYAXIAL 3Z
7125525	MAS PLIF SCREW, 5.5X25MM POTYAXIAL 3Z
7125530	MAS PLIF SCREW, 5.5X30MM POLYAXIAL 3Z
7125535	MAS PLIF SCREW. 5.5X35MM POLYAXIAL 3Z



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
7125540	MAS PLIF SCREW, 5.5X40MM POLYAXIAL 3Z
7242020	MAS PLIF FIXATION ROD, 20MM PRE-BENT
7242025	MAS PLIF FIXATION ROD, 25MM PRE-BENT
7242030	MAS PLIF FIXATION ROD, 30MM PRE-BENT
7242035	MAS PLIF FIXATION ROD, 35MM PRE-BENT
7242040	MAS PLIF FIXATION ROD, 40MM PRE-BENT
7242045	MAS PLIF FIXATION ROD, 45MM PRE-BENT
7242050	MAS PLIF FIXATION ROD. 50MM PRE-BENT
7242055	MAS PLIF FIXATION ROD. 55MM PRE-BENT
7242060	MAS PLIF FIXATION ROD, 60MM PRE-BENT
7242065	MAS PLIF FIXATION ROD, 65MM PRE-BENT
7242070	MAS PLIF FIXATION ROD. 70MM PRE-BENT
7242075	MAS PLIF FIXATION ROD, 75MM PRE-BENT
7242080	MAS PLIF FIXATION ROD, 80MM PRE-BENT
7242085	MAS PLIF FIXATION ROD, 85MM PRE-BENT
7242090	MAS PLIF FIXATION ROD, 90MM PRE-BENT
7242095	MAS PLIF FIXATION ROD, 95MM PRE-BENT
7242100	MAS PLIF FIXATION ROD, 100MM PRE-BENT
7242110	MAS PLIF FIXATION ROD, 110MM PRE-BENT
7242120	MAS PLIF FIXATION ROD, 120MM PRE-BENT
7291020	MAS PLIF ROD, 20MM
7292020	MAS PLIF ROD, 20MM PRE-BENT
7292130	MAS PLIF ROD, 130MM PRE-BENT
7292140	MAS PLIF ROD, 140MM PRE-BENT
7292150	MAS PLIF ROD, 150MM PRE-BENT
7292300	MAS PLIF ROD, 300MM PRE-BENT
7884525	MAS PLIF SCREW, 4.5 X 25MM POLYAXIAL 3Z
7884530	MAS PLIF SCREW, 4.5 X 30MM POLYAXIAL 3Z
7884535	MAS PLIF SCREW. 4.5 X 35MM POLYAXIAL 3Z
7884540	MAS PLIF SCREW, 4.5 X 40MM POLYAXIAL 3Z
7885025	MAS PLIF SCREW, 5.0 X 25MM POLYAXIAL 3Z
7885030	MAS PLIF SCREW, 5.0 X 30MM POLYAXIAL 3Z
7885035	MAS PLIF SCREW, 5.0 X 35MM POLYAXIAL 3Z
7885040	MAS PLIF SCREW. 5.0 X 40MM POLYAXIAL 3Z
7885525	MAS PLIF SCREW. 5.5 X 25MM POLYAXIAL 3Z
7885530	MAS PLIF SCREW. 5.5 X 30MM POLYAXIAL 3Z
7885535	MAS PLIF SCREW. 5.5 X 35MM POLYAXIAL 3Z
7885540	MAS PLIF SCREW, 5.5 X 40MM POLYAXIAL 3Z
8800000	PRECEPT SET SCREW



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8804525	PRECEPT SCREW. 4.5X25MM POLYAXIAL
8804530	PRECEPT SCREW. 4.5X30MM POLYAXIAL
8804535	PRECEPT SCREW, 4.5X35MM POLYAXIAL
8804540	PRECEPT SCREW, 4.5X40MM POLYAXIAL
8804545	PRECEPT SCREW, 4.5X45MM POLYAXIAL
8805530	PRECEPT SCREW, 5.5X30MM POLYAXIAL
8805535	PRECEPT SCREW, 5.5X35MM POLYAXIAL
8805540	PRECEPT SCREW, 5.5X40MM POLYAXIAL
8805545	PRECEPT SCREW, 5.5X45MM POLYAXIAL
8805550	PRECEPT SCREW, 5.5X50MM POLYAXIAL
8805555	PRECEPT SCREW, 5.5X55MM POLYAXIAL
8806530	PRECEPT SCREW, 6.5X30MM POLYAXIAL
8806535	PRECEPT SCREW, 6.5X35MM POLYAXIAL
8806540	PRECEPT SCREW, 6.5X40MM POLYAXIAL
8806545	PRECEPT SCREW, 6.5X45MM POLYAXIAL
8806550	PRECEPT SCREW, 6.5X50MM POLYAXIAL
8806555	PRECEPT SCREW, 6.5X55MM POLYAXIAL
8807040	PRECEPT SCREW, 7.0X40MM POLYAXIAL
8807045	PRECEPT SCREW, 7.0X45MM POLYAXIAL
8807500	PRECEPT SCREW, 7.5X100MM POLYAXIAL
8807530	PRECEPT SCREW, 7.5X30MM POLYAXIAL
8807535	PRECEPT SCREW. 7.5X35MM POLYAXIAL
8807540	PRECEPT SCREW, 7.5X40MM POLYAXIAL
8807545	PRECEPT SCREW, 7.5X45MM POLYAXIAL
8807550	PRECEPT SCREW, 7.5X50MM POLYAXIAL
8807555	PRECEPT SCREW, 7.5X55MM POLYAXIAL
8807560	PRECEPT SCREW, 7.5X60MM POLYAXIAL
8807570	PRECEPT SCREW, 7.5X70MM POLYAXIAL
8807580	PRECEPT SCREW, 7.5X80MM POLYAXIAL
8807590	PRECEPT SCREW, 7.5X90MM POLYAXIAL
8808500	PRECEPT SCREW, 8.5X100MM POLYAXIAL
8808535	PRECEPT SCREW, 8.5X35MM POLYAXIAL
8808540	PRECEPT SCREW, 8.5X40MM POLYAXIAL
8808545	PRECEPT SCREW, 8.5X45MM POLYAXIAL
8808550	PRECEPT SCREW, 8.5X50MM POLYAXIAL
8808555	PRECEPT SCREW, 8.5X55MM POLYAXIAL
8808560	PRECEPT SCREW, 8.5X60MM POLYAXIAL
8808570	PRECEPT SCREW. 8.5X70MM POLYAXIAL
8808580	PRECEPT SCREW, 8.5X80MM POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8808590	PRECEPT SCREW, 8.5X90MM POLYAXIAL
8815550	PRECEPT SCREW. 5.5X50MM REDUCTION
8815555	PRECEPT SCREW 5.5X55MM REDUCTION
8816535	PRECEPT SCREW 6.5X35MM REDUCTION
8816545	PRECEPT SCREW 6.5X45MM REDUCTION
8824525	PRECEPT REDUCTION SCREW 4.5X25MM
8824530	PRECEPT REDUCTION SCREW 4.5X30MM
8824535	PRECEPT REDUCTION SCREW 4.5X35MM
8824540	PRECEPT REDUCTION SCREW 4.5X40MM
8824545	PRECEPT REDUCTION SCREW. 4.5X45MM
8825530	PRECEPT REDUCTION SCREW 5.5X30MM
8825535	PRECEPT REDUCTION SCREW 5.5X35MM
8825540	PRECEPT REDUCTION SCREW. 5.5X40MM
8825545	PRECEPT REDUCTION SCREW. 5.5X45MM
8825550	PRECEPT REDUCTION SCREW. 5.5X50MM
8825555	PRECEPT REDUCTION SCREW 5.5X55MM
8826530	PRECEPT REDUCTION SCREW 6.5X30MM
8826535	PRECEPT REDUCTION SCREW 6.5X35MM
8826540	PRECEPT REDUCTION SCREW, 6.5X40MM
8826545	PRECEPT REDUCTION SCREW. 6.5X45MM
8826550	PRECEPT REDUCTION SCREW 6.5X50MM
8826555	PRECEPT REDUCTION SCREW 6.5X55MM
8827530	PRECEPT REDUCTION SCREW 7.5X30MM
8827535	PRECEPT REDUCTION SCREW 7.5X35MM
8827540	PRECEPT REDUCTION SCREW 7.5X40MM
8827545	PRECEPT REDUCTION SCREW 7.5X45MM
8827550	PRECEPT REDUCTION SCREW 7.5X50MM
8827555	PRECEPT REDUCTION SCREW 7.5X55MM
8828540	PRECEPT REDUCTION SCREW, 8.5X40MM
8828545	PRECEPT REDUCTION SCREW, 8.5X45MM
8828550	PRECEPT REDUCTION SCREW, 8.5X50MM
8828555	PRECEPT REDUCTION SCREW, 8.5X55MM
8828560	PRECEPT REDUCTION SCREW, 8.5X60MM
8841020	PRECEPT ROD, 20MM
8841025	PRECEPT ROD, 25MM
8841030	PRECEPT ROD. 30MM
8841035	PRECEPT ROD, 35MM
8841040	PRECEPT ROD, 40MM
8841045	PRECEPT ROD, 45MM



República de Colombia  
Ministerio de Salud y Protección Social  
Instituto Nacional de Vigilancia de Medicamentos y Alimentos – INVIMA

**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**  
**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8841050	PRECEPT ROD, 50MM
8841055	PRECEPT ROD, 55MM
8841060	PRECEPT ROD, 60MM
8841065	PRECEPT ROD, 65MM
8841075	PRECEPT ROD, 75MM
8841085	PRECEPT ROD, 85MM
8841095	PRECEPT ROD, 95MM
8842120	PRECEPT ROD, 120MM STRAIGHT
8842300	PRECEPT ROD, 300MM STRAIGHT
8842500	PRECEPT ROD, 500MM STRAIGHT
8842600	PRECEPT ROD, 600MM STRAIGHT
8843500	PRECEPT TI ROD, 500MM PREBENT KYPHOTIC
8845120	PRECEPT COCR ROD, 120MM STRAIGHT
8845150	PRECEPT COCR ROD, 150MM STRAIGHT
8845300	PRECEPT COCR ROD, 300MM STRAIGHT
8845500	PRECEPT COCR ROD, 500MM STRAIGHT
8845600	PRECEPT COCR ROD, 600MM STRAIGHT
8846500	PRECEPT COCR ROD, 500MM PREBENT KYPHOTIC
8847020	PRECEPT TI ROD, 20MM PREBENT LORDOTIC
8847025	PRECEPT TI ROD, 25MM PREBENT LORDOTIC
8847030	PRECEPT TI ROD, 30MM PREBENT LORDOTIC
8847035	PRECEPT TI ROD, 35MM PREBENT LORDOTIC
8847040	PRECEPT TI ROD, 40MM PREBENT LORDOTIC
8847045	PRECEPT TI ROD, 45MM PREBENT LORDOTIC
8847050	PRECEPT TI ROD. 50MM PREBENT LORDOTIC
8847055	PRECEPT TI ROD, 55MM PREBENT LORDOTIC
8847060	PRECEPT TI ROD, 60MM PREBENT LORDOTIC
8847065	PRECEPT TI ROD, 65MM PREBENT LORDOTIC
8847070	PRECEPT TI ROD, 70MM PREBENT LORDOTIC
8847075	PRECEPT TI ROD, 75MM PREBENT LORDOTIC
8847080	PRECEPT TI ROD. 80MM PREBENT LORDOTIC
8847085	PRECEPT TI ROD, 85MM PREBENT LORDOTIC
8847090	PRECEPT TI ROD. 90MM PREBENT LORDOTIC
8847095	PRECEPT TI ROD, 95MM PREBENT LORDOTIC
8847100	PRECEPT TI ROD, 100MM PREBENT LORDOTIC
8847110	PRECEPT TI ROD, 110MM PREBENT LORDOTIC
8847120	PRECEPT TI ROD, 120MM PREBENT LORDOTIC
8847130	PRECEPT TI ROD, 130MM PREBENT LORDOTIC
8847140	PRECEPT TI ROD, 140MM PREBENT LORDOTIC



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8847150	PRECEPT TI ROD. 150MM PREBENT LORDOTIC
8847160	PRECEPT TI ROD, 160MM PREBENT LORDOTIC
8847170	PRECEPT TI ROD. 170MM PREBENT LORDOTIC
8847180	PRECEPT TI ROD. 180MM PREBENT LORDOTIC
8847190	PRECEPT TI ROD, 190MM PREBENT LORDOTIC
8847200	PRECEPT TI ROD, 200MM PREBENT LORDOTIC
8848035	PRECEPT TI ROD, 35MM STRAIGHT
8848045	PRECEPT TI ROD, 45MM STRAIGHT
8848050	PRECEPT TI ROD. 50MM STRAIGHT
8848055	PRECEPT TI ROD. 55MM STRAIGHT
8848060	PRECEPT TI ROD, 60MM STRAIGHT
8848070	PRECEPT TI ROD, 70MM STRAIGHT
8848075	PRECEPT TI ROD, 75MM STRAIGHT
8848080	PRECEPT TI ROD, 80MM STRAIGHT
8848090	PRECEPT TI ROD, 90MM STRAIGHT
8848100	PRECEPT TI ROD, 100MM STRAIGHT
8848110	PRECEPT TI ROD, 110MM STRAIGHT
8848120	PRECEPT TI ROD, 120MM STRAIGHT
8848130	PRECEPT TI ROD, 130MM STRAIGHT
8848140	PRECEPT TI ROD, 140MM STRAIGHT
8848150	PRECEPT TI ROD, 150MM STRAIGHT
8848160	PRECEPT TI ROD, 160MM STRAIGHT
8848170	PRECEPT TI ROD, 170MM STRAIGHT
8848180	PRECEPT TI ROD, 180MM STRAIGHT
8848190	PRECEPT TI ROD, 190MM STRAIGHT
8848200	PRECEPT TI ROD, 200MM STRAIGHT
8848210	PRECEPT TI ROD, 210MM STRAIGHT
8848220	PRECEPT TI ROD, 220MM STRAIGHT
8848230	PRECEPT TI ROD, 230MM STRAIGHT
8848240	PRECEPT TI ROD, 240MM STRAIGHT
8848250	PRECEPT TI ROD, 250MM STRAIGHT
8848300	PRECEPT TI ROD, 300MM STRAIGHT
8848500	PRECEPT TI ROD. 500MM STRAIGHT
8848600	PRECEPT TI ROD, 600MM STRAIGHT
8884530	PRECEPT SCREW, 4.5X30MM FIXED
8884535	PRECEPT SCREW, 4.5X35MM FIXED
8884540	PRECEPT SCREW, 4.5X40MM FIXED
8884545	PRECEPT SCREW, 4.5X45MM FIXED
8885535	PRECEPT SCREW, 5.5X35MM FIXED



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8885540	PRECEPT SCREW. 5.5X40MM FIXED
8885545	PRECEPT SCREW, 5.5X45MM FIXED
8885550	PRECEPT SCREW. 5.5X50MM FIXED
8885555	PRECEPT SCREW, 5.5X55MM FIXED
8886535	PRECEPT SCREW, 6.5X35MM FIXED
8886540	PRECEPT SCREW, 6.5X40MM FIXED
8886545	PRECEPT SCREW, 6.5X45MM FIXED
8886550	PRECEPT SCREW, 6.5X50MM FIXED
8886555	PRECEPT SCREW, 6.5X55MM FIXED
8887540	PRECEPT SCREW. 7.5X40MM FIXED
8887545	PRECEPT SCREW, 7.5X45MM FIXED
8887550	PRECEPT SCREW, 7.5X50MM FIXED
8887555	PRECEPT SCREW, 7.5X55MM FIXED
8887560	PRECEPT SCREW, 7.5X60MM FIXED
8888550	PRECEPT SCREW, 8.5X50MM FIXED
8888555	PRECEPT SCREW, 8.5X55MM FIXED
8888560	PRECEPT SCREW. 8.5X60MM FIXED
8890000	PRECEPT TULIP, MODULAR
8894525	PRECEPT SCREW SHANK, 4.5X25MM MODULAR
8894530	PRECEPT SCREW SHANK, 4.5X30MM MODULAR
8894535	PRECEPT SCREW SHANK, 4.5X35MM MODULAR
8894540	PRECEPT SCREW SHANK, 4.5X40MM MODULAR
8894545	PRECEPT SCREW SHANK, 4.5X45MM MODULAR
8894550	PRECEPT SCREW SHANK, 4.5X50MM MODULAR
8895530	PRECEPT SCREW SHANK, 5.5X30MM MODULAR
8895535	PRECEPT SCREW SHANK. 5.5X35MM MODULAR
8895540	PRECEPT SCREW SHANK, 5.5X40MM MODULAR
8895545	PRECEPT SCREW SHANK. 5.5X45MM MODULAR
8895550	PRECEPT SCREW SHANK, 5.5X50MM MODULAR
8895555	PRECEPT SCREW SHANK, 5.5X55MM MODULAR
8895560	PRECEPT SCREW SHANK, 5.5X60MM MODULAR
8896040	PRECEPT SCREW SHANK, 6.0X40MM MODULAR
8896045	PRECEPT SCREW SHANK. 6.0X45MM MODULAR
8896050	PRECEPT SCREW SHANK, 6.0X50MM MODULAR
8896055	PRECEPT SCREW SHANK, 6.0X55MM MODULAR
8896530	PRECEPT SCREW SHANK, 6.5X30MM MODULAR
8896535	PRECEPT SCREW SHANK, 6.5X35MM MODULAR
8896540	PRECEPT SCREW SHANK, 6.5X40MM MODULAR
8886545	PRECEPT SCREW SHANK, 6.5X45MM MODULAR



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
8896550	PRECEPT SCREW SHANK, 6.5X50MM MODULAR
8896555	PRECEPT SCREW SHANK, 6.5X55MM MODULAR
8896560	PRECEPT SCREW SHANK. 6.5X60MM MODULAR
8897045	PRECEPT SCREW SHANK. 7.0X45MM MODULAR
8897050	PRECEPT SCREW SHANK, 7.0X50MM MODULAR
8897055	PRECEPT SCREW SHANK, 7.0X55MM MODULAR
8897500	PRECEPT SCREW SHANK, 7.5X100MM MODULAR
8897510	PRECEPT SCREW SHANK, 7.5X110MM MODULAR
8897520	PRECEPT SCREW SHANK. 7.5X120MM MODULAR
8897530	PRECEPT SCREW SHANK, 7.5X30MM MODULAR
8897535	PRECEPT SCREW SHANK. 7.5X35MM MODULAR
8897540	PRECEPT SCREW SHANK, 7.5X40MM MODULAR
8897545	PRECEPT SCREW SHANK. 7.5X45MM MODULAR
8897550	PRECEPT SCREW SHANK, 7.5X50MM MODULAR
8897555	PRECEPT SCREW SHANK, 7.5X55MM MODULAR
8897560	PRECEPT SCREW SHANK, 7.5X60MM MODULAR
8897570	PRECEPT SCREW SHANK, 7.5X70MM MODULAR
8897580	PRECEPT SCREW SHANK, 7.5X80MM MODULAR
8897590	PRECEPT SCREW SHANK, 7.5X90MM MODULAR
8898500	PRECEPT SCREW SHANK, 8.5X100MM MODULAR
8898510	PRECEPT SCREW SHANK. 8.5X110MM MODULAR
8898520	PRECEPT SCREW SHANK 8.5X120MM MODULAR
8898530	PRECEPT SCREW SHANK 8.5X30MM MODULAR
8898535	PRECEPT SCREW SHANK 8.5X35MM MODULAR
8898540	PRECEPT SCREW SHANK 8.5X40MM MODULAR
8898545	PRECEPT SCREW SHANK 8.5X45MM MODULAR
8898550	PRECEPT SCREW SHANK 8.5X50MM MODULAR
8898555	PRECEPT SCREW SHANK 8.5X55MM MODULAR
8898570	PRECEPT SCREW SHANK 8.5X70MM MODULAR
8898580	PRECEPT SCREW SHANK, 8.5X80MM MODULAR
8898590	PRECEPT SCREW SHANK, 8.5X90MM MODULAR
8945535	PRECEPT SCREW 5.5X35MM SOLID UNIPLANAR
7070000006	MAS PLIF TULIP MODULAR
7070003006	MAS PLIF FIXATION LOCK SCREW
7094525006	MAS PLIF SHANK 4.5X25MM 3Z
7094530006	MAS PLIF SHANK 4.5X30MM 3Z
7094535006	MAS PLIF SHANK 4.5X35MM 3Z
7094540006	MAS PLIF SHANK 4.5X40MM 3Z
7095025006	MAS PLIF SHANK. 5.0X25MM 3Z



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
7095030006	MAS PLIF SHANK 5.0X30MM 3Z
7095035006	MAS PLIF SHANK, 5.0X35MM 3Z
7095040006	MAS PLIF SHANK, 5.0X40MM 3Z
7095525006	MAS PLIF SHANK, 5.5X25MM 3Z
7095530006	MAS PLIF SHANK, 5.5X30MM 3Z
7095535006	MAS PLIF SHANK. 5.5X35MM 3Z
7095540006	MAS PLIF SHANK, 5.5X40MM 3Z
7095545006	MAS PLIF SHANK. 5.5X45MM 3Z
7095550006	MAS PLIF SHANK, 5.5X50MM 3Z
7096530006	MAS PLIF SHANK, 6.5X30MM 3Z
7096535006	MAS PLIF SHANK, 6.5X35MM 3Z
7096540006	MAS PLIF SHANK, 6.5X40MM 3Z
7096545006	MAS PLIF SHANK, 6.5X45MM 3Z
7096550006	MAS PLIF SHANK, 6.5X50MM 3Z
7096555006	MAS PLIF SHANK, 6.5X55MM 3Z
7097530006	MAS PLIF SHANK, 7.5X30MM 3Z
7097535006	MAS PLIF SHANK, 7.5X35MM 3Z
7097540006	MAS PLIF SHANK, 7.5X40MM 3Z
7097545006	MAS PLIF SHANK, 7.5X45MM 3Z
7097550006	MAS PLIF SHANK. 7.5X50MM 3Z
7097555006	MAS PLIF SHANK. 7.5X55MM 3Z
7097560006	MAS PLIF SHANK. 7.5X60MM 3Z
7098535006	MAS PLIF SHANK, 8.5X35MM 3Z
7098540006	MAS PLIF SHANK, 8.5X40MM 3Z
7098545006	MAS PLIF SHANK, 8.5X45MM 3Z
7098550006	MAS PLIF SHANK, 8.5X50MM 3Z
7098555006	MAS PLIF SHANK, 8.5X55MM 3Z
7098560006	MAS PLIF SHANK, 8.5X60MM 3Z
7125535006	MAS PLIF SCREW, 5.5X35MM POLYAXIAL 3Z
7125540006	MAS PLIF SCREW, 5.5X40MM POLYAXIAL 3Z
7125545006	MAS PLIF SCREW, 5.5X45MM POLYAXIAL 3Z
7126535006	MAS PLIF SCREW, 6.5X35MM POLYAXIAL 3Z
7126540006	MAS PLIF SCREW, 6.5X40MM POLYAXIAL 3Z
7126545006	MAS PLIF SCREW, 6.5X45MM POLYAXIAL 3Z
7126550006	MAS PLIF SCREW, 6.5X50MM POLYAXIAL 3Z
8800000006	PRECEPT SET SCREW
8800535006	PRECEPT SCREW, 10.5X35MM POLYAXIAL C.R.
8800540006	PRECEPT SCREW, 10.5X40MM POLYAXIAL C.R.
8800545006	PRECEPT SCREW, 10.5X45MM POLYAXIAL C.R.



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8800550006	PRECEPT SCREW. 10.5X50MM POLYAXIAL C.R.
8804525006	PRECEPT SCREW, 4.5X25MM POLYAXIAL
8804530006	PRECEPT SCREW, 4.5X30MM POLYAXIAL
8804535006	PRECEPT SCREW, 4.5X35MM POLYAXIAL
8804540006	PRECEPT SCREW. 4.5X40MM POLYAXIAL
8804545006	PRECEPT SCREW, 4.5X45MM POLYAXIAL
8805525006	PRECEPT SCREW, 5.5X25MM POLYAXIAL
8805530006	PRECEPT SCREW, 5.5X30MM POLYAXIAL
8805535006	PRECEPT SCREW, 5.5X35MM POLYAXIAL
8805540006	PRECEPT SCREW, 5.5X40MM POLYAXIAL
8805545006	PRECEPT SCREW, 5.5X45MM POLYAXIAL
8805550006	PRECEPT SCREW, 5.5X50MM POLYAXIAL
8805555006	PRECEPT SCREW, 5.5X55MM POLYAXIAL
8806500006	PRECEPT SCREW, 6.5X100MM POLYAXIAL
8806525006	PRECEPT SCREW, 6.5X25MM POLYAXIAL
8806530006	PRECEPT SCREW, 6.5X30MM POLYAXIAL
8806535006	PRECEPT SCREW, 6.5X35MM POLYAXIAL
8806540006	PRECEPT SCREW, 6.5X40MM POLYAXIAL
8806545006	PRECEPT SCREW, 6.5X45MM POLYAXIAL
8806550006	PRECEPT SCREW, 6.5X50MM POLYAXIAL
8806555006	PRECEPT SCREW, 6.5X55MM POLYAXIAL
8806560006	PRECEPT SCREW, 6.5X60MM POLYAXIAL
8806570006	PRECEPT SCREW, 6.5X70MM POLYAXIAL
8806580006	PRECEPT SCREW, 6.5X80MM POLYAXIAL
8806590006	PRECEPT SCREW. 6.5X90MM POLYAXIAL
8807500006	PRECEPT SCREW, 7.5X100MM POLYAXIAL C.R.
8807530006	PRECEPT SCREW, 7.5X30MM POLYAXIAL C.R.
8807535006	PRECEPT SCREW, 7.5X35MM POLYAXIAL C.R.
8807540006	PRECEPT SCREW, 7.5X40MM POLYAXIAL C.R.
8807545006	PRECEPT SCREW, 7.5X45MM POLYAXIAL C.R.
8807550006	PRECEPT SCREW, 7.5X50MM POLYAXIAL C.R.
8807555006	PRECEPT SCREW, 7.5X55MM POLYAXIAL C.R.
8807560006	PRECEPT SCREW, 7.5X60MM POLYAXIAL C.R.
8807570006	PRECEPT SCREW, 7.5X70MM POLYAXIAL C.R.
8807580006	PRECEPT SCREW, 7.5X80MM POLYAXIAL C.R.
8807590006	PRECEPT SCREW, 7.5X90MM POLYAXIAL C.R.
8808500006	PRECEPT SCREW, 8.5X100MM POLYAXIAL C.R.
8808540006	PRECEPT SCREW, 8.5X40MM POLYAXIAL C.R.
8808545006	PRECEPT SCREW, 8.5X45MM POLYAXIAL C.R.



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8808550006	PRECEPT SCREW, 8.5X50MM POLYAXIAL C.R.
8808555006	PRECEPT SCREW, 8.5X55MM POLYAXIAL C.R.
8808560006	PRECEPT SCREW, 8.5X60MM POLYAXIAL C.R.
8808570006	PRECEPT SCREW, 8.5X70MM POLYAXIAL C.R.
8808580006	PRECEPT SCREW, 8.5X80MM POLYAXIAL C.R.
8808590006	PRECEPT SCREW, 8.5X90MM POLYAXIAL C.R.
8809500006	PRECEPT SCREW, 9.5X100MM POLYAXIAL C.R.
8809535006	PRECEPT SCREW, 9.5X35MM POLYAXIAL C.R.
8809540006	PRECEPT SCREW, 9.5X40MM POLYAXIAL C.R.
8809545006	PRECEPT SCREW, 9.5X45MM POLYAXIAL C.R.
8809550006	PRECEPT SCREW, 9.5X50MM POLYAXIAL C.R.
8809590006	PRECEPT SCREW, 9.5X90MM POLYAXIAL C.R.
8824525006	PRECEPT REDUCTION SCREW, 4.5X25MM
8824530006	PRECEPT REDUCTION SCREW, 4.5X30MM
8824535006	PRECEPT REDUCTION SCREW, 4.5X35MM
8824540006	PRECEPT REDUCTION SCREW, 4.5X40MM
8824545006	PRECEPT REDUCTION SCREW, 4.5X45MM
8824550006	PRECEPT REDUCTION SCREW, 4.5X50MM
8824555006	PRECEPT REDUCTION SCREW, 4.5X55MM
8824560006	PRECEPT REDUCTION SCREW, 4.5X60MM
8825025006	PRECEPT REDUCTION SCREW, 5.0X25MM
8825030006	PRECEPT REDUCTION SCREW, 5.0X30MM
8825035006	PRECEPT REDUCTION SCREW, 5.0X35MM
8825040006	PRECEPT REDUCTION SCREW, 5.0X40MM
8825045006	PRECEPT REDUCTION SCREW, 5.0X45MM
8825050006	PRECEPT REDUCTION SCREW, 5.0X50MM
8825055006	PRECEPT REDUCTION SCREW, 5.0X55MM
8825060006	PRECEPT REDUCTION SCREW, 5.0X60MM
8825525006	PRECEPT REDUCTION SCREW, 5.5X25MM
8825530006	PRECEPT REDUCTION SCREW, 5.5X30MM
8825535006	PRECEPT REDUCTION SCREW, 5.5X35MM
8825540006	PRECEPT REDUCTION SCREW, 5.5X40MM
8825545006	PRECEPT REDUCTION SCREW, 5.5X45MM
8825550006	PRECEPT REDUCTION SCREW, 5.5X50MM
8825555006	PRECEPT REDUCTION SCREW, 5.5X55MM
8825560006	PRECEPT REDUCTION SCREW, 5.5X60MM
8826025006	PRECEPT REDUCTION SCREW, 6.0X25MM
8826030006	PRECEPT REDUCTION SCREW, 6.0X30MM
8826035006	PRECEPT REDUCTION SCREW, 6.0X35MM



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8826040006	PRECEPT REDUCTION SCREW, 6.0X40MM
8826045006	PRECEPT REDUCTION SCREW, 6.0X45MM
8826050006	PRECEPT REDUCTION SCREW, 6.0X50MM
8826055006	PRECEPT REDUCTION SCREW, 6.0X55MM
8826060006	PRECEPT REDUCTION SCREW, 6.0X60MM
8826525006	PRECEPT REDUCTION SCREW, 6.5X25MM
8826530006	PRECEPT REDUCTION SCREW. 6.5X30MM
8826535006	PRECEPT REDUCTION SCREW, 6.5X35MM
8826540006	PRECEPT REDUCTION SCREW, 6.5X40MM
8826545006	PRECEPT REDUCTION SCREW, 6.5X45MM
8826550006	PRECEPT REDUCTION SCREW. 6.5X50MM
8826555006	PRECEPT REDUCTION SCREW, 6.5X55MM
8826560006	PRECEPT REDUCTION SCREW, 6.5X60MM
8826570006	PRECEPT REDUCTION SCREW, 6.5X70MM
8826580006	PRECEPT REDUCTION SCREW, 6.5X80MM
8827500006	PRECEPT SCREW, 7.5X100MM REDUCTION C.R.
8827530006	PRECEPT SCREW, 7.5X30MM REDUCTION C.R.
8827535006	PRECEPT SCREW, 7.5X35MM REDUCTION C.R.
8827540006	PRECEPT SCREW. 7.5X40MM REDUCTION C.R.
8827545006	PRECEPT SCREW, 7.5X45MM REDUCTION C.R.
8827550006	PRECEPT SCREW, 7.5X50MM REDUCTION C.R.
8827555006	PRECEPT SCREW, 7.5X55MM REDUCTION C.R.
8827560006	PRECEPT SCREW, 7.5X60MM REDUCTION C.R.
8827570006	PRECEPT SCREW, 7.5X70MM REDUCTION C.R.
8827580006	PRECEPT SCREW, 7.5X80MM REDUCTION C.R.
8827590006	PRECEPT SCREW, 7.5X90MM REDUCTION C.R.
8828500006	PRECEPT SCREW, 8.5X100MM REDUCTION C.R.
8828535006	PRECEPT SCREW, 8.5X35MM REDUCTION C.R.
8828540006	PRECEPT SCREW, 8.5X40MM REDUCTION C.R.
8828545006	PRECEPT SCREW, 8.5X45MM REDUCTION C.R.
8828550006	PRECEPT SCREW, 8.5X50MM REDUCTION C.R.
8828555006	PRECEPT SCREW, 8.5X55MM REDUCTION C.R.
8828560006	PRECEPT SCREW, 8.5X60MM REDUCTION C.R.
8828570006	PRECEPT SCREW, 8.5X70MM REDUCTION C.R.
8828580006	PRECEPT SCREW, 8.5X80MM REDUCTION C.R.
8828590006	PRECEPT SCREW, 8.5X90MM REDUCTION C.R.
8847020006	PRECEPT TI ROD, 20MM PREBENT LORDOTIC
8847025006	PRECEPT TI ROD, 25MM PREBENT LORDOTIC
8847030006	PRECEPT TI ROD 30MM PREBENT LORDOTIC



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8847035006	PRECEPT TI. ROD, 35MM PREBENT LORDOTIC
8847040006	PRECEPT TI ROD, 40MM PREBENT LORDOTIC
8847045006	PRECEPT TI ROD, 45MM PREBENT LORDOTIC
8847050006	PRECEPT TI ROD, 50MM PREBENT LORDOTIC
8847055006	PRECEPT TI ROD 55MM PREBENT LORDOTIC
8847060006	PRECEPT TI ROD 60MM PREBENT LORDOTIC
8847065006	PRECEPT TI ROD 65MM PREBENT LORDOTIC
8847070006	PRECEPT TI ROD 70MM PREBENT LORDOTIC
8847075006	PRECEPT TI ROD 75MM PREBENT LORDOTIC
8847080006	PRECEPT TI ROD 80MM PREBENT LORDOTIC
8847085006	PRECEPT TI ROD 85MM PREBENT LORDOTIC
8847090006	PRECEPT TI ROD 90MM PREBENT LORDOTIC
8847095006	PRECEPT TI ROD 95MM PREBENT LORDOTIC
8847100006	PRECEPT TI ROD 100MM PREBENT LORDOTIC
8847110006	PRECEPT TI ROD, 110MM PREBENT LORDOTIC
8847120006	PRECEPT TI ROD, 120MM PREBENT LORDOTIC
8847130006	PRECEPT TI ROD, 130MM PREBENT LORDOTIC
8847140006	PRECEPT TI ROD, 140MM PREBENT LORDOTIC
8847150006	PRECEPT TI ROD, 150MM PREBENT LORDOTIC
8847160006	PRECEPT TI ROD, 160MM PREBENT LORDOTIC
8848040006	PRECEPT TI ROD, 40MM STRAIGHT
8848045006	PRECEPT TI ROD, 45MM STRAIGHT
8848050006	PRECEPT TI ROD, 50MM STRAIGHT
8848060006	PRECEPT TI ROD. 60MM STRAIGHT
8848070006	PRECEPT TI ROD, 70MM STRAIGHT
8848080006	PRECEPT TI ROD, 80MM STRAIGHT
8848090006	PRECEPT TI ROD, 90MM STRAIGHT
8848100006	PRECEPT TI ROD, 100MM STRAIGHT
8848110006	PRECEPT TI ROD, 110MM STRAIGHT
8848120006	PRECEPT TI ROD, 120MM STRAIGHT
8848130006	PRECEPT TI ROD, 130MM STRAIGHT
8848140006	PRECEPT TI ROD, 140MM STRAIGHT
8848150006	PRECEPT TI ROD, 150MM STRAIGHT
8848160006	PRECEPT TI ROD, 160MM STRAIGHT
8848170006	PRECEPT TI ROD, 170MM STRAIGHT
8848180006	PRECEPT TI ROD, 180MM STRAIGHT
8848190006	PRECEPT TI ROD, 190MM STRAIGHT
8848200006	PRECEPT TI ROD, 200MM STRAIGHT
8848210006	PRECEPT TI ROD, 210MM STRAIGHT



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8848220006	PRECEPT TI ROD, 220MM STRAIGHT
8848230006	PRECEPT TI ROD, 230MM STRAIGHT
8848240006	PRECEPT TI ROD, 240MM STRAIGHT
8848250006	PRECEPT TI ROD, 250MM STRAIGHT
8848300006	PRECEPT TI ROD, 300MM STRAIGHT
8848500006	PRECEPT TI ROD, 500MM STRAIGHT
8875530994	PRECEPT SCREW, 5.5X30MM FENESTRATED
8875535994	PRECEPT SCREW. 5.5X35MM FENESTRATED
8875540994	PRECEPT SCREW 5.5X 40MM FENESTRATED
8875545994	PRECEPT SCREW 5.5X45MM FENESTRATED
8875550994	PRECEPT SCREW, 5.5X 50MM FENESTRATED
8875555994	PRECEPT SCREW, 5.5X55MM FENESTRATED
8875560994	PRECEPT SCREW, 5.5X60MM FENESTRATED
8876530994	PRECEPT SCREW. 6.5X30MM FENESTRATED
8876535994	PRECEPT SCREW, 6.5X35MM FENESTRATED
8876540994	PRECEPT SCREW, 6.5X40MM FENESTRATED
8876545994	PRECEPT SCREW, 6.5X45MM FENESTRATED
8876550994	PRECEPT SCREW, 6.5X50MM FENESTRATED
8876555994	PRECEPT SCREW, 6.5X55MM FENESTRATED
8876560994	PRECEPT SCREW, 6.5X60MM FENESTRATED
8877535994	PRECEPT SCREW. 7.5X35MM FENESTRATED
8877540994	PRECEPT SCREW, 7.5X40MM FENESTRATED
8877545994	PRECEPT SCREW, 7.5X45MM FENESTRATED
8877550994	PRECEPT SCREW, 7.5X50MM FENESTRATED
8877555994	PRECEPT SCREW, 7.5X55MM FENESTRATED
8877560994	PRECEPT SCREW, 7.5X60MM FENESTRATED
8878500994	PRECEPT SCREW, 8.5X100MM FENESTRATED
8878510994	PRECEPT SCREW, 8.5X110MM FENESTRATED
8878520994	PRECEPT SCREW. 8.5X120MM FENESTRATED
8878535994	PRECEPT SCREW, 8.5X35MM FENESTRATED
8878540994	PRECEPT SCREW, 8.5X40MM FENESTRATED
8878545994	PRECEPT SCREW. 8.5X45MM FENESTRATED
8878550994	PRECEPT SCREW. 8.5X50MM FENESTRATED
8878555994	PRECEPT SCREW. 8.5X55MM FENESTRATED
8878560994	PRECEPT SCREW. 8.5X60MM FENESTRATED
8878580994	PRECEPT SCREW. 8.5X80MM FENESTRATED
8878590994	PRECEPT SCREW. 8.5X90MM FENESTRATED
8879500994	PRECEPT SCREW. 9.5X100MM FENESTRATED
8879510994	PRECEPT SCREW. 9.5X110MM FENESTRATED



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8879520994	PRECEPT SCREW, 9.5X120MM FENESTRATED
8879580994	PRECEPT SCREW. 9.5X80MM FENESTRATED
8879590994	PRECEPT SCREW, 9.5X90MM FENESTRATED
7070000P	MAS PLIF TULIP, MODULAR
7070003P	MAS PLIF FIXATION LOCK SCREW
7094525P	MAS PLIF SHANK, 4.5X25MM 3Z
7094530P	MAS PLIF SHANK, 4.5X30MM 3Z
7094535P	MAS PLIF SHANK, 4.5X35MM 3Z
7094540P	MAS PLIF SHANK, 4.5X40MM 3Z
7095025P	MAS PLIF SHANK, 5.0X25MM 3Z
7095030P	MAS PLIF SHANK, 5.0X30MM 3Z
7095035P	MAS PLIF SHANK, 5.0X35MM 3Z
7095040P	MAS PLIF SHANK, 5.0X40MM 3Z
7095525P	MAS PLIF SHANK, 5.5X25MM 3Z
7095530P	MAS PLIF SHANK, 5.5X30MM 3Z
7095535P	MAS PLIF SHANK, 5.5X35MM 3Z
7095540P	MAS PLIF SHANK, 5,5X40MM 3Z
7096530P	MAS PLIF SHANK, 6.5X30MM 3Z
7096535P	MAS PLIF SHANK, 6.5X35MM 3Z
7096540P	MAS PLIF SHANK, 6.5X40MM 3Z
7096545P	MAS PLIF SHANK, 6.5X45MM 3Z
7124525P	MAS PLIF SCREW, 4.5X25MM POLYAXIAL 3Z
7124530P	MAS PLIF SCREW, 4.5X30MM POLYAXIAL 3Z
7124535P	MAS PLIF SCREW, 4.5X35MM POLYAXIAL 3Z
7124540P	MAS PLIF SCREW, 4.5X40MM POLYAXIAL 3Z
7125025P	MAS PLIF SCREW, 5.0X25MM POLYAXIAL 3Z
7125030P	MAS PLIF SCREW, 5.0X30MM POLYAXIAL 3Z
7125035P	MAS PLIF SCREW, 5.0X35MM POLYAXIAL 3Z
7125040P	MAS PLIF SCREW, 5.0X40MM POLYAXIAL 3Z
7125525P	MAS PLIF SCREW, 5.5X25MM POLYAXIAL 3Z
7125530P	MAS PLIF SCREW, 5.5X30MM POLYAXIAL 3Z
7125535P	MAS PLIF SCREW. 5.5X35MM POLYAXIAL 3Z
7125540P	MAS PLIF SCREW, 5.5X40MM POLYAXIAL 3Z
8800000P	PRECEPT SET SCREW
8800500A	PRECEPT SCREW, 10.5X100MM POLYAXIAL C.R
8804525A	PRECEPT SCREW, 4.5X25MM POLYAXIAL
8804525AP	PRECEPT SCREW. 4.5X25MM POLYAXIAL S R.
8804530A	PRECEPT SCREW, 4.5X30MM POLYAXIAL
8804530AP	PRECEPT SCREW, 4.5X30MM POLYAXIAL S.R.



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8804535A	PRECEPT SCREW, 4.5X35MM POLYAXIAL
8804535AP	PRECEPT SCREW, 4.5X35MM POLYAXIAL S.R.
8804540A	PRECEPT SCREW, 4.5X40MM POLYAXIAL
8804540AP	PRECEPT SCREW, 4.5X40MM POLYAXIAL S.R.
8804545A	PRECEPT SCREW, 4.5X45MM POLYAXIAL
8804545AP	PRECEPT SCREW, 4.5X45MM POLYAXIAL S.R.
8804550A	PRECEPT SCREW, 4.5X50MM POLYAXIAL S.R.
8804555A	PRECEPT SCREW, 4.5X55MM POLYAXIAL S.R.
8804560A	PRECEPT SCREW. 4.5X60MM POLYAXIAL S.R.
8805530A	PRECEPT SCEW, 5.5X30MM POLYAXIAL
8805530AP	PRECEPT SCREW, 5.5X30MM POLYAXIAL S.R.
8805535A	PRECEPT SCREW, 5.5X35MM POLYAXIAL
8805535AP	PRECEPT SCREW, 5,5X35MM POLYAXIAL S.R.
8805540A	PRECEPT SCREW, 5.5X40MM POLYAXIAL
8805540AP	PRECEPT SCREW. 5.5X40MM POLYAXIAL S.R.
8805545A	PRECEPT SCREW, 5.5X45MM POLYAXIAL
8805545AP	PRECEPT SCREW. 5.5X45MM POLYAXIAL S.R.
8805550A	PRECEPT SCREW, 5.5X50MM POLYAXIAL
8805550AP	PRECEPT SCREW, 5.5X50MM POLYAXIAL S.R.
8805555A	PRECEPT SCREW, 5.5X55MM POLYAXIAL
8805555AP	PRECEPT SCREW, 5.5X55MM POLYAXIAL S.R.
8805560A	PRECEPT SCREW, 5.5X60MM POLYAXIAL S.R.
8805560AP	PRECEPT SCREW, 5.5X60MM POLYAXIAL S.R.
8806530A	PRECEPT SCREW, 6.5X30MM POLYAXIAL
8806530AP	PRECEPT SCREW. 6.5X30MM POLYAXIAL S.R
8806535A	PRECEPT SCREW 6.5X35MM POLYAXIAL
8806535AP	PRECEPT SCREW 6.5X35MM POLYAXIAL S R
8806540A	PRECEPT SCREW 6.5X40MM POLYAXIAL
8806540AP	PRECEPT SCREW 6.5X40MM POLYAXIAL S R
8806545A	PRECEPT SCREW 6.5X45MM POLYAXIAL
8806545AP	PRECEPT SCREW 6.5X45MM POLYAXIAL S R
8806550A	PRECEPT SCREW 6.5X50MM POLYAXIAL
8806550AP	PRECEPT SCREW, 6.5X50MM POLYAXIAL S.R
8806555A	PRECEPT SCREW 6.5X55MM POLYAXIAL
8806555AP	PRECEPT SCREW, 6.5X55MM POLYAXIAL S.R
8806560A	PRECEPT SCREW, 6.5X60MM POLYAXIAL S R
8806560AP	PRECEPT SCREW, 6.5X60MM POLYAXIAL S.R
8806565A	PRECEPT SCREW 6.5X65MM POLYAXIAL S.R
8807500A	PRECEPT SCREW 7.5X100MM POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8807500AP	PRECEPT SCREW. 7.5X100MM POLYAXIAL C.R
8807530A	PRECEPT SCREW 7.5X30MM POLYAXIAL
8807530AP	PRECEPT SCREW. 7.5X30MM POLYAXIAL C R
8807535A	PRECEPT SCREW 7.5X35MM POLYAXIAL
8807535AP	PRECEPT SCREW. 7.5X35MM POLYAXIAL C R
8807540A	PRECEPT SCREW 7.5X40MM POLYAXIAL
8807540AP	PRECEPT SCREW 7.5X40MM POLYAXIAL C R
8807545A	PRECEPT SCREW 7.5X45MM POLYAXIAL
8807545AP	PRECEPT SCREW 7.5X45MM POLYAXIAL C R
8807550A	PRECEPT SCREW 7.5X50MM POLYAXIAL
8807550AP	PRECEPT SCREW, 7.5X50MM POLYAXIAL C R
8807555A	PRECEPT SCREW 7.5X55MM POLYAXIAL
8807555AP	PRECEPT SCREW, 7.5X55MM POLYAXIAL C Q
8807560A	PRECEPT SCREW 7.5X60MM POLYAXIAL
8807560AP	PRECEPT SCREW. 7.5X60MM POLYAXIAL C R
8807565A	PRECEPT SCREW 7.5X65MM POLYAXIAL C R
8807570A	PRECEPT SCREW 7.5X70MM POLYAXIAL
8807570AP	PRECEPT SCREW, 7.5X70MM POLYAXIAL C.R.
8807580A	PRECEPT SCREW, 7.5X80MM POLYAXIAL
8807580AP	PRECEPT SCREW, 7.5X80MM POLYAXIAL C.R.
8807590A	PRECEPT SCREW, 7.5X90MM POLYAXIAL
8807590AP	PRECEPT SCREW, 7.5X90MM POLYAXIAL C.R.
8808500A	PRECEPT SCREW, 8.5X100MM POLYAXIAL
8808500AP	PRECEPT SCREW, 8.5X100MM POLYAXIAL C.R.
8808530A	PRECEPT SCREW, 8.5X30MM POLYAXIAL
8808530AP	PRECEPT SCREW, 8.5X30MM POLYAXIAL C.R.
8808535A	PRECEPT SCREW, 8.5X35MM POLYAXIAL
8808535AP	PRECEPT SCREW, 8.5X35MM POLYAXIAL C.R.
8808540A	PRECEPT SCREW, 8.5X40MM POLYAXIAL
8808540AP	PRECEPT SCREW, 8.5X40MM POLYAXIAL C.R.
8808545A	PRECEPT SCREW, 8.5X45MM POLYAXIAL
8808545AP	PRECEPT SCREW, 8.5X45MM POLYAXIAL C.R.
8808550A	PRECEPT SCREW, 8.5X50MM POLYAXIAL
8808550AP	PRECEPT SCREW, 8.5X50MM POLYAXIAL C.R.
8808555A	PRECEPT SCREW, 8.5X55MM POLYAXIAL
8808555AP	PRECEPT SCREW, 8.5X55MM POLYAXIAL C.R.
8808560A	PRECEPT SCREW, 8.5X60MM POLYAXIAL
8808560AP	PRECEPT SCREW, 8.5X60MM POLYAXIAL C.R.
8808565A	PRECEPT SCREW, 8.5X65MM POLYAXIAL C.R.



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8808570A	PRECEPT SCREW, 8.5X70MM POLYAXIAL
8808570AP	PRECEPT SCREW, 8.5X70MM POLYAXIAL C.R.
8808580A	PRECEPT SCREW, 8.5X80MM POLYAXIAL
8808580AP	PRECEPT SCREW, 8.5X80MM POLYAXIAL C.R.
8808590A	PRECEPT SCREW, 8.5X90MM POLYAXIAL
8808590AP	PRECEPT SCREW, 8.5X90MM POLYAXIAL C.R.
8809500A	PRECEPT SCREW. 9.5X100MM POLYAXIAL C.R.
8809500AP	PRECEPT SCREW, 9.5X100MM POLYAXIAL C.R.
8809540A	PRECEPT SCREW, 9.5X40MM POLYAXIAL C.R.
8809545A	PRECEPT SCREW, 9.5X45MM POLYAXIAL C.R.
8809550A	PRECEPT SCREW, 9.5X50MM POLYAXIAL C.R.
8809555A	PRECEPT SCREW, 9.5X55MM POLYAXIAL C.R.
8809565A	PRECEPT SCREW, 9.5X65MM POLYAXIAL C.R.
8809570A	PRECEPT SCREW, 9.5X70MM POLYAXIAL C.R.
8809570AP	PRECEPT SCREW, 9.5X70MM POLYAXIAL C,R.
8809580A	PRECEPT SCREW, 9.5X80MM POLYAXIAL C.R.
8809580AP	PRECEPT SCREW, 9.5X80MM POLYAXIAL C.R.
8809590A	PRECEPT SCREW, 9.5X90MM POLYAXIAL C.R.
8809590AP	PRECEPT SCREW, 9.5X90MM POLYAXIAL C.R.
8834525A	PRECEPT SCREW, 4.5X25MM POLYAXIAL TL
8834530A	PRECEPT SCREW, 4.5X30MM POLYAXIAL TL
8834535A	PRECEPT SCREW, 4.5X35MM POLYAXIAL TL
8834540A	PRECEPT SCREW, 4.5X40MM POLYAXIAL TL
8834545A	PRECEPT SCREW, 4.5X45MM POLYAXIAL TL
8835530A	PRECEPT SCREW, 5.5X30MM POLYAXIAL TL
8835535A	PRECEPT SCREW, 5.5X35MM POLYAXIAL TL
8835540A	PRECEPT SCREW, 5.5X40MM POLYAXIAL TL
8835545A	PRECEPT SCREW, 5.5X45MM POLYAXIAL TL
8835550A	PRECEPT SCREW, 5.5X50MM POLYAXIAL TL
8835555A	PRECEPT SCREW, 5.5X55MM POLYAXIAL TL
8836530A	PRECEPT SCREW, 6.5X30MM POLYAXIAL TL
8836535A	PRECEPT SCREW, 6.5X35MM POLYAXIAL TL
8836540A	PRECEPT SCREW, 6.5X40MM POLYAXIAL TL
8836545A	PRECEPT SCREW, 6.5X45MM POLYAXIAL TL
8836550A	PRECEPT SCREW, 6.5X50MM POLYAXIAL TL
8836555A	PRECEPT SCREW, 6.5X55MM POLYAXIAL TL
8837530A	PRECEPT SCREW, 7.5X30MM POLYAXIAL TL
8837535A	PRECEPT SCREW, 7.5X35MM POLYAXIAL TL
8837540A	PRECEPT SCREW, 7.5X40MM POLYAXIAL TL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8837545A	PRECEPT SCREW, 7.5X45MM POLYAXIAL TL
8837550A	PRECEPT SCREW, 7.5X50MM POLYAXIAL TL
8837555A	PRECEPT SCREW, 7.5X55MM POLYAXIAL TL
8845300P	PRECEPT COCR ROD, 300MM STRAIGHT
8845500P	PRECEPT COCR ROD, 500MM STRAIGHT
8847020P	PRECEPT TI ROD, 20MM PREBENT
8847025P	PRECEPT TI ROD, 25MM PREBENT
8847030P	PRECEPT TI ROD, 30MM PREBENT
8847035P	PRECEPT TI ROD, 35MM PREBENT
8847040P	PRECEPT TI ROD, 40MM PREBENT
8847045P	PRECEPT TI ROD, 45MM PREBENT
8847050P	PRECEPT TI ROD, 50MM PREBENT
8847055P	PRECEPT TI ROD, 55MM PREBENT
8847060P	PRECEPT TI ROD, 60MM PREBENT
8847065P	PRECEPT TI ROD. 65MM PREBENT
8847070P	PRECEPT TI ROD, 70MM PREBENT
8847075P	PRECEPT TI ROD, 75MM PREBENT
8847080P	PRECEPT TI ROD. 80MM PREBENT
8847085P	PRECEPT TI ROD. 85MM PREBENT
8847090P	PRECEPT TI ROD, 90MM PREBENT
8847095P	PRECEPT TI ROD, 95MM PREBENT
8847100P	PRECEPT TI ROD, 100MM PREBENT
8847110P	PRECEPT TI ROD, 110MM PREBENT
8847120P	PRECEPT TI ROD, 120MM PREBENT
8847130P	PRECEPT TI ROD, 130MM PREBENT
8847140P	PRECEPT TI ROD, 140MM PREBENT
8847150P	PRECEPT TI ROD, 150MM PREBENT
8847160P	PRECEPT TI ROD, 160MM PREBENT
8848160P	PRECEPT TI ROD, 160MM STRAIGHT
8848170P	PRECEPT TI ROD, 170MM STRAIGHT
8848180P	PRECEPT TI ROD, 180MM STRAIGHT
8848190P	PRECEPT TI ROD, 190MM STRAIGHT
8848200P	PRECEPT TI ROD, 200MM STRAIGHT
8848210P	PRECEPT TI ROD. 210MM STRAIGHT
8848220P	PRECEPT TI ROD, 220MM STRAIGHT
8848230P	PRECEPT TI ROD, 230MM STRAIGHT
8848240P	PRECEPT TI ROD. 240MM STRAIGHT
8848250P	PRECEPT TI ROD, 250MM STRAIGHT
8848300P	PRECEPT TI ROD, 300MM STRAIGHT



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8848500P	PRECEPT TI ROD, 500MM STRAIGHT
8890000P	PRECEPT TULIP, MODULAR
8895530P	PRECEPT SCREW SHANK, 5.5X30MM MODULAR
8895535P	PRECEPT SCREW SHANK, 5.5X35MM MODULAR
8895540P	PRECEPT SCREW SHANK, 5.5X40MM MODULAR
8895545P	PRECEPT SCREW SHANK. 5.5X45MM MODULAR
8895550P	PRECEPT SCREW SHANK. 5.5X50MM MODULAR
8895555P	PRECEPT SCREW SHANK. 5.5X55MM MODULAR
8896535P	PRECEPT SCREW SHANK, 6.5X35MM MODULAR
8896540P	PRECEPT SCEW SHANK, 6.5X40MM MODULAR
8896545P	PRECEPT SCREW SHANK, 6.5X45MM MODULAR
8896550P	PRECEPT SCREW SHANK, 6.5X50MM MODULAR
8886555P	PRECEPT SCREW SHANK, 6.5X55MM MODULAR
8897530P	PRECEPT SCREW SHANK, 7.5X30MM MODULAR
8897535P	PRECEPT SCREW SHANK, 7.5X35MM MODULAR
8897540P	PRECEPT SCREW SHANK, 7.5X40MM MODULAR
8897545P	PRECEPT SCREW SHANK, 7.5X45MM MODULAR
8897550P	PRECEPT SCREW SHANK, 7.5X50MM MODULAR
8897555P	PRECEPT SCREW SHANK, 7.5X55MM MODULAR
1635999	RELIN LOCK SCREW, 5.5MM OPEN TULIP
10010006	GSB TULIP, MOD
10014525	GSB SCREW, 4.5X25MM 2C POLYAXIAL
10014530	GSB SCREW, 4.5X30MM 2C POLYAXIAL
10014535	GSB SCREW, 4.5X35MM 2C POLYAXIAL
10014540	GSB SCREW, 4.5X40MM 2C POLYAXIAL
10015040	GSB SCREW, 5.0X40MM 2C POLYAXIAL
10015530	GSB SCREW, 5.5X30MM 2C POLYAXIAL
10015535	GSB SCREW, 5.5X35MM 2C POLYAXIAL
10015540	GSB SCREW, 5.5X40MM 2C POLYAXIAL
10015545	GSB SCREW, 5.5X45MM 2C POLYAXIAL
10015550	GSB SCREW, 5.5X50MM 2C POLYAXIAL
10016530	GSB SCREW, 6.5X30MM 2C POLYAXIAL
10016535	GSB SCREW, 6.5X35MM 2C POLYAXIAL
10016540	GSB SCREW, 6.5X40MM 2C POLYAXIAL
10016545	GSB SCREW, 6.5X45MM 2C POLYAXIAL
10016550	GSB SCREW, 6.5X50MM 2C POLYAXIAL
10016555	GSB SCREW, 6.5X55MM 2C POLYAXIAL
10016560	GSB SCREW, 6.5X60MM 2C POLYAXIAL
10016570	GSB SCREW, 6.5X70MM 2C POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
10017535	GSB SCREW, 7.5X35MM 2C POLYAXIAL
10017540	GSB SCREW, 7.5X40MM 2C POLYAXIAL
10017545	GSB SCREW, 7.5X45MM 2C POLYAXIAL
10017550	GSB SCREW, 7.5X50MM 2C POLYAXIAL
10017555	GSB SCREW, 7.5X55MM 2C POLYAXIAL
10017560	GSB SCREW, 7.5X60MM 2C POLYAXIAL
10017570	GSB SCREW, 7.5X70MM 2C POLYAXIAL
10017580	GSB SCREW, 7.5X80MM 2C POLYAXIAL
10017590	GSB SCREW, 7.5X90MM 2C POLYAXIAL
10018510	GSB SCREW, 8.5X100MM 2C POLYAXIAL
10018535	GSB SCREW, 8.5X35MM 2C POLYAXIAL
10018540	GSB SCREW, 8.5X40MM 2C POLYAXIAL
10018545	GSB SCREW, 8.5X45MM 2C POLYAXIAL
10018550	GSB SCREW, 8.5X50MM 2C POLYAXIAL
10018555	GSB SCREW, 8.5X55MM 2C POLYAXIAL
10018560	GSB SCREW, 8.5X60MM 2C POLYAXIAL
10018570	GSB SCREW, 8.5X70MM 2C POLYAXIAL
10018580	GSB SCREW, 8.5X80MM 2C POLYAXIAL
10018590	GSB SCREW, 8.5X90MM 2C POLYAXIAL
10025530	GSB SCREW, 5.5X30MM 2C MAS RED
10025535	GSB SCREW, 5.5X35MM 2C MAS RED
10025540	GSB SCREW, 5.5X40MM 2C MAS RED
10025545	GSB SCREW, 5.5X45MM 2C MAS RED
10025550	GSB SCREW, 5.5X50MM 2C MAS RED
10025555	GSB SCREW, 5.5X55MM 2C MAS RED
10026530	GSB SCREW, 6.5X30MM 2C MAS RED
10026535	GSB SCREW, 6.5X35MM 2C MAS RED
10026540	GSB SCREW, 6.5X40MM 2C MAS RED
10026545	GSB SCREW, 6.5X45MM 2C MAS RED
10026550	GSB SCREW, 6.5X50MM 2C MAS RED
10026555	GSB SCREW, 6.5X55MM 2C MAS RED
10027530	GSB SCREW, 7.5X30MM 2C MAS RED
10027535	GSB SCREW, 7.5X35MM 2C MAS RED
10027540	GSB SCREW, 7.5X40MM 2C MAS RED
10027545	GSB SCREW, 7.5X45MM 2C MAS RED
10027550	GSB SCREW, 7.5X50MM 2C MAS RED
10027555	GSB SCREW, 7.5X55MM 2C MAS RED
10028545	GSB SCREW, 8.5X45MM 2C MAS RED
10028550	GSB SCREW, 8.5X50MM 2C MAS RED



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
10050320	RELINE-O X-CONN, 20MM 5.0 ROD FIXED
10050322	RELINE-O X-CONN, 22.5MM 5.0 ROD FIXED
10050325	RELINE-O X-CONN, 25MM 5.0 ROD FIXED
10050327	RELINE-O X-CONN, 27.5MM 5.0 ROD FIXED
10050330	RELINE-O X-CONN, 30MM 5.0 ROD FIXED
10050332	RELINE-O X-CONN, 32.5MM 5.0 ROD FIXED
10050335	RELINE-O X-CONN, 35MM 5.0 ROD FIXED
10050337	RELINE-O X-CONN, 37.5MM 5.0 ROD FIXED
10050340	RELINE-O X-CONN, 40MM 5.0 ROD FIXED
10050430	RELINE-O X-CONN, 30-35MM 5.0 LP ADJ
10050435	RELINE-O X-CONN, 35-42MM 5.0 LP ADJ
10050440	RELINE-O X-CONN, 40-50MM 5.0 LP ADJ
10050445	RELINE-O X-CONN, 45-65MM 5.0 LP ADJ
10055020	GSB X-LINK, 20MM FIXED 5.5 ROD
10055022	GSB X-LINK, 22.5MM FIXED 5.5 ROD
10055025	GSB X-LINK, 25MM FIXED 5.5 ROD
10055027	GSB X-LINK, 27.5MM FIXED 5.5 ROD
10055030	GSB X-LINK, 30MM FIXED 5.5 ROD
10055032	GSB X-LINK, 32.5MM FIXED 5.5 ROD
10055035	GSB X-LINK, 35MM FIXED 5.5 ROD
10055037	GSB X-LINK, 37.5MM FIXED 5.5 ROD
10055040	GSB X-LINK, 40MM FIXED 5.5 ROD
10055045	GSB X-LINK, 45MM FIXED 5.5 ROD
10055050	GSB X-LINK, 50MM FIXED 5.5 ROD
10055055	GSB X-LINK, 55MM FIXED 5.5 ROD
10055060	GSB X-LINK, 60MM FIXED 5.5 ROD
10055145	GSB X-LINK, 45-50MM ADJ. 5.5 ROD
10055150	GSB X-LINK, 50-60MM ADJ. 5.5 ROD
10055160	GSB X-LINK, 60-70MM ADJ. 5.5 ROD
10055225	GSB X-LINK, 25-30MM ADJ. 5.5 ROD LP
10055230	GSB X-LINK, 30-35MM ADJ. 5.5 ROD LP
10055235	GSB X-LINK, 35-42MM ADJ. 5.5 ROD LP
10055240	GSB X-LINK, 40-50MM ADJ. 5.5 ROD LP
10055245	GSB X-LINK, 45-65MM ADJ. 5.5 ROD LP
10055320	RELINE-O X-CONN, 20MM 5.5 ROD FIXED
10055322	RELINE-O X-CONN, 22.5MM 5.5 ROD FIXED
10055325	RELINE-O X-CONN, 25MM 5.5 ROD FIXED
10055327	RELINE-O X-CONN, 27.5MM 5.5 ROD FIXED
10055330	RELINE-O X-CONN, 30MM 5.5 ROD FIXED

**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022****Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
10055332	RELINE-O X-CONN, 32.5MM 5.5 ROD FIXED
10055335	RELINE-O X-CONN, 35MM 5.5 ROD FIXED
10055337	RELINE-O X-CONN, 37.5MM 5.5 ROD FIXED
10055340	RELINE-O X-CONN, 40MM 5.5 ROD FIXED
10055345	RELINE-O X-CONN, 45MM 5.5 ROD FIXED
10055350	RELINE-O X-CONN, 50MM 5.5 ROD FIXED
10055355	RELINE-O X-CONN, 55MM 5.5 ROD FIXED
10055430	RELINE-O X-CONN, 30-35MM 5.5 LP ADJ
10055435	RELINE-O X-CONN, 35-42MM 5.5 LP ADJ
10055440	RELINE-O X-CONN, 40-50MM 5.5 LP ADJ
10055445	RELINE-O X-CONN, 45-65MM 5.5 LP ADJ
10055455	RELINE-O X-CONN, 55-80MM 5.5 LP ADJ
10056430	RELINE-O X-CON, 30-35MM 5.5/6.0 LP ADJ
10056435	RELINE-O X-CON, 35-42MM 5.5/6.0 LP ADJ
10056440	RELINE-O X-CON, 40-50MM 5.5/6.0 LP ADJ
10056445	RELINE-O X-CON, 45-65MM 5.5/6.0 LP ADJ
10056455	RELINE-O X-CON, 55-80MM 5.5/6.0 LP ADJ
10060020	GSB X-LINK, 20MM FIXED 6.0 ROD
10060022	GSB X-LINK, 22.5MM FIXED 6.0 ROD
10060025	GSB X-LINK, 25MM FIXED 6.0 ROD
10060027	GSB X-LINK, 27.5MM FIXED 6.0 ROD
10060030	GSB X-LINK, 30MM FIXED 6.0 ROD
10060032	GSB X-LINK, 32.5MM FIXED 6.0 ROD
10060035	GSB X-LINK, 35MM FIXED 6.0 ROD
10060037	GSB X-LINK, 37.5MM FIXED 6.0 ROD
10060040	GSB X-LINK, 40MM FIXED 6.0 ROD
10060045	GSB X-LINK, 45MM FIXED 6.0 ROD
10060050	GSB X-LINK, 50MM FIXED 6.0 ROD
10060055	GSB X-LINK, 55MM FIXED 6.0 ROD
10060060	GSB X-LINK, 60MM FIXED 6.0 ROD
10060320	RELINE-O X-CONN, 20MM 6.0 ROD FIXED
10060322	RELINE-O X-CONN, 22.5MM 6.0 ROD FIXED
10060325	RELINE-O X-CONN, 25MM 6.0 ROD FIXED
10060327	RELINE-O X-CONN, 27.5MM 6.0 ROD FIXED
10060330	RELINE-O X-CONN, 30MM 6.0 ROD FIXED
10060332	RELINE-O X-CONN, 32.5MM 6.0 ROD FIXED
10060335	RELINE-O X-CONN, 35MM 6.0 ROD FIXED
10060337	RELINE-O X-CONN, 37.5MM 6.0 ROD FIXED
10060340	RELINE-O X-CONN, 40MM 6.0 ROD FIXED



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
10060345	RELIN-O X-CONN, 45MM 6.0 ROD FIXED
10060430	RELIN-O X-CONN, 30-35MM 6.0 LP ADJ
10060435	RELIN-O X-CONN, 35-42MM 6.0 LP ADJ
10060440	RELIN-O X-CONN, 40-50MM 6.0 LP ADJ
10060445	RELIN-O X-CONN, 45-65MM 6.0 LP ADJ
10060455	RELIN-O X-CONN, 55-80MM 6.0 LP ADJ
10103516	RELIN HS SCREW, 3.5 X 16MM
10103520	RELIN HS SCREW, 3.5 X 20MM
10103524	RELIN HS SCREW, 3.5 X 24MM
10104016	RELIN HS SCREW, 4.0 X 16MM
10104020	RELIN HS SCREW, 4.0 X 20MM
10104024	RELIN HS SCREW, 4.0 X 24MM
10145311	RSS CROSS CONN, 11MM FIXED 4.5 ROD
10145313	RSS CROSS CONN, 13MM FIXED 4.5 ROD
10145315	RSS CROSS CONN, 15MM FIXED 4.5 ROD
10145317	RSS CROSS CONN, 17.5MM FIXED 4.5 ROD
10145320	RSS CROSS CONN, 20MM FIXED 4.5 ROD
10145322	RSS CROSS CONN, 22.5MM FIXED 4.5 ROD
10145325	RSS CROSS CONN, 25MM FIXED 4.5 ROD
10145327	RSS CROSS CONN, 27.5MM FIXED 4.5 ROD
10145330	RSS CROSS CONN, 30MM FIXED 4.5 ROD
10145332	RSS CROSS CONN, 32.5MM FIXED 4.5 ROD
10145335	RSS CROSS CONN, 35MM FIXED 4.5 ROD
10145337	RSS CROSS CONN, 37.5MM FIXED 4.5 ROD
10145340	RSS CROSS CONN, 40MM FIXED 4.5 ROD
10145342	RSS CROSS CONN, 42.5MM FIXED 4.5 ROD
10145345	RSS CROSS CONN, 45MM FIXED 4.5 ROD
10145347	RSS CROSS CONN, 47.5MM FIXED 4.5 ROD
10145350	RSS CROSS CONN, 50MM FIXED 4.5 ROD
10145352	RSS CROSS CONN, 52.5MM FIXED 4.5 ROD
10145355	RSS CROSS CONN, 55MM FIXED 4.5 ROD
10145357	RSS CROSS CONN, 57.5MM FIXED 4.5 ROD
10145360	RSS CROSS CONN, 60MM FIXED 4.5 ROD
10145362	RSS CROSS CONN, 62.5MM FIXED 4.5 ROD
10145365	RSS CROSS CONN, 65MM FIXED 4.5 ROD
10145367	RSS CROSS CONN, 67.5MM FIXED 4.5 ROD
10145370	RSS CROSS CONN, 70MM FIXED 4.5 ROD
10145372	RSS CROSS CONN, 72.5MM FIXED 4.5 ROD
10145375	RSS CROSS CONN, 75MM FIXED 4.5 ROD



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
10145377	RSS CROSS CONN, 77.5MM FIXED 4.5 ROD
10145380	RSS CROSS CONN, 80MM FIXED 4.5 ROD
10145422	RSS CROSS CONN, 22-25MM LP ADJ 4.5 ROD
10145425	RSS CROSS CONN, 25-30MM LP ADJ 4.5 ROD
10145430	RSS CROSS CONN, 30-40MM LP ADJ 4.5 ROD
10145440	RSS CROSS CONN, 40-60MM LP ADJ 4.5 ROD
10145460	RSS CROSS CONN, 60-100MM LP ADJ 4.5 ROD
10150311	RSS CROSS CONN, 11MM FIXED 5.0 ROD
10150313	RSS CROSS CONN, 13MM FIXED 5.0 ROD
10150315	RSS CROSS CONN, 15MM FIXED 5.0 ROD
10150317	RSS CROSS CONN, 17.5MM FIXED 5.0 ROD
10150320	RSS CROSS CONN, 20MM FIXED 5.0 ROD
10150322	RSS CROSS CONN, 22.5MM FIXED 5.0 ROD
10150325	RSS CROSS CONN, 25MM FIXED 5.0 ROD
10150327	RSS CROSS CONN, 27.5MM FIXED 5.0 ROD
10150330	RSS CROSS CONN, 30MM FIXED 5.0 ROD
10150332	RSS CROSS CONN, 32.5MM FIXED 5.0 ROD
10150335	RSS CROSS CONN, 35MM FIXED 5.0 ROD
10150422	RSS CROSS CONN, 22-25MM LP ADJ 5.0 ROD
10150425	RSS CROSS CONN, 25-30MM LP ADJ 5.0 ROD
10150430	RSS X-CONNECTOR, 30-40MM 5.0 ROD LP ADJ.
10150440	RSS X-CONNECTOR, 40-60MM 5.0 ROD LP ADJ.
10150460	RSS CROSS CONN, 60-100MM LP ADJ 5.0 ROD
10155025	GSB TI ROD, 5.5X25MM LORDOTIC PIVOT MAS
10155030	GSB TI ROD, 5.5X30MM LORDOTIC PIVOT MAS
10155035	GSB TI ROD, 5.5X35MM LORDOTIC PIVOT MAS
10155040	GSB TI ROD, 5.5X40MM LORDOTIC PIVOT MAS
10155045	GSB TI ROD, 5.5X45MM LORDOTIC PIVOT MAS
10155050	GSB TI ROD, 5.5X50MM LORDOTIC PIVOT MAS
10155055	GSB TI ROD, 5.5X55MM LORDOTIC PIVOT MAS
10155060	GSB TI ROD, 5.5X60MM LORDOTIC PIVOT MAS
10250108	RSS ILIAC SADDLE, 5.0MM STRAIGHT Ø8
10250110	RSS ILIAC SADDLE, 5.0MM STRAIGHT Ø10
10250112	RSS ILIAC SADDLE, 5.0MM STRAIGHT Ø12
10250114	RSS ILIAC SADDLE, 5.0MM STRAIGHT Ø14
10250116	RSS ILIAC SADDLE, 5.0MM STRAIGHT Ø16
10250118	RSS ILIAC SADDLE, 5.0MM STRAIGHT Ø18
10250120	RSS ILIAC SADDLE, 5.0MM STRAIGHT Ø20
10250208	RSS ILIAC SADDLE, 5.0MM ANTERIOR ANG Ø8



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
10250210	RSS ILIAC SADDLE, 5.0MM ANTERIOR ANG Ø10
10250212	RSS ILIAC SADDLE, 5.0MM ANTERIOR ANG Ø12
10250214	RSS ILIAC SADDLE, 5.0MM ANTERIOR ANG Ø14
10250216	RSS ILIAC SADDLE, 5.0MM ANTERIOR ANG Ø16
10250218	RSS ILIAC SADDLE, 5.0MM ANTERIOR ANG Ø18
10250220	RSS ILIAC SADDLE, 5.0MM ANTERIOR ANG Ø20
10250308	RSS ILIAC SADDLE, 5.0MM ANG LEFT Ø8
10250310	RSS ILIAC SADDLE, 5.0MM ANG LEFT Ø10
10250312	RSS ILIAC SADDLE, 5.0MM ANG LEFT Ø12
10250314	RSS ILIAC SADDLE, 5.0MM ANG LEFT Ø14
10250316	RSS ILIAC SADDLE, 5.0MM ANG LEFT Ø16
10250318	RSS ILIAC SADDLE, 5.0MM ANG LEFT Ø18
10250320	RSS ILIAC SADDLE, 5.0MM ANG LEFT Ø20
10250408	RSS ILIAC SADDLE, 5.0MM ANG RIGHT Ø8
10250410	RSS ILIAC SADDLE, 5.0MM ANG RIGHT Ø10
10250412	RSS ILIAC SADDLE, 5.0MM ANG RIGHT Ø12
10250414	RSS ILIAC SADDLE, 5.0MM ANG RIGHT Ø14
10250416	RSS ILIAC SADDLE, 5.0MM ANG RIGHT Ø16
10250418	RSS ILIAC SADDLE, 5.0MM ANG RIGHT Ø18
10250420	RSS ILIAC SADDLE, 5.0MM ANG RIGHT Ø20
10255108	RSS ILIAC SADDLE, 5.5MM STRAIGHT Ø8
10255110	RSS ILIAC SADDLE, 5.5MM STRAIGHT Ø10
10255112	RSS ILIAC SADDLE, 5.5MM STRAIGHT Ø12
10255114	RSS ILIAC SADDLE, 5.5MM STRAIGHT Ø14
10255116	RSS ILIAC SADDLE, 5.5MM STRAIGHT Ø16
10255118	RSS ILIAC SADDLE, 5.5MM STRAIGHT Ø18
10255120	RSS ILIAC SADDLE, 5.5MM STRAIGHT Ø20
10255208	RSS ILIAC SADDLE, 5.5MM ANT ANG Ø8
10255210	RSS ILIAC SADDLE, 5.5MM ANT ANG Ø10
10255212	RSS ILIAC SADDLE, 5.5MM ANT ANG Ø12
10255214	RSS ILIAC SADDLE, 5.5MM ANT ANG Ø14
10255216	RSS ILIAC SADDLE, 5.5MM ANT ANG Ø16
10255218	RSS ILIAC SADDLE, 5.5MM ANT ANG Ø18
10255220	RSS ILIAC SADDLE, 5.5MM ANT ANG Ø20
10255308	RSS ILIAC SADDLE, 5.5MM ANG LEFT Ø8
10255310	RSS ILIAC SADDLE, 5.5MM ANG LEFT Ø10
10255312	RSS ILIAC SADDLE, 5.5MM ANG LEFT Ø12
10255314	RSS ILIAC SADDLE, 5.5MM ANG LEFT Ø14
10255316	RSS ILIAC SADDLE, 5.5MM ANG LEFT Ø16



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
10255318	RSS ILIAC SADDLE, 5.5MM ANG LEFT Ø18
10255320	RSS ILIAC SADDLE, 5.5MM ANG LEFT Ø20
10255408	RSS ILIAC SADDLE, 5.5MM ANG RIGHT Ø8
10255410	RSS ILIAC SADDLE, 5.5MM ANG RIGHT Ø10
10255412	RSS ILIAC SADDLE, 5.5MM ANG RIGHT Ø12
10255414	RSS ILIAC SADDLE, 5.5MM ANG RIGHT Ø14
10255416	RSS ILIAC SADDLE, 5.5MM ANG RIGHT Ø16
10255418	RSS ILIAC SADDLE, 5.5MM ANG RIGHT Ø18
10255420	RSS ILIAC SADDLE, 5.5MM ANG RIGHT Ø20
10260108	RSS ILIAC SADDLE, 6.0MM STRAIGHT Ø8
10260110	RSS ILIAC SADDLE, 6.0MM STRAIGHT Ø10
10260112	RSS ILIAC SADDLE, 6.0MM STRAIGHT Ø12
10260114	RSS ILIAC SADDLE, 6.0MM STRAIGHT Ø14
10260116	RSS ILIAC SADDLE, 6.0MM STRAIGHT Ø16
10260118	RSS ILIAC SADDLE, 6.0MM STRAIGHT Ø18
10260120	RSS ILIAC SADDLE, 6.0MM STRAIGHT Ø20
10260208	RSS ILIAC SADDLE, 6.0MM ANT ANG Ø8
10260210	RSS ILIAC SADDLE, 6.0MM ANT ANG Ø10
10260212	RSS ILIAC SADDLE, 6.0MM ANT ANG Ø12
10260214	RSS ILIAC SADDLE, 6.0MM ANT ANG Ø14
10260216	RSS ILIAC SADDLE, 6.0MM ANT ANG Ø16
10260218	RSS ILIAC SADDLE, 6.0MM ANT ANG Ø18
10260220	RSS ILIAC SADDLE, 6.0MM ANT ANG Ø20
10260308	RSS ILIAC SADDLE, 6.0MM ANG LEFT Ø8
10260310	RSS ILIAC SADDLE, 6.0MM ANG LEFT Ø10
10260312	RSS ILIAC SADDLE, 6.0MM ANG LEFT Ø12
10260314	RSS ILIAC SADDLE, 6.0MM ANG LEFT Ø14
10260316	RSS ILIAC SADDLE, 6.0MM ANG LEFT Ø16
10260318	RSS ILIAC SADDLE, 6.0MM ANG LEFT Ø18
10260320	RSS ILIAC SADDLE, 6.0MM ANG LEFT Ø20
10260408	RSS ILIAC SADDLE, 6.0MM ANG RIGHT Ø8
10260410	RSS ILIAC SADDLE, 6.0MM ANG RIGHT Ø10
10260412	RSS ILIAC SADDLE, 6.0MM ANG RIGHT Ø12
10260414	RSS ILIAC SADDLE, 6.0MM ANG RIGHT Ø14
10260416	RSS ILIAC SADDLE, 6.0MM ANG RIGHT Ø16
10260418	RSS ILIAC SADDLE, 6.0MM ANG RIGHT Ø18
10260420	RSS ILIAC SADDLE, 6.0MM ANG RIGHT Ø20
11014025	GSB SCREW, 4.0X25MM 2S POLYAXIAL
11014030	GSB SCREW, 4.0X30MM 2S POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
11014035	GSB SCREW, 4.0X35MM 2S POLYAXIAL
11014040	GSB SCREW, 4.0X40MM 2S POLYAXIAL
11014525	GSB SCREW, 4.5X25MM 2S POLYAXIAL
11014530	GSB SCREW, 4.5X30MM 2S POLYAXIAL
11014535	GSB SCREW, 4.5X35MM 2S POLYAXIAL
11014540	GSB SCREW, 4.5X40MM 2S POLYAXIAL
11014545	GSB SCREW, 4.5X45MM 2S POLYAXIAL
11015025	GSB SCREW, 5.0X25MM 2S POLYAXIAL
11015030	GSB SCREW, 5.0X30MM 2S POLYAXIAL
11015035	GSB SCREW, 5.0X35MM 2S POLYAXIAL
11015040	GSB SCREW, 5.0X40MM 2S POLYAXIAL
11015045	GSB SCREW, 5.0X45MM 2S POLYAXIAL
11015525	GSB SCREW, 5.5X25MM 2S POLYAXIAL
11015530	GSB SCREW, 5.5X30MM 2S POLYAXIAL
11015535	GSB SCREW, 5.5X35MM 2S POLYAXIAL
11015540	GSB SCREW, 5.5X40MM 2S POLYAXIAL
11015545	GSB SCREW, 5.5X45MM 2S POLYAXIAL
11015550	GSB SCREW, 5.5X50MM 2S POLYAXIAL
11015555	GSB SCREW, 5.5X55MM 2S POLYAXIAL
11016530	GSB SCREW, 6.5X30MM 2S POLYAXIAL
11016535	GSB SCREW, 6.5X35MM 2S POLYAXIAL
11016540	GSB SCREW, 6.5X40MM 2S POLYAXIAL
11016545	GSB SCREW, 6.5X45MM 2S POLYAXIAL
11016550	GSB SCREW, 6.5X50MM 2S POLYAXIAL
11016555	GSB SCREW, 6.5X55MM 2S POLYAXIAL
11016560	GSB SCREW, 6.5X60MM 2S POLYAXIAL
11016570	GSB SCREW, 6.5X70MM 2S POLYAXIAL
11017530	GSB SCREW, 7.5X30MM 2S POLYAXIAL
11017535	GSB SCREW, 7.5X35MM 2S POLYAXIAL
11017540	GSB SCREW, 7.5X40MM 2S POLYAXIAL
11017545	GSB SCREW, 7.5X45MM 2S POLYAXIAL
11017550	GSB SCREW, 7.5X50MM 2S POLYAXIAL
11017555	GSB SCREW, 7.5X55MM 2S POLYAXIAL
11017560	GSB SCREW, 7.5X60MM 2S POLYAXIAL
11017570	GSB SCREW, 7.5X70MM 2S POLYAXIAL
11017580	GSB SCREW, 7.5X80MM 2S POLYAXIAL
11017590	GSB SCREW, 7.5X90MM 2S POLYAXIAL
11018510	GSB SCREW, 8.5X100MM 2S POLYAXIAL
11018540	GSB SCREW, 8.5X40MM 2S POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
11018545	GSB SCREW, 8.5X45MM 2S POLYAXIAL
11018550	GSB SCREW, 8.5X50MM 2S POLYAXIAL
11018555	GSB SCREW, 8.5X55MM 2S POLYAXIAL
11018560	GSB SCREW, 8.5X60MM 2S POLYAXIAL
11018570	GSB SCREW, 8.5X70MM 2S POLYAXIAL
11018580	GSB SCREW, 8.5X80MM 2S POLYAXIAL
11018590	GSB SCREW, 8.5X90MM 2S POLYAXIAL
11019510	GSB SCREW, 9.5X100MM 2S POLYAXIAL
11019560	GSB SCREW, 9.5X60MM 2S POLYAXIAL
11019570	GSB SCREW, 9.5X70MM 2S POLYAXIAL
11019580	GSB SCREW, 9.5X80MM 2S POLYAXIAL
11019590	GSB SCREW, 9.5X90MM 2S POLYAXIAL
11024035	GSB SCREW, 4.0X35MM 2S REDUCTION
11024530	GSB SCREW, 4.5X30MM 2S REDUCTION
11024535	GSB SCREW, 4.5X35MM 2S REDUCTION
11024540	GSB SCREW, 4.5X40MM 2S REDUCTION
11024545	GSB SCREW, 4.5X45MM 2S REDUCTION
11025530	GSB SCREW, 5.5X30MM 2S REDUCTION
11025535	GSB SCREW, 5.5X35MM 2S REDUCTION
11025540	GSB SCREW, 5.5X40MM 2S REDUCTION
11025545	GSB SCREW, 5.5X45MM 2S REDUCTION
11025550	GSB SCREW, 5.5X50MM 2S REDUCTION
11025555	GSB SCREW, 5.5X55MM 2S REDUCTION
11026535	GSB SCREW, 6.5X35MM 2S REDUCTION
11026540	GSB SCREW, 6.5X40MM 2S REDUCTION
11026545	GSB SCREW, 6.5X45MM 2S REDUCTION
11026550	GSB SCREW, 6.5X50MM 2S REDUCTION
11026555	GSB SCREW, 6.5X55MM 2S REDUCTION
11027535	GSB SCREW, 7.5X35MM 2S REDUCTION
11027540	GSB SCREW, 7.5X40MM 2S REDUCTION
11027545	GSB SCREW, 7.5X45MM 2S REDUCTION
11027550	GSB SCREW, 7.5X50MM 2S REDUCTION
11027555	GSB SCREW, 7.5X55MM 2S REDUCTION
11028535	GSB SCREW, 8.5X35MM 2S REDUCTION
11028540	GSB SCREW, 8.5X40MM 2S REDUCTION
11028545	GSB SCREW, 8.5X45MM 2S REDUCTION
11028550	GSB SCREW, 8.5X50MM 2S REDUCTION
11028555	GSB SCREW, 8.5X55MM 2S REDUCTION
11036560	GSB SCREW, 6.5X60MM 2S ILIAC CLOSED



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
11036570	GSB SCREW, 6.5X70MM 2S ILIAC CLOSED
11037560	GSB SCREW, 7.5X60MM 2S ILIAC CLOSED
11037570	GSB SCREW, 7.5X70MM 2S ILIAC CLOSED
11037580	GSB SCREW, 7.5X80MM 2S ILIAC CLOSED
11037590	GSB SCREW, 7.5X90MM 2S ILIAC CLOSED
11038510	GSB SCREW, 8.5X100MM 2S ILIAC CLOSED
11038560	GSB SCREW, 8.5X60MM 2S ILIAC CLOSED
11038570	GSB SCREW, 8.5X70MM 2S ILIAC CLOSED
11038580	GSB SCREW, 8.5X80MM 2S ILIAC CLOSED
11038590	GSB SCREW, 8.5X90MM 2S ILIAC CLOSED
11039510	GSB SCREW, 9.5X100MM 2S ILIAC CLOSED
11039560	GSB SCREW, 9.5X60MM 2S ILIAC CLOSED
11039570	GSB SCREW, 9.5X70MM 2S ILIAC CLOSED
11039580	GSB SCREW, 9.5X80MM 2S ILIAC CLOSED
11039590	GSB SCREW, 9.5X90MM 2S ILIAC CLOSED
11044530	GSB SCREW, 4.5X30MM 2S UNIPLANAR
11044535	GSB SCREW, 4.5X35MM 2S UNIPLANAR
11044540	GSB SCREW, 4.5X40MM 2S UNIPLANAR
11045030	GSB SCREW, 5.0X30MM 2S UNIPLANAR
11045035	GSB SCREW, 5.0X35MM 2S UNIPLANAR
11045040	GSB SCREW, 5.0X40MM 2S UNIPLANAR
11045045	GSB SCREW, 5.0X45MM 2S UNIPLANAR
11045530	GSB SCREW, 5.5X30MM 2S UNIPLANAR
11045535	GSB SCREW, 5.5X35MM 2S UNIPLANAR
11045540	GSB SCREW, 5.5X40MM 2S UNIPLANAR
11045545	GSB SCREW, 5.5X45MM 2S UNIPLANAR
11045550	GSB SCREW, 5.5X50MM 2S UNIPLANAR
11045555	GSB SCREW, 5.5X55MM 2S UNIPLANAR
11046530	GSB SCREW, 6.5X30MM 2S UNIPLANAR
11046535	GSB SCREW, 6.5X35MM 2S UNIPLANAR
11046540	GSB SCREW, 6.5X40MM 2S UNIPLANAR
11046545	GSB SCREW, 6.5X45MM 2S UNIPLANAR
11046550	GSB SCREW, 6.5X50MM 2S UNIPLANAR
11046555	GSB SCREW, 6.5X55MM 2S UNIPLANAR
11084530	GSB SCREW, 4.5X30MM 2S MONOAXIAL
11084535	GSB SCREW, 4.5X35MM 2S MONOAXIAL
11084540	GSB SCREW, 4.5X40MM 2S MONOAXIAL
11085530	GSB SCREW, 5.5X30MM 2S MONOAXIAL
11085535	GSB SCREW, 5.5X35MM 2S MONOAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
11085540	GSB SCREW, 5.5X40MM 2S MONOAXIAL
11085545	GSB SCREW, 5.5X45MM 2S MONOAXIAL
11085550	GSB SCREW, 5.5X50MM 2S MONOAXIAL
11086535	GSB SCREW, 6.5X35MM 2S MONOAXIAL
11086540	GSB SCREW, 6.5X40MM 2S MONOAXIAL
11086545	GSB SCREW, 6.5X45MM 2S MONOAXIAL
11086550	GSB SCREW, 6.5X50MM 2S MONOAXIAL
11087540	GSB SCREW, 7.5X40MM 2S MONOAXIAL
11087545	GSB SCREW, 7.5X45MM 2S MONOAXIAL
11087550	GSB SCREW, 7.5X50MM 2S MONOAXIAL
11087555	GSB SCREW, 7.5X55MM 2S MONOAXIAL
11095530	GSB SCREW, 5.5X30MM 4S POLYAXIAL
11095535	GSB SCREW, 5.5X35MM 4S POLYAXIAL
11095540	GSB SCREW, 5.5X40MM 4S POLYAXIAL
11095545	GSB SCREW, 5.5X45MM 4S POLYAXIAL
11095550	GSB SCREW, 5.5X50MM 4S POLYAXIAL
11096535	GSB SCREW, 6.5X35MM 4S POLYAXIAL
11096540	GSB SCREW, 6.5X40MM 4S POLYAXIAL
11096545	GSB SCREW, 6.5X45MM 4S POLYAXIAL
11096550	GSB SCREW, 6.5X50MM 4S POLYAXIAL
11096555	GSB SCREW, 6.5X55MM 4S POLYAXIAL
11097535	GSB SCREW, 7.5X35MM 4S POLYAXIAL
11097540	GSB SCREW, 7.5X40MM 4S POLYAXIAL
11097545	GSB SCREW, 7.5X45MM 4S POLYAXIAL
11097550	GSB SCREW, 7.5X50MM 4S POLYAXIAL
11097555	GSB SCREW, 7.5X55MM 4S POLYAXIAL
11100020	GSB CONN, 20MM 5.0/6.0 ROD OPEN OFFSET
11100030	GSB CONN, 30MM 5.0/6.0 ROD OPEN OFFSET
11100060	GSB CONN, 60MM 5.0/6.0 ROD OPEN OFFSET
11100120	GSB CONN, 20MM 5.0/6.0 ROD CLOSED OFFSET
11100130	GSB CONN, 30MM 5.0/6.0 ROD CLOSED OFFSET
11100160	GSB CONN, 60MM 5.0/6.0 ROD CLOSED OFFSET
11101255	RELIN-O CONN, 5-6/5-6MM O-O MED
11101260	RELIN-O CONN, 5-6/6.35MM O-O MED
11101355	RELIN-O CONN, 5-6/5-6MM O-O LG
11102255	RELIN-O CONN, 5-6/5-6MM O-H MED
11102260	RELIN-O CONN, 5-6/6.35MM O-H MED
11102355	RELIN-O CONN, 5-6/5-6MM O-H LG
11103055	RELIN-O CONN, 5-6/5-6MM ROTATING O-O



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
11103060	RELINE-O CONN, 5-6/6.35MM ROTATING O-O
11104055	RELINE-O CONN, 5-6/5-6MM O-H ROTATING
11104060	RELINE-O CONN, 5-6/6.35MM O-H ROTATING
11106155	RELINE-O CONN, 5-6/5-6MM 2H PARALLEL SM
11106160	RELINE-O CONN, 6.0/6.35MM 2H PARALLEL SM
11106255	RELINE-O CONN, 5-6/5-6MM 2H PARALLEL MED
11106260	RELINE-O CONN, 6.0/6.35MM 2H PARALLEL MED
11106355	RELINE-O CONN, 5-6/5-6MM 2H PARALLEL LG
11107155	RELINE-O CONN, 5-6/5-6MM 4H PARALLEL SM
11107160	RELINE-O CONN, 6.0/6.35MM 4H PARALLEL SM
11107255	RELINE-O CONN, 5-6/5-6MM 4H PARALLEL MED
11107260	RELINE-O CONN, 6.0/6.35MM 4H PARALLEL MED
11108055	RELINE-O CONN, 5-6/5-6MM 2H IN-LINE
11108060	RELINE-O CONN, 5-6/6.35MM 2H IN-LINE
11109055	RELINE-O CONN, 5-6/5-6MM 4H IN-LINE
11109060	RELINE-O CONN, 6.0/6.35MM 4H IN-LINE
11111000	GSB LOCK SCREW, CLOSED TULIP
11111550	GSB LOCK SCREW, 5.5MM OPEN TULIP
11111600	GSB LOCK SCREW, 6.0MM OPEN TULIP
11155025	GSB TI ROD, 5.5X25MM LORDOTIC MAS
11155030	GSB TI ROD, 5.5X30MM LORDOTIC MAS
11155035	GSB TI ROD, 5.5X35MM LORDOTIC MAS
11155040	GSB TI ROD, 5.5X40MM LORDOTIC MAS
11155045	GSB TI ROD, 5.5X45MM LORDOTIC MAS
11155050	GSB TI ROD, 5.5X50MM LORDOTIC MAS
11155055	GSB TI ROD, 5.5X55MM LORDOTIC MAS
11155060	GSB TI ROD, 5.5X60MM LORDOTIC MAS
11155065	GSB TI ROD, 5.5X65MM LORDOTIC MAS
11155070	GSB TI ROD, 5.5X70MM LORDOTIC MAS
11155075	GSB TI ROD, 5.5X75MM LORDOTIC MAS
11155080	GSB TI ROD, 5.5X80MM LORDOTIC MAS
11155085	GSB TI ROD, 5.5X85MM LORDOTIC MAS
11155090	GSB TI ROD, 5.5X90MM LORDOTIC MAS
11155095	GSB TI ROD, 5.5X95MM LORDOTIC MAS
11155100	GSB TI ROD, 5.5X100MM LORDOTIC MAS
11155110	GSB TI ROD, 5.5X110MM LORDOTIC MAS
11155120	GSB TI ROD, 5.5X120MM LORDOTIC MAS
11155130	GSB TI ROD, 5.5X130MM LORDOTIC MAS
11155140	GSB TI ROD, 5.5X140MM LORDOTIC MAS



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
11155150	GSB TI ROD, 5.5X150MM LORDOTIC MAS
11155160	GSB TI ROD, 5.5X160MM LORDOTIC MAS
11200020	RELINE-O CONNECTOR, 20MM OPEN OFFSET
11200030	RELINE-O CONNECTOR, 30MM OPEN OFFSET
11200060	RELINE-O CONNECTOR, 60MM OPEN OFFSET
11200120	RELINE-O CONNECTOR, 20MM CLOSED OFFSET
11200130	RELINE-O CONNECTOR, 30MM CLOSED OFFSET
11200160	RELINE-O CONNECTOR, 60MM CLOSED OFFSET
11200220	RELINE-O CONNECTOR, 200MM OPEN OFFSET
11200410	RELINE-O CONN, 10MM OFFSET WITH TAB
11200420	RELINE-O CONN, 20MM OFFSET WITH TAB
11200430	RELINE-O CONN, 30MM OFFSET WITH TAB
11200440	RELINE-O CONN, 40MM OFFSET WITH TAB
11200460	RELINE-O CONN, 60MM OFFSET WITH TAB
11355025	RELINE MAS TI ROD, 5.5X25MM LORDOTIC
11355030	RELINE MAS TI ROD, 5.5X30MM LORDOTIC
11355035	RELINE MAS TI ROD, 5.5X35MM LORDOTIC
11355040	RELINE MAS TI ROD, 5.5X40MM LORDOTIC
11355045	RELINE MAS TI ROD, 5.5X45MM LORDOTIC
11355050	RELINE MAS TI ROD, 5.5X50MM LORDOTIC
11355055	RELINE MAS TI ROD, 5.5X55MM LORDOTIC
11355060	RELINE MAS TI ROD, 5.5X60MM LORDOTIC
11355065	RELINE MAS TI ROD, 5.5X65MM LORDOTIC
11355070	RELINE MAS TI ROD, 5.5X70MM LORDOTIC
11355075	RELINE MAS TI ROD, 5.5X75MM LORDOTIC
11355080	RELINE MAS TI ROD, 5.5X80MM LORDOTIC
11355085	RELINE MAS TI ROD, 5.5X85MM LORDOTIC
11355090	RELINE MAS TI ROD, 5.5X90MM LORDOTIC
11355095	RELINE MAS TI ROD, 5.5X95MM LORDOTIC
11355100	RELINE MAS TI ROD, 5.5X100MM LORDOTIC
11355105	RELINE MAS TI ROD, 5.5X105MM LORDOTIC
11355110	RELINE MAS TI ROD, 5.5X110MM LORDOTIC
11355120	RELINE MAS TI ROD, 5.5X120MM LORDOTIC
11355130	RELINE MAS TI ROD, 5.5X130MM LORDOTIC
11355140	RELINE MAS TI ROD, 5.5X140MM LORDOTIC
11355150	RELINE MAS TI ROD, 5.5X150MM LORDOTIC
11355160	RELINE MAS TI ROD, 5.5X160MM LORDOTIC
11360025	RELINE MAS TI ROD, 6.0X25MM LORDOTIC
11360030	RELINE MAS TI ROD, 6.0X30MM LORDOTIC



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
11360035	RELINE MAS TI ROD, 6.0X35MM LORDOTIC
11360040	RELINE MAS TI ROD, 6.0X40MM LORDOTIC
11360045	RELINE MAS TI ROD, 6.0X45MM LORDOTIC
11360050	RELINE MAS TI ROD, 6.0X50MM LORDOTIC
11360055	RELINE MAS TI ROD, 6.0X55MM LORDOTIC
11360060	RELINE MAS TI ROD, 6.0X60MM LORDOTIC
11360065	RELINE MAS TI ROD, 6.0X65MM LORDOTIC
11360070	RELINE MAS TI ROD, 6.0X70MM LORDOTIC
11360075	RELINE MAS TI ROD, 6.0X75MM LORDOTIC
11360080	RELINE MAS TI ROD, 6.0X80MM LORDOTIC
11360085	RELINE MAS TI ROD, 6.0X85MM LORDOTIC
11360090	RELINE MAS TI ROD, 6.0X90MM LORDOTIC
11360095	RELINE MAS TI ROD, 6.0X95MM LORDOTIC
11360100	RELINE MAS TI ROD, 6.0X100MM LORDOTIC
11360110	RELINE MAS TI ROD, 6.0X110MM LORDOTIC
11360120	RELINE MAS TI ROD, 6.0X120MM LORDOTIC
11360140	RELINE MAS TI ROD, 6.0X140MM LORDOTIC
11360160	RELINE MAS TI ROD, 6.0X160MM LORDOTIC
11455025	RELINE MAS COCR ROD, 5.5X25MM LORDOTIC
11455030	RELINE MAS COCR ROD, 5.5X30MM LORDOTIC
11455035	RELINE MAS COCR ROD, 5.5X35MM LORDOTIC
11455040	RELINE MAS COCR ROD, 5.5X40MM LORDOTIC
11455045	RELINE MAS COCR ROD, 5.5X45MM LORDOTIC
11455050	RELINE MAS COCR ROD, 5.5X50MM LORDOTIC
11455055	RELINE MAS COCR ROD, 5.5X55MM LORDOTIC
11455060	RELINE MAS COCR ROD, 5.5X60MM LORDOTIC
11455065	RELINE MAS COCR ROD, 5.5X65MM LORDOTIC
11455070	RELINE MAS COCR ROD, 5.5X70MM LORDOTIC
11455080	RELINE MAS COCR ROD, 5.5X80MM LORDOTIC
11455090	RELINE MAS COCR ROD, 5.5X90MM LORDOTIC
11455095	RELINE MAS COCR ROD, 5.5X95MM LORDOTIC
11455100	RELINE MAS COCR ROD, 5.5X100MM LORDOTIC
11455110	RELINE MAS COCR ROD, 5.5X110MM LORDOTIC
11455120	RELINE MAS COCR ROD, 5.5X120MM LORDOTIC
11455130	RELINE MAS COCR ROD, 5.5X130MM LORDOTIC
11455140	RELINE MAS COCR ROD, 5.5X140MM LORDOTIC
11455150	RELINE MAS COCR ROD, 5.5X150MM LORDOTIC
11455160	RELINE MAS COCR ROD, 5.5X160MM LORDOTIC
11560160	RELINE-O CONN, OFFSET 6.0MM CLOSED 60MM



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
12100150	RSS CONN, 4.5/5.0-4.5/5.0MM ROD O-O
12100160	RSS CONN, 4.5/5.0-5.0/6.0MM ROD O-O
12100260	RSS CONN, 4.5/5.0-5.0/6.0MM ROD O-H
12100550	RSS CONN, 4.5/5.0-4.5/5.0MM ROD O-C
12100650	RSS CONN,4.5/5.0-4.5/5.0MM ROD 2H PAR SM
12100660	RSS CONN,4.5/5.0-5.0/6.0MM ROD 2H PAR SM
12100760	RSS CONN,4.5/5.0-5.0/6.0MM ROD 4H PAR SM
12100850	RSS CONN,4.5/5.0-4.5/5.0MM ROD 2H INLINE
12100950	RSS CONN,4.5/5.0-4.5/5.0MM ROD 4H INLINE
12101260	RSS CONNECTOR, 60MM OFFSET OPEN
12101360	RSS CONNECTOR, 60MM OFFSET CLOSED
12101450	RSS CONNECTOR, 50MM IN-LINE EXTENDED
12101470	RSS CONNECTOR, 70MM IN-LINE EXTENDED
12101490	RSS CONNECTOR, 90MM IN-LINE EXTENDED
12155300	GSB TI ROD, 5.5X300MM STRAIGHT MAS
12155500	GSB TI ROD, 5.5X500MM STRAIGHT MAS
12255150	GSB COCR ROD, 5.5X150MM STRAIGHT MAS
12255300	GSB COCR ROD, 5.5X300MM STRAIGHT MAS
12255500	GSB COCR ROD, 5.5X500MM STRAIGHT MAS
12355040	RELIN MAS TI ROD, 5.5X40MM STRAIGHT
12355045	RELIN MAS TI ROD, 5.5X45MM STRAIGHT
12355050	RELIN MAS TI ROD, 5.5X50MM STRAIGHT
12355055	RELIN MAS TI ROD, 5.5X55MM STRAIGHT
12355060	RELIN MAS TI ROD, 5.5X60MM STRAIGHT
12355065	RELIN MAS TI ROD, 5.5X65MM STRAIGHT
12355070	RELIN MAS TI ROD, 5.5X70MM STRAIGHT
12355075	RELIN MAS TI ROD, 5.5X75MM STRAIGHT
12355080	RELIN MAS TI ROD, 5.5X80MM STRAIGHT
12355085	RELIN MAS TI ROD, 5.5X85MM STRAIGHT
12355090	RELIN MAS TI ROD, 5.5X90MM STRAIGHT
12355095	RELIN MAS TI ROD, 5.5X95MM STRAIGHT
12355100	RELIN MAS TI ROD, 5.5X100MM STRAIGHT
12355110	RELIN MAS TI ROD, 5.5X110MM STRAIGHT
12355120	RELIN MAS TI ROD, 5.5X120MM STRAIGHT
12355130	RELIN MAS TI ROD, 5.5X130MM STRAIGHT
12355140	RELIN MAS TI ROD, 5.5X140MM STRAIGHT
12355150	RELIN MAS TI ROD, 5.5X150MM STRAIGHT
12355160	RELIN MAS TI ROD, 5.5X160MM STRAIGHT
12355170	RELIN MAS TI ROD, 5.5X170MM STRAIGHT



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
12355180	RELIN MAS TI ROD, 5.5X180MM STRAIGHT
12355190	RELIN MAS TI ROD, 5.5X190MM STRAIGHT
12355200	RELIN MAS TI ROD, 5.5X200MM STRAIGHT
12355210	RELIN MAS TI ROD, 5.5X210MM STRAIGHT
12355220	RELIN MAS TI ROD, 5.5X220MM STRAIGHT
12355230	RELIN MAS TI ROD, 5.5X230MM STRAIGHT
12355240	RELIN MAS TI ROD, 5.5X240MM STRAIGHT
12355250	RELIN MAS TI ROD, 5.5X250MM STRAIGHT
12355300	RELIN MAS TI ROD, 5.5X300MM STRAIGHT
12355500	RELIN MAS TI ROD, 5.5X500MM STRAIGHT
12360150	RELIN MAS TI ROD, 6.0X150MM STRAIGHT
12360200	RELIN MAS TI ROD, 6.0X200MM STRAIGHT
12360300	RELIN MAS TI ROD, 6.0X300MM STRAIGHT
12360500	RELIN MAS TI ROD, 6.0X500MM STRAIGHT
12455070	RELIN MAS COCR ROD, 5.5X70MM STRAIGHT
12455075	RELIN MAS COCR ROD, 5.5X75MM STRAIGHT
12455080	RELIN MAS COCR ROD, 5.5X80MM STRAIGHT
12455140	RELIN MAS COCR ROD, 5.5X140MM STRAIGHT
12455150	RELIN MAS COCR ROD, 5.5X150MM STRAIGHT
12455200	RELIN MAS COCR ROD, 5.5X200MM STRAIGHT
12455300	RELIN MAS COCR ROD, 5.5X300MM STRAIGHT
12455500	RELIN MAS COCR ROD, 5.5X500MM STRAIGHT
12455600	RELIN MAS COCR ROD, 5.5X600MM STRAIGHT
12460500	RELIN MAS COCR ROD, 6.0X500MM STRAIGHT
13000000	RELIN LOCK SCREW, CLOSED TULIP
13010510	RELIN-O SCREW, 10.5X100MM 2S POLY ILIAC
13010511	RELIN-O SCREW, 10.5X110MM 2S POLY ILIAC
13010512	RELIN-O SCREW, 10.5X120MM 2S POLY ILIAC
13010535	RELIN-O SCREW, 10.5X35MM 2S POLYAXIAL
13010540	RELIN-O SCREW, 10.5X40MM 2S POLYAXIAL
13010545	RELIN-O SCREW, 10.5X45MM 2S POLYAXIAL
13010550	RELIN-O SCREW, 10.5X50MM 2S POLYAXIAL
13010555	RELIN-O SCREW, 10.5X55MM 2S POLYAXIAL
13010560	RELIN-O SCREW, 10.5X60MM 2S POLY ILIAC
13010570	RELIN-O SCREW, 10.5X70MM 2S POLY ILIAC
13010580	RELIN-O SCREW, 10.5X80MM 2S POLY ILIAC
13010590	RELIN-O SCREW, 10.5X90MM 2S POLY ILIAC
13011510	RELIN-O SCREW, 11.5X100MM 2S POLY ILIAC
13011545	RELIN-O SCREW, 11.5X45MM 2S POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13011550	RELIN-O SCREW, 11.5X50MM 2S POLYAXIAL
13011570	RELIN-O SCREW, 11.5X70MM 2S POLY ILIAC
13011580	RELIN-O SCREW, 11.5X80MM 2S POLY ILIAC
13011590	RELIN-O SCREW, 11.5X90MM 2S POLY ILIAC
13014020	RELIN-O SCREW, 4.0X20MM 2S POLYAXIAL
13014025	RELIN-O SCREW, 4.0X25MM 2S POLYAXIAL
13014030	RELIN-O SCREW, 4.0X30MM 2S POLYAXIAL
13014035	RELIN-O SCREW, 4.0X35MM 2S POLYAXIAL
13014040	RELIN-O SCREW, 4.0X40MM 2S POLYAXIAL
13014045	RELIN-O SCREW, 4.0X45MM 2S POLYAXIAL
13014050	RELIN-O SCREW, 4.0X50MM 2S POLYAXIAL
13014520	RELIN-O SCREW, 4.5X20MM 2S POLYAXIAL
13014525	RELIN-O SCREW, 4.5X25MM 2S POLYAXIAL
13014530	RELIN-O SCREW, 4.5X30MM 2S POLYAXIAL
13014535	RELIN-O SCREW, 4.5X35MM 2S POLYAXIAL
13014540	RELIN-O SCREW, 4.5X40MM 2S POLYAXIAL
13014545	RELIN-O SCREW, 4.5X45MM 2S POLYAXIAL
13014550	RELIN-O SCREW, 4.5X50MM 2S POLYAXIAL
13014555	RELIN-O SCREW, 4.5X55MM 2S POLYAXIAL
13015020	RELIN-O SCREW, 5.0X20MM 2S POLYAXIAL
13015025	RELIN-O SCREW, 5.0X25MM 2S POLYAXIAL
13015030	RELIN-O SCREW, 5.0X30MM 2S POLYAXIAL
13015035	RELIN-O SCREW, 5.0X35MM 2S POLYAXIAL
13015040	RELIN-O SCREW, 5.0X40MM 2S POLYAXIAL
13015045	RELIN-O SCREW, 5.0X45MM 2S POLYAXIAL
13015050	RELIN-O SCREW, 5.0X50MM 2S POLYAXIAL
13015055	RELIN-O SCREW, 5.0X55MM 2S POLYAXIAL
13015060	RELIN-O SCREW, 5.0X60MM 2S POLY ILIAC
13015520	RELIN-O SCREW, 5.5X20MM 2S POLYAXIAL
13015525	RELIN-O SCREW, 5.5X25MM 2S POLYAXIAL
13015530	RELIN-O SCREW, 5.5X30MM 2S POLYAXIAL
13015535	RELIN-O SCREW, 5.5X35MM 2S POLYAXIAL
13015540	RELIN-O SCREW, 5.5X40MM 2S POLYAXIAL
13015545	RELIN-O SCREW, 5.5X45MM 2S POLYAXIAL
13015550	RELIN-O SCREW, 5.5X50MM 2S POLYAXIAL
13015555	RELIN-O SCREW, 5.5X55MM 2S POLYAXIAL
13015560	RELIN-O SCREW, 5.5X60MM 2S POLY ILIAC
13015565	RELIN-O SCREW, 5.5X65MM 2S POLY ILIAC
13016030	RELIN-O SCREW, 6.0X30MM 2S POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13016035	RELIN-O SCREW, 6.0X35MM 2S POLYAXIAL
13016040	RELIN-O SCREW, 6.0X40MM 2S POLYAXIAL
13016045	RELIN-O SCREW, 6.0X45MM 2S POLYAXIAL
13016050	RELIN-O SCREW, 6.0X50MM 2S POLYAXIAL
13016055	RELIN-O SCREW, 6.0X55MM 2S POLYAXIAL
13016060	RELIN-O SCREW, 6.0X60MM 2S POLY ILIAC
13016510	RELIN-O SCREW, 6.5X100MM 2S POLY ILIAC
13016520	RELIN-O SCREW, 6.5X20MM 2S POLYAXIAL
13016525	RELIN-O SCREW, 6.5X25MM 2S POLYAXIAL
13016530	RELIN-O SCREW, 6.5X30MM 2S POLYAXIAL
13016535	RELIN-O SCREW, 6.5X35MM 2S POLYAXIAL
13016540	RELIN-O SCREW, 6.5X40MM 2S POLYAXIAL
13016545	RELIN-O SCREW, 6.5X45MM 2S POLYAXIAL
13016550	RELIN-O SCREW, 6.5X50MM 2S POLYAXIAL
13016555	RELIN-O SCREW, 6.5X55MM 2S POLYAXIAL
13016560	RELIN-O SCREW, 6.5X60MM 2S POLY ILIAC
13016565	RELIN-O SCREW, 6.5X65MM 2S POLY ILIAC
13016570	RELIN-O SCREW, 6.5X70MM 2S POLY ILIAC
13016580	RELIN-O SCREW, 6.5X80MM 2S POLY ILIAC
13016590	RELIN-O SCREW, 6.5X90MM 2S POLY ILIAC
13017030	RELIN-O SCREW, 7.0X30MM 2S POLYAXIAL
13017035	RELIN-O SCREW, 7.0X35MM 2S POLYAXIAL
13017040	RELIN-O SCREW, 7.0X40MM 2S POLYAXIAL
13017045	RELIN-O SCREW, 7.0X45MM 2S POLYAXIAL
13017050	RELIN-O SCREW, 7.0X50MM 2S POLYAXIAL
13017055	RELIN-O SCREW, 7.0X55MM 2S POLYAXIAL
13017060	RELIN-O SCREW, 7.0X60MM 2S POLY ILIAC
13017510	RELIN-O SCREW, 7.5X100MM 2S POLY ILIAC
13017525	RELIN-O SCREW, 7.5X25MM 2S POLYAXIAL
13017530	RELIN-O SCREW, 7.5X30MM 2S POLYAXIAL
13017535	RELIN-O SCREW, 7.5X35MM 2S POLYAXIAL
13017540	RELIN-O SCREW, 7.5X40MM 2S POLYAXIAL
13017545	RELIN-O SCREW, 7.5X45MM 2S POLYAXIAL
13017550	RELIN-O SCREW, 7.5X50MM 2S POLYAXIAL
13017555	RELIN-O SCREW, 7.5X55MM 2S POLYAXIAL
13017560	RELIN-O SCREW, 7.5X60MM 2S POLY ILIAC
13017565	RELIN-O SCREW, 7.5X65MM 2S POLY ILIAC
13017570	RELIN-O SCREW, 7.5X70MM 2S POLY ILIAC
13017580	RELIN-O SCREW, 7.5X80MM 2S POLY ILIAC



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13017590	RELIN-O SCREW, 7.5X90MM 2S POLY ILIAC
13018510	RELIN-O SCREW, 8.5X100MM 2S POLY ILIAC
13018511	RELIN-O SCREW, 8.5X110MM 2S POLY ILIAC
13018512	RELIN-O SCREW, 8.5X120MM 2S POLY ILIAC
13018530	RELIN-O SCREW, 8.5X30MM 2S POLYAXIAL
13018535	RELIN-O SCREW, 8.5X35MM 2S POLYAXIAL
13018540	RELIN-O SCREW, 8.5X40MM 2S POLYAXIAL
13018545	RELIN-O SCREW, 8.5X45MM 2S POLYAXIAL
13018550	RELIN-O SCREW, 8.5X50MM 2S POLYAXIAL
13018555	RELIN-O SCREW, 8.5X55MM 2S POLYAXIAL
13018560	RELIN-O SCREW, 8.5X60MM 2S POLY ILIAC
13018565	RELIN-O SCREW, 8.5X65MM 2S POLY ILIAC
13018570	RELIN-O SCREW, 8.5X70MM 2S POLY ILIAC
13018580	RELIN-O SCREW, 8.5X80MM 2S POLY ILIAC
13018590	RELIN-O SCREW, 8.5X90MM 2S POLY ILIAC
13019510	RELIN-O SCREW, 9.5X100MM 2S POLY ILIAC
13019511	RELIN-O SCREW, 9.5X110MM 2S POLY ILIAC
13019512	RELIN-O SCREW, 9.5X120MM 2S POLY ILIAC
13019530	RELIN-O SCREW, 9.5X30MM 2S POLYAXIAL
13019535	RELIN-O SCREW, 9.5X35MM 2S POLYAXIAL
13019540	RELIN-O SCREW, 9.5X40MM 2S POLYAXIAL
13019545	RELIN-O SCREW, 9.5X45MM 2S POLYAXIAL
13019550	RELIN-O SCREW, 9.5X50MM 2S POLYAXIAL
13019555	RELIN-O SCREW, 9.5X55MM 2S POLYAXIAL
13019560	RELIN-O SCREW, 9.5X60MM 2S POLY ILIAC
13019565	RELIN-O SCREW, 9.5X65MM 2S POLY ILIAC
13019570	RELIN-O SCREW, 9.5X70MM 2S POLY ILIAC
13019580	RELIN-O SCREW, 9.5X80MM 2S POLY ILIAC
13019590	RELIN-O SCREW, 9.5X90MM 2S POLY ILIAC
13024030	RELIN-O SCREW, 4.0X30MM 2S REDUCTION
13024035	RELIN-O SCREW, 4.0X35MM 2S REDUCTION
13024040	RELIN-O SCREW, 4.0X40MM 2S REDUCTION
13024525	RELIN-O SCREW, 4.5X25MM 2S REDUCTION
13024530	RELIN-O SCREW, 4.5X30MM 2S REDUCTION
13024535	RELIN-O SCREW, 4.5X35MM 2S REDUCTION
13024540	RELIN-O SCREW, 4.5X40MM 2S REDUCTION
13024545	RELIN-O SCREW, 4.5X45MM 2S REDUCTION
13024550	RELIN-O SCREW, 4.5X50MM 2S REDUCTION
13024555	RELIN-O SCREW, 4.5X55MM 2S REDUCTION



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13025030	RELIN-O SCREW, 5.0X30MM 2S REDUCTION
13025035	RELIN-O SCREW, 5.0X35MM 2S REDUCTION
13025040	RELIN-O SCREW, 5.0X40MM 2S REDUCTION
13025525	RELIN-O SCREW, 5.5X25MM 2S REDUCTION
13025530	RELIN-O SCREW, 5.5X30MM 2S REDUCTION
13025535	RELIN-O SCREW, 5.5X35MM 2S REDUCTION
13025540	RELIN-O SCREW, 5.5X40MM 2S REDUCTION
13025545	RELIN-O SCREW, 5.5X45MM 2S REDUCTION
13025550	RELIN-O SCREW, 5.5X50MM 2S REDUCTION
13025555	RELIN-O SCREW, 5.5X55MM 2S REDUCTION
13026530	RELIN-O SCREW, 6.5X30MM 2S REDUCTION
13026535	RELIN-O SCREW, 6.5X35MM 2S REDUCTION
13026540	RELIN-O SCREW, 6.5X40MM 2S REDUCTION
13026545	RELIN-O SCREW, 6.5X45MM 2S REDUCTION
13026550	RELIN-O SCREW, 6.5X50MM 2S REDUCTION
13026555	RELIN-O SCREW, 6.5X55MM 2S REDUCTION
13026560	RELIN-O SCREW, 6.5X60MM 2S REDUCTION
13027530	RELIN-O SCREW, 7.5X30MM 2S REDUCTION
13027535	RELIN-O SCREW, 7.5X35MM 2S REDUCTION
13027540	RELIN-O SCREW, 7.5X40MM 2S REDUCTION
13027545	RELIN-O SCREW, 7.5X45MM 2S REDUCTION
13027550	RELIN-O SCREW, 7.5X50MM 2S REDUCTION
13027555	RELIN-O SCREW, 7.5X55MM 2S REDUCTION
13027560	RELIN-O SCREW, 7.5X60MM 2S REDUCTION
13028512	RELIN-O SCREW, 8.5X120MM 2S REDUCTION
13028535	RELIN-O SCREW, 8.5X35MM 2S REDUCTION
13028540	RELIN-O SCREW, 8.5X40MM 2S REDUCTION
13028545	RELIN-O SCREW, 8.5X45MM 2S REDUCTION
13028550	RELIN-O SCREW, 8.5X50MM 2S REDUCTION
13028555	RELIN-O SCREW, 8.5X55MM 2S REDUCTION
13028560	RELIN-O SCREW, 8.5X60MM 2S REDUCTION
13036510	RELIN-O SCREW, 6.5X100MM 2S ILIAC CLSD
13036560	RELIN-O SCREW, 6.5X60MM 2S ILIAC CLOSED
13036570	RELIN-O SCREW, 6.5X70MM 2S ILIAC CLOSED
13036580	RELIN-O SCREW, 6.5X80MM 2S ILIAC CLOSED
13036590	RELIN-O SCREW, 6.5X90MM 2S ILIAC CLOSED
13037510	RELIN-O SCREW, 7.5X100MM 2S ILIAC CLSD
13037560	RELIN-O SCREW, 7.5X60MM 2S ILIAC CLOSED
13037570	RELIN-O SCREW, 7.5X70MM 2S ILIAC CLOSED



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13037580	RELIN-O SCREW, 7.5X80MM 2S ILIAC CLOSED
13037590	RELIN-O SCREW, 7.5X90MM 2S ILIAC CLOSED
13038510	RELIN-O SCREW, 8.5X100MM 2S ILIAC CLSD
13038560	RELIN-O SCREW, 8.5X60MM 2S ILIAC CLOSED
13038570	RELIN-O SCREW, 8.5X70MM 2S ILIAC CLOSED
13038580	RELIN-O SCREW, 8.5X80MM 2S ILIAC CLOSED
13038590	RELIN-O SCREW, 8.5X90MM 2S ILIAC CLOSED
13039510	RELIN-O SCREW, 9.5X100MM 2S ILIAC CLSD
13039560	RELIN-O SCREW, 9.5X60MM 2S ILIAC CLOSED
13039570	RELIN-O SCREW, 9.5X70MM 2S ILIAC CLOSED
13039580	RELIN-O SCREW, 9.5X80MM 2S ILIAC CLOSED
13039590	RELIN-O SCREW, 9.5X90MM 2S ILIAC CLOSED
13044025	RELIN-O SCREW, 4.0X25MM 2S UNIPLANAR
13044030	RELIN-O SCREW, 4.0X30MM 2S UNIPLANAR
13044035	RELIN-O SCREW, 4.0X35MM 2S UNIPLANAR
13044040	RELIN-O SCREW, 4.0X40MM 2S UNIPLANAR
13044520	RELIN-O SCREW, 4.5X20MM 2S UNIPLANAR
13044525	RELIN-O SCREW, 4.5X25MM 2S UNIPLANAR
13044530	RELIN-O SCREW, 4.5X30MM 2S UNIPLANAR
13044535	RELIN-O SCREW, 4.5X35MM 2S UNIPLANAR
13044540	RELIN-O SCREW, 4.5X40MM 2S UNIPLANAR
13044545	RELIN-O SCREW, 4.5X45MM 2S UNIPLANAR
13044550	RELIN-O SCREW, 4.5X50MM 2S UNIPLANAR
13045020	RELIN-O SCREW, 5.0X20MM 2S UNIPLANAR
13045025	RELIN-O SCREW, 5.0X25MM 2S UNIPLANAR
13045030	RELIN-O SCREW, 5.0X30MM 2S UNIPLANAR
13045035	RELIN-O SCREW, 5.0X35MM 2S UNIPLANAR
13045040	RELIN-O SCREW, 5.0X40MM 2S UNIPLANAR
13045045	RELIN-O SCREW, 5.0X45MM 2S UNIPLANAR
13045520	RELIN-O SCREW, 5.5X20MM 2S UNIPLANAR
13045525	RELIN-O SCREW, 5.5X25MM 2S UNIPLANAR
13045530	RELIN-O SCREW, 5.5X30MM 2S UNIPLANAR
13045535	RELIN-O SCREW, 5.5X35MM 2S UNIPLANAR
13045540	RELIN-O SCREW, 5.5X40MM 2S UNIPLANAR
13045545	RELIN-O SCREW, 5.5X45MM 2S UNIPLANAR
13045550	RELIN-O SCREW, 5.5X50MM 2S UNIPLANAR
13045555	RELIN-O SCREW, 5.5X55MM 2S UNIPLANAR
13046025	RELIN-O SCREW, 6.0X25MM 2S UNIPLANAR
13046030	RELIN-O SCREW, 6.0X30MM 2S UNIPLANAR



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13046035	RELIN-O SCREW, 6.0X35MM 2S UNIPLANAR
13046040	RELIN-O SCREW, 6.0X40MM 2S UNIPLANAR
13046045	RELIN-O SCREW, 6.0X45MM 2S UNIPLANAR
13046050	RELIN-O SCREW, 6.0X50MM 2S UNIPLANAR
13046525	RELIN-O SCREW, 6.5X25MM 2S UNIPLANAR
13046530	RELIN-O SCREW, 6.5X30MM 2S UNIPLANAR
13046535	RELIN-O SCREW, 6.5X35MM 2S UNIPLANAR
13046540	RELIN-O SCREW, 6.5X40MM 2S UNIPLANAR
13046545	RELIN-O SCREW, 6.5X45MM 2S UNIPLANAR
13046550	RELIN-O SCREW, 6.5X50MM 2S UNIPLANAR
13046555	RELIN-O SCREW, 6.5X55MM 2S UNIPLANAR
13047030	RELIN-O SCREW, 7.0X30MM 2S UNIPLANAR
13047035	RELIN-O SCREW, 7.0X35MM 2S UNIPLANAR
13047040	RELIN-O SCREW, 7.0X40MM 2S UNIPLANAR
13047045	RELIN-O SCREW, 7.0X45MM 2S UNIPLANAR
13047050	RELIN-O SCREW, 7.0X50MM 2S UNIPLANAR
13047530	RELIN-O SCREW, 7.5X30MM 2S UNIPLANAR
13047535	RELIN-O SCREW, 7.5X35MM 2S UNIPLANAR
13047540	RELIN-O SCREW, 7.5X40MM 2S UNIPLANAR
13047545	RELIN-O SCREW, 7.5X45MM 2S UNIPLANAR
13047550	RELIN-O SCREW, 7.5X50MM 2S UNIPLANAR
13047555	RELIN-O SCREW, 7.5X55MM 2S UNIPLANAR
13054525	RELIN-O SCREW, 4.5X25MM 2S RED UNI
13054530	RELIN-O SCREW, 4.5X30MM 2S RED UNI
13054535	RELIN-O SCREW, 4.5X35MM 2S RED UNI
13054540	RELIN-O SCREW, 4.5X40MM 2S RED UNI
13054545	RELIN-O SCREW, 4.5X45MM 2S RED UNI
13054550	RELIN-O SCREW, 4.5X50MM 2S RED UNI
13055025	RELIN-O SCREW, 5.0X25MM 2S RED UNI
13055030	RELIN-O SCREW, 5.0X30MM 2S RED UNI
13055035	RELIN-O SCREW, 5.0X35MM 2S RED UNI
13055040	RELIN-O SCREW, 5.0X40MM 2S RED UNI
13055045	RELIN-O SCREW, 5.0X45MM 2S RED UNI
13055525	RELIN-O SCREW, 5.5X25MM 2S RED UNI
13055530	RELIN-O SCREW, 5.5X30MM 2S RED UNI
13055535	RELIN-O SCREW, 5.5X35MM 2S REDUNI
13055540	RELIN-O SCREW, 5.5X40MM 2S REDUNI
13055545	RELIN-O SCREW, 5.5X45MM 2S REDUNI
13055550	RELIN-O SCREW, 5.5X50MM 2S REDUNI



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13055555	RELIN-O SCREW, 5.5X55MM 2S REDUNI
13056530	RELIN-O SCREW, 6.5X30MM 2S REDUNI
13056535	RELIN-O SCREW, 6.5X35MM 2S REDUNI
13056540	RELIN-O SCREW, 6.5X40MM 2S REDUNI
13056545	RELIN-O SCREW, 6.5X45MM 2S REDUNI
13056550	RELIN-O SCREW, 6.5X50MM 2S REDUNI
13056555	RELIN-O SCREW, 6.5X55MM 2S REDUNI
13064525	RELIN-O SCREW, 4.5X25MM 2S PLS
13064530	RELIN-O SCREW, 4.5X30MM 2S PLS
13064535	RELIN-O SCREW, 4.5X35MM 2S PLS
13064540	RELIN-O SCREW, 4.5X40MM 2S PLS
13064545	RELIN-O SCREW, 4.5X45MM 2S PLS
13065530	RELIN-O SCREW, 5.5X30MM 2S PLS
13065535	RELIN-O SCREW, 5.5X35MM 2S PLS
13065540	RELIN-O SCREW, 5.5X40MM 2S PLS
13065545	RELIN-O SCREW, 5.5X45MM 2S PLS
13065550	RELIN-O SCREW, 5.5X50MM 2S PLS
13066530	RELIN-O SCREW, 6.5X30MM 2S PLS
13066535	RELIN-O SCREW, 6.5X35MM 2S PLS
13066540	RELIN-O SCREW, 6.5X40MM 2S PLS
13066545	RELIN-O SCREW, 6.5X45MM 2S PLS
13066550	RELIN-O SCREW, 6.5X50MM 2S PLS
13067535	RELIN-O SCREW, 7.5X35MM 2S PLS
13067540	RELIN-O SCREW, 7.5X40MM 2S PLS
13067545	RELIN-O SCREW, 7.5X45MM 2S PLS
13067550	RELIN-O SCREW, 7.5X50MM 2S PLS
13067555	RELIN-O SCREW, 7.5X55MM 2S PLS
13084025	RELIN-O SCREW, 4.0X25MM 2S MONOAXIAL
13084030	RELIN-O SCREW, 4.0X30MM 2S MONOAXIAL
13084035	RELIN-O SCREW, 4.0X35MM 2S MONOAXIAL
13084040	RELIN-O SCREW, 4.0X40MM 2S MONOAXIAL
13084045	RELIN-O SCREW, 4.0X45MM 2S MONOAXIAL
13084530	RELIN-O SCREW, 4.5X30MM 2S MONOAXIAL
13084535	RELIN-O SCREW, 4.5X35MM 2S MONOAXIAL
13084540	RELIN-O SCREW, 4.5X40MM 2S MONOAXIAL
13084545	RELIN-O SCREW, 4.5X45MM 2S MONOAXIAL
13085530	RELIN-O SCREW, 5.5X30MM 2S MONOAXIAL
13085535	RELIN-O SCREW, 5.5X35MM 2S MONOAXIAL
13085540	RELIN-O SCREW, 5.5X40MM 2S MONOAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13085545	RELIN-O SCREW, 5.5X45MM 2S MONOAXIAL
13085550	RELIN-O SCREW, 5.5X50MM 2S MONOAXIAL
13085555	RELIN-O SCREW, 5.5X55MM 2S MONOAXIAL
13085560	RELIN-O SCREW, 5.5X60MM 2S MONOAXIAL
13086530	RELIN-O SCREW, 6.5X30MM 2S MONOAXIAL
13086535	RELIN-O SCREW, 6.5X35MM 2S MONOAXIAL
13086540	RELIN-O SCREW, 6.5X40MM 2S MONOAXIAL
13086545	RELIN-O SCREW, 6.5X45MM 2S MONOAXIAL
13086550	RELIN-O SCREW, 6.5X50MM 2S MONOAXIAL
13086555	RELIN-O SCREW, 6.5X55MM 2S MONOAXIAL
13086560	RELIN-O SCREW, 6.5X60MM 2S MONOAXIAL
13087535	RELIN-O SCREW, 7.5X35MM 2S MONOAXIAL
13087540	RELIN-O SCREW, 7.5X40MM 2S MONOAXIAL
13087545	RELIN-O SCREW, 7.5X45MM 2S MONOAXIAL
13087550	RELIN-O SCREW, 7.5X50MM 2S MONOAXIAL
13087555	RELIN-O SCREW, 7.5X55MM 2S MONOAXIAL
13087560	RELIN-O SCREW, 7.5X60MM 2S MONOAXIAL
13088535	RELIN-O SCREW, 8.5X35MM 2S MONOAXIAL
13088540	RELIN-O SCREW, 8.5X40MM 2S MONOAXIAL
13088545	RELIN-O SCREW, 8.5X45MM 2S MONOAXIAL
13088550	RELIN-O SCREW, 8.5X50MM 2S MONOAXIAL
13088555	RELIN-O SCREW, 8.5X55MM 2S MONOAXIAL
13100000	RSS LOCK SCREW, CLOSED TULIP
13184025	RELIN-O MOD SCREW SHANK, 4.0X25MM 2S
13184045	RELIN-O MOD SCREW SHANK, 4.0X45MM 2S
13184540	RELIN-O MOD SCREW SHANK, 4.5X40MM 2S
13184545	RELIN-O MOD SCREW SHANK, 4.5X45MM 2S
13185535	RELIN-O MOD SCREW SHANK, 5.5X35MM 2S
13185540	RELIN-O MOD SCREW SHANK, 5.5X40MM 2S
13185545	RELIN-O MOD SCREW SHANK, 5.5X45MM 2S
13185550	RELIN-O MOD SCREW SHANK, 5.5X50MM 2S
13185555	RELIN-O MOD SCREW SHANK, 5.5X55MM 2S
13186535	RELIN-O MOD SCREW SHANK, 6.5X35MM 2S
13186540	RELIN-O MOD SCREW SHANK, 6.5X40MM 2S
13186545	RELIN-O MOD SCREW SHANK, 6.5X45MM 2S
13186550	RELIN-O MOD SCREW SHANK, 6.5X50MM 2S
13186555	RELIN-O MOD SCREW SHANK, 6.5X55MM 2S
13186560	RELIN-O MOD SCREW SHANK, 6.5X60MM 2S
13187535	RELIN-O MOD SCREW SHANK, 7.5X35MM 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13187540	RELIN-O MOD SCREW SHANK, 7.5X40MM 2S
13187545	RELIN-O MOD SCREW SHANK, 7.5X45MM 2S
13187550	RELIN-O MOD SCREW SHANK, 7.5X50MM 2S
13187555	RELIN-O MOD SCREW SHANK, 7.5X55MM 2S
13187560	RELIN-O MOD SCREW SHANK, 7.5X60MM 2S
13188535	RELIN-O MOD SCREW SHANK, 8.5X35MM 2S
13188540	RELIN-O MOD SCREW SHANK, 8.5X40MM 2S
13188545	RELIN-O MOD SCREW SHANK, 8.5X45MM 2S
13188550	RELIN-O MOD SCREW SHANK, 8.5X50MM 2S
13188555	RELIN-O MOD SCREW SHANK, 8.5X55MM 2S
13215535	RELIN-O SCREW, 5.5X35MM DBL FEN POLY
13215540	RELIN-O SCREW, 5.5X40MM DBL FEN POLY
13215545	RELIN-O SCREW, 5.5X45MM DBL FEN POLY
13215550	RELIN-O SCREW, 5.5X50MM DBL FEN POLY
13215555	RELIN-O SCREW, 5.5X55MM DBL FEN POLY
13215560	RELIN-O SCREW, 5.5X60MM DBL FEN POLY
13216535	RELIN-O SCREW, 6.5X35MM DBL FEN POLY
13216540	RELIN-O SCREW, 6.5X40MM DBL FEN POLY
13216545	RELIN-O SCREW, 6.5X45MM DBL FEN POLY
13216550	RELIN-O SCREW, 6.5X50MM DBL FEN POLY
13216555	RELIN-O SCREW, 6.5X55MM DBL FEN POLY
13216560	RELIN-O SCREW, 6.5X60MM DBL FEN POLY
13217511	RELIN-O SCREW, 7.5X110MM DBL FEN POLY
13217535	RELIN-O SCREW, 7.5X35MM DBL FEN POLY
13217540	RELIN-O SCREW, 7.5X40MM DBL FEN POLY
13217545	RELIN-O SCREW, 7.5X45MM DBL FEN POLY
13217550	RELIN-O SCREW, 7.5X50MM DBL FEN POLY
13217555	RELIN-O SCREW, 7.5X55MM DBL FEN POLY
13217560	RELIN-O SCREW, 7.5X60MM DBL FEN POLY
13218535	RELIN-O SCREW, 8.5X35MM DBL FEN POLY
13218540	RELIN-O SCREW, 8.5X40MM DBL FEN POLY
13218545	RELIN-O SCREW, 8.5X45MM DBL FEN POLY
13218550	RELIN-O SCREW, 8.5X50MM DBL FEN POLY
13218555	RELIN-O SCREW, 8.5X55MM DBL FEN POLY
13218560	RELIN-O SCREW, 8.5X60MM DBL FEN POLY
13218570	RELIN-O SCREW, 8.5X70MM DBL FEN POLY
13218580	RELIN-O SCREW, 8.5X80MM DBL FEN POLY
13218590	RELIN-O SCREW, 8.5X90MM DBL FEN POLY
13219090	RELIN-O SCREW, 9.0X90MM DBL FEN POLY



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13219570	RELIN-O SCREW, 9.5X70MM DBL FEN POLY
13219580	RELIN-O SCREW, 9.5X80MM DBL FEN POLY
13219590	RELIN-O SCREW, 9.5X90MM DBL FEN POLY
13450000	RSS LOCK SCREW, 4.5/5.0MM
13550000	RELIN LOCK SCREW, 5.5MM OPEN TULIP
13600000	RELIN LOCK SCREW, 6.0MM OPENTULIP
14155035	GSB TI ROD, 5.5X35MM LORDOTIC
14155040	GSB TI ROD, 5.5X40MM LORDOTIC
14155045	GSB TI ROD, 5.5X45MM LORDOTIC
14155050	GSB TI ROD, 5.5X50MM LORDOTIC
14155055	GSB TI ROD, 5.5X55MM LORDOTIC
14155060	GSB TI ROD, 5.5X60MM LORDOTIC
14155065	GSB TI ROD, 5.5X65MM LORDOTIC
14155070	GSB TI ROD, 5.5X70MM LORDOTIC
14155075	GSB TI ROD, 5.5X75MM LORDOTIC
14155080	GSB TI ROD, 5.5X80MM LORDOTIC
14155090	GSB TI ROD, 5.5X90MM LORDOTIC
14155100	GSB TI ROD, 5.5X100MM LORDOTIC
14155110	GSB TI ROD, 5.5X110MM LORDOTIC
14155120	GSB TI ROD, 5.5X120MM LORDOTIC
14210530	RELIN SCREW, 10.5X30MM 2FS POLYAXIAL
14210535	RELIN SCREW, 10.5X35MM 2FS POLYAXIAL
14210540	RELIN SCREW, 10.5X40MM 2FS POLYAXIAL
14210545	RELIN SCREW, 10.5X45MM 2FS POLYAXIAL
14210550	RELIN SCREW, 10.5X50MM 2FS POLYAXIAL
14210555	RELIN SCREW, 10.5X55MM 2FS POLYAXIAL
14210560	RELIN SCREW, 10.5X60MM 2FS POLYAXIAL
14211030	RELIN SCREW, 10.0X30MM 2FS POLYAXIAL
14211035	RELIN SCREW, 10.0X35MM 2FS POLYAXIAL
14211040	RELIN SCREW, 10.0X40MM 2FS POLYAXIAL
14211045	RELIN SCREW, 10.0X45MM 2FS POLYAXIAL
14211050	RELIN SCREW, 10.0X50MM 2FS POLYAXIAL
14211055	RELIN SCREW, 10.0X55MM 2FS POLYAXIAL
14211060	RELIN SCREW, 10.0X60MM 2FS POLYAXIAL
14214530	RELIN SCREW, 4.5X30MM 2FS POLYAXIAL
14214535	RELIN SCREW, 4.5X35MM 2FS POLYAXIAL
14214540	RELIN SCREW, 4.5X40MM 2FS POLYAXIAL
14214545	RELIN SCREW, 4.5X45MM 2FS POLYAXIAL
14214550	RELIN SCREW, 4.5X50MM 2FS POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14214555	RELINER SCREW, 4.5X55MM 2FS POLYAXIAL
14214560	RELINER SCREW, 4.5X60MM 2FS POLYAXIAL
14215030	RELINER SCREW, 5.0X30MM 2FS POLYAXIAL
14215035	RELINER SCREW, 5.0X35MM 2FS POLYAXIAL
14215040	RELINER SCREW, 5.0X40MM 2FS POLYAXIAL
14215045	RELINER SCREW, 5.0X45MM 2FS POLYAXIAL
14215050	RELINER SCREW, 5.0X50MM 2FS POLYAXIAL
14215055	RELINER SCREW, 5.0X55MM 2FS POLYAXIAL
14215060	RELINER SCREW, 5.0X60MM 2FS POLYAXIAL
14215530	RELINER SCREW, 5.5X30MM 2FS POLYAXIAL
14215535	RELINER SCREW, 5.5X35MM 2FS POLYAXIAL
14215540	RELINER SCREW, 5.5X40MM 2FS POLYAXIAL
14215545	RELINER SCREW, 5.5X45MM 2FS POLYAXIAL
14215550	RELINER SCREW, 5.5X50MM 2FS POLYAXIAL
14215555	RELINER SCREW, 5.5X55MM 2FS POLYAXIAL
14215560	RELINER SCREW, 5.5X60MM 2FS POLYAXIAL
14216030	RELINER SCREW, 6.0X30MM 2FS POLYAXIAL
14216035	RELINER SCREW, 6.0X35MM 2FS POLYAXIAL
14216040	RELINER SCREW, 6.0X40MM 2FS POLYAXIAL
14216045	RELINER SCREW, 6.0X45MM 2FS POLYAXIAL
14216050	RELINER SCREW, 6.0X50MM 2FS POLYAXIAL
14216055	RELINER SCREW, 6.0X55MM 2FS POLYAXIAL
14216060	RELINER SCREW, 6.0X60MM 2FS POLYAXIAL
14216530	RELINER SCREW, 6.5X30MM 2FS POLYAXIAL
14216535	RELINER SCREW, 6.5X35MM 2FS POLYAXIAL
14216540	RELINER SCREW, 6.5X40MM 2FS POLYAXIAL
14216545	RELINER SCREW, 6.5X45MM 2FS POLYAXIAL
14216550	RELINER SCREW, 6.5X50MM 2FS POLYAXIAL
14216555	RELINER SCREW, 6.5X55MM 2FS POLYAXIAL
14216560	RELINER SCREW, 6.5X60MM 2FS POLYAXIAL
14217030	RELINER SCREW, 7.0X30MM 2FS POLYAXIAL
14217035	RELINER SCREW, 7.0X35MM 2FS POLYAXIAL
14217040	RELINER SCREW, 7.0X40MM 2FS POLYAXIAL
14217045	RELINER SCREW, 7.0X45MM 2FS POLYAXIAL
14217050	RELINER SCREW, 7.0X50MM 2FS POLYAXIAL
14217055	RELINER SCREW, 7.0X55MM 2FS POLYAXIAL
14217060	RELINER SCREW, 7.0X60MM 2FS POLYAXIAL
14217530	RELINER SCREW, 7.5X30MM 2FS POLYAXIAL
14217535	RELINER SCREW, 7.5X35MM 2FS POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14217540	RELINER SCREW, 7.5X40MM 2FS POLYAXIAL
14217545	RELINER SCREW, 7.5X45MM 2FS POLYAXIAL
14217550	RELINER SCREW, 7.5X50MM 2FS POLYAXIAL
14217555	RELINER SCREW, 7.5X55MM 2FS POLYAXIAL
14217560	RELINER SCREW, 7.5X60MM 2FS POLYAXIAL
14218030	RELINER SCREW, 8.0X30MM 2FS POLYAXIAL
14218035	RELINER SCREW, 8.0X35MM 2FS POLYAXIAL
14218040	RELINER SCREW, 8.0X40MM 2FS POLYAXIAL
14218045	RELINER SCREW, 8.0X45MM 2FS POLYAXIAL
14218050	RELINER SCREW, 8.0X50MM 2FS POLYAXIAL
14218055	RELINER SCREW, 8.0X55MM 2FS POLYAXIAL
14218060	RELINER SCREW, 8.0X60MM 2FS POLYAXIAL
14218530	RELINER SCREW, 8.5X30MM 2FS POLYAXIAL
14218535	RELINER SCREW, 8.5X35MM 2FS POLYAXIAL
14218540	RELINER SCREW, 8.5X40MM 2FS POLYAXIAL
14218545	RELINER SCREW, 8.5X45MM 2FS POLYAXIAL
14218550	RELINER SCREW, 8.5X50MM 2FS POLYAXIAL
14218555	RELINER SCREW, 8.5X55MM 2FS POLYAXIAL
14218560	RELINER SCREW, 8.5X60MM 2FS POLYAXIAL
14219030	RELINER SCREW, 9.0X30MM 2FS POLYAXIAL
14219035	RELINER SCREW, 9.0X35MM 2FS POLYAXIAL
14219040	RELINER SCREW, 9.0X40MM 2FS POLYAXIAL
14219045	RELINER SCREW, 9.0X45MM 2FS POLYAXIAL
14219050	RELINER SCREW, 9.0X50MM 2FS POLYAXIAL
14219055	RELINER SCREW, 9.0X55MM 2FS POLYAXIAL
14219060	RELINER SCREW, 9.0X60MM 2FS POLYAXIAL
14219530	RELINER SCREW, 9.5X30MM 2FS POLYAXIAL
14219535	RELINER SCREW, 9.5X35MM 2FS POLYAXIAL
14219540	RELINER SCREW, 9.5X40MM 2FS POLYAXIAL
14219545	RELINER SCREW, 9.5X45MM 2FS POLYAXIAL
14219550	RELINER SCREW, 9.5X50MM 2FS POLYAXIAL
14219555	RELINER SCREW, 9.5X55MM 2FS POLYAXIAL
14219560	RELINER SCREW, 9.5X60MM 2FS POLYAXIAL
14255080	GSB COCR ROD, 5.5X80MM LORDOTIC
14355025	RELINER-O TI ROD, 5.5X25MM LORDOTIC
14355030	RELINER-O TI ROD, 5.5X30MM LORDOTIC
14355035	RELINER-O TI ROD, 5.5X35MM LORDOTIC
14355040	RELINER-O TI ROD, 5.5X40MM LORDOTIC
14355045	RELINER-O TI ROD, 5.5X45MM LORDOTIC



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14355050	RELIN-O TI ROD, 5.5X50MM LORDOTIC
14355055	RELIN-O TI ROD, 5.5X55MM LORDOTIC
14355060	RELIN-O TI ROD, 5.5X60MM LORDOTIC
14355065	RELIN-O TI ROD, 5.5X65MM LORDOTIC
14355070	RELIN-O TI ROD, 5.5X70MM LORDOTIC
14355075	RELIN-O TI ROD, 5.5X75MM LORDOTIC
14355080	RELIN-O TI ROD, 5.5X80MM LORDOTIC
14355085	RELIN-O TI ROD, 5.5X85MM LORDOTIC
14355090	RELIN-O TI ROD, 5.5X90MM LORDOTIC
14355095	RELIN-O TI ROD, 5.5X95MM LORDOTIC
14355100	RELIN-O TI ROD, 5.5X100MM LORDOTIC
14355110	RELIN-O TI ROD, 5.5X110MM LORDOTIC
14355120	RELIN-O TI ROD, 5.5X120MM LORDOTIC
14355140	RELIN-O TI ROD, 5.5X140MM LORDOTIC
14355150	RELIN-O TI ROD, 5.5X150MM LORDOTIC
14355160	RELIN-O TI ROD, 5.5X160MM LORDOTIC
14355180	RELIN-O TI ROD, 5.5X180MM LORDOTIC
14355200	RELIN-O TI ROD, 5.5X200MM LORDOTIC
14360030	RELIN-O TI ROD, 6.0X30MM LORDOTIC
14360035	RELIN-O TI ROD, 6.0X35MM LORDOTIC
14360040	RELIN-O TI ROD, 6.0X40MM LORDOTIC
14360045	RELIN-O TI ROD, 6.0X45MM LORDOTIC
14360050	RELIN-O TI ROD, 6.0X50MM LORDOTIC
14360055	RELIN-O TI ROD, 6.0X55MM LORDOTIC
14360060	RELIN-O TI ROD, 6.0X60MM LORDOTIC
14360065	RELIN-O TI ROD, 6.0X65MM LORDOTIC
14360070	RELIN-O TI ROD, 6.0X70MM LORDOTIC
14360075	RELIN-O TI ROD, 6.0X75MM LORDOTIC
14360080	RELIN-O TI ROD, 6.0X80MM LORDOTIC
14360085	RELIN-O TI ROD, 6.0X85MM LORDOTIC
14360090	RELIN-O TI ROD, 6.0X90MM LORDOTIC
14360095	RELIN-O TI ROD, 6.0X95MM LORDOTIC
14360100	RELIN-O TI ROD, 6.0X100MM LORDOTIC
14360110	RELIN-O TI ROD, 6.0X110MM LORDOTIC
14360120	RELIN-O TI ROD, 6.0X120MM LORDOTIC
14450025	RELIN-O COCR ROD, 5.0X25MM LORDOTIC
14450030	RELIN-O COCR ROD, 5.0X30MM LORDOTIC
14450035	RELIN-O COCR ROD, 5.0X35MM LORDOTIC
14450040	RELIN-O COCR ROD, 5.0X40MM LORDOTIC



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14450045	RELINE-O COCR ROD, 5.0X45MM LORDOTIC
14450050	RELINE-O COCR ROD, 5.0X50MM LORDOTIC
14450055	RELINE-O COCR ROD, 5.0X55MM LORDOTIC
14450060	RELINE-O COCR ROD, 5.0X60MM LORDOTIC
14450065	RELINE-O COCR ROD, 5.0X65MM LORDOTIC
14450070	RELINE-O COCR ROD, 5.0X70MM LORDOTIC
14450075	RELINE-O COCR ROD, 5.0X75MM LORDOTIC
14450080	RELINE-O COCR ROD, 5.0X80MM LORDOTIC
14450085	RELINE-O COCR ROD, 5.0X85MM LORDOTIC
14450090	RELINE-O COCR ROD, 5.0X90MM LORDOTIC
14450095	RELINE-O COCR ROD, 5.0X95MM LORDOTIC
14450100	RELINE-O COCR ROD, 5.0X100MM LORDOTIC
14450110	RELINE-O COCR ROD, 5.0X110MM LORDOTIC
14450120	RELINE-O COCR ROD, 5.0X120MM LORDOTIC
14455030	RELINE-O COCR ROD, 5.5X30MM LORDOTIC
14455035	RELINE-O COCR ROD, 5.5X35MM LORDOTIC
14455040	RELINE-O COCR ROD, 5.5X40MM LORDOTIC
14455045	RELINE-O COCR ROD, 5.5X45MM LORDOTIC
14455050	RELINE-O COCR ROD, 5.5X50MM LORDOTIC
14455055	RELINE-O COCR ROD, 5.5X55MM LORDOTIC
14455060	RELINE-O COCR ROD, 5.5X60MM LORDOTIC
14455065	RELINE-O COCR ROD, 5.5X65MM LORDOTIC
14455070	RELINE-O COCR ROD, 5.5X70MM LORDOTIC
14455075	RELINE-O COCR ROD, 5.5X75MM LORDOTIC
14455080	RELINE-O COCR ROD, 5.5X80MM LORDOTIC
14455085	RELINE-O COCR ROD, 5.5X85MM LORDOTIC
14455090	RELINE-O COCR ROD, 5.5X90MM LORDOTIC
14455095	RELINE-O COCR ROD, 5.5X95MM LORDOTIC
14455100	RELINE-O COCR ROD, 5.5X100MM LORDOTIC
14455110	RELINE-O COCR ROD, 5.5X110MM LORDOTIC
14455120	RELINE-O COCR ROD, 5.5X120MM LORDOTIC
14455200	RELINE-O COCR ROD, 5.5X200MM LORDOTIC
14460030	RELINE-O COCR ROD, 6.0X30MM LORDOTIC
14460035	RELINE-O COCR ROD, 6.0X35MM LORDOTIC
14460040	RELINE-O COCR ROD, 6.0X40MM LORDOTIC
14460045	RELINE-O COCR ROD, 6.0X45MM LORDOTIC
14460050	RELINE-O COCR ROD, 6.0X50MM LORDOTIC
14460055	RELINE-O COCR ROD, 6.0X55MM LORDOTIC
14460060	RELINE-O COCR ROD, 6.0X60MM LORDOTIC



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
14460065	RELIN-O COCR ROD, 6.0X65MM LORDOTIC
14460070	RELIN-O COCR ROD, 6.0X70MM LORDOTIC
14460075	RELIN-O COCR ROD, 6.0X75MM LORDOTIC
14460080	RELIN-O COCR ROD, 6.0X80MM LORDOTIC
14460085	RELIN-O COCR ROD, 6.0X85MM LORDOTIC
14460090	RELIN-O COCR ROD, 6.0X90MM LORDOTIC
14460095	RELIN-O COCR ROD, 6.0X95MM LORDOTIC
14460100	RELIN-O COCR ROD, 6.0X100MM LORDOTIC
14460110	RELIN-O COCR ROD, 6.0X110MM LORDOTIC
14460120	RELIN-O COCR ROD, 6.0X120MM LORDOTIC
14460200	RELIN-O COCR ROD, 6.0X200MM LORDOTIC
14514017	RSS SCREW, 4.0X17.5MM POLYAXIAL 2S
14514020	RSS SCREW, 4.0X20MM POLYAXIAL 2S
14514022	RSS SCREW, 4.0X22.5MM POLYAXIAL 2S
14514025	RSS SCREW, 4.0X25MM POLYAXIAL 2S
14514027	RSS SCREW, 4.0X27.5MM POLYAXIAL 2S
14514030	RSS SCREW, 4.0X30MM POLYAXIAL 2S
14514035	RSS SCREW, 4.0X35MM POLYAXIAL 2S
14514040	RSS SCREW, 4.0X40MM POLYAXIAL 2S
14514517	RSS SCREW, 4.5X17.5MM POLYAXIAL 2S
14514520	RSS SCREW, 4.5X20MM POLYAXIAL 2S
14514522	RSS SCREW, 4.5X22.5MM POLYAXIAL 2S
14514525	RSS SCREW, 4.5X25MM POLYAXIAL 2S
14514527	RSS SCREW, 4.5X27.5MM POLYAXIAL 2S
14514530	RSS SCREW, 4.5X30MM POLYAXIAL 2S
14514535	RSS SCREW, 4.5X35MM POLYAXIAL 2S
14514540	RSS SCREW, 4.5X40MM POLYAXIAL 2S
14514545	RSS SCREW, 4.5X45MM POLYAXIAL 2S
14515017	RSS SCREW, 5.0X17.5MM POLYAXIAL 2S
14515020	RSS SCREW, 5.0X20MM POLYAXIAL 2S
14515022	RSS SCREW, 5.0X22.5MM POLYAXIAL 2S
14515025	RSS SCREW, 5.0X25MM POLYAXIAL 2S
14515027	RSS SCREW, 5.0X27.5MM POLYAXIAL 2S
14515030	RSS SCREW, 5.0X30MM POLYAXIAL 2S
14515035	RSS SCREW, 5.0X35MM POLYAXIAL 2S
14515040	RSS SCREW, 5.0X40MM POLYAXIAL 2S
14515045	RSS SCREW, 5.0X45MM POLYAXIAL 2S
14515517	RSS SCREW, 5.5X17.5MM POLYAXIAL 2S
14515520	RSS SCREW, 5.5X20MM POLYAXIAL 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14515522	RSS SCREW, 5.5X22.5MM POLYAXIAL 2S
14515525	RSS SCREW, 5.5X25MM POLYAXIAL 2S
14515527	RSS SCREW, 5.5X27.5MM POLYAXIAL 2S
14515530	RSS SCREW, 5.5X30MM POLYAXIAL 2S
14515535	RSS SCREW, 5.5X35MM POLYAXIAL 2S
14515540	RSS SCREW, 5.5X40MM POLYAXIAL 2S
14515545	RSS SCREW, 5.5X45MM POLYAXIAL 2S
14515550	RSS SCREW, 5.5X50MM POLYAXIAL 2S
14515555	RSS SCREW, 5.5X55MM POLYAXIAL 2S
14515560	RSS SCREW, 5.5X60MM POLYAXIAL 2S
14515565	RSS SCREW, 5.5X65MM POLYAXIAL 2S
14515570	RSS SCREW, 5.5X70MM POLYAXIAL 2S
14515575	RSS SCREW, 5.5X75MM POLYAXIAL 2S
14515580	RSS SCREW, 5.5X80MM POLYAXIAL 2S
14516020	RSS SCREW, 6.0X20MM POLYAXIAL 2S
14516025	RSS SCREW, 6.0X25MM POLYAXIAL 2S
14516030	RSS SCREW, 6.0X30MM POLYAXIAL 2S
14516035	RSS SCREW, 6.0X35MM POLYAXIAL 2S
14516040	RSS SCREW, 6.0X40MM POLYAXIAL 2S
14516045	RSS SCREW, 6.0X45MM POLYAXIAL 2S
14516050	RSS SCREW, 6.0X50MM POLYAXIAL 2S
14516055	RSS SCREW, 6.0X55MM POLYAXIAL 2S
14516060	RSS SCREW, 6.0X60MM POLYAXIAL 2S
14516520	RSS SCREW, 6.5X20MM POLYAXIAL 2S
14516525	RSS SCREW, 6.5X25MM POLYAXIAL 2S
14516530	RSS SCREW, 6.5X30MM POLYAXIAL 2S
14516535	RSS SCREW, 6.5X35MM POLYAXIAL 2S
14516540	RSS SCREW, 6.5X40MM POLYAXIAL 2S
14516545	RSS SCREW, 6.5X45MM POLYAXIAL 2S
14516550	RSS SCREW, 6.5X50MM POLYAXIAL 2S
14516555	RSS SCREW, 6.5X55MM POLYAXIAL 2S
14516560	RSS SCREW, 6.5X60MM POLYAXIAL 2S
14516565	RSS SCREW, 6.5X65MM POLYAXIAL 2S
14516570	RSS SCREW, 6.5X70MM POLYAXIAL 2S
14516575	RSS SCREW, 6.5X75MM POLYAXIAL 2S
14516580	RSS SCREW, 6.5X80MM POLYAXIAL 2S
14516590	RSS SCREW, 6.5X90MM POLYAXIAL 2S
14517530	RSS SCREW, 7.5X30MM POLYAXIAL 2S
14517535	RSS SCREW, 7.5X35MM POLYAXIAL 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14517540	RSS SCREW, 7.5X40MM POLYAXIAL 2S
14517545	RSS SCREW, 7.5X45MM POLYAXIAL 2S
14517550	RSS SCREW, 7.5X50MM POLYAXIAL 2S
14517555	RSS SCREW, 7.5X55MM POLYAXIAL 2S
14517560	RSS SCREW, 7.5X60MM POLYAXIAL 2S
14517565	RSS SCREW, 7.5X65MM POLYAXIAL 2S
14517570	RSS SCREW, 7.5X70MM POLYAXIAL 2S
14517575	RSS SCREW, 7.5X75MM POLYAXIAL 2S
14517580	RSS SCREW, 7.5X80MM POLYAXIAL 2S
14517590	RSS SCREW, 7.5X90MM POLYAXIAL 2S
14518560	RSS SCREW, 8.5X60MM POLYAXIAL 2S
14518565	RSS SCREW, 8.5X65MM POLYAXIAL 2S
14518570	RSS SCREW, 8.5X70MM POLYAXIAL 2S
14518580	RSS SCREW, 8.5X80MM POLYAXIAL 2S
14518590	RSS SCREW, 8.5X90MM POLYAXIAL 2S
14524020	RSS SCREW, 4.0X20MM REDUCTION 2S
14524025	RSS SCREW, 4.0X25MM REDUCTION 2S
14524030	RSS SCREW, 4.0X30MM REDUCTION 2S
14524035	RSS SCREW, 4.0X35MM REDUCTION 2S
14524520	RSS SCREW, 4.5X20MM REDUCTION 2S
14524525	RSS SCREW, 4.5X25MM REDUCTION 2S
14524530	RSS SCREW, 4.5X30MM REDUCTION 2S
14524535	RSS SCREW, 4.5X35MM REDUCTION 2S
14524540	RSS SCREW, 4.5X40MM REDUCTION 2S
14525020	RSS SCREW, 5.0X20MM REDUCTION 2S
14525025	RSS SCREW, 5.0X25MM REDUCTION 2S
14525027	RSS SCREW, 5.0X27.5MM REDUCTION 2S
14525030	RSS SCREW, 5.0X30MM REDUCTION 2S
14525035	RSS SCREW, 5.0X35MM REDUCTION 2S
14525040	RSS SCREW, 5.0X40MM REDUCTION 2S
14525045	RSS SCREW, 5.0X45MM REDUCTION 2S
14525520	RSS SCREW, 5.5X20MM REDUCTION 2S
14525525	RSS SCREW, 5.5X25MM REDUCTION 2S
14525530	RSS SCREW, 5.5X30MM REDUCTION 2S
14525535	RSS SCREW, 5.5X35MM REDUCTION 2S
14525540	RSS SCREW, 5.5X40MM REDUCTION 2S
14525545	RSS SCREW, 5.5X45MM REDUCTION 2S
14526020	RSS SCREW, 6.0X20MM REDUCTION 2S
14526025	RSS SCREW, 6.0X25MM REDUCTION 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14526027	RSS SCREW, 6.0X27.5MM REDUCTION 2S
14526030	RSS SCREW, 6.0X30MM REDUCTION 2S
14526035	RSS SCREW, 6.0X35MM REDUCTION 2S
14526040	RSS SCREW, 6.0X40MM REDUCTION 2S
14526045	RSS SCREW, 6.0X45MM REDUCTION 2S
14526525	RSS SCREW, 6.5X25MM REDUCTION 2S
14526530	RSS SCREW, 6.5X30MM REDUCTION 2S
14526535	RSS SCREW, 6.5X35MM REDUCTION 2S
14526540	RSS SCREW, 6.5X40MM REDUCTION 2S
14526545	RSS SCREW, 6.5X45MM REDUCTION 2S
14544025	RSS SCREW, 4.0X25MM UNIPLANAR 2S
14544030	RSS SCREW, 4.0X30MM UNIPLANAR 2S
14544035	RSS SCREW, 4.0X35MM UNIPLANAR 2S
14544520	RSS SCREW, 4.5X20MM UNIPLANAR 2S
14544525	RSS SCREW, 4.5X25MM UNIPLANAR 2S
14544530	RSS SCREW, 4.5X30MM UNIPLANAR 2S
14544535	RSS SCREW, 4.5X35MM UNIPLANAR 2S
14544540	RSS SCREW, 4.5X40MM UNIPLANAR 2S
14545020	RSS SCREW, 5.0X20MM UNIPLANAR 2S
14545025	RSS SCREW, 5.0X25MM UNIPLANAR 2S
14545030	RSS SCREW, 5.0X30MM UNIPLANAR 2S
14545035	RSS SCREW, 5.0X35MM UNIPLANAR 2S
14545040	RSS SCREW, 5.0X40MM UNIPLANAR 2S
14545045	RSS SCREW, 5.0X45MM UNIPLANAR 2S
14545520	RSS SCREW, 5.5X20MM UNIPLANAR 2S
14545525	RSS SCREW, 5.5X25MM UNIPLANAR 2S
14545530	RSS SCREW, 5.5X30MM UNIPLANAR 2S
14545535	RSS SCREW, 5.5X35MM UNIPLANAR 2S
14545540	RSS SCREW, 5.5X40MM UNIPLANAR 2S
14545545	RSS SCREW, 5.5X45MM UNIPLANAR 2S
14546025	RSS SCREW, 6.0X25MM UNIPLANAR 2S
14546030	RSS SCREW, 6.0X30MM UNIPLANAR 2S
14546035	RSS SCREW, 6.0X35MM UNIPLANAR 2S
14546040	RSS SCREW, 6.0X40MM UNIPLANAR 2S
14546045	RSS SCREW, 6.0X45MM UNIPLANAR 2S
14546525	RSS SCREW, 6.5X25MM UNIPLANAR 2S
14546530	RSS SCREW, 6.5X30MM UNIPLANAR 2S
14546535	RSS SCREW, 6.5X35MM UNIPLANAR 2S
14546540	RSS SCREW, 6.5X40MM UNIPLANAR 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14546545	RSS SCREW, 6.5X45MM UNIPLANAR 2S
14554010	RSS SCREW, 4.0X100MM UNI-REDUCTION 2S
14554011	RSS SCREW, 4.0X110MM UNI-REDUCTION 2S
14554012	RSS SCREW, 4.0X120MM UNI-REDUCTION 2S
14554017	RSS SCREW, 4.0X17.5MM UNI-REDUCTION 2S
14554020	RSS SCREW, 4.0X20MM UNI-REDUCTION 2S
14554022	RSS SCREW, 4.0X22.5MM UNI-REDUCTION 2S
14554025	RSS SCREW, 4.0X25MM UNI-REDUCTION 2S
14554027	RSS SCREW, 4.0X27.5MM UNI-REDUCTION 2S
14554030	RSS SCREW, 4.0X30MM UNI-REDUCTION 2S
14554035	RSS SCREW, 4.0X35MM UNI-REDUCTION 2S
14554040	RSS SCREW, 4.0X40MM UNI-REDUCTION 2S
14554045	RSS SCREW, 4.0X45MM UNI-REDUCTION 2S
14554050	RSS SCREW, 4.0X50MM UNI-REDUCTION 2S
14554055	RSS SCREW, 4.0X55MM UNI-REDUCTION 2S
14554060	RSS SCREW, 4.0X60MM UNI-REDUCTION 2S
14554065	RSS SCREW, 4.0X65MM UNI-REDUCTION 2S
14554070	RSS SCREW, 4.0X70MM UNI-REDUCTION 2S
14554075	RSS SCREW, 4.0X75MM UNI-REDUCTION 2S
14554080	RSS SCREW, 4.0X80MM UNI-REDUCTION 2S
14554085	RSS SCREW, 4.0X85MM UNI-REDUCTION 2S
14554090	RSS SCREW, 4.0X90MM UNI-REDUCTION 2S
14554095	RSS SCREW, 4.0X95MM UNI-REDUCTION 2S
14554510	RSS SCREW, 4.5X100MM UNI-REDUCTION 2S
14554511	RSS SCREW, 4.5X110MM UNI-REDUCTION 2S
14554512	RSS SCREW, 4.5X120MM UNI-REDUCTION 2S
14554517	RSS SCREW, 4.5X17.5MM UNI-REDUCTION 2S
14554520	RSS SCREW, 4.5X20MM UNI-REDUCTION 2S
14554522	RSS SCREW, 4.5X22.5MM UNI-REDUCTION 2S
14554525	RSS SCREW, 4.5X25MM UNI-REDUCTION 2S
14554527	RSS SCREW, 4.5X27.5MM UNI-REDUCTION 2S
14554530	RSS SCREW, 4.5X30MM UNI-REDUCTION 2S
14554535	RSS SCREW, 4.5X35MM UNI-REDUCTION 2S
14554540	RSS SCREW, 4.5X40MM UNI-REDUCTION 2S
14554545	RSS SCREW, 4.5X45MM UNI-REDUCTION 2S
14554550	RSS SCREW, 4.5X50MM UNI-REDUCTION 2S
14554555	RSS SCREW, 4.5X55MM UNI-REDUCTION 2S
14554560	RSS SCREW, 4.5X60MM UNI-REDUCTION 2S
14554565	RSS SCREW, 4.5X65MM UNI-REDUCTION 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
14554570	RSS SCREW, 4.5X70MM UNI-REDUCTION 2S
14554575	RSS SCREW, 4.5X75MM UNI-REDUCTION 2S
14554580	RSS SCREW, 4.5X80MM UNI-REDUCTION 2S
14554585	RSS SCREW, 4.5X85MM UNI-REDUCTION 2S
14554590	RSS SCREW, 4.5X90MM UNI-REDUCTION 2S
14554595	RSS SCREW, 4.5X95MM UNI-REDUCTION 2S
14555010	RSS SCREW, 5.0X100MM UNI-REDUCTION 2S
14555011	RSS SCREW, 5.0X110MM UNI-REDUCTION 2S
14555012	RSS SCREW, 5.0X120MM UNI-REDUCTION 2S
14555017	RSS SCREW, 5.0X17.5MM UNI-REDUCTION 2S
14555020	RSS SCREW, 5.0X20MM UNI-REDUCTION 2S
14555022	RSS SCREW, 5.0X22.5MM UNI-REDUCTION 2S
14555025	RSS SCREW, 5.0X25MM UNI-REDUCTION 2S
14555027	RSS SCREW, 5.0X27.5MM UNI-REDUCTION 2S
14555030	RSS SCREW, 5.0X30MM UNI-REDUCTION 2S
14555035	RSS SCREW, 5.0X35MM UNI-REDUCTION 2S
14555040	RSS SCREW, 5.0X40MM UNI-REDUCTION 2S
14555045	RSS SCREW, 5.0X45MM UNI-REDUCTION 2S
14555050	RSS SCREW, 5.0X50MM UNI-REDUCTION 2S
14555055	RSS SCREW, 5.0X55MM UNI-REDUCTION 2S
14555060	RSS SCREW, 5.0X60MM UNI-REDUCTION 2S
14555065	RSS SCREW, 5.0X65MM UNI-REDUCTION 2S
14555070	RSS SCREW, 5.0X70MM UNI-REDUCTION 2S
14555075	RSS SCREW, 5.0X75MM UNI-REDUCTION 2S
14555080	RSS SCREW, 5.0X80MM UNI-REDUCTION 2S
14555085	RSS SCREW, 5.0X85MM UNI-REDUCTION 2S
14555090	RSS SCREW, 5.0X90MM UNI-REDUCTION 2S
14555095	RSS SCREW, 5.0X95MM UNI-REDUCTION 2S
14555510	RSS SCREW, 5.5X100MM UNI-REDUCTION 2S
14555511	RSS SCREW, 5.5X110MM UNI-REDUCTION 2S
14555512	RSS SCREW, 5.5X120MM UNI-REDUCTION 2S
14555517	RSS SCREW, 5.5X17.5MM UNI-REDUCTION 2S
14555520	RSS SCREW, 5.5X20MM UNI-REDUCTION 2S
14555522	RSS SCREW, 5.5X22.5MM UNI-REDUCTION 2S
14555525	RSS SCREW, 5.5X25MM UNI-REDUCTION 2S
14555527	RSS SCREW, 5.5X27.5MM UNI-REDUCTION 2S
14555530	RSS SCREW, 5.5X30MM UNI-REDUCTION 2S
14555535	RSS SCREW, 5.5X35MM UNI-REDUCTION 2S
14555540	RSS SCREW, 5.5X40MM UNI-REDUCTION 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14555545	RSS SCREW, 5.5X45MM UNI-REDUCTION 2S
14555550	RSS SCREW, 5.5X50MM UNI-REDUCTION 2S
14555555	RSS SCREW, 5.5X55MM UNI-REDUCTION 2S
14555560	RSS SCREW, 5.5X60MM UNI-REDUCTION 2S
14555565	RSS SCREW, 5.5X65MM UNI-REDUCTION 2S
14555570	RSS SCREW, 5.5X70MM UNI-REDUCTION 2S
14555575	RSS SCREW, 5.5X75MM UNI-REDUCTION 2S
14555580	RSS SCREW, 5.5X80MM UNI-REDUCTION 2S
14555585	RSS SCREW, 5.5X85MM UNI-REDUCTION 2S
14555590	RSS SCREW, 5.5X90MM UNI-REDUCTION 2S
14555595	RSS SCREW, 5.5X95MM UNI-REDUCTION 2S
14556010	RSS SCREW, 6.0X100MM UNI-REDUCTION 2S
14556011	RSS SCREW, 6.0X110MM UNI-REDUCTION 2S
14556012	RSS SCREW, 6.0X120MM UNI-REDUCTION 2S
14556017	RSS SCREW, 6.0X17.5MM UNI-REDUCTION 2S
14556020	RSS SCREW, 6.0X20MM UNI-REDUCTION 2S
14556022	RSS SCREW, 6.0X22.5MM UNI-REDUCTION 2S
14556025	RSS SCREW, 6.0X25MM UNI-REDUCTION 2S
14556027	RSS SCREW, 6.0X27.5MM UNI-REDUCTION 2S
14556030	RSS SCREW, 6.0X30MM UNI-REDUCTION 2S
14556035	RSS SCREW, 6.0X35MM UNI-REDUCTION 2S
14556040	RSS SCREW, 6.0X40MM UNI-REDUCTION 2S
14556045	RSS SCREW, 6.0X45MM UNI-REDUCTION 2S
14556050	RSS SCREW, 6.0X50MM UNI-REDUCTION 2S
14556055	RSS SCREW, 6.0X55MM UNI-REDUCTION 2S
14556060	RSS SCREW, 6.0X60MM UNI-REDUCTION 2S
14556065	RSS SCREW, 6.0X65MM UNI-REDUCTION 2S
14556070	RSS SCREW, 6.0X70MM UNI-REDUCTION 2S
14556075	RSS SCREW, 6.0X75MM UNI-REDUCTION 2S
14556080	RSS SCREW, 6.0X80MM UNI-REDUCTION 2S
14556085	RSS SCREW, 6.0X85MM UNI-REDUCTION 2S
14556090	RSS SCREW, 6.0X90MM UNI-REDUCTION 2S
14556095	RSS SCREW, 6.0X95MM UNI-REDUCTION 2S
14556510	RSS SCREW, 6.5X100MM UNI-REDUCTION 2S
14556511	RSS SCREW, 6.5X110MM UNI-REDUCTION 2S
14556512	RSS SCREW, 6.5X120MM UNI-REDUCTION 2S
14556517	RSS SCREW, 6.5X17.5MM UNI-REDUCTION 2S
14556520	RSS SCREW, 6.5X20MM UNI-REDUCTION 2S
14556522	RSS SCREW, 6.5X22.5MM UNI-REDUCTION 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14556525	RSS SCREW, 6.5X25MM UNI-REDUCTION 2S
14556527	RSS SCREW, 6.5X27.5MM UNI-REDUCTION 2S
14556530	RSS SCREW, 6.5X30MM UNI-REDUCTION 2S
14556535	RSS SCREW, 6.5X35MM UNI-REDUCTION 2S
14556540	RSS SCREW, 6.5X40MM UNI-REDUCTION 2S
14556545	RSS SCREW, 6.5X45MM UNI-REDUCTION 2S
14556550	RSS SCREW, 6.5X50MM UNI-REDUCTION 2S
14556555	RSS SCREW, 6.5X55MM UNI-REDUCTION 2S
14556560	RSS SCREW, 6.5X60MM UNI-REDUCTION 2S
14556565	RSS SCREW, 6.5X65MM UNI-REDUCTION 2S
14556570	RSS SCREW, 6.5X70MM UNI-REDUCTION 2S
14556575	RSS SCREW, 6.5X75MM UNI-REDUCTION 2S
14556580	RSS SCREW, 6.5X80MM UNI-REDUCTION 2S
14556585	RSS SCREW, 6.5X85MM UNI-REDUCTION 2S
14556590	RSS SCREW, 6.5X90MM UNI-REDUCTION 2S
14556595	RSS SCREW, 6.5X95MM UNI-REDUCTION 2S
14584017	RSS SCREW, 4.0X17.5MM MONOAXIAL 2S
14584020	RSS SCREW, 4.0X20MM MONOAXIAL 2S
14584022	RSS SCREW, 4.0X22.5MM MONOAXIAL 2S
14584025	RSS SCREW, 4.0X25MM MONOAXIAL 2S
14584027	RSS SCREW, 4.0X27.5MM MONOAXIAL 2S
14584030	RSS SCREW, 4.0X30MM MONOAXIAL 2S
14584035	RSS SCREW, 4.0X35MM MONOAXIAL 2S
14584517	RSS SCREW, 4.5X17.5MM MONOAXIAL 2S
14584520	RSS SCREW, 4.5X20MM MONOAXIAL 2S
14584522	RSS SCREW, 4.5X22.5MM MONOAXIAL 2S
14584525	RSS SCREW, 4.5X25MM MONOAXIAL 2S
14584527	RSS SCREW, 4.5X27.5MM MONOAXIAL 2S
14584530	RSS SCREW, 4.5X30MM MONOAXIAL 2S
14584535	RSS SCREW, 4.5X35MM MONOAXIAL 2S
14584540	RSS SCREW, 4.5X40MM MONOAXIAL 2S
14585020	RSS SCREW, 5.0X20MM MONOAXIAL 2S
14585022	RSS SCREW, 5.0X22.5MM MONOAXIAL 2S
14585025	RSS SCREW, 5.0X25MM MONOAXIAL 2S
14585027	RSS SCREW, 5.0X27.5MM MONOAXIAL 2S
14585030	RSS SCREW, 5.0X30MM MONOAXIAL 2S
14585035	RSS SCREW, 5.0X35MM MONOAXIAL 2S
14585040	RSS SCREW, 5.0X40MM MONOAXIAL 2S
14585045	RSS SCREW, 5.0X45MM MONOAXIAL 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14585520	RSS SCREW, 5.5X20MM MONOAXIAL 2S
14585522	RSS SCREW, 5.5X22.5MM MONOAXIAL 2S
14585525	RSS SCREW, 5.5X25MM MONOAXIAL 2S
14585527	RSS SCREW, 5.5X27.5MM MONOAXIAL 2S
14585530	RSS SCREW, 5.5X30MM MONOAXIAL 2S
14585535	RSS SCREW, 5.5X35MM MONOAXIAL 2S
14585540	RSS SCREW, 5.5X40MM MONOAXIAL 2S
14585545	RSS SCREW, 5.5X45MM MONOAXIAL 2S
14586025	RSS SCREW, 6.0X25MM MONOAXIAL 2S
14586030	RSS SCREW, 6.0X30MM MONOAXIAL 2S
14586035	RSS SCREW, 6.0X35MM MONOAXIAL 2S
14586040	RSS SCREW, 6.0X40MM MONOAXIAL 2S
14586045	RSS SCREW, 6.0X45MM MONOAXIAL 2S
14586530	RSS SCREW, 6.5X30MM MONOAXIAL 2S
14586535	RSS SCREW, 6.5X35MM MONOAXIAL 2S
14586540	RSS SCREW, 6.5X40MM MONOAXIAL 2S
14586545	RSS SCREW, 6.5X45MM MONOAXIAL 2S
14611111	RMM MOD TULIP
14615030	RMM SCREW, 5.0X30MM POLYAXIAL 4S
14615035	RMM SCREW, 5.0X35MM POLYAXIAL 4S
14615040	RMM SCREW, 5.0X40MM POLYAXIAL 4S
14615045	RMM SCREW, 5.0X45MM POLYAXIAL 4S
14615050	RMM SCREW, 5.0X50MM POLYAXIAL 4S
14615530	RMM SCREW, 5.5X30MM POLYAXIAL 4S
14615535	RMM SCREW, 5.5X35MM POLYAXIAL 4S
14615540	RMM SCREW, 5.5X40MM POLYAXIAL 4S
14615545	RMM SCREW, 5.5X45MM POLYAXIAL 4S
14615550	RMM SCREW, 5.5X50MM POLYAXIAL 4S
14616530	RMM SCREW, 6.5X30MM POLYAXIAL 4S
14616535	RMM SCREW, 6.5X35MM POLYAXIAL 4S
14616540	RMM SCREW, 6.5X40MM POLYAXIAL 4S
14616545	RMM SCREW, 6.5X45MM POLYAXIAL 4S
14616550	RMM SCREW, 6.5X50MM POLYAXIAL 4S
14617530	RMM SCREW, 7.5X30MM POLYAXIAL 4S
14617535	RMM SCREW, 7.5X35MM POLYAXIAL 4S
14617540	RMM SCREW, 7.5X40MM POLYAXIAL 4S
14617545	RMM SCREW, 7.5X45MM POLYAXIAL 4S
14617550	RMM SCREW, 7.5X50MM POLYAXIAL 4S
14618540	RMM SCREW, 8.5X40MM POLYAXIAL 4S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
15055200	GSB CPTI ROD, 5.5X200MM STRAIGHT
15055300	GSB CPTI ROD, 5.5X300MM STRAIGHT
15055400	GSB CPTI ROD, 5.5X400MM STRAIGHT
15060300	GSB CPTI ROD, 6.0X300MM STRAIGHT
15060400	GSB CPTI ROD, 6.0X400MM STRAIGHT
15060500	GSB CPTI ROD, 6.0X500MM STRAIGHT
15105300	RELIN-O TI ROD, 300MM 5.5-5.0 TAPERED
15105400	RELIN-O TI ROD, 400MM 5.5-5.0 TAPERED
15105500	RELIN-O TI ROD, 500MM 5.5-5.0 TAPERED
15155200	GSB TI ROD, 5.5X200MM STRAIGHT
15155300	GSB TI ROD, 5.5X300MM STRAIGHT
15155500	GSB TI ROD, 5.5X500MM STRAIGHT
15160500	GSB TI ROD, 6.0X500MM STRAIGHT
15165500	RELIN-O TI ROD, 500MM 6.0-5.0 TAPERED
15205300	RELIN-O COCR ROD, 300MM 5.5-5.0 TAPERED
15205400	RELIN-O COCR ROD, 400MM 5.5-5.0 TAPERED
15205500	RELIN-O COCR ROD, 500MM 5.5-5.0 TAPERED
15205600	RELIN-O COCR ROD, 600MM 5.5-5.0 TAPERED
15255300	GSB COCR ROD, 5.5X300MM STRAIGHT
15255500	GSB COCR ROD, 5.5X500MM STRAIGHT
15260500	GSB COCR ROD, 6.0X500MM STRAIGHT
15265500	RELIN-O COCR ROD, 500MM 6.0-5.0 TAPERED
15265600	RELIN-O COCR ROD, 600MM 6.0-5.0 TAPERED
15350300	RELIN-O TI ROD, 5.0X300MM STRAIGHT
15350500	RELIN-O TI ROD, 5.0X500MM STRAIGHT
15355035	RELIN-O TI ROD, 5.5X35MM STRAIGHT
15355040	RELIN-O TI ROD, 5.5X40MM STRAIGHT
15355045	RELIN-O TI ROD, 5.5X45MM STRAIGHT
15355050	RELIN-O TI ROD, 5.5X50MM STRAIGHT
15355060	RELIN-O TI ROD, 5.5X60MM STRAIGHT
15355070	RELIN-O TI ROD, 5.5X70MM STRAIGHT
15355075	RELIN-O TI ROD, 5.5X75MM STRAIGHT
15355080	RELIN-O TI ROD, 5.5X80MM STRAIGHT
15355085	RELIN-O TI ROD, 5.5X85MM STRAIGHT
15355090	RELIN-O TI ROD, 5.5X90MM STRAIGHT
15355095	RELIN-O TI ROD, 5.5X95MM STRAIGHT
15355100	RELIN-O TI ROD, 5.5X100MM STRAIGHT
15355110	RELIN-O TI ROD, 5.5X110MM STRAIGHT
15355120	RELIN-O TI ROD, 5.5X120MM STRAIGHT



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
15355130	RELIN-O TI ROD, 5.5X130MM STRAIGHT
15355140	RELIN-O TI ROD, 5.5X140MM STRAIGHT
15355150	RELIN-O TI ROD, 5.5X150MM STRAIGHT
15355160	RELIN-O TI ROD, 5.5X160MM STRAIGHT
15355170	RELIN-O TI ROD, 5.5X170MM STRAIGHT
15355180	RELIN-O TI ROD, 5.5X180MM STRAIGHT
15355190	RELIN-O TI ROD, 5.5X190MM STRAIGHT
15355200	RELIN-O TI ROD, 5.5X200MM STRAIGHT
15355300	RELIN-O TI ROD, 5.5X300MM STRAIGHT
15355400	RELIN-O TI ROD, 5.5X400MM STRAIGHT
15355500	RELIN-O TI ROD, 5.5X500MM STRAIGHT
15355600	RELIN-O TI ROD, 5.5X600MM STRAIGHT
15360035	RELIN-O TI ROD, 6.0X35MM STRAIGHT
15360040	RELIN-O TI ROD, 6.0X40MM STRAIGHT
15360045	RELIN-O TI ROD, 6.0X45MM STRAIGHT
15360050	RELIN-O TI ROD, 6.0X50MM STRAIGHT
15360060	RELIN-O TI ROD, 6.0X60MM STRAIGHT
15360070	RELIN-O TI ROD, 6.0X70MM STRAIGHT
15360080	RELIN-O TI ROD, 6.0X80MM STRAIGHT
15360090	RELIN-O TI ROD, 6.0X90MM STRAIGHT
15360100	RELIN-O TI ROD, 6.0X100MM STRAIGHT
15360200	RELIN-O TI ROD, 6.0X200MM STRAIGHT
15360300	RELIN-O TI ROD, 6.0X300MM STRAIGHT
15360500	RELIN-O TI ROD, 6.0X500MM STRAIGHT
15360600	RELIN-O TI ROD, 6.0X600MM STRAIGHT
15450200	RELIN-O COCR ROD, 5.0X200MM STRAIGHT
15450300	RELIN-O COCR ROD, 5.0X300MM STRAIGHT
15450500	RELIN-O COCR ROD, 5.0X500MM STRAIGHT
15455300	RELIN-O COCR ROD, 5.5X300MM STRAIGHT
15455400	RELIN-O COCR ROD, 5.5X400MM STRAIGHT
15455500	RELIN-O COCR ROD, 5.5X500MM STRAIGHT
15455600	RELIN-O COCR ROD, 5.5X600MM STRAIGHT
15460300	RELIN-O COCR ROD, 6.0X300MM STRAIGHT
15460500	RELIN-O COCR ROD, 6.0X500MM STRAIGHT
15460600	RELIN-O COCR ROD, 6.0X600MM STRAIGHT
16010510	RELIN MAS SCREW, 10.5X100MM 2C POLYILIAC
16010511	RELIN MAS SCREW, 10.5X110MM 2C POLYILIAC
16010512	RELIN MAS SCREW, 10.5X120MM 2C POLYILIAC
16010580	RELIN MAS SCREW, 10.5X80MM 2C POLYILIAC



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
16010590	RELIN MAS SCREW,10.5X90MM 2C POLYILIAC
16011510	RELIN MAS SCREW,11.5X100MM 2C POLYILIAC
16011590	RELIN MAS SCREW,11.5X90MM 2C POLYILIAC
16012510	RELIN MAS SCREW,12.5X100MM 2C POLYILIAC
16012511	RELIN MAS SCREW,12.5X110MM 2C POLYILIAC
16012512	RELIN MAS SCREW,12.5X120MM 2C POLYILIAC
16012550	RELIN MAS SCREW, 12.5X50MM 2CPOLYAXIAL
16012560	RELIN MAS SCREW,12.5X60MM 2C POLYILIAC
16012570	RELIN MAS SCREW,12.5X70MM 2C POLYILIAC
16012580	RELIN MAS SCREW,12.5X80MM 2C POLYILIAC
16012590	RELIN MAS SCREW,12.5X90MM 2C POLYILIAC
16014525	RELIN MAS SCREW, 4.5X25MM 2C POLYAXIAL
16014530	RELIN MAS SCREW, 4.5X30MM 2C POLYAXIAL
16014535	RELIN MAS SCREW, 4.5X35MM 2C POLYAXIAL
16014540	RELIN MAS SCREW, 4.5X40MM 2C POLYAXIAL
16014545	RELIN MAS SCREW, 4.5X45MM 2C POLYAXIAL
16014550	RELIN MAS SCREW, 4.5X50MM 2C POLYAXIAL
16014555	RELIN MAS SCREW, 4.5X55MM 2C POLYAXIAL
16015030	RELIN MAS SCREW, 5.0X30MM 2C POLYAXIAL
16015035	RELIN MAS SCREW, 5.0X35MM 2C POLYAXIAL
16015040	RELIN MAS SCREW, 5.0X40MM 2C POLYAXIAL
16015045	RELIN MAS SCREW, 5.0X45MM 2C POLYAXIAL
16015050	RELIN MAS SCREW, 5.0X50MM 2C POLYAXIAL
16015055	RELIN MAS SCREW, 5.0X55MM 2C POLYAXIAL
16015525	RELIN MAS SCREW, 5.5X25MM 2C POLYAXIAL
16015530	RELIN MAS SCREW, 5.5X30MM 2C POLYAXIAL
16015535	RELIN MAS SCREW, 5.5X35MM 2C POLYAXIAL
16015540	RELIN MAS SCREW, 5.5X40MM 2C POLYAXIAL
16015545	RELIN MAS SCREW, 5.5X45MM 2C POLYAXIAL
16015550	RELIN MAS SCREW, 5.5X50MM 2C POLYAXIAL
16015555	RELIN MAS SCREW, 5.5X55MM 2C POLYAXIAL
16015560	RELIN MAS SCREW, 5.5X60MM 2C POLYILIAC
16015565	RELIN MAS SCREW, 5.5X65MM 2C POLYILIAC
16016030	RELIN MAS SCREW, 6.0X30MM 2C POLYAXIAL
16016035	RELIN MAS SCREW, 6.0X35MM 2C POLYAXIAL
16016040	RELIN MAS SCREW, 6.0X40MM 2C POLYAXIAL
16016045	RELIN MAS SCREW, 6.0X45MM 2C POLYAXIAL
16016050	RELIN MAS SCREW, 6.0X50MM 2C POLYAXIAL
16016055	RELIN MAS SCREW, 6.0X55MM 2C POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
16016510	RELIN MAS SCREW,6.5X100MM 2C POLYILIAC
16016525	RELIN MAS SCREW, 6.5X25MM 2C POLYAXIAL
16016530	RELIN MAS SCREW, 6.5X30MM 2C POLYAXIAL
16016535	RELIN MAS SCREW, 6.5X35MM 2C POLYAXIAL
16016540	RELIN MAS SCREW, 6.5X40MM 2C POLYAXIAL
16016545	RELIN MAS SCREW, 6.5X45MM 2C POLYAXIAL
16016550	RELIN MAS SCREW, 6.5X50MM 2C POLYAXIAL
16016555	RELIN MAS SCREW, 6.5X55MM 2C POLYAXIAL
16016560	RELIN MAS SCREW, 6.5X60MM 2C POLYILIAC
16016565	RELIN MAS SCREW, 6.5X65MM 2C POLYILIAC
16016570	RELIN MAS SCREW, 6.5X70MM 2C POLYILIAC
16017510	RELIN MAS SCREW,7.5X100MM 2C POLYILIAC
16017525	RELIN MAS SCREW, 7.5X25MM 2C POLYAXIAL
16017530	RELIN MAS SCREW, 7.5X30MM 2C POLYAXIAL
16017535	RELIN MAS SCREW, 7.5X35MM 2C POLYAXIAL
16017540	RELIN MAS SCREW, 7.5X40MM 2C POLYAXIAL
16017545	RELIN MAS SCREW, 7.5X45MM 2C POLYAXIAL
16017550	RELIN MAS SCREW, 7.5X50MM 2C POLYAXIAL
16017555	RELIN MAS SCREW, 7.5X55MM 2C POLYAXIAL
16017560	RELIN MAS SCREW, 7.5X60MM 2C POLYILIAC
16017565	RELIN MAS SCREW, 7.5X65MM 2C POLYILIAC
16017570	RELIN MAS SCREW, 7.5X70MM 2C POLYILIAC
16017580	RELIN MAS SCREW, 7.5X80MM 2C POLYILIAC
16017590	RELIN MAS SCREW, 7.5X90MM 2C POLYILIAC
16018510	RELIN MAS SCREW,8.5X100MM 2C POLYILIAC
16018530	RELIN MAS SCREW, 8.5X30MM 2C POLYAXIAL
16018535	RELIN MAS SCREW, 8.5X35MM 2C POLYAXIAL
16018540	RELIN MAS SCREW, 8.5X40MM 2C POLYAXIAL
16018545	RELIN MAS SCREW, 8.5X45MM 2C POLYAXIAL
16018550	RELIN MAS SCREW, 8.5X50MM 2C POLYAXIAL
16018555	RELIN MAS SCREW, 8.5X55MM 2C POLYAXIAL
16018560	RELIN MAS SCREW, 8.5X60MM 2C POLYILIAC
16018565	RELIN MAS SCREW, 8.5X65MM 2C POLYILIAC
16018570	RELIN MAS SCREW, 8.5X70MM 2C POLYILIAC
16018580	RELIN MAS SCREW, 8.5X80MM 2C POLYILIAC
16018590	RELIN MAS SCREW, 8.5X90MM 2C POLYILIAC
16019510	RELIN MAS SCREW,9.5X100MM 2C POLYILIAC
16019512	RELIN MAS SCREW,9.5X120MM 2C POLYILIAC
16019535	RELIN MAS SCREW, 9.5X35MM 2C POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
16019540	RELIN MAS SCREW, 9.5X40MM 2C POLYAXIAL
16019545	RELIN MAS SCREW, 9.5X45MM 2C POLYAXIAL
16019550	RELIN MAS SCREW, 9.5X50MM 2C POLYAXIAL
16019555	RELIN MAS SCREW, 9.5X55MM 2C POLYAXIAL
16019560	RELIN MAS SCREW, 9.5X60MM 2C POLYILIAC
16019570	RELIN MAS SCREW, 9.5X70MM 2C POLYILIAC
16019580	RELIN MAS SCREW, 9.5X80MM 2C POLYILIAC
16019590	RELIN MAS SCREW, 9.5X90MM 2C POLYILIAC
16024525	RELIN MAS RED SCREW, 4.5X25MM POLY 2C
16024530	RELIN MAS RED SCREW, 4.5X30MM POLY 2C
16024535	RELIN MAS RED SCREW, 4.5X35MM POLY 2C
16024540	RELIN MAS RED SCREW, 4.5X40MM POLY 2C
16024545	RELIN MAS RED SCREW, 4.5X45MM POLY 2C
16024550	RELIN MAS RED SCREW, 4.5X50MM POLY 2C
16024555	RELIN MAS RED SCREW, 4.5X55MM POLY 2C
16025530	RELIN MAS RED SCREW, 5.5X30MM POLY 2C
16025535	RELIN MAS RED SCREW, 5.5X35MM POLY 2C
16025540	RELIN MAS RED SCREW, 5.5X40MM POLY 2C
16025545	RELIN MAS RED SCREW, 5.5X45MM POLY 2C
16025550	RELIN MAS RED SCREW, 5.5X50MM POLY 2C
16025555	RELIN MAS RED SCREW, 5.5X55MM POLY 2C
16026040	RELIN MAS RED SCREW, 6.0X40MM POLY 2C
16026530	RELIN MAS RED SCREW, 6.5X30MM POLY 2C
16026535	RELIN MAS RED SCREW, 6.5X35MM POLY 2C
16026540	RELIN MAS RED SCREW, 6.5X40MM POLY 2C
16026545	RELIN MAS RED SCREW, 6.5X45MM POLY 2C
16026550	RELIN MAS RED SCREW, 6.5X50MM POLY 2C
16026555	RELIN MAS RED SCREW, 6.5X55MM POLY 2C
16026560	RELIN MAS-R SCREW,6.5X60MM POLYILIAC 2C
16027510	RELIN MAS-R SCREW,7.5X100 POLYILIAC 2C
16027535	RELIN MAS RED SCREW, 7.5X35MM POLY 2C
16027540	RELIN MAS RED SCREW, 7.5X40MM POLY 2C
16027545	RELIN MAS RED SCREW, 7.5X45MM POLY 2C
16027550	RELIN MAS RED SCREW, 7.5X50MM POLY 2C
16027555	RELIN MAS RED SCREW, 7.5X55MM POLY 2C
16027560	RELIN MAS-R SCREW,7.5X60MM POLYILIAC 2C
16028510	RELIN MAS-R SCREW,8.5X100 POLYILIAC 2C
16028535	RELIN MAS RED SCREW, 8.5X35MM POLY 2C
16028540	RELIN MAS RED SCREW, 8.5X40MM POLY 2C



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
16028545	RELIN MAS RED SCREW, 8.5X45MM POLY 2C
16028550	RELIN MAS RED SCREW, 8.5X50MM POLY 2C
16028555	RELIN MAS RED SCREW, 8.5X55MM POLY 2C
16028560	RELIN MAS-R SCREW,8.5X60MM POLYILIAC 2C
16028570	RELIN MAS-R SCREW,8.5X70MM POLYILIAC 2C
16028580	RELIN MAS-R SCREW,8.5X80MM POLYILIAC 2C
16028590	RELIN MAS-R SCREW,8.5X90MM POLYILIAC 2C
16029510	RELIN MAS-R SCREW,9.5X100 POLYILIAC 2C
16029511	RELIN MAS-R SCREW,9.5X110 POLYILIAC 2C
16029535	RELIN MAS RED SCREW, 9.5X35MM POLY 2C
16029540	RELIN MAS RED SCREW, 9.5X40MM POLY 2C
16029545	RELIN MAS RED SCREW, 9.5X45MM POLY 2C
16029550	RELIN MAS RED SCREW, 9.5X50MM POLY 2C
16029555	RELIN MAS RED SCREW, 9.5X55MM POLY 2C
16029560	RELIN MAS-R SCREW,9.5X60MM POLYILIAC 2C
16029570	RELIN MAS-R SCREW,9.5X70MM POLYILIAC 2C
16029580	RELIN MAS-R SCREW,9.5X80MM POLYILIAC 2C
16029590	RELIN MAS-R SCREW,9.5X90MM POLYILIAC 2C
16084525	RELIN MAS TRAUMA SCREW, 4.5X25MM CANN
16084530	RELIN MAS TRAUMA SCREW, 4.5X30MM CANN
16084535	RELIN MAS TRAUMA SCREW, 4.5X35MM CANN
16084540	RELIN MAS TRAUMA SCREW, 4.5X40MM CANN
16084545	RELIN MAS TRAUMA SCREW, 4.5X45MM CANN
16084550	RELIN MAS TRAUMA SCREW, 4.5X50MM CANN
16084555	RELIN MAS TRAUMA SCREW, 4.5X55MM CANN
16085525	RELIN MAS TRAUMA SCREW, 5.5X25MM CANN
16085530	RELIN MAS TRAUMA SCREW, 5.5X30MM CANN
16085535	RELIN MAS TRAUMA SCREW, 5.5X35MM CANN
16085540	RELIN MAS TRAUMA SCREW, 5.5X40MM CANN
16085545	RELIN MAS TRAUMA SCREW, 5.5X45MM CANN
16085550	RELIN MAS TRAUMA SCREW, 5.5X50MM CANN
16085555	RELIN MAS TRAUMA SCREW, 5.5X55MM CANN
16085560	RELIN MAS TRAUMA SCREW, 5.5X60MM CANN
16086525	RELIN MAS TRAUMA SCREW, 6.5X25MM CANN
16086530	RELIN MAS TRAUMA SCREW, 6.5X30MM CANN
16086535	RELIN MAS TRAUMA SCREW, 6.5X35MM CANN
16086555	RELIN MAS TRAUMA SCREW, 6.5X55MM CANN
16086560	RELIN MAS TRAUMA SCREW, 6.5X60MM CANN
16087525	RELIN MAS TRAUMA SCREW, 7.5X25MM CANN



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
16087530	RELINE MAS TRAUMA SCREW, 7.5X30MM CANN
16087535	RELINE MAS TRAUMA SCREW, 7.5X35MM CANN
16087540	RELINE MAS TRAUMA SCREW, 7.5X40MM CANN
16087545	RELINE MAS TRAUMA SCREW, 7.5X45MM CANN
16087550	RELINE MAS TRAUMA SCREW, 7.5X50MM CANN
16087555	RELINE MAS TRAUMA SCREW, 7.5X55MM CANN
16087560	RELINE MAS TRAUMA SCREW, 7.5X60MM CANN
16088540	RELINE MAS TRAUMA SCREW, 8.5X40MM CANN
16088545	RELINE MAS TRAUMA SCREW, 8.5X45MM CANN
16088550	RELINE MAS TRAUMA SCREW, 8.5X50MM CANN
16088555	RELINE MAS TRAUMA SCREW, 8.5X55MM CANN
16155500	GSB TI ROD, 5.5X500MM TL
16171111	RELINE TULIP, MOD
16184530	RELINE MAS MOD SCREW SHANK, 4.5X30MM 2C
16184535	RELINE MAS MOD SCREW SHANK, 4.5X35MM 2C
16184540	RELINE MAS MOD SCREW SHANK, 4.5X40MM 2C
16184545	RELINE MAS MOD SCREW SHANK, 4.5X45MM 2C
16184550	RELINE MAS MOD SCREW SHANK, 4.5X50MM 2C
16185530	RELINE MAS MOD SCREW SHANK, 5.5X30MM 2C
16185535	RELINE MAS MOD SCREW SHANK, 5.5X35MM 2C
16185540	RELINE MAS MOD SCREW SHANK, 5.5X40MM 2C
16185545	RELINE MAS MOD SCREW SHANK, 5.5X45MM 2C
16185550	RELINE MAS MOD SCREW SHANK, 5.5X50MM 2C
16185555	RELINE MAS MOD SCREW SHANK, 5.5X55MM 2C
16185560	RELINE MAS MOD SCREW SHANK, 5.5X60MM 2C
16185565	RELINE MAS MOD SCREW SHANK, 5.5X65MM 2C
16186530	RELINE MAS MOD SCREW SHANK, 6.5X30MM 2C
16186535	RELINE MAS MOD SCREW SHANK, 6.5X35MM 2C
16186540	RELINE MAS MOD SCREW SHANK, 6.5X40MM 2C
16186545	RELINE MAS MOD SCREW SHANK, 6.5X45MM 2C
16186550	RELINE MAS MOD SCREW SHANK, 6.5X50MM 2C
16186555	RELINE MAS MOD SCREW SHANK, 6.5X55MM 2C
16186560	RELINE MAS MOD SCREW SHANK, 6.5X60MM 2C
16186565	RELINE MAS MOD SCREW SHANK, 6.5X65MM 2C
16187530	RELINE MAS MOD SCREW SHANK, 7.5X30MM 2C
16187535	RELINE MAS MOD SCREW SHANK, 7.5X35MM 2C
16187540	RELINE MAS MOD SCREW SHANK, 7.5X40MM 2C
16187545	RELINE MAS MOD SCREW SHANK, 7.5X45MM 2C
16187550	RELINE MAS MOD SCREW SHANK, 7.5X50MM 2C



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
16187555	RELIN MAS MOD SCREW SHANK, 7.5X55MM 2C
16187560	RELIN MAS MOD SCREW SHANK, 7.5X60MM 2C
16187565	RELIN MAS MOD SCREW SHANK, 7.5X65MM 2C
16187570	RELIN MAS MOD SCREW SHANK, 7.5X70MM 2C
16187580	RELIN MAS MOD SCREW SHANK, 7.5X80MM 2C
16187590	RELIN MAS MOD SCREW SHANK, 7.5X90MM 2C
16188530	RELIN MAS MOD SCREW SHANK, 8.5X30MM 2C
16188535	RELIN MAS MOD SCREW SHANK, 8.5X35MM 2C
16188540	RELIN MAS MOD SCREW SHANK, 8.5X40MM 2C
16188545	RELIN MAS MOD SCREW SHANK, 8.5X45MM 2C
16188550	RELIN MAS MOD SCREW SHANK, 8.5X50MM 2C
16188555	RELIN MAS MOD SCREW SHANK, 8.5X55MM 2C
16188560	RELIN MAS MOD SCREW SHANK, 8.5X60MM 2C
16188565	RELIN MAS MOD SCREW SHANK, 8.5X65MM 2C
16188570	RELIN MAS MOD SCREW SHANK, 8.5X70MM 2C
16188580	RELIN MAS MOD SCREW SHANK, 8.5X80MM 2C
16188590	RELIN MAS MOD SCREW SHANK, 8.5X90MM 2C
16189535	RELIN MAS MOD SCREW SHANK, 9.5X35MM 2C
16189540	RELIN MAS MOD SCREW SHANK, 9.5X40MM 2C
16189545	RELIN MAS MOD SCREW SHANK, 9.5X45MM 2C
16189550	RELIN MAS MOD SCREW SHANK, 9.5X50MM 2C
16191111	RELIN MAS MOD, REDUCTION EXTENSION
16210510	RELIN SCREW, 10.5X100MM 2FC POLYAXIAL
16210530	RELIN SCREW, 10.5X30MM 2FC POLYAXIAL
16210535	RELIN SCREW, 10.5X35MM 2FC POLYAXIAL
16210540	RELIN SCREW, 10.5X40MM 2FC POLYAXIAL
16210545	RELIN SCREW, 10.5X45MM 2FC POLYAXIAL
16210550	RELIN SCREW, 10.5X50MM 2FC POLYAXIAL
16210555	RELIN SCREW, 10.5X55MM 2FC POLYAXIAL
16210560	RELIN SCREW, 10.5X60MM 2FC POLYAXIAL
16210565	RELIN SCREW, 10.5X65MM 2FC POLYAXIAL
16210570	RELIN SCREW, 10.5X70MM 2FC POLYAXIAL
16210575	RELIN SCREW, 10.5X75MM 2FC POLYAXIAL
16210580	RELIN SCREW, 10.5X80MM 2FC POLYAXIAL
16210585	RELIN SCREW, 10.5X85MM 2FC POLYAXIAL
16210590	RELIN SCREW, 10.5X90MM 2FC POLYAXIAL
16210595	RELIN SCREW, 10.5X95MM 2FC POLYAXIAL
16211010	RELIN SCREW, 10.0X100MM 2FC POLYAXIAL
16211030	RELIN SCREW, 10.0X30MM 2FC POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
16211035	RELINER SCREW, 10.0X35MM 2FC POLYAXIAL
16211040	RELINER SCREW, 10.0X40MM 2FC POLYAXIAL
16211045	RELINER SCREW, 10.0X45MM 2FC POLYAXIAL
16211050	RELINER SCREW, 10.0X50MM 2FC POLYAXIAL
16211055	RELINER SCREW, 10.0X55MM 2FC POLYAXIAL
16211060	RELINER SCREW, 10.0X60MM 2FC POLYAXIAL
16211065	RELINER SCREW, 10.0X65MM 2FC POLYAXIAL
16211070	RELINER SCREW, 10.0X70MM 2FC POLYAXIAL
16211075	RELINER SCREW, 10.0X75MM 2FC POLYAXIAL
16211080	RELINER SCREW, 10.0X80MM 2FC POLYAXIAL
16211085	RELINER SCREW, 10.0X85MM 2FC POLYAXIAL
16211090	RELINER SCREW, 10.0X90MM 2FC POLYAXIAL
16211095	RELINER SCREW, 10.0X95MM 2FC POLYAXIAL
16214510	RELINER SCREW, 4.5X100MM 2FC POLYAXIAL
16214530	RELINER SCREW, 4.5X30MM 2FC POLYAXIAL
16214535	RELINER SCREW, 4.5X35MM 2FC POLYAXIAL
16214540	RELINER SCREW, 4.5X40MM 2FC POLYAXIAL
16214545	RELINER SCREW, 4.5X45MM 2FC POLYAXIAL
16214550	RELINER SCREW, 4.5X50MM 2FC POLYAXIAL
16214555	RELINER SCREW, 4.5X55MM 2FC POLYAXIAL
16214560	RELINER SCREW, 4.5X60MM 2FC POLYAXIAL
16214565	RELINER SCREW, 4.5X65MM 2FC POLYAXIAL
16214570	RELINER SCREW, 4.5X70MM 2FC POLYAXIAL
16214575	RELINER SCREW, 4.5X75MM 2FC POLYAXIAL
16214580	RELINER SCREW, 4.5X80MM 2FC POLYAXIAL
16214585	RELINER SCREW, 4.5X85MM 2FC POLYAXIAL
16214590	RELINER SCREW, 4.5X90MM 2FC POLYAXIAL
16214595	RELINER SCREW, 4.5X95MM 2FC POLYAXIAL
16215010	RELINER SCREW, 5.0X100MM 2FC POLYAXIAL
16215030	RELINER SCREW, 5.0X30MM 2FC POLYAXIAL
16215035	RELINER SCREW, 5.0X35MM 2FC POLYAXIAL
16215040	RELINER SCREW, 5.0X40MM 2FC POLYAXIAL
16215045	RELINER SCREW, 5.0X45MM 2FC POLYAXIAL
16215050	RELINER SCREW, 5.0X50MM 2FC POLYAXIAL
16215055	RELINER SCREW, 5.0X55MM 2FC POLYAXIAL
16215060	RELINER SCREW, 5.0X60MM 2FC POLYAXIAL
16215065	RELINER SCREW, 5.0X65MM 2FC POLYAXIAL
16215070	RELINER SCREW, 5.0X70MM 2FC POLYAXIAL
16215075	RELINER SCREW, 5.0X75MM 2FC POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
16215080	RELINER SCREW, 5.0X80MM 2FC POLYAXIAL
16215085	RELINER SCREW, 5.0X85MM 2FC POLYAXIAL
16215090	RELINER SCREW, 5.0X90MM 2FC POLYAXIAL
16215095	RELINER SCREW, 5.0X95MM 2FC POLYAXIAL
16215510	RELINER SCREW, 5.5X100MM 2FC POLYAXIAL
16215530	RELINER SCREW, 5.5X30MM 2FC POLYAXIAL
16215535	RELINER SCREW, 5.5X35MM 2FC POLYAXIAL
16215540	RELINER SCREW, 5.5X40MM 2FC POLYAXIAL
16215545	RELINER SCREW, 5.5X45MM 2FC POLYAXIAL
16215550	RELINER SCREW, 5.5X50MM 2FC POLYAXIAL
16215555	RELINER SCREW, 5.5X55MM 2FC POLYAXIAL
16215560	RELINER SCREW, 5.5X60MM 2FC POLYAXIAL
16215565	RELINER SCREW, 5.5X65MM 2FC POLYAXIAL
16215570	RELINER SCREW, 5.5X70MM 2FC POLYAXIAL
16215575	RELINER SCREW, 5.5X75MM 2FC POLYAXIAL
16215580	RELINER SCREW, 5.5X80MM 2FC POLYAXIAL
16215585	RELINER SCREW, 5.5X85MM 2FC POLYAXIAL
16215590	RELINER SCREW, 5.5X90MM 2FC POLYAXIAL
16215595	RELINER SCREW, 5.5X95MM 2FC POLYAXIAL
16216010	RELINER SCREW, 6.0X100MM 2FC POLYAXIAL
16216030	RELINER SCREW, 6.0X30MM 2FC POLYAXIAL
16216035	RELINER SCREW, 6.0X35MM 2FC POLYAXIAL
16216040	RELINER SCREW, 6.0X40MM 2FC POLYAXIAL
16216045	RELINER SCREW, 6.0X45MM 2FC POLYAXIAL
16216050	RELINER SCREW, 6.0X50MM 2FC POLYAXIAL
16216055	RELINER SCREW, 6.0X55MM 2FC POLYAXIAL
16216060	RELINER SCREW, 6.0X60MM 2FC POLYAXIAL
16216065	RELINER SCREW, 6.0X65MM 2FC POLYAXIAL
16216070	RELINER SCREW, 6.0X70MM 2FC POLYAXIAL
16216075	RELINER SCREW, 6.0X75MM 2FC POLYAXIAL
16216080	RELINER SCREW, 6.0X80MM 2FC POLYAXIAL
16216085	RELINER SCREW, 6.0X85MM 2FC POLYAXIAL
16216090	RELINER SCREW, 6.0X90MM 2FC POLYAXIAL
16216095	RELINER SCREW, 6.0X95MM 2FC POLYAXIAL
16216510	RELINER SCREW, 6.5X100MM 2FC POLYAXIAL
16216530	RELINER SCREW, 6.5X30MM 2FC POLYAXIAL
16216535	RELINER SCREW, 6.5X35MM 2FC POLYAXIAL
16216540	RELINER SCREW, 6.5X40MM 2FC POLYAXIAL
16216545	RELINER SCREW, 6.5X45MM 2FC POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
16216550	RELINER SCREW, 6.5X50MM 2FC POLYAXIAL
16216555	RELINER SCREW, 6.5X55MM 2FC POLYAXIAL
16216560	RELINER SCREW, 6.5X60MM 2FC POLYAXIAL
16216565	RELINER SCREW, 6.5X65MM 2FC POLYAXIAL
16216570	RELINER SCREW, 6.5X70MM 2FC POLYAXIAL
16216575	RELINER SCREW, 6.5X75MM 2FC POLYAXIAL
16216580	RELINER SCREW, 6.5X80MM 2FC POLYAXIAL
16216585	RELINER SCREW, 6.5X85MM 2FC POLYAXIAL
16216590	RELINER SCREW, 6.5X90MM 2FC POLYAXIAL
16216595	RELINER SCREW, 6.5X95MM 2FC POLYAXIAL
16217010	RELINER SCREW, 7.0X100MM 2FC POLYAXIAL
16217030	RELINER SCREW, 7.0X30MM 2FC POLYAXIAL
16217035	RELINER SCREW, 7.0X35MM 2FC POLYAXIAL
16217040	RELINER SCREW, 7.0X40MM 2FC POLYAXIAL
16217045	RELINER SCREW, 7.0X45MM 2FC POLYAXIAL
16217050	RELINER SCREW, 7.0X50MM 2FC POLYAXIAL
16217055	RELINER SCREW, 7.0X55MM 2FC POLYAXIAL
16217060	RELINER SCREW, 7.0X60MM 2FC POLYAXIAL
16217065	RELINER SCREW, 7.0X65MM 2FC POLYAXIAL
16217070	RELINER SCREW, 7.0X70MM 2FC POLYAXIAL
16217075	RELINER SCREW, 7.0X75MM 2FC POLYAXIAL
16217080	RELINER SCREW, 7.0X80MM 2FC POLYAXIAL
16217085	RELINER SCREW, 7.0X85MM 2FC POLYAXIAL
16217090	RELINER SCREW, 7.0X90MM 2FC POLYAXIAL
16217095	RELINER SCREW, 7.0X95MM 2FC POLYAXIAL
16217510	RELINER SCREW, 7.5X100MM 2FC POLYAXIAL
16217530	RELINER SCREW, 7.5X30MM 2FC POLYAXIAL
16217535	RELINER SCREW, 7.5X35MM 2FC POLYAXIAL
16217540	RELINER SCREW, 7.5X40MM 2FC POLYAXIAL
16217545	RELINER SCREW, 7.5X45MM 2FC POLYAXIAL
16217550	RELINER SCREW, 7.5X50MM 2FC POLYAXIAL
16217555	RELINER SCREW, 7.5X55MM 2FC POLYAXIAL
16217560	RELINER SCREW, 7.5X60MM 2FC POLYAXIAL
16217565	RELINER SCREW, 7.5X65MM 2FC POLYAXIAL
16217570	RELINER SCREW, 7.5X70MM 2FC POLYAXIAL
16217575	RELINER SCREW, 7.5X75MM 2FC POLYAXIAL
16217580	RELINER SCREW, 7.5X80MM 2FC POLYAXIAL
16217585	RELINER SCREW, 7.5X85MM 2FC POLYAXIAL
16217590	RELINER SCREW, 7.5X90MM 2FC POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
16217595	RELINER SCREW, 7.5X95MM 2FC POLYAXIAL
16218010	RELINER SCREW, 8.0X100MM 2FC POLYAXIAL
16218030	RELINER SCREW, 8.0X30MM 2FC POLYAXIAL
16218035	RELINER SCREW, 8.0X35MM 2FC POLYAXIAL
16218040	RELINER SCREW, 8.0X40MM 2FC POLYAXIAL
16218045	RELINER SCREW, 8.0X45MM 2FC POLYAXIAL
16218050	RELINER SCREW, 8.0X50MM 2FC POLYAXIAL
16218055	RELINER SCREW, 8.0X55MM 2FC POLYAXIAL
16218060	RELINER SCREW, 8.0X60MM 2FC POLYAXIAL
16218065	RELINER SCREW, 8.0X65MM 2FC POLYAXIAL
16218070	RELINER SCREW, 8.0X70MM 2FC POLYAXIAL
16218075	RELINER SCREW, 8.0X75MM 2FC POLYAXIAL
16218080	RELINER SCREW, 8.0X80MM 2FC POLYAXIAL
16218085	RELINER SCREW, 8.0X85MM 2FC POLYAXIAL
16218090	RELINER SCREW, 8.0X90MM 2FC POLYAXIAL
16218095	RELINER SCREW, 8.0X95MM 2FC POLYAXIAL
16218510	RELINER SCREW, 8.5X100MM 2FC POLYAXIAL
16218530	RELINER SCREW, 8.5X30MM 2FC POLYAXIAL
16218535	RELINER SCREW, 8.5X35MM 2FC POLYAXIAL
16218540	RELINER SCREW, 8.5X40MM 2FC POLYAXIAL
16218545	RELINER SCREW, 8.5X45MM 2FC POLYAXIAL
16218550	RELINER SCREW, 8.5X50MM 2FC POLYAXIAL
16218555	RELINER SCREW, 8.5X55MM 2FC POLYAXIAL
16218560	RELINER SCREW, 8.5X60MM 2FC POLYAXIAL
16218565	RELINER SCREW, 8.5X65MM 2FC POLYAXIAL
16218570	RELINER SCREW, 8.5X70MM 2FC POLYAXIAL
16218575	RELINER SCREW, 8.5X75MM 2FC POLYAXIAL
16218580	RELINER SCREW, 8.5X80MM 2FC POLYAXIAL
16218585	RELINER SCREW, 8.5X85MM 2FC POLYAXIAL
16218590	RELINER SCREW, 8.5X90MM 2FC POLYAXIAL
16218595	RELINER SCREW, 8.5X95MM 2FC POLYAXIAL
16219010	RELINER SCREW, 9.0X100MM 2FC POLYAXIAL
16219030	RELINER SCREW, 9.0X30MM 2FC POLYAXIAL
16219035	RELINER SCREW, 9.0X35MM 2FC POLYAXIAL
16219040	RELINER SCREW, 9.0X40MM 2FC POLYAXIAL
16219045	RELINER SCREW, 9.0X45MM 2FC POLYAXIAL
16219050	RELINER SCREW, 9.0X50MM 2FC POLYAXIAL
16219055	RELINER SCREW, 9.0X55MM 2FC POLYAXIAL
16219060	RELINER SCREW, 9.0X60MM 2FC POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
16219065	RELINER SCREW, 9.0X65MM 2FC POLYAXIAL
16219070	RELINER SCREW, 9.0X70MM 2FC POLYAXIAL
16219075	RELINER SCREW, 9.0X75MM 2FC POLYAXIAL
16219080	RELINER SCREW, 9.0X80MM 2FC POLYAXIAL
16219085	RELINER SCREW, 9.0X85MM 2FC POLYAXIAL
16219090	RELINER SCREW, 9.0X90MM 2FC POLYAXIAL
16219095	RELINER SCREW, 9.0X95MM 2FC POLYAXIAL
16219510	RELINER SCREW, 9.5X100MM 2FC POLYAXIAL
16219530	RELINER SCREW, 9.5X30MM 2FC POLYAXIAL
16219535	RELINER SCREW, 9.5X35MM 2FC POLYAXIAL
16219540	RELINER SCREW, 9.5X40MM 2FC POLYAXIAL
16219545	RELINER SCREW, 9.5X45MM 2FC POLYAXIAL
16219550	RELINER SCREW, 9.5X50MM 2FC POLYAXIAL
16219555	RELINER SCREW, 9.5X55MM 2FC POLYAXIAL
16219560	RELINER SCREW, 9.5X60MM 2FC POLYAXIAL
16219565	RELINER SCREW, 9.5X65MM 2FC POLYAXIAL
16219570	RELINER SCREW, 9.5X70MM 2FC POLYAXIAL
16219575	RELINER SCREW, 9.5X75MM 2FC POLYAXIAL
16219580	RELINER SCREW, 9.5X80MM 2FC POLYAXIAL
16219585	RELINER SCREW, 9.5X85MM 2FC POLYAXIAL
16219590	RELINER SCREW, 9.5X90MM 2FC POLYAXIAL
16219595	RELINER SCREW, 9.5X95MM 2FC POLYAXIAL
16255500	GSB COCR ROD, 5.5X500MM TL
16355400	RELINER-O TI ROD 5.5X400MM TL
16355500	RELINER TI ROD, 5.5X500MM TL
16360500	RELINER TI ROD, 6.0X500MM TL
16455400	RELINER-O COCR ROD 5.5X400MM TL
16455500	RELINER-O COCR ROD, 5.5X500MM TL
16460400	RELINER-O COCR ROD 6.0X400MM TL
16460500	RELINER-O COCR ROD, 6.0X500MM TL
17155001	RELINER-O TI ROD, 5.5MM LATERAL OFF SMALL
17160001	RELINER-O TI ROD, 6.0MM LATERAL OFF SMALL
17255001	RELINER-O COCR ROD, 5.5MM LATERAL OFF SM
17260001	RELINER-O COCR ROD, 6.0MM LATERAL OFF SM
17355300	RELINER-O TI ROD, 5.5X300MM J-SHAPE
17355400	RELINER-O TI ROD, 5.5X400MM J-SHAPE
17355500	RELINER-O TI ROD, 5.5X500MM J-SHAPE
17355600	RELINER-O TI ROD, 5.5X600MM J-SHAPE
17360300	RELINER-O TI ROD, 6.0X300MM J-SHAPE



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
17360400	RELIN-O TI ROD, 6.0X400MM J-SHAPE
17360500	RELIN-O TI ROD, 6.0X500MM J-SHAPE
17360600	RELIN-O TI ROD, 6.0X600MM J-SHAPE
17455300	RELIN-O COCR ROD, 5.5X300MM J-SHAPE
17455400	RELIN-O COCR ROD, 5.5X400MM J-SHAPE
17455500	RELIN-O COCR ROD, 5.5X500MM J-SHAPE
17455600	RELIN-O COCR ROD, 5.5X600MM J-SHAPE
17460300	RELIN-O COCR ROD, 6.0X300MM J-SHAPE
17460400	RELIN-O COCR ROD, 6.0X400MM J-SHAPE
17460500	RELIN-O COCR ROD, 6.0X500MM J-SHAPE
17460600	RELIN-O COCR ROD, 6.0X600MM J-SHAPE
18155002	RELIN-O TI ROD, 5.5MM LATERAL OFF MED
18160002	RELIN-O TI ROD, 6.0MM LATERAL OFF MED
18255002	RELIN-O COCR ROD, 5.5MM LATERAL OFF MED
18260002	RELIN-O COCR ROD, 6.0MM LATERAL OFF MED
19155003	RELIN-O TI ROD, 5.5MM LATERAL OFF LG
19160003	RELIN-O TI ROD, 6.0MM LATERAL OFF LG
19255003	RELIN-O COCR ROD, 5.5MM LATERAL OFF LG
19260003	RELIN-O COCR ROD, 6.0MM LATERAL OFF LG
20345015	RSS TI ROD, 4.5X15MM STRAIGHT
20345020	RSS TI ROD, 4.5X20MM STRAIGHT
20345025	RSS TI ROD, 4.5X25MM STRAIGHT
20345030	RSS TI ROD, 4.5X30MM STRAIGHT
20345035	RSS TI ROD, 4.5X35MM STRAIGHT
20345040	RSS TI ROD, 4.5X40MM STRAIGHT
20345045	RSS TI ROD, 4.5X45MM STRAIGHT
20345050	RSS TI ROD, 4.5X50MM STRAIGHT
20345055	RSS TI ROD, 4.5X55MM STRAIGHT
20345060	RSS TI ROD, 4.5X60MM STRAIGHT
20345065	RSS TI ROD, 4.5X65MM STRAIGHT
20345070	RSS TI ROD, 4.5X70MM STRAIGHT
20345075	RSS TI ROD, 4.5X75MM STRAIGHT
20345080	RSS TI ROD, 4.5X80MM STRAIGHT
20345085	RSS TI ROD, 4.5X85MM STRAIGHT
20345090	RSS TI ROD, 4.5X90MM STRAIGHT
20345095	RSS TI ROD, 4.5X95MM STRAIGHT
20345100	RSS TI ROD, 4.5X100MM STRAIGHT
20345110	RSS TI ROD, 4.5X110MM STRAIGHT
20345120	RSS TI ROD, 4.5X120MM STRAIGHT



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
20345130	RSS TI ROD, 4.5X130MM STRAIGHT
20345140	RSS TI ROD, 4.5X140MM STRAIGHT
20345150	RSS TI ROD, 4.5X150MM STRAIGHT
20345160	RSS TI ROD, 4.5X160MM STRAIGHT
20345170	RSS TI ROD, 4.5X170MM STRAIGHT
20345180	RSS TI ROD, 4.5X180MM STRAIGHT
20345190	RSS TI ROD, 4.5X190MM STRAIGHT
20345200	RSS TI ROD, 4.5X200MM STRAIGHT
20345300	RSS TI ROD, 4.5X300MM STRAIGHT
20345400	RSS TI ROD, 4.5X400MM STRAIGHT
20345500	RSS TI ROD, 4.5X500MM STRAIGHT
20345600	RSS TI ROD, 4.5X600MM STRAIGHT
20445300	RSS COCR ROD, 4.5X300MM STRAIGHT
20445500	RSS COCR ROD, 4.5X500MM STRAIGHT
13016540P3	RELIN-O SCREW, 6.5X40MM 2S POLYAXIAL
13016545P3	RELIN-O SCREW, 6.5X45MM 2S POLYAXIAL
13016550P3	RELIN-O SCREW, 6.5X50MM 2S POLYAXIAL
10050320P	RELIN-O X-CONN, 20MM 5.0 ROD FIXED
10050322P	RELIN-O X-CONN, 22.5MM 5.0 ROD FIXED
10050327P	RELIN-O X-CONN, 27.5MM 5.0 ROD FIXED
10050330P	RELIN-O X-CONN, 30MM 5.0 ROD FIXED
10050332P	RELIN-O X-CONN, 32.5MM 5.0 ROD FIXED
10050335P	RELIN-O X-CONN, 35MM 5.0 ROD FIXED
10050430P	RELIN-O X-CONN, 30-35MM 5.0 LP ADJ
10050435P	RELIN-O X-CONN, 35-42MM 5.0 LP ADJ
10050440P	RELIN-O X-CONN, 40-50MM 5.0 LP ADJ
10050445P	RELIN-O X-CONN, 45-65MM 5.0 LP ADJ
10055320P	RELIN-O X-CONN, 20MM 5.5 ROD FIXED
10055322P	RELIN-O X-CONN, 22.5MM 5.5 ROD FIXED
10055325P	RELIN-O X-CONN, 25MM 5.5 ROD FIXED
10055327P	RELIN-O X-CONN, 27.5MM 5.5 ROD FIXED
10055330P	RELIN-O X-CONN, 30MM 5.5 ROD FIXED
10055332P	RELIN-O X-CONN, 32.5MM 5.5 ROD FIXED
10055335P	RELIN-O X-CONN, 35MM 5.5 ROD FIXED
10055340P	RELIN-O X-CONN, 40MM 5.5 ROD FIXED
10055345P	RELIN-O X-CONN, 45MM 5.5 ROD FIXED
10055350P	RELIN-O X-CONN, 50MM 5.5 ROD FIXED
10055355P	RELIN-O X-CONN, 55MM 5.5 ROD FIXED
10055430P	RELIN-O X-CONN, 30-35MM 5.5 LP ADJ



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
10055435P	RELIN-O X-CONN, 35-42MM 5.5 LP ADJ
10055440P	RELIN-O X-CONN, 40-50MM 5.5 LP ADJ
10055445P	RELIN-O X-CONN, 45-65MM 5.5 LP ADJ
10055455P	RELIN-O X-CONN, 55-80MM 5.5 LP ADJ
10056430P	RELIN-O X-CON, 30-35MM 5.5/6.0 LP ADJ
10056435P	RELIN-O X-CON, 35-42MM 5.5/6.0 LP ADJ
10056440P	RELIN-O X-CON, 40-50MM 5.5/6.0 LP ADJ
10056445P	RELIN-O X-CON, 45-65MM 5.5/6.0 LP ADJ
10056455P	RELIN-O X-CON, 55-80MM 5.5/6.0 LP ADJ
10060320P	RELIN-O X-CONN, 20MM 6.0 ROD FIXED
10060322P	RELIN-O X-CONN, 22.5MM 6.0 ROD FIXED
10060325P	RELIN-O X-CONN, 25MM 6.0 ROD FIXED
10060327P	RELIN-O X-CONN, 27.5MM 6.0 ROD FIXED
10060330P	RELIN-O X-CONN, 30MM 6.0 ROD FIXED
10060332P	RELIN-O X-CONN, 32.5MM 6.0 ROD FIXED
10060335P	RELIN-O X-CONN, 35MM 6.0 ROD FIXED
10060340P	RELIN-O X-CONN, 40MM 6.0 ROD FIXED
10060430P	RELIN-O X-CONN, 30-35MM 6.0 LP ADJ
10060435P	RELIN-O X-CONN, 35-42MM 6.0 LP ADJ
10060440P	RELIN-O X-CONN, 40-50MM 6.0 LP ADJ
10060445P	RELIN-O X-CONN, 45-65MM 6.0 LP ADJ
10060455P	RELIN-O X-CONN, 55-80MM 6.0 LP ADJ
10145311P	RSS CROSS CONN, 11MM FIXED 4.5 ROD
10145313P	RSS CROSS CONN, 13MM FIXED 4.5 ROD
10145315P	RSS CROSS CONN, 15MM FIXED 4.5 ROD
10145317P	RSS CROSS CONN, 17.5MM FIXED 4.5 ROD
10145320P	RSS CROSS CONN, 20MM FIXED 4.5 ROD
10145322P	RSS CROSS CONN, 22.5MM FIXED 4.5 ROD
10145325P	RSS CROSS CONN, 25MM FIXED 4.5 ROD
10145327P	RSS CROSS CONN, 27.5MM FIXED 4.5 ROD
10145330P	RSS CROSS CONN, 30MM FIXED 4.5 ROD
10145332P	RSS CROSS CONN, 32.5MM FIXED 4.5 ROD
10145335P	RSS CROSS CONN, 35MM FIXED 4.5 ROD
10145425P	RSS CROSS CONN, 25-30MM LP ADJ4.5 ROD
10145430P	RSS CROSS CONN, 30-40MM LP ADJ4.5 ROD
10145440P	RSS CROSS CONN, 40-60MM LP ADJ4.5 ROD
10150311P	RSS CROSS CONN, 11MM FIXED 5.0 ROD
10150313P	RSS CROSS CONN, 13MM FIXED 5.0 ROD
10150315P	RSS CROSS CONN, 15MM FIXED 5.0 ROD



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
10150317P	RSS CROSS CONN, 17.5MM FIXED 5.0 ROD
10150320P	RSS CROSS CONN, 20MM FIXED 5.0 ROD
10150322P	RSS CROSS CONN, 22.5MM FIXED 5.0 ROD
10150325P	RSS CROSS CONN, 25MM FIXED 5.0 ROD
10150327P	RSS CROSS CONN, 27.5MM FIXED 5.0 ROD
10150330P	RSS CROSS CONN, 30MM FIXED 5.0 ROD
10150332P	RSS CROSS CONN, 32.5MM FIXED 5.0 ROD
10150335P	RSS CROSS CONN, 35MM FIXED 5.0 ROD
10150425P	RSS CROSS CONN, 25-30MM LP ADJ 5.0 ROD
10150430P	RSS X-CONNECTOR, 30-40MM 5.0 ROD LPADJ.
10150440P	RSS X-CONNECTOR, 40-60MM 5.0 ROD LPADJ.
10250112P	RSS ILIAC SADDLE, 5.0MM STRAIGHT Ø12
10250114P	RSS ILIAC SADDLE, 5.0MM STRAIGHT Ø14
10255110P	RSS ILIAC SADDLE, 5.5MM STRAIGHT Ø10
10255112P	RSS ILIAC SADDLE, 5.5MM STRAIGHT Ø12
10255114P	RSS ILIAC SADDLE, 5.5MM STRAIGHT Ø14
10255116P	RSS ILIAC SADDLE, 5.5MM STRAIGHT Ø16
10255210P	RSS ILIAC SADDLE, 5.5MM ANT ANG Ø10
10255212P	RSS ILIAC SADDLE, 5.5MM ANT ANG Ø12
10255214P	RSS ILIAC SADDLE, 5.5MM ANT ANG Ø14
10255216P	RSS ILIAC SADDLE, 5.5MM ANT ANG Ø16
11106155P	RELIN-O CONN, 5-6/5-6MM 2H PARALLEL SM
11106160P	RELIN-O CONN, 6.0/6.35MM 2H PARALLEL SM
11106255P	RELIN-O CONN, 5-6/5-6MM 2H PARALLEL MED
11106260P	RELIN-O CONN,6.0/6.35MM 2H PARALLEL MED
11107155P	RELIN-O CONN, 5-6/5-6MM 4H PARALLEL SM
11107255P	RELIN-O CONN, 5-6/5-6MM 4H PARALLEL MED
11108055P	RELIN-O CONN, 5-6/5-6MM 2H IN-LINE
11108060P	RELIN-O CONN, 5-6/6.35MM 2H IN-LINE
11109055P	RELIN-O CONN, 5-6/5-6MM 4H IN-LINE
11109060P	RELIN-O CONN, 6.0/6.35MM 4H IN-LINE
11355025P	RELIN MAS TI ROD, 5.5X25MM LORDOTIC
11355030P	RELIN MAS TI ROD, 5.5X30MM LORDOTIC
11355035P	RELIN MAS TI ROD, 5.5X35MM LORDOTIC
11355040P	RELIN MAS TI ROD, 5.5X40MM LORDOTIC
11355045P	RELIN MAS TI ROD, 5.5X45MM LORDOTIC
11355050P	RELIN MAS TI ROD, 5.5X50MM LORDOTIC
11355055P	RELIN MAS TI ROD, 5.5X55MM LORDOTIC
11355060P	RELIN MAS TI ROD, 5.5X60MM LORDOTIC



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
11355065P	RELINE MAS TI ROD, 5.5X65MM LORDOTIC
11355070P	RELINE MAS TI ROD, 5.5X70MM LORDOTIC
11355075P	RELINE MAS TI ROD, 5.5X75MM LORDOTIC
11355080P	RELINE MAS TI ROD, 5.5X80MM LORDOTIC
11355085P	RELINE MAS TI ROD, 5.5X85MM LORDOTIC
11355090P	RELINE MAS TI ROD, 5.5X90MM LORDOTIC
11355095P	RELINE MAS TI ROD, 5.5X95MM LORDOTIC
11355100P	RELINE MAS TI ROD, 5.5X100MM LORDOTIC
11355105P	RELINE MAS TI ROD, 5.5X105MM LORDOTIC
11355110P	RELINE MAS TI ROD, 5.5X110MM LORDOTIC
11355120P	RELINE MAS TI ROD, 5.5X120MM LORDOTIC
11355130P	RELINE MAS TI ROD, 5.5X130MM LORDOTIC
11355140P	RELINE MAS TI ROD, 5.5X140MM LORDOTIC
11355150P	RELINE MAS TI ROD, 5.5X150MM LORDOTIC
11355160P	RELINE MAS TI ROD, 5.5X160MM LORDOTIC
11360035P	RELINE MAS TI ROD, 6.0X35MM LORDOTIC
11360040P	RELINE MAS TI ROD, 6.0X40MM LORDOTIC
11360045P	RELINE MAS TI ROD, 6.0X45MM LORDOTIC
11360050P	RELINE MAS TI ROD, 6.0X50MM LORDOTIC
11360055P	RELINE MAS TI ROD, 6.0X55MM LORDOTIC
11360060P	RELINE MAS TI ROD, 6.0X60MM LORDOTIC
11360065P	RELINE MAS TI ROD, 6.0X65MM LORDOTIC
11360070P	RELINE MAS TI ROD, 6.0X70MM LORDOTIC
11360075P	RELINE MAS TI ROD, 6.0X75MM LORDOTIC
11360080P	RELINE MAS TI ROD, 6.0X80MM LORDOTIC
11360085P	RELINE MAS TI ROD, 6.0X85MM LORDOTIC
11360090P	RELINE MAS TI ROD, 6.0X90MM LORDOTIC
11360100P	RELINE MAS TI ROD, 6.0X100MM LORDOTIC
11360110P	RELINE MAS TI ROD, 6.0X110MM LORDOTIC
11360120P	RELINE MAS TI ROD, 6.0X120MM LORDOTIC
11360140P	RELINE MAS TI ROD, 6.0X140MM LORDOTIC
11360160P	RELINE MAS TI ROD, 6.0X160MM LORDOTIC
11455025P	RELINE MAS COCR ROD, 5.5X25MM LORDOTIC
11455030P	RELINE MAS COCR ROD, 5.5X30MM LORDOTIC
11455035P	RELINE MAS COCR ROD, 5.5X35MM LORDOTIC
11455040P	RELINE MAS COCR ROD, 5.5X40MM LORDOTIC
11455045P	RELINE MAS COCR ROD, 5.5X45MM LORDOTIC
11455050P	RELINE MAS COCR ROD, 5.5X50MM LORDOTIC
11455055P	RELINE MAS COCR ROD, 5.5X55MM LORDOTIC



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
11455060P	RELINE MAS COCR ROD, 5.5X60MM LORDOTIC
11455065P	RELINE MAS COCR ROD, 5.5X65MM LORDOTIC
11455070P	RELINE MAS COCR ROD, 5.5X70MM LORDOTIC
11455080P	RELINE MAS COCR ROD, 5.5X80MM LORDOTIC
11455090P	RELINE MAS COCR ROD, 5.5X90MM LORDOTIC
11455100P	RELINE MAS COCR ROD, 5.5X100MM LORDOTIC
11455110P	RELINE MAS COCR ROD, 5.5X110MM LORDOTIC
11455120P	RELINE MAS COCR ROD, 5.5X120MM LORDOTIC
11455130P	RELINE MAS COCR ROD, 5.5X130MM LORDOTIC
11455140P	RELINE MAS COCR ROD, 5.5X140MM LORDOTIC
11455150P	RELINE MAS COCR ROD, 5.5X150MM LORDOTIC
11455160P	RELINE MAS COCR ROD, 5.5X160MM LORDOTIC
12100650P	RSS CONN,4.5/5.0-4.5/5.0MM ROD 2H PAR SM
12100660P	RSS CONN,4.5/5.0-5.0/6.0MM ROD 2H PAR SM
12100760P	RSS CONN,4.5/5.0-5.0/6.0MM ROD 4H PAR SM
12100850P	RSS CONN,4.5/5.0-4.5/5.0MM ROD 2H INLINE
12100950P	RSS CONN,4.5/5.0-4.5/5.0MM ROD 4H INLINE
12101450P	RSS CONNECTOR, 50MM IN-LINE EXTENDED
12101470P	RSS CONNECTOR, 70MM IN-LINE EXTENDED
12101490P	RSS CONNECTOR, 90MM IN-LINE EXTENDED
12355040P	RELINE MAS TI ROD, 5.5X40MM STRAIGHT
12355045P	RELINE MAS TI ROD, 5.5X45MM STRAIGHT
12355050P	RELINE MAS TI ROD, 5.5X50MM STRAIGHT
12355055P	RELINE MAS TI ROD, 5.5X55MM STRAIGHT
12355060P	RELINE MAS TI ROD, 5.5X60MM STRAIGHT
12355065P	RELINE MAS TI ROD, 5.5X65MM STRAIGHT
12355070P	RELINE MAS TI ROD, 5.5X70MM STRAIGHT
12355075P	RELINE MAS TI ROD, 5.5X75MM STRAIGHT
12355080P	RELINE MAS TI ROD, 5.5X80MM STRAIGHT
12355085P	RELINE MAS TI ROD, 5.5X85MM STRAIGHT
12355090P	RELINE MAS TI ROD, 5.5X90MM STRAIGHT
12355095P	RELINE MAS TI ROD, 5.5X95MM STRAIGHT
12355100P	RELINE MAS TI ROD, 5.5X100MM STRAIGHT
12355110P	RELINE MAS TI ROD, 5.5X110MM STRAIGHT
12355120P	RELINE MAS TI ROD, 5.5X120MM STRAIGHT
12355130P	RELINE MAS TI ROD, 5.5X130MM STRAIGHT
12355140P	RELINE MAS TI ROD, 5.5X140MM STRAIGHT
12355150P	RELINE MAS TI ROD, 5.5X150MM STRAIGHT
12355160P	RELINE MAS TI ROD, 5.5X160MM STRAIGHT



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
12355170P	RELIN MAS TI ROD, 5.5X170MM STRAIGHT
12355180P	RELIN MAS TI ROD, 5.5X180MM STRAIGHT
12355190P	RELIN MAS TI ROD, 5.5X190MM STRAIGHT
12355200P	RELIN MAS TI ROD, 5.5X200MM STRAIGHT
12355210P	RELIN MAS TI ROD, 5.5X210MM STRAIGHT
12355220P	RELIN MAS TI ROD, 5.5X220MM STRAIGHT
12355230P	RELIN MAS TI ROD, 5.5X230MM STRAIGHT
12355240P	RELIN MAS TI ROD, 5.5X240MM STRAIGHT
12355250P	RELIN MAS TI ROD, 5.5X250MM STRAIGHT
12355300P	RELIN MAS TI ROD, 5.5X300MM STRAIGHT
12355500P	RELIN MAS TI ROD, 5.5X500MM STRAIGHT
12360150P	RELIN MAS TI ROD, 6.0X150MM STRAIGHT
12360200P	RELIN MAS TI ROD, 6.0X200MM STRAIGHT
12360300P	RELIN MAS TI ROD, 6.0X300MM STRAIGHT
12360500P	RELIN MAS TI ROD, 6.0X500MM STRAIGHT
12455150P	RELIN MAS COCR ROD, 5.5X150MM STRAIGHT
12455200P	RELIN MAS COCR ROD, 5.5X200MM STRAIGHT
12455300P	RELIN MAS COCR ROD, 5.5X300MM STRAIGHT
12455500P	RELIN MAS COCR ROD, 5.5X500MM STRAIGHT
12460500P	RELIN MAS COCR ROD, 6.0X500MM STRAIGHT
13000000P	RELIN LOCK SCREW, CLOSED TULIP
13100000P	RSS LOCK SCREW, CLOSED TULIP
13184540P	RELIN-O MOD SCREW SHANK, 4.5X40MM 2S
13184545P	RELIN-O MOD SCREW SHANK, 4.5X45MM 2S
13185535P	RELIN-O MOD SCREW SHANK, 5.5X35MM 2S
13185540P	RELIN-O MOD SCREW SHANK, 5.5X40MM 2S
13185545P	RELIN-O MOD SCREW SHANK, 5.5X45MM 2S
13185550P	RELIN-O MOD SCREW SHANK, 5.5X50MM 2S
13185555P	RELIN-O MOD SCREW SHANK, 5.5X55MM 2S
13186535P	RELIN-O MOD SCREW SHANK, 6.5X35MM 2S
13186540P	RELIN-O MOD SCREW SHANK, 6.5X40MM 2S
13186545P	RELIN-O MOD SCREW SHANK, 6.5X45MM 2S
13186550P	RELIN-O MOD SCREW SHANK, 6.5X50MM 2S
13186555P	RELIN-O MOD SCREW SHANK, 6.5X55MM 2S
13186560P	RELIN-O MOD SCREW SHANK, 6.5X60MM 2S
13187535P	RELIN-O MOD SCREW SHANK, 7.5X35MM 2S
13187540P	RELIN-O MOD SCREW SHANK, 7.5X40MM 2S
13187545P	RELIN-O MOD SCREW SHANK, 7.5X45MM 2S
13187550P	RELIN-O MOD SCREW SHANK, 7.5X50MM 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
13187555P	RELIN-O MOD SCREW SHANK, 7.5X55MM 2S
13187560P	RELIN-O MOD SCREW SHANK, 7.5X60MM 2S
13188535P	RELIN-O MOD SCREW SHANK, 8.5X35MM 2S
13188540P	RELIN-O MOD SCREW SHANK, 8.5X40MM 2S
13188545P	RELIN-O MOD SCREW SHANK, 8.5X45MM 2S
13188550P	RELIN-O MOD SCREW SHANK, 8.5X50MM 2S
13188555P	RELIN-O MOD SCREW SHANK, 8.5X55MM 2S
13450000P	RSS LOCK SCREW, 4.5-5.0MM
13550000P	RELIN LOCK SCREW, 5.5MM OPEN TULIP
13600000P	RELIN LOCK SCREW, 6.0MM OPEN TULIP
14355025P	RELIN-O TI ROD, 5.5X25MM LORDOTIC
14355030P	RELIN-O TI ROD, 5.5X30MM LORDOTIC
14355035P	RELIN-O TI ROD, 5.5X35MM LORDOTIC
14355040P	RELIN-O TI ROD, 5.5X40MM LORDOTIC
14355045P	RELIN-O TI ROD, 5.5X45MM LORDOTIC
14355050P	RELIN-O TI ROD, 5.5X50MM LORDOTIC
14355055P	RELIN-O TI ROD, 5.5X55MM LORDOTIC
14355060P	RELIN-O TI ROD, 5.5X60MM LORDOTIC
14355065P	RELIN-O TI ROD, 5.5X65MM LORDOTIC
14355070P	RELIN-O TI ROD, 5.5X70MM LORDOTIC
14355075P	RELIN-O TI ROD, 5.5X75MM LORDOTIC
14355080P	RELIN-O TI ROD, 5.5X80MM LORDOTIC
14355085P	RELIN-O TI ROD, 5.5X85MM LORDOTIC
14355090P	RELIN-O TI ROD, 5.5X90MM LORDOTIC
14355095P	RELIN-O TI ROD, 5.5X95MM LORDOTIC
14355100P	RELIN-O TI ROD, 5.5X100MM LORDOTIC
14355110P	RELIN-O TI ROD, 5.5X110MM LORDOTIC
14355120P	RELIN-O TI ROD, 5.5X120MM LORDOTIC
14355140P	RELIN-O TI ROD, 5.5X140MM LORDOTIC
14355150P	RELIN-O TI ROD, 5.5X150MM LORDOTIC
14355160P	RELIN-O TI ROD, 5.5X160MM LORDOTIC
14355180P	RELIN-O TI ROD, 5.5X180MM LORDOTIC
14355200P	RELIN-O TI ROD, 5.5X200MM LORDOTIC
14360030P	RELIN-O TI ROD, 6.0X30MM LORDOTIC
14360035P	RELIN-O TI ROD, 6.0X35MM LORDOTIC
14360040P	RELIN-O TI ROD, 6.0X40MM LORDOTIC
14360045P	RELIN-O TI ROD, 6.0X45MM LORDOTIC
14360050P	RELIN-O TI ROD, 6.0X50MM LORDOTIC
14360055P	RELIN-O TI ROD, 6.0X55MM LORDOTIC



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14360060P	RELIN-O TI ROD, 6.0X60MM LORDOTIC
14360065P	RELIN-O TI ROD, 6.0X65MM LORDOTIC
14360070P	RELIN-O TI ROD, 6.0X70MM LORDOTIC
14360075P	RELIN-O TI ROD, 6.0X75MM LORDOTIC
14360080P	RELIN-O TI ROD, 6.0X80MM LORDOTIC
14360090P	RELIN-O TI ROD, 6.0X90MM LORDOTIC
14360100P	RELIN-O TI ROD, 6.0X100MM LORDOTIC
14360110P	RELIN-O TI ROD, 6.0X110MM LORDOTIC
14360120P	RELIN-O TI ROD, 6.0X120MM LORDOTIC
14450025P	RELIN-O COCR ROD, 5.0X25MM LORDOTIC
14450030P	RELIN-O COCR ROD, 5.0X30MM LORDOTIC
14450035P	RELIN-O COCR ROD, 5.0X35MM LORDOTIC
14450040P	RELIN-O COCR ROD, 5.0X40MM LORDOTIC
14450045P	RELIN-O COCR ROD, 5.0X45MM LORDOTIC
14450050P	RELIN-O COCR ROD, 5.0X50MM LORDOTIC
14450055P	RELIN-O COCR ROD, 5.0X55MM LORDOTIC
14450060P	RELIN-O COCR ROD, 5.0X60MM LORDOTIC
14450065P	RELIN-O COCR ROD, 5.0X65MM LORDOTIC
14450070P	RELIN-O COCR ROD, 5.0X70MM LORDOTIC
14450075P	RELIN-O COCR ROD, 5.0X75MM LORDOTIC
14450080P	RELIN-O COCR ROD, 5.0X80MM LORDOTIC
14450085P	RELIN-O COCR ROD, 5.0X85MM LORDOTIC
14450090P	RELIN-O COCR ROD, 5.0X90MM LORDOTIC
14450095P	RELIN-O COCR ROD, 5.0X95MM LORDOTIC
14450100P	RELIN-O COCR ROD, 5.0X100MM LORDOTIC
14450110P	RELIN-O COCR ROD, 5.0X110MM LORDOTIC
14455030P	RELIN-O COCR ROD, 5.5X30MM LORDOTIC
14455035P	RELIN-O COCR ROD, 5.5X35MM LORDOTIC
14455040P	RELIN-O COCR ROD, 5.5X40MM LORDOTIC
14455045P	RELIN-O COCR ROD, 5.5X45MM LORDOTIC
14455050P	RELIN-O COCR ROD, 5.5X50MM LORDOTIC
14455055P	RELIN-O COCR ROD, 5.5X55MM LORDOTIC
14455060P	RELIN-O COCR ROD, 5.5X60MM LORDOTIC
14455065P	RELIN-O COCR ROD, 5.5X65MM LORDOTIC
14455070P	RELIN-O COCR ROD, 5.5X70MM LORDOTIC
14455075P	RELIN-O COCR ROD, 5.5X75MM LORDOTIC
14455080P	RELIN-O COCR ROD, 5.5X80MM LORDOTIC
14455085P	RELIN-O COCR ROD, 5.5X85MM LORDOTIC
14455090P	RELIN-O COCR ROD, 5.5X90MM LORDOTIC



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
14455095P	RELIN-O COCR ROD, 5.5X95MM LORDOTIC
14455100P	RELIN-O COCR ROD, 5.5X100MM LORDOTIC
14455110P	RELIN-O COCR ROD, 5.5X110MM LORDOTIC
14455120P	RELIN-O COCR ROD, 5.5X120MM LORDOTIC
14455200P	RELIN-O COCR ROD, 5.5X200MM LORDOTIC
14460030P	RELIN-O COCR ROD, 6.0X30MM LORDOTIC
14460035P	RELIN-O COCR ROD, 6.0X35MM LORDOTIC
14460040P	RELIN-O COCR ROD, 6.0X40MM LORDOTIC
14460045P	RELIN-O COCR ROD, 6.0X45MM LORDOTIC
14460050P	RELIN-O COCR ROD, 6.0X50MM LORDOTIC
14460055P	RELIN-O COCR ROD, 6.0X55MM LORDOTIC
14460060P	RELIN-O COCR ROD, 6.0X60MM LORDOTIC
14460065P	RELIN-O COCR ROD, 6.0X65MM LORDOTIC
14460070P	RELIN-O COCR ROD, 6.0X70MM LORDOTIC
14460075P	RELIN-O COCR ROD, 6.0X75MM LORDOTIC
14460080P	RELIN-O COCR ROD, 6.0X80MM LORDOTIC
14460085P	RELIN-O COCR ROD, 6.0X85MM LORDOTIC
14460090P	RELIN-O COCR ROD, 6.0X90MM LORDOTIC
14460095P	RELIN-O COCR ROD, 6.0X95MM LORDOTIC
14460100P	RELIN-O COCR ROD, 6.0X100MM LORDOTIC
14460110P	RELIN-O COCR ROD, 6.0X110MM LORDOTIC
14460120P	RELIN-O COCR ROD, 6.0X120MM LORDOTIC
14460200P	RELIN-O COCR ROD, 6.0X200MM LORDOTIC
15105300P	RELIN-O TI ROD, 300MM 5.5-5.0 TAPERED
15105500P	RELIN-O TI ROD, 500MM 5.5-5.0 TAPERED
15165500P	RELIN-O TI ROD, 500MM 6.0-5.0 TAPERED
15205300P	RELIN-O COCR ROD, 300MM 5.5-5.0 TAPERED
15205500P	RELIN-O COCR ROD, 500MM 5.5-5.0 TAPERED
15265500P	RELIN-O COCR ROD, 500MM 6.0-5.0 TAPERED
15350300P	RELIN-O TI ROD, 5.0X300MM STRAIGHT
15350500P	RELIN-O TI ROD, 5.0X500MM STRAIGHT
15355035P	RELIN-O TI ROD, 5.5X35MM STRAIGHT
15355040P	RELIN-O TI ROD, 5.5X40MM STRAIGHT
15355045P	RELIN-O TI ROD, 5.5X45MM STRAIGHT
15355050P	RELIN-O TI ROD, 5.5X50MM STRAIGHT
15355060P	RELIN-O TI ROD, 5.5X60MM STRAIGHT
15355070P	RELIN-O TI ROD, 5.5X70MM STRAIGHT
15355080P	RELIN-O TI ROD, 5.5X80MM STRAIGHT
15355085P	RELIN-O TI ROD, 5.5X85MM STRAIGHT



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
15355090P	RELIN-O TI ROD, 5.5X90MM STRAIGHT
15355100P	RELIN-O TI ROD, 5.5X100MM STRAIGHT
15355110P	RELIN-O TI ROD, 5.5X110MM STRAIGHT
15355120P	RELIN-O TI ROD, 5.5X120MM STRAIGHT
15355130P	RELIN-O TI ROD, 5.5X130MM STRAIGHT
15355140P	RELIN-O TI ROD, 5.5X140MM STRAIGHT
15355150P	RELIN-O TI ROD, 5.5X150MM STRAIGHT
15355160P	RELIN-O TI ROD, 5.5X160MM STRAIGHT
15355170P	RELIN-O TI ROD, 5.5X170MM STRAIGHT
15355180P	RELIN-O TI ROD, 5.5X180MM STRAIGHT
15355190P	RELIN-O TI ROD, 5.5X190MM STRAIGHT
15355200P	RELIN-O TI ROD, 5.5X200MM STRAIGHT
15355300P	RELIN-O TI ROD, 5.5X300MM STRAIGHT
15355500P	RELIN-O TI ROD, 5.5X500MM STRAIGHT
15360035P	RELIN-O TI ROD, 6.0X35MM STRAIGHT
15360040P	RELIN-O TI ROD, 6.0X40MM STRAIGHT
15360045P	RELIN-O TI ROD, 6.0X45MM STRAIGHT
15360050P	RELIN-O TI ROD, 6.0X50MM STRAIGHT
15360060P	RELIN-O TI ROD, 6.0X60MM STRAIGHT
15360070P	RELIN-O TI ROD, 6.0X70MM STRAIGHT
15360080P	RELIN-O TI ROD, 6.0X80MM STRAIGHT
15360090P	RELIN-O TI ROD, 6.0X90MM STRAIGHT
15360100P	RELIN-O TI ROD, 6.0X100MM STRAIGHT
15360200P	RELIN-O TI ROD, 6.0X200MM STRAIGHT
15360300P	RELIN-O TI ROD, 6.0X300MM STRAIGHT
15360500P	RELIN-O TI ROD, 6.0X500MM STRAIGHT
15450200P	RELIN-O COCR ROD, 5.0X200MM STRAIGHT
15450300P	RELIN-O COCR ROD, 5.0X300MM STRAIGHT
15450500P	RELIN-O COCR ROD, 5.0X500MM STRAIGHT
15455300P	RELIN-O COCR ROD, 5.5X300MM STRAIGHT
15455500P	RELIN-O COCR ROD, 5.5X500MM STRAIGHT
15460300P	RELIN-O COCR ROD, 6.0X300MM STRAIGHT
15460500P	RELIN-O COCR ROD, 6.0X500MM STRAIGHT
16184530P	RELIN MAS MOD SCREW SHANK, 4.5X30MM 2C
16184535P	RELIN MAS MOD SCREW SHANK, 4.5X35MM 2C
16184540P	RELIN MAS MOD SCREW SHANK, 4.5X40MM 2C
16184545P	RELIN MAS MOD SCREW SHANK, 4.5X45MM 2C
16184550P	RELIN MAS MOD SCREW SHANK, 4.5X50MM 2C
16185530P	RELIN MAS MOD SCREW SHANK, 5.5X30MM 2C



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
16185535P	RELIN MAS MOD SCREW SHANK, 5.5X35MM 2C
16185540P	RELIN MAS MOD SCREW SHANK, 5.5X40MM 2C
16185545P	RELIN MAS MOD SCREW SHANK, 5.5X45MM 2C
16185550P	RELIN MAS MOD SCREW SHANK, 5.5X50MM 2C
16185555P	RELIN MAS MOD SCREW SHANK, 5.5X55MM 2C
16185560P	RELIN MAS MOD SCREW SHANK, 5.5X60MM 2C
16185565P	RELIN MAS MOD SCREW SHANK, 5.5X65MM 2C
16186530P	RELIN MAS MOD SCREW SHANK, 6.5X30MM 2C
16186535P	RELIN MAS MOD SCREW SHANK, 6.5X35MM 2C
16186540P	RELIN MAS MOD SCREW SHANK, 6.5X40MM 2C
16186545P	RELIN MAS MOD SCREW SHANK, 6.5X45MM 2C
16186550P	RELIN MAS MOD SCREW SHANK, 6.5X50MM 2C
16186555P	RELIN MAS MOD SCREW SHANK, 6.5X55MM 2C
16186560P	RELIN MAS MOD SCREW SHANK, 6.5X60MM 2C
16186565P	RELIN MAS MOD SCREW SHANK, 6.5X65MM 2C
16187530P	RELIN MAS MOD SCREW SHANK, 7.5X30MM 2C
16187535P	RELIN MAS MOD SCREW SHANK, 7.5X35MM 2C
16187540P	RELIN MAS MOD SCREW SHANK, 7.5X40MM 2C
16187545P	RELIN MAS MOD SCREW SHANK, 7.5X45MM 2C
16187550P	RELIN MAS MOD SCREW SHANK, 7.5X50MM 2C
16187555P	RELIN MAS MOD SCREW SHANK, 7.5X55MM 2C
16187560P	RELIN MAS MOD SCREW SHANK, 7.5X60MM 2C
16187565P	RELIN MAS MOD SCREW SHANK, 7.5X65MM 2C
16187580P	RELIN MAS MOD SCREW SHANK, 7.5X80MM 2C
16188535P	RELIN MAS MOD SCREW SHANK, 8.5X35MM 2C
16188540P	RELIN MAS MOD SCREW SHANK, 8.5X40MM 2C
16188545P	RELIN MAS MOD SCREW SHANK, 8.5X45MM 2C
16188550P	RELIN MAS MOD SCREW SHANK, 8.5X50MM 2C
16188555P	RELIN MAS MOD SCREW SHANK, 8.5X55MM 2C
16188560P	RELIN MAS MOD SCREW SHANK, 8.5X60MM 2C
16188570P	RELIN MAS MOD SCREW SHANK, 8.5X70MM 2C
16188580P	RELIN MAS MOD SCREW SHANK, 8.5X80MM 2C
16188590P	RELIN MAS MOD SCREW SHANK, 8.5X90MM 2C
16189540P	RELIN MAS MOD SCREW SHANK, 9.5X40MM 2C
16189545P	RELIN MAS MOD SCREW SHANK, 9.5X45MM 2C
16189550P	RELIN MAS MOD SCREW SHANK, 9.5X50MM 2C
16355400P	RELIN-O TI ROD 5.5X400MM TL
16355500P	RELIN TI ROD, 5.5X500MM TL
16360400P	RELIN-O TI ROD 6.0X400MM TL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
16360500P	RELIN TI ROD, 6.0X500MM TL
16455400P	RELIN-O COCR ROD 5.5X400MM TL
16455500P	RELIN-O COCR ROD, 5.5X500MM TL
16460400P	RELIN-O COCR ROD 6.0X400MM TL
16460500P	RELIN-O COCR ROD, 6.0X500MM TL
18155002P	RELIN-O TI ROD, 5.5MM LATERAL OFF MED
18160002P	RELIN-O TI ROD, 6.0MM LATERAL OFF MED
18255002P	RELIN-O COCR ROD, 5.5MM LATERAL OFF MED
18260002P	RELIN-O COCR ROD, 6.0MM LATERAL OFF MED
20345300P	RSS TI ROD, 4.5X300MM STRAIGHT
20345500P	RSS TI ROD, 4.5X500MM STRAIGHT
20445300P	RSS COCR ROD, 4.5X300MM STRAIGHT
20445500P	RSS COCR ROD, 4.5X500MM STRAIGHT
13104025P	RELIN-O SCREW+LS, 4.0X25MM TRAC POLY
13104030P	RELIN-O SCREW+LS, 4.0X30MM TRAC POLY
13104035P	RELIN-O SCREW+LS, 4.0X35MM TRAC POLY
13104040P	RELIN-O SCREW+LS, 4.0X40MM TRAC POLY
13104045P	RELIN-O SCREW+LS, 4.0X45MM TRAC POLY
13104525P	RELIN-O SCREW+LS, 4.5X25MM TRAC POLY
13104530P	RELIN-O SCREW+LS, 4.5X30MM TRAC POLY
13104535P	RELIN-O SCREW+LS, 4.5X35MM TRAC POLY
13104540P	RELIN-O SCREW+LS, 4.5X40MM TRAC POLY
13104545P	RELIN-O SCREW+LS, 4.5X45MM TRAC POLY
13104550P	RELIN-O SCREW+LS, 4.5X50MM TRAC POLY
13104555P	RELIN-O SCREW+LS, 4.5X55MM TRAC POLY
13105025P	RELIN-O SCREW+LS, 5.0X25MM TRAC POLY
13105030P	RELIN-O SCREW+LS, 5.0X30MM TRAC POLY
13105035P	RELIN-O SCREW+LS, 5.0X35MM TRAC POLY
13105040P	RELIN-O SCREW+LS, 5.0X40MM TRAC POLY
13105045P	RELIN-O SCREW+LS, 5.0X45MM TRAC POLY
13105050P	RELIN-O SCREW+LS, 5.0X50MM TRAC POLY
13105055P	RELIN-O SCREW+LS, 5.0X55MM TRAC POLY
13105525P	RELIN-O SCREW+LS, 5.5X25MM TRAC POLY
13105530P	RELIN-O SCREW+LS, 5.5X30MM TRAC POLY
13105535P	RELIN-O SCREW+LS, 5.5X35MM TRAC POLY
13105540P	RELIN-O SCREW+LS, 5.5X40MM TRAC POLY
13105545P	RELIN-O SCREW+LS, 5.5X45MM TRAC POLY
13105550P	RELIN-O SCREW+LS, 5.5X50MM TRAC POLY
13105555P	RELIN-O SCREW+LS, 5.5X55MM TRAC POLY



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13105560P	RELIN-O SCREW+LS, 5.5X60MM TRAC POLY
13106530P	RELIN-O SCREW+LS, 6.5X30MM TRAC POLY
13106535P	RELIN-O SCREW+LS, 6.5X35MM TRAC POLY
13106540P	RELIN-O SCREW+LS, 6.5X40MM TRAC POLY
13106545P	RELIN-O SCREW+LS, 6.5X45MM TRAC POLY
13106550P	RELIN-O SCREW+LS, 6.5X50MM TRAC POLY
13106555P	RELIN-O SCREW+LS, 6.5X55MM TRAC POLY
13106560P	RELIN-O SCREW+LS, 6.5X60MM TRAC POLY
13107530P	RELIN-O SCREW+LS, 7.5X30MM TRAC POLY
13107535P	RELIN-O SCREW+LS, 7.5X35MM TRAC POLY
13107540P	RELIN-O SCREW+LS, 7.5X40MM TRAC POLY
13107545P	RELIN-O SCREW+LS, 7.5X45MM TRAC POLY
13107550P	RELIN-O SCREW+LS, 7.5X50MM TRAC POLY
13107555P	RELIN-O SCREW+LS, 7.5X55MM TRAC POLY
13107560P	RELIN-O SCREW+LS, 7.5X60MM TRAC POLY
13108530P	RELIN-O SCREW+LS, 8.5X30MM TRAC POLY
13108535P	RELIN-O SCREW+LS, 8.5X35MM TRAC POLY
13108540P	RELIN-O SCREW+LS, 8.5X40MM TRAC POLY
13108545P	RELIN-O SCREW+LS, 8.5X45MM TRAC POLY
13108550P	RELIN-O SCREW+LS, 8.5X50MM TRAC POLY
13108555P	RELIN-O SCREW+LS, 8.5X55MM TRAC POLY
13109540P	RELIN-O SCREW+LS, 9.5X40MM TRAC POLY
13109545P	RELIN-O SCREW+LS, 9.5X45MM TRAC POLY
13109550P	RELIN-O SCREW+LS, 9.5X50MM TRAC POLY
13191110P	RELIN-O MOD+LS, REDUCTION EXT
13225535P	RELIN-O SCREW+LS, 5.5X35MM 2 FENPOLY
13225540P	RELIN-O SCREW+LS, 5.5X40MM 2 FENPOLY
13225545P	RELIN-O SCREW+LS, 5.5X45MM 2 FENPOLY
13225550P	RELIN-O SCREW+LS, 5.5X50MM 2 FENPOLY
13225555P	RELIN-O SCREW+LS, 5.5X55MM 2 FENPOLY
13225560P	RELIN-O SCREW+LS, 5.5X60MM 2 FENPOLY
13226535P	RELIN-O SCREW+LS, 6.5X35MM 2 FENPOLY
13226540P	RELIN-O SCREW+LS, 6.5X40MM 2 FENPOLY
13226545P	RELIN-O SCREW+LS, 6.5X45MM 2 FENPOLY
13226550P	RELIN-O SCREW+LS, 6.5X50MM 2 FENPOLY
13226555P	RELIN-O SCREW+LS, 6.5X55MM 2 FENPOLY
13226560P	RELIN-O SCREW+LS, 6.5X60MM 2 FENPOLY
13227535P	RELIN-O SCREW+LS, 7.5X35MM 2 FENPOLY
13227540P	RELIN-O SCREW+LS, 7.5X40MM 2 FENPOLY



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13227545P	RELIN-O SCREW+LS, 7.5X45MM 2 FENPOLY
13227550P	RELIN-O SCREW+LS, 7.5X50MM 2 FENPOLY
13227555P	RELIN-O SCREW+LS, 7.5X55MM 2 FENPOLY
13227560P	RELIN-O SCREW+LS, 7.5X60MM 2 FENPOLY
13228535P	RELIN-O SCREW+LS, 8.5X35MM 2 FENPOLY
13228540P	RELIN-O SCREW+LS, 8.5X40MM 2 FENPOLY
13228545P	RELIN-O SCREW+LS, 8.5X45MM 2 FENPOLY
13228550P	RELIN-O SCREW+LS, 8.5X50MM 2 FENPOLY
13228555P	RELIN-O SCREW+LS, 8.5X55MM 2 FENPOLY
13228560P	RELIN-O SCREW+LS, 8.5X60MM 2 FENPOLY
13228580P	RELIN-O SCREW+LS, 8.5X80MM 2 FENPOLY
13310510P	RELIN-O SCREW+LS, 10.5X100MM 2S P.I.
13310511P	RELIN-O SCREW+LS, 10.5X110MM 2S P.I.
13310512P	RELIN-O SCREW+LS, 10.5X120MM 2S P.I.
13310535P	RELIN-O SCREW+LS, 10.5X35MM 2S POLY
13310540P	RELIN-O SCREW+LS, 10.5X40MM 2S POLY
13310545P	RELIN-O SCREW+LS, 10.5X45MM 2S POLY
13310550P	RELIN-O SCREW+LS, 10.5X50MM 2S POLY
13310555P	RELIN-O SCREW+LS, 10.5X55MM 2S POLY
13310560P	RELIN-O SCREW+LS, 10.5X60MM 2S P.I.
13310570P	RELIN-O SCREW+LS, 10.5X70MM 2S P.I.
13310580P	RELIN-O SCREW+LS, 10.5X80MM 2S P.I.
13310590P	RELIN-O SCREW+LS, 10.5X90MM 2S P.I.
13311510P	RELIN-O SCREW+LS, 11.5X100MM 2S P.I.
13311580P	RELIN-O SCREW+LS, 11.5X80MM 2S P.I.
13311590P	RELIN-O SCREW+LS, 11.5X90MM 2S P.I.
13314020P	RELIN-O SCREW+LS, 4.0X20MM 2S POLY
13314025P	RELIN-O SCREW+LS, 4.0X25MM 2S POLY
13314030P	RELIN-O SCREW+LS, 4.0X30MM 2S POLY
13314035P	RELIN-O SCREW+LS, 4.0X35MM 2S POLY
13314040P	RELIN-O SCREW+LS, 4.0X40MM 2S POLY
13314045P	RELIN-O SCREW+LS, 4.0X45MM 2S POLY
13314520P	RELIN-O SCREW+LS, 4.5X20MM 2S POLY
13314525P	RELIN-O SCREW+LS, 4.5X25MM 2S POLY
13314530P	RELIN-O SCREW+LS, 4.5X30MM 2S POLY
13314535P	RELIN-O SCREW+LS, 4.5X35MM 2S POLY
13314540P	RELIN-O SCREW+LS, 4.5X40MM 2S POLY
13314545P	RELIN-O SCREW+LS, 4.5X45MM 2S POLY
13314550P	RELIN-O SCREW+LS, 4.5X50MM 2S POLY



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13314555P	RELIN-O SCREW+LS, 4.5X55MM 2S POLY
13315020P	RELIN-O SCREW+LS, 5.0X20MM 2S POLY
13315025P	RELIN-O SCREW+LS, 5.0X25MM 2S POLY
13315030P	RELIN-O SCREW+LS, 5.0X30MM 2S POLY
13315035P	RELIN-O SCREW+LS, 5.0X35MM 2S POLY
13315040P	RELIN-O SCREW+LS, 5.0X40MM 2S POLY
13315045P	RELIN-O SCREW+LS, 5.0X45MM 2S POLY
13315050P	RELIN-O SCREW+LS, 5.0X50MM 2S POLY
13315055P	RELIN-O SCREW+LS, 5.0X55MM 2S POLY
13315520P	RELIN-O SCREW+LS, 5.5X20MM 2S POLY
13315525P	RELIN-O SCREW+LS, 5.5X25MM 2S POLY
13315530P	RELIN-O SCREW+LS, 5.5X30MM 2S POLY
13315535P	RELIN-O SCREW+LS, 5.5X35MM 2S POLY
13315540P	RELIN-O SCREW+LS, 5.5X40MM 2S POLY
13315545P	RELIN-O SCREW+LS, 5.5X45MM 2S POLY
13315550P	RELIN-O SCREW+LS, 5.5X50MM 2S POLY
13315555P	RELIN-O SCREW+LS, 5.5X55MM 2S POLY
13315560P	RELIN-O SCREW+LS, 5.5X60MM 2S P.I.
13316030P	RELIN-O SCREW+LS, 6.0X30MM 2S POLY
13316035P	RELIN-O SCREW+LS, 6.0X35MM 2S POLY
13316040P	RELIN-O SCREW+LS, 6.0X40MM 2S POLY
13316045P	RELIN-O SCREW+LS, 6.0X45MM 2S POLY
13316050P	RELIN-O SCREW+LS, 6.0X50MM 2S POLY
13316055P	RELIN-O SCREW+LS, 6.0X55MM 2S POLY
13316060P	RELIN-O SCREW+LS, 6.0X60MM 2S P.I.
13316525P	RELIN-O SCREW+LS, 6.5X25MM 2S POLY
13316530P	RELIN-O SCREW+LS, 6.5X30MM 2S POLY
13316535P	RELIN-O SCREW+LS, 6.5X35MM 2S POLY
13316540P	RELIN-O SCREW+LS, 6.5X40MM 2S POLY
13316545P	RELIN-O SCREW+LS, 6.5X45MM 2S POLY
13316550P	RELIN-O SCREW+LS, 6.5X50MM 2S POLY
13316555P	RELIN-O SCREW+LS, 6.5X55MM 2S POLY
13316560P	RELIN-O SCREW+LS, 6.5X60MM 2S P.I.
13316565P	RELIN-O SCREW+LS, 6.5X65MM 2S P.I.
13316570P	RELIN-O SCREW+LS, 6.5X70MM 2S P.I.
13316580P	RELIN-O SCREW+LS, 6.5X80MM 2S P.I.
13316590P	RELIN-O SCREW+LS, 6.5X90MM 2S P.I.
13317030P	RELIN-O SCREW+LS, 7.0X30MM 2S POLY
13317035P	RELIN-O SCREW+LS, 7.0X35MM 2S POLY



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13317040P	RELIN-O SCREW+LS, 7.0X40MM 2S POLY
13317045P	RELIN-O SCREW+LS, 7.0X45MM 2S POLY
13317050P	RELIN-O SCREW+LS, 7.0X50MM 2S POLY
13317055P	RELIN-O SCREW+LS, 7.0X55MM 2S POLY
13317060P	RELIN-O SCREW+LS, 7.0X60MM 2S P.I.
13317510P	RELIN-O SCREW+LS, 7.5X100MM 2S P.I.
13317525P	RELIN-O SCREW+LS, 7.5X25MM 2S POLY
13317530P	RELIN-O SCREW+LS, 7.5X30MM 2S POLY
13317535P	RELIN-O SCREW+LS, 7.5X35MM 2S POLY
13317540P	RELIN-O SCREW+LS, 7.5X40MM 2S POLY
13317545P	RELIN-O SCREW+LS, 7.5X45MM 2S POLY
13317550P	RELIN-O SCREW+LS, 7.5X50MM 2S POLY
13317555P	RELIN-O SCREW+LS, 7.5X55MM 2S POLY
13317560P	RELIN-O SCREW+LS, 7.5X60MM 2S P.I.
13317565P	RELIN-O SCREW+LS, 7.5X65MM 2S P.I.
13317570P	RELIN-O SCREW+LS, 7.5X70MM 2S P.I.
13317580P	RELIN-O SCREW+LS, 7.5X80MM 2S P.I.
13317590P	RELIN-O SCREW+LS, 7.5X90MM 2S P.I.
13318510P	RELIN-O SCREW+LS, 8.5X100MM 2S P.I.
13318511P	RELIN-O SCREW+LS, 8.5X110MM 2S P.I.
13318512P	RELIN-O SCREW+LS, 8.5X120MM 2S P.I.
13318530P	RELIN-O SCREW+LS, 8.5X30MM 2S POLY
13318535P	RELIN-O SCREW+LS, 8.5X35MM 2S POLY
13318540P	RELIN-O SCREW+LS, 8.5X40MM 2S POLY
13318545P	RELIN-O SCREW+LS, 8.5X45MM 2S POLY
13318550P	RELIN-O SCREW+LS, 8.5X50MM 2S POLY
13318555P	RELIN-O SCREW+LS, 8.5X55MM 2S POLY
13318560P	RELIN-O SCREW+LS, 8.5X60MM 2S P.I.
13318565P	RELIN-O SCREW+LS, 8.5X65MM 2S P.I.
13318570P	RELIN-O SCREW+LS, 8.5X70MM 2S P.I.
13318580P	RELIN-O SCREW+LS, 8.5X80MM 2S P.I.
13318590P	RELIN-O SCREW+LS, 8.5X90MM 2S P.I.
13319510P	RELIN-O SCREW+LS, 9.5X100MM 2S P.I.
13319511P	RELIN-O SCREW+LS, 9.5X110MM 2S P.I.
13319512P	RELIN-O SCREW+LS, 9.5X120MM 2S P.I.
13319530P	RELIN-O SCREW+LS, 9.5X30MM 2S POLY
13319535P	RELIN-O SCREW+LS, 9.5X35MM 2S POLY
13319540P	RELIN-O SCREW+LS, 9.5X40MM 2S POLY
13319545P	RELIN-O SCREW+LS, 9.5X45MM 2S POLY



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13319550P	RELIN-O SCREW+LS, 9.5X50MM 2S POLY
13319555P	RELIN-O SCREW+LS, 9.5X55MM 2S POLY
13319560P	RELIN-O SCREW+LS, 9.5X60MM 2S P.I.
13319565P	RELIN-O SCREW+LS, 9.5X65MM 2S P.I.
13319570P	RELIN-O SCREW+LS, 9.5X70MM 2S P.I.
13319580P	RELIN-O SCREW+LS, 9.5X80MM 2S P.I.
13319590P	RELIN-O SCREW+LS, 9.5X90MM 2S P.I.
13324030P	RELIN-O SCREW+LS, 4.0X30MM 2S RED
13324035P	RELIN-O SCREW+LS, 4.0X35MM 2S RED
13324040P	RELIN-O SCREW+LS, 4.0X40MM 2S RED
13324525P	RELIN-O SCREW+LS, 4.5X25MM 2S RED
13324530P	RELIN-O SCREW+LS, 4.5X30MM 2S RED
13324535P	RELIN-O SCREW+LS, 4.5X35MM 2S RED
13324540P	RELIN-O SCREW+LS, 4.5X40MM 2S RED
13324545P	RELIN-O SCREW+LS, 4.5X45MM 2S RED
13325030P	RELIN-O SCREW+LS, 5.0X30MM 2S RED
13325035P	RELIN-O SCREW+LS, 5.0X35MM 2S RED
13325040P	RELIN-O SCREW+LS, 5.0X40MM 2S RED
13325525P	RELIN-O SCREW+LS, 5.5X25MM 2S RED
13325530P	RELIN-O SCREW+LS, 5.5X30MM 2S RED
13325535P	RELIN-O SCREW+LS, 5.5X35MM 2S RED
13325540P	RELIN-O SCREW+LS, 5.5X40MM 2S RED
13325545P	RELIN-O SCREW+LS, 5.5X45MM 2S RED
13325550P	RELIN-O SCREW+LS, 5.5X50MM 2S RED
13325555P	RELIN-O SCREW+LS, 5.5X55MM 2S RED
13326530P	RELIN-O SCREW+LS, 6.5X30MM 2S RED
13326535P	RELIN-O SCREW+LS, 6.5X35MM 2S RED
13326540P	RELIN-O SCREW+LS, 6.5X40MM 2S RED
13326545P	RELIN-O SCREW+LS, 6.5X45MM 2S RED
13326550P	RELIN-O SCREW+LS, 6.5X50MM 2S RED
13326555P	RELIN-O SCREW+LS, 6.5X55MM 2S RED
13327535P	RELIN-O SCREW+LS, 7.5X35MM 2S RED
13327540P	RELIN-O SCREW+LS, 7.5X40MM 2S RED
13327545P	RELIN-O SCREW+LS, 7.5X45MM 2S RED
13327550P	RELIN-O SCREW+LS, 7.5X50MM 2S RED
13327555P	RELIN-O SCREW+LS, 7.5X55MM 2S RED
13327560P	RELIN-O SCREW+LS, 7.5X60MM 2S RED
13328535P	RELIN-O SCREW+LS, 8.5X35MM 2S RED
13328540P	RELIN-O SCREW+LS, 8.5X40MM 2S RED



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
13328545P	RELIN-O SCREW+LS, 8.5X45MM 2S RED
13328550P	RELIN-O SCREW+LS, 8.5X50MM 2S RED
13328555P	RELIN-O SCREW+LS, 8.5X55MM 2S RED
13337510P	RELIN-O SCREW+LS, 7.5X100MM 2S I.C.
13337560P	RELIN-O SCREW+LS, 7.5X60MM 2S I.C.
13337570P	RELIN-O SCREW+LS, 7.5X70MM 2S I.C.
13337580P	RELIN-O SCREW+LS, 7.5X80MM 2S I.C.
13337590P	RELIN-O SCREW+LS, 7.5X90MM 2S I.C.
13338510P	RELIN-O SCREW+LS, 8.5X100MM 2S I.C.
13338560P	RELIN-O SCREW+LS, 8.5X60MM 2S I.C.
13338570P	RELIN-O SCREW+LS, 8.5X70MM 2S I.C.
13338580P	RELIN-O SCREW+LS, 8.5X80MM 2S I.C.
13338590P	RELIN-O SCREW+LS, 8.5X90MM 2S I.C.
13339510P	RELIN-O SCREW+LS, 9.5X100MM 2S I.C.
13339560P	RELIN-O SCREW+LS, 9.5X60MM 2S I.C.
13339570P	RELIN-O SCREW+LS, 9.5X70MM 2S I.C.
13339580P	RELIN-O SCREW+LS, 9.5X80MM 2S I.C.
13339590P	RELIN-O SCREW+LS, 9.5X90MM 2S I.C.
13344025P	RELIN-O SCREW+LS, 4.0X25MM 2S UNI
13344030P	RELIN-O SCREW+LS, 4.0X30MM 2S UNI
13344035P	RELIN-O SCREW+LS, 4.0X35MM 2S UNI
13344040P	RELIN-O SCREW+LS, 4.0X40MM 2S UNI
13344520P	RELIN-O SCREW+LS, 4.5X20MM 2S UNI
13344525P	RELIN-O SCREW+LS, 4.5X25MM 2S UNI
13344530P	RELIN-O SCREW+LS, 4.5X30MM 2S UNI
13344535P	RELIN-O SCREW+LS, 4.5X35MM 2S UNI
13344540P	RELIN-O SCREW+LS, 4.5X40MM 2S UNI
13344545P	RELIN-O SCREW+LS, 4.5X45MM 2S UNI
13345025P	RELIN-O SCREW+LS, 5.0X25MM 2S UNI
13345030P	RELIN-O SCREW+LS, 5.0X30MM 2S UNI
13345035P	RELIN-O SCREW+LS, 5.0X35MM 2S UNI
13345040P	RELIN-O SCREW+LS, 5.0X40MM 2S UNI
13345045P	RELIN-O SCREW+LS, 5.0X45MM 2S UNI
13345525P	RELIN-O SCREW+LS, 5.5X25MM 2S UNI
13345530P	RELIN-O SCREW+LS, 5.5X30MM 2S UNI
13345535P	RELIN-O SCREW+LS, 5.5X35MM 2S UNI
13345540P	RELIN-O SCREW+LS, 5.5X40MM 2S UNI
13345545P	RELIN-O SCREW+LS, 5.5X45MM 2S UNI
13345550P	RELIN-O SCREW+LS, 5.5X50MM 2S UNI



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
13345555P	RELIN-O SCREW+LS, 5.5X55MM 2S UNI
13346030P	RELIN-O SCREW+LS, 6.0X30MM 2S UNI
13346035P	RELIN-O SCREW+LS, 6.0X35MM 2S UNI
13346040P	RELIN-O SCREW+LS, 6.0X40MM 2S UNI
13346045P	RELIN-O SCREW+LS, 6.0X45MM 2S UNI
13346525P	RELIN-O SCREW+LS, 6.5X25MM 2S UNI
13346530P	RELIN-O SCREW+LS, 6.5X30MM 2S UNI
13346535P	RELIN-O SCREW+LS, 6.5X35MM 2S UNI
13346540P	RELIN-O SCREW+LS, 6.5X40MM 2S UNI
13346545P	RELIN-O SCREW+LS, 6.5X45MM 2S UNI
13346550P	RELIN-O SCREW+LS, 6.5X50MM 2S UNI
13346555P	RELIN-O SCREW+LS, 6.5X55MM 2S UNI
13347035P	RELIN-O SCREW+LS, 7.0X35MM 2S UNI
13347040P	RELIN-O SCREW+LS, 7.0X40MM 2S UNI
13347045P	RELIN-O SCREW+LS, 7.0X45MM 2S UNI
13347530P	RELIN-O SCREW+LS, 7.5X30MM 2S UNI
13347535P	RELIN-O SCREW+LS, 7.5X35MM 2S UNI
13347540P	RELIN-O SCREW+LS, 7.5X40MM 2S UNI
13347545P	RELIN-O SCREW+LS, 7.5X45MM 2S UNI
13347550P	RELIN-O SCREW+LS, 7.5X50MM 2S UNI
13347555P	RELIN-O SCREW+LS, 7.5X55MM 2S UNI
13354525P	RELIN-O SCREW+LS, 4.5X25MM 2S REDUNI
13354530P	RELIN-O SCREW+LS, 4.5X30MM 2S REDUNI
13354535P	RELIN-O SCREW+LS, 4.5X35MM 2S REDUNI
13354540P	RELIN-O SCREW+LS, 4.5X40MM 2S REDUNI
13354545P	RELIN-O SCREW+LS, 4.5X45MM 2S REDUNI
13355025P	RELIN-O SCREW+LS, 5.0X25MM 2S REDUNI
13355030P	RELIN-O SCREW+LS, 5.0X30MM 2S REDUNI
13355035P	RELIN-O SCREW+LS, 5.0X35MM 2S REDUNI
13355040P	RELIN-O SCREW+LS, 5.0X40MM 2S REDUNI
13355045P	RELIN-O SCREW+LS, 5.0X45MM 2S REDUNI
13355525P	RELIN-O SCREW+LS, 5.5X25MM 2S REDUNI
13355530P	RELIN-O SCREW+LS, 5.5X30MM 2S REDUNI
13355535P	RELIN-O SCREW+LS, 5.5X35MM 2S REDUNI
13355540P	RELIN-O SCREW+LS, 5.5X40MM 2S REDUNI
13355545P	RELIN-O SCREW+LS, 5.5X45MM 2S REDUNI
13355550P	RELIN-O SCREW+LS, 5.5X50MM 2S REDUNI
13355555P	RELIN-O SCREW+LS, 5.5X55MM 2S REDUNI
13356530P	RELIN-O SCREW+LS, 6.5X30MM 2S REDUNI



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
13356535P	RELIN-O SCREW+LS, 6.5X35MM 2S REDUNI
13356540P	RELIN-O SCREW+LS, 6.5X40MM 2S REDUNI
13356545P	RELIN-O SCREW+LS, 6.5X45MM 2S REDUNI
13356550P	RELIN-O SCREW+LS, 6.5X50MM 2S REDUNI
13356555P	RELIN-O SCREW+LS, 6.5X55MM 2S REDUNI
13364530P	RELIN-O SCREW+LS, 4.5X30MM 2S PLS
13364535P	RELIN-O SCREW+LS, 4.5X35MM 2S PLS
13364540P	RELIN-O SCREW+LS, 4.5X40MM 2S PLS
13364545P	RELIN-O SCREW+LS, 4.5X45MM 2S PLS
13365530P	RELIN-O SCREW+LS, 5.5X30MM 2S PLS
13365535P	RELIN-O SCREW+LS, 5.5X35MM 2S PLS
13365540P	RELIN-O SCREW+LS, 5.5X40MM 2S PLS
13365545P	RELIN-O SCREW+LS, 5.5X45MM 2S PLS
13365550P	RELIN-O SCREW+LS, 5.5X50MM 2S PLS
13366530P	RELIN-O SCREW+LS, 6.5X30MM 2S PLS
13366535P	RELIN-O SCREW+LS, 6.5X35MM 2S PLS
13366540P	RELIN-O SCREW+LS, 6.5X40MM 2S PLS
13366545P	RELIN-O SCREW+LS, 6.5X45MM 2S PLS
13366550P	RELIN-O SCREW+LS, 6.5X50MM 2S PLS
13367535P	RELIN-O SCREW+LS, 7.5X35MM 2S PLS
13367540P	RELIN-O SCREW+LS, 7.5X40MM 2S PLS
13367545P	RELIN-O SCREW+LS, 7.5X45MM 2S PLS
13367550P	RELIN-O SCREW+LS, 7.5X50MM 2S PLS
13384025P	RELIN-O SCREW+LS, 4.0X25MM 2S MONO
13384030P	RELIN-O SCREW+LS, 4.0X30MM 2S MONO
13384035P	RELIN-O SCREW+LS, 4.0X35MM 2S MONO
13384040P	RELIN-O SCREW+LS, 4.0X40MM 2S MONO
13384530P	RELIN-O SCREW+LS, 4.5X30MM 2S MONO
13384535P	RELIN-O SCREW+LS, 4.5X35MM 2S MONO
13384540P	RELIN-O SCREW+LS, 4.5X40MM 2S MONO
13384545P	RELIN-O SCREW+LS, 4.5X45MM 2S MONO
13385030P	RELIN-O SCREW+LS, 5.0X30MM 2S MONO
13385035P	RELIN-O SCREW+LS, 5.0X35MM 2S MONO
13385040P	RELIN-O SCREW+LS, 5.0X40MM 2S MONO
13385530P	RELIN-O SCREW+LS, 5.5X30MM 2S MONO
13385535P	RELIN-O SCREW+LS, 5.5X35MM 2S MONO
13385540P	RELIN-O SCREW+LS, 5.5X40MM 2S MONO
13385545P	RELIN-O SCREW+LS, 5.5X45MM 2S MONO
13385550P	RELIN-O SCREW+LS, 5.5X50MM 2S MONO



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
13385555P	RELIN-O SCREW+LS, 5.5X55MM 2S MONO
13386530P	RELIN-O SCREW+LS, 6.5X30MM 2S MONO
13386535P	RELIN-O SCREW+LS, 6.5X35MM 2S MONO
13386540P	RELIN-O SCREW+LS, 6.5X40MM 2S MONO
13386545P	RELIN-O SCREW+LS, 6.5X45MM 2S MONO
13386550P	RELIN-O SCREW+LS, 6.5X50MM 2S MONO
13386555P	RELIN-O SCREW+LS, 6.5X55MM 2S MONO
13387530P	RELIN-O SCREW+LS, 7.5X30MM 2S MONO
13387535P	RELIN-O SCREW+LS, 7.5X35MM 2S MONO
13387540P	RELIN-O SCREW+LS, 7.5X40MM 2S MONO
13387545P	RELIN-O SCREW+LS, 7.5X45MM 2S MONO
13387550P	RELIN-O SCREW+LS, 7.5X50MM 2S MONO
13387555P	RELIN-O SCREW+LS, 7.5X55MM 2S MONO
13388540P	RELIN-O SCREW+LS, 8.5X40MM 2S MONO
13450001P	RSS LOCK SCREW, 4.5-5.0MM/4PK
14114535P	RELIN SCREW+LS, 4.5X35MM 2FS POLY
14114540P	RELIN SCREW+LS, 4.5X40MM 2FS POLY
14114545P	RELIN SCREW+LS, 4.5X45MM 2FS POLY
14115535P	RELIN SCREW+LS, 5.5X35MM 2FS POLY
14115540P	RELIN SCREW+LS, 5.5X40MM 2FS POLY
14115545P	RELIN SCREW+LS, 5.5X45MM 2FS POLY
14115550P	RELIN SCREW+LS, 5.5X50MM 2FS POLY
14116535P	RELIN SCREW+LS, 6.5X35MM 2FS POLY
14116540P	RELIN SCREW+LS, 6.5X40MM 2FS POLY
14116545P	RELIN SCREW+LS, 6.5X45MM 2FS POLY
14116550P	RELIN SCREW+LS, 6.5X50MM 2FS POLY
14116555P	RELIN SCREW+LS, 6.5X55MM 2FS POLY
14117540P	RELIN SCREW+LS, 7.5X40MM 2FS POLY
14117545P	RELIN SCREW+LS, 7.5X45MM 2FS POLY
14117550P	RELIN SCREW+LS, 7.5X50MM 2FS POLY
14117555P	RELIN SCREW+LS, 7.5X55MM 2FS POLY
14118540P	RELIN SCREW+LS, 8.5X40MM 2FS POLY
14118545P	RELIN SCREW+LS, 8.5X45MM 2FS POLY
14118550P	RELIN SCREW+LS, 8.5X50MM 2FS POLY
14118555P	RELIN SCREW+LS, 8.5X55MM 2FS POLY
14414017P	RSS SCREW+LS, 4.0X17.5MM POLY 2S
14414020P	RSS SCREW+LS, 4.0X20MM POLY 2S
14414022P	RSS SCREW+LS, 4.0X22.5MM POLY 2S
14414025P	RSS SCREW+LS, 4.0X25MM POLY 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14414027P	RSS SCREW+LS, 4.0X27.5MM POLY 2S
14414030P	RSS SCREW+LS, 4.0X30MM POLY 2S
14414035P	RSS SCREW+LS, 4.0X35MM POLY 2S
14414040P	RSS SCREW+LS, 4.0X40MM POLY 2S
14414417P	RSS SCREW+LS, 4.5X17.5MM POLY 2S
14414420P	RSS SCREW+LS, 4.5X20MM POLY 2S
14414422P	RSS SCREW+LS, 4.5X22.5MM POLY 2S
14414425P	RSS SCREW+LS, 4.5X25MM POLY 2S
14414427P	RSS SCREW+LS, 4.5X27.5MM POLY 2S
14414430P	RSS SCREW+LS, 4.5X30MM POLY 2S
14414435P	RSS SCREW+LS, 4.5X35MM POLY 2S
14414440P	RSS SCREW+LS, 4.5X40MM POLY 2S
14415017P	RSS SCREW+LS, 5.0X17.5MM POLY 2S
14415020P	RSS SCREW+LS, 5.0X20MM POLY 2S
14415022P	RSS SCREW+LS, 5.0X22.5MM POLY 2S
14415025P	RSS SCREW+LS, 5.0X25MM POLY 2S
14415027P	RSS SCREW+LS, 5.0X27.5MM POLY 2S
14415030P	RSS SCREW+LS, 5.0X30MM POLY 2S
14415035P	RSS SCREW+LS, 5.0X35MM POLY 2S
14415040P	RSS SCREW+LS, 5.0X40MM POLY 2S
14415520P	RSS SCREW+LS, 5.5X20MM POLY 2S
14415522P	RSS SCREW+LS, 5.5X22.5MM POLY 2S
14415525P	RSS SCREW+LS, 5.5X25MM POLY 2S
14415527P	RSS SCREW+LS, 5.5X27.5MM POLY 2S
14415530P	RSS SCREW+LS, 5.5X30MM POLY 2S
14415535P	RSS SCREW+LS, 5.5X35MM POLY 2S
14415540P	RSS SCREW+LS, 5.5X40MM POLY 2S
14415545P	RSS SCREW+LS, 5.5X45MM POLY 2S
14415550P	RSS SCREW+LS, 5.5X50MM POLY 2S
14415555P	RSS SCREW+LS, 5.5X55MM POLY 2S
14415560P	RSS SCREW+LS, 5.5X60MM POLY 2S
14415565P	RSS SCREW+LS, 5.5X65MM POLY 2S
14415570P	RSS SCREW+LS, 5.5X70MM POLY 2S
14415575P	RSS SCREW+LS, 5.5X75MM POLY 2S
14415580P	RSS SCREW+LS, 5.5X80MM POLY 2S
14416025P	RSS SCREW+LS, 6.0X25MM POLY 2S
14416030P	RSS SCREW+LS, 6.0X30MM POLY 2S
14416035P	RSS SCREW+LS, 6.0X35MM POLY 2S
14416040P	RSS SCREW+LS, 6.0X40MM POLY 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14416045P	RSS SCREW+LS, 6.0X45MM POLY 2S
14416050P	RSS SCREW+LS, 6.0X50MM POLY 2S
14416525P	RSS SCREW+LS, 6.5X25MM POLY 2S
14416530P	RSS SCREW+LS, 6.5X30MM POLY 2S
14416535P	RSS SCREW+LS, 6.5X35MM POLY 2S
14416540P	RSS SCREW+LS, 6.5X40MM POLY 2S
14416545P	RSS SCREW+LS, 6.5X45MM POLY 2S
14416550P	RSS SCREW+LS, 6.5X50MM POLY 2S
14416555P	RSS SCREW+LS, 6.5X55MM POLY 2S
14416560P	RSS SCREW+LS, 6.5X60MM POLY 2S
14416565P	RSS SCREW+LS, 6.5X65MM POLY 2S
14416570P	RSS SCREW+LS, 6.5X70MM POLY 2S
14416575P	RSS SCREW+LS, 6.5X75MM POLY 2S
14416580P	RSS SCREW+LS, 6.5X80MM POLY 2S
14416590P	RSS SCREW+LS, 6.5X90MM POLY 2S
14417530P	RSS SCREW+LS, 7.5X30MM POLY 2S
14417535P	RSS SCREW+LS, 7.5X35MM POLY 2S
14417540P	RSS SCREW+LS, 7.5X40MM POLY 2S
14417545P	RSS SCREW+LS, 7.5X45MM POLY 2S
14417550P	RSS SCREW+LS, 7.5X50MM POLY 2S
14417555P	RSS SCREW+LS, 7.5X55MM POLY 2S
14417560P	RSS SCREW+LS, 7.5X60MM POLY 2S
14417565P	RSS SCREW+LS, 7.5X65MM POLY 2S
14417570P	RSS SCREW+LS, 7.5X70MM POLY 2S
14417575P	RSS SCREW+LS, 7.5X75MM POLY 2S
14417580P	RSS SCREW+LS, 7.5X80MM POLY 2S
14417590P	RSS SCREW+LS, 7.5X90MM POLY 2S
14418560P	RSS SCREW+LS, 8.5X60MM POLY 2S
14418565P	RSS SCREW+LS, 8.5X65MM POLY 2S
14418570P	RSS SCREW+LS, 8.5X70MM POLY 2S
14418580P	RSS SCREW+LS, 8.5X80MM POLY 2S
14418590P	RSS SCREW+LS, 8.5X90MM POLY 2S
14424017P	RSS SCREW+LS, 4.0X17.5MM RED 2S
14424020P	RSS SCREW+LS, 4.0X20MM RED 2S
14424025P	RSS SCREW+LS, 4.0X25MM RED 2S
14424030P	RSS SCREW+LS, 4.0X30MM RED 2S
14424035P	RSS SCREW+LS, 4.0X35MM RED 2S
14424517P	RSS SCREW+LS, 4.5X17.5MM RED 2S
14424520P	RSS SCREW+LS, 4.5X20MM RED 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14424525P	RSS SCREW+LS, 4.5X25MM RED 2S
14424530P	RSS SCREW+LS, 4.5X30MM RED 2S
14424535P	RSS SCREW+LS, 4.5X35MM RED 2S
14424540P	RSS SCREW+LS, 4.5X40MM RED 2S
14425017P	RSS SCREW+LS, 5.0X17.5MM RED 2S
14425020P	RSS SCREW+LS, 5.0X20MM RED 2S
14425025P	RSS SCREW+LS, 5.0X25MM RED 2S
14425030P	RSS SCREW+LS, 5.0X30MM RED 2S
14425035P	RSS SCREW+LS, 5.0X35MM RED 2S
14425040P	RSS SCREW+LS, 5.0X40MM RED 2S
14425045P	RSS SCREW+LS, 5.0X45MM RED 2S
14426030P	RSS SCREW+LS, 6.0X30MM RED 2S
14426035P	RSS SCREW+LS, 6.0X35MM RED 2S
14426040P	RSS SCREW+LS, 6.0X40MM RED 2S
14426045P	RSS SCREW+LS, 6.0X45MM RED 2S
14426530P	RSS SCREW+LS, 6.5X30MM RED 2S
14426535P	RSS SCREW+LS, 6.5X35MM RED 2S
14426540P	RSS SCREW+LS, 6.5X40MM RED 2S
14426545P	RSS SCREW+LS, 6.5X45MM RED 2S
14444025P	RSS SCREW+LS, 4.0X25MM UNI 2S
14444030P	RSS SCREW+LS, 4.0X30MM UNI 2S
14444035P	RSS SCREW+LS, 4.0X35MM UNI 2S
14444517P	RSS SCREW+LS, 4.5X17.5MM UNI 2S
14444525P	RSS SCREW+LS, 4.5X25MM UNI 2S
14444530P	RSS SCREW+LS, 4.5X30MM UNI 2S
14444535P	RSS SCREW+LS, 4.5X35MM UNI 2S
14444540P	RSS SCREW+LS, 4.5X40MM UNI 2S
14445025P	RSS SCREW+LS, 5.0X25MM UNI 2S
14445030P	RSS SCREW+LS, 5.0X30MM UNI 2S
14445035P	RSS SCREW+LS, 5.0X35MM UNI 2S
14445040P	RSS SCREW+LS, 5.0X40MM UNI 2S
14445045P	RSS SCREW+LS, 5.0X45MM UNI 2S
14445525P	RSS SCREW+LS, 5.5X25MM UNI 2S
14445530P	RSS SCREW+LS, 5.5X30MM UNI 2S
14445535P	RSS SCREW+LS, 5.5X35MM UNI 2S
14445540P	RSS SCREW+LS, 5.5X40MM UNI 2S
14445545P	RSS SCREW+LS, 5.5X45MM UNI 2S
14446030P	RSS SCREW+LS, 6.0X30MM UNI 2S
14446035P	RSS SCREW+LS, 6.0X35MM UNI 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
14446040P	RSS SCREW+LS, 6.0X40MM UNI 2S
14446045P	RSS SCREW+LS, 6.0X45MM UNI 2S
14446530P	RSS SCREW+LS, 6.5X30MM UNI 2S
14446535P	RSS SCREW+LS, 6.5X35MM UNI 2S
14446540P	RSS SCREW+LS, 6.5X40MM UNI 2S
14446545P	RSS SCREW+LS, 6.5X45MM UNI 2S
14455530P	RSS SCREW+LS, 5.5X30MM UNI-RED 2S
14455535P	RSS SCREW+LS, 5.5X35MM UNI-RED 2S
14455540P	RSS SCREW+LS, 5.5X40MM UNI-RED 2S
14456030P	RSS SCREW+LS, 6.0X30MM UNI-RED 2S
14456035P	RSS SCREW+LS, 6.0X35MM UNI-RED 2S
14456040P	RSS SCREW+LS, 6.0X40MM UNI-RED 2S
14456045P	RSS SCREW+LS, 6.0X45MM UNI-RED 2S
14456050P	RSS SCREW+LS, 6.0X50MM UNI-RED 2S
14484017P	RSS SCREW+LS, 4.0X17.5MM MONO 2S
14484020P	RSS SCREW+LS, 4.0X20MM MONO 2S
14484025P	RSS SCREW+LS, 4.0X25MM MONO 2S
14484030P	RSS SCREW+LS, 4.0X30MM MONO 2S
14484035P	RSS SCREW+LS, 4.0X35MM MONO 2S
14484517P	RSS SCREW+LS, 4.5X17.5MM MONO 2S
14484520P	RSS SCREW+LS, 4.5X20MM MONO 2S
14484525P	RSS SCREW+LS, 4.5X25MM MONO 2S
14484530P	RSS SCREW+LS, 4.5X30MM MONO 2S
14484535P	RSS SCREW+LS, 4.5X35MM MONO 2S
14484540P	RSS SCREW+LS, 4.5X40MM MONO 2S
14485017P	RSS SCREW+LS, 5.0X17.5MM MONO 2S
14485020P	RSS SCREW+LS, 5.0X20MM MONO 2S
14485025P	RSS SCREW+LS, 5.0X25MM MONO 2S
14485030P	RSS SCREW+LS, 5.0X30MM MONO 2S
14485035P	RSS SCREW+LS, 5.0X35MM MONO 2S
14485040P	RSS SCREW+LS, 5.0X40MM MONO 2S
14485045P	RSS SCREW+LS, 5.0X45MM MONO 2S
14485520P	RSS SCREW+LS, 5.5X20MM MONO 2S
14485527P	RSS SCREW+LS, 5.5X27.5MM MONO 2S
14486030P	RSS SCREW+LS, 6.0X30MM MONO 2S
14486035P	RSS SCREW+LS, 6.0X35MM MONO 2S
14486040P	RSS SCREW+LS, 6.0X40MM MONO 2S
14486045P	RSS SCREW+LS, 6.0X45MM MONO 2S
14486530P	RSS SCREW+LS, 6.5X30MM MONO 2S



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
14486535P	RSS SCREW+LS, 6.5X35MM MONO 2S
14486540P	RSS SCREW+LS, 6.5X40MM MONO 2S
14486545P	RSS SCREW+LS, 6.5X45MM MONO 2S
16171110P	RELINE TULIP+LS, MOD
16291111P	RELINE MAS MOD+LS, REDUCTION EXT
16314535P	RELINE SCREW+LS, 4.5X35MM 2FC POLY
16314540P	RELINE SCREW+LS, 4.5X40MM 2FC POLY
16314545P	RELINE SCREW+LS, 4.5X45MM 2FC POLY
16315535P	RELINE SCREW+LS, 5.5X35MM 2FC POLY
16315540P	RELINE SCREW+LS, 5.5X40MM 2FC POLY
16315545P	RELINE SCREW+LS, 5.5X45MM 2FC POLY
16315550P	RELINE SCREW+LS, 5.5X50MM 2FC POLY
16315555P	RELINE SCREW+LS, 5.5X55MM 2FC POLY
16315560P	RELINE SCREW+LS, 5.5X60MM 2FC POLY
16316535P	RELINE SCREW+LS, 6.5X35MM 2FC POLY
16316540P	RELINE SCREW+LS, 6.5X40MM 2FC POLY
16316545P	RELINE SCREW+LS, 6.5X45MM 2FC POLY
16316550P	RELINE SCREW+LS, 6.5X50MM 2FC POLY
16316555P	RELINE SCREW+LS, 6.5X55MM 2FC POLY
16316560P	RELINE SCREW+LS, 6.5X60MM 2FC POLY
16317535P	RELINE SCREW+LS, 7.5X35MM 2FC POLY
16317540P	RELINE SCREW+LS, 7.5X40MM 2FC POLY
16317545P	RELINE SCREW+LS, 7.5X45MM 2FC POLY
16317550P	RELINE SCREW+LS, 7.5X50MM 2FC POLY
16317555P	RELINE SCREW+LS, 7.5X55MM 2FC POLY
16317560P	RELINE SCREW+LS, 7.5X60MM 2FC POLY
16318535P	RELINE SCREW+LS, 8.5X35MM 2FC POLY
16318540P	RELINE SCREW+LS, 8.5X40MM 2FC POLY
16318545P	RELINE SCREW+LS, 8.5X45MM 2FC POLY
16318550P	RELINE SCREW+LS, 8.5X50MM 2FC POLY
16318555P	RELINE SCREW+LS, 8.5X55MM 2FC POLY
16318560P	RELINE SCREW+LS, 8.5X60MM 2FC POLY
16324525P	RELINE MAS RED SCREW+LS, 4.5X25POLY2C
16324530P	RELINE MAS RED SCREW+LS, 4.5X30POLY2C
16324535P	RELINE MAS RED SCREW+LS, 4.5X35POLY2C
16324540P	RELINE MAS RED SCREW+LS, 4.5X40POLY2C
16324545P	RELINE MAS RED SCREW+LS, 4.5X45POLY2C
16325530P	RELINE MAS RED SCREW+LS, 5.5X30POLY2C
16325535P	RELINE MAS RED SCREW+LS, 5.5X35POLY2C



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
16325540P	RELIN MAS RED SCREW+LS, 5.5X40POLY2C
16325545P	RELIN MAS RED SCREW+LS, 5.5X45POLY2C
16325550P	RELIN MAS RED SCREW+LS, 5.5X50POLY2C
16325555P	RELIN MAS RED SCREW+LS, 5.5X55POLY2C
16326530P	RELIN MAS RED SCREW+LS, 6.5X30POLY2C
16326535P	RELIN MAS RED SCREW+LS, 6.5X35POLY2C
16326540P	RELIN MAS RED SCREW+LS, 6.5X40POLY2C
16326545P	RELIN MAS RED SCREW+LS, 6.5X45POLY2C
16326550P	RELIN MAS RED SCREW+LS, 6.5X50POLY2C
16326555P	RELIN MAS RED SCREW+LS, 6.5X55POLY2C
16326560P	RELIN MAS-R SCREW+LS, 6.5X60MM P.I.2C
16327535P	RELIN MAS RED SCREW+LS, 7.5X35POLY2C
16327540P	RELIN MAS RED SCREW+LS, 7.5X40POLY2C
16327545P	RELIN MAS RED SCREW+LS, 7.5X45POLY2C
16327550P	RELIN MAS RED SCREW+LS, 7.5X50POLY2C
16327555P	RELIN MAS RED SCREW+LS, 7.5X55POLY2C
16327560P	RELIN MAS-R SCREW+LS, 7.5X60MM P.I.2C
16328510P	RELIN MAS-R SCREW+LS, 8.5X100 P.I.2C
16328535P	RELIN MAS RED SCREW+LS, 8.5X35POLY2C
16328540P	RELIN MAS RED SCREW+LS, 8.5X40POLY2C
16328545P	RELIN MAS RED SCREW+LS, 8.5X45POLY2C
16328550P	RELIN MAS RED SCREW+LS, 8.5X50POLY2C
16328555P	RELIN MAS RED SCREW+LS, 8.5X55POLY2C
16328560P	RELIN MAS-R SCREW+LS, 8.5X60MM P.I.2C
16328570P	RELIN MAS-R SCREW+LS, 8.5X70MM P.I.2C
16328580P	RELIN MAS-R SCREW+LS, 8.5X80MM P.I.2C
16328590P	RELIN MAS-R SCREW+LS, 8.5X90MM P.I.2C
16329510P	RELIN MAS-R SCREW+LS, 9.5X100 P.I.2C
16329511P	RELIN MAS-R SCREW+LS, 9.5X110 P.I.2C
16329545P	RELIN MAS RED SCREW+LS, 9.5X45POLY2C
16329570P	RELIN MAS-R SCREW+LS, 9.5X70MM P.I.2C
16329580P	RELIN MAS-R SCREW+LS, 9.5X80MM P.I.2C
16329590P	RELIN MAS-R SCREW+LS, 9.5X90MM P.I.2C
16384530P	RELIN MAS-T SCREW+LS, 4.5X30MM CANN
16384535P	RELIN MAS-T SCREW+LS, 4.5X35MM CANN
16384540P	RELIN MAS-T SCREW+LS, 4.5X40MM CANN
16384545P	RELIN MAS-T SCREW+LS, 4.5X45MM CANN
16384550P	RELIN MAS-T SCREW+LS, 4.5X50MM CANN
16384555P	RELIN MAS-T SCREW+LS, 4.5X55MM CANN



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
16385525P	RELIN MAS-T SCREW+LS, 5.5X25MM CANN
16385530P	RELIN MAS-T SCREW+LS, 5.5X30MM CANN
16385535P	RELIN MAS-T SCREW+LS, 5.5X35MM CANN
16385540P	RELIN MAS-T SCREW+LS, 5.5X40MM CANN
16385545P	RELIN MAS-T SCREW+LS, 5.5X45MM CANN
16385550P	RELIN MAS-T SCREW+LS, 5.5X50MM CANN
16385555P	RELIN MAS-T SCREW+LS, 5.5X55MM CANN
16385560P	RELIN MAS-T SCREW+LS, 5.5X60MM CANN
16386525P	RELIN MAS-T SCREW+LS, 6.5X25MM CANN
16386535P	RELIN MAS-T SCREW+LS, 6.5X35MM CANN
16386540P	RELIN MAS-T SCREW+LS, 6.5X40MM CANN
16386545P	RELIN MAS-T SCREW+LS, 6.5X45MM CANN
16386550P	RELIN MAS-T SCREW+LS, 6.5X50MM CANN
16386555P	RELIN MAS-T SCREW+LS, 6.5X55MM CANN
16386560P	RELIN MAS-T SCREW+LS, 6.5X60MM CANN
16387535P	RELIN MAS-T SCREW+LS, 7.5X35MM CANN
16387540P	RELIN MAS-T SCREW+LS, 7.5X40MM CANN
16387545P	RELIN MAS-T SCREW+LS, 7.5X45MM CANN
16387550P	RELIN MAS-T SCREW+LS, 7.5X50MM CANN
16387555P	RELIN MAS-T SCREW+LS, 7.5X55MM CANN
16387560P	RELIN MAS-T SCREW+LS, 7.5X60MM CANN
16410510P	RELIN MAS SCREW+LS, 10.5X100MM 2CP.L.
16410590P	RELIN MAS SCREW+LS, 10.5X90MM 2CP.I.
16414525P	RELIN MAS SCREW+LS, 4.5X25MM 2CPOLY
16414530P	RELIN MAS SCREW+LS, 4.5X30MM 2CPOLY
16414535P	RELIN MAS SCREW+LS, 4.5X35MM 2CPOLY
16414540P	RELIN MAS SCREW+LS, 4.5X40MM 2CPOLY
16414545P	RELIN MAS SCREW+LS, 4.5X45MM 2CPOLY
16414550P	RELIN MAS SCREW+LS, 4.5X50MM 2CPOLY
16415030P	RELIN MAS SCREW+LS, 5.0X30MM 2CPOLY
16415035P	RELIN MAS SCREW+LS, 5.0X35MM 2CPOLY
16415040P	RELIN MAS SCREW+LS, 5.0X40MM 2CPOLY
16415045P	RELIN MAS SCREW+LS, 5.0X45MM 2CPOLY
16415050P	RELIN MAS SCREW+LS, 5.0X50MM 2CPOLY
16415055P	RELIN MAS SCREW+LS, 5.0X55MM 2CPOLY
16415525P	RELIN MAS SCREW+LS, 5.5X25MM 2CPOLY
16415530P	RELIN MAS SCREW+LS, 5.5X30MM 2CPOLY
16415535P	RELIN MAS SCREW+LS, 5.5X35MM 2CPOLY
16415540P	RELIN MAS SCREW+LS, 5.5X40MM 2CPOLY



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
16415545P	RELIN MAS SCREW+LS, 5.5X45MM 2CPOLY
16415550P	RELIN MAS SCREW+LS, 5.5X50MM 2CPOLY
16415555P	RELIN MAS SCREW+LS, 5.5X55MM 2CPOLY
16415560P	RELIN MAS SCREW+LS, 5.5X60MM 2CP.I.
16415565P	RELIN MAS SCREW+LS, 5.5X65MM 2CP.I.
16416040P	RELIN MAS SCREW+LS, 6.0X40MM 2CPOLY
16416045P	RELIN MAS SCREW+LS, 6.0X45MM 2CPOLY
16416050P	RELIN MAS SCREW+LS, 6.0X50MM 2CPOLY
16416525P	RELIN MAS SCREW+LS, 6.5X25MM 2CPOLY
16416530P	RELIN MAS SCREW+LS, 6.5X30MM 2CPOLY
16416535P	RELIN MAS SCREW+LS, 6.5X35MM 2CPOLY
16416540P	RELIN MAS SCREW+LS, 6.5X40MM 2CPOLY
16416545P	RELIN MAS SCREW+LS, 6.5X45MM 2CPOLY
16416550P	RELIN MAS SCREW+LS, 6.5X50MM 2CPOLY
16416555P	RELIN MAS SCREW+LS, 6.5X55MM 2CPOLY
16416560P	RELIN MAS SCREW+LS, 6.5X60MM 2CP.I.
16416565P	RELIN MAS SCREW+LS, 6.5X65MM 2CP.I.
16417510P	RELIN MAS SCREW+LS, 7.5X100MM 2CP.I.
16417525P	RELIN MAS SCREW+LS, 7.5X25MM 2CPOLY
16417530P	RELIN MAS SCREW+LS, 7.5X30MM 2CPOLY
16417535P	RELIN MAS SCREW+LS, 7.5X35MM 2CPOLY
16417540P	RELIN MAS SCREW+LS, 7.5X40MM 2CPOLY
16417545P	RELIN MAS SCREW+LS, 7.5X45MM 2CPOLY
16417550P	RELIN MAS SCREW+LS, 7.5X50MM 2CPOLY
16417555P	RELIN MAS SCREW+LS, 7.5X55MM 2CPOLY
16417560P	RELIN MAS SCREW+LS, 7.5X60MM 2CP.I.
16417565P	RELIN MAS SCREW+LS, 7.5X65MM 2CP.I.
16417570P	RELIN MAS SCREW+LS, 7.5X70MM 2CP.I.
16417580P	RELIN MAS SCREW+LS, 7.5X80MM 2CP.I.
16417590P	RELIN MAS SCREW+LS, 7.5X90MM 2CP.I.
16418510P	RELIN MAS SCREW+LS, 8.5X100MM 2CP.I.
16418530P	RELIN MAS SCREW+LS, 8.5X30MM 2CPOLY
16418535P	RELIN MAS SCREW+LS, 8.5X35MM 2CPOLY
16418540P	RELIN MAS SCREW+LS, 8.5X40MM 2CPOLY
16418545P	RELIN MAS SCREW+LS, 8.5X45MM 2CPOLY
16418550P	RELIN MAS SCREW+LS, 8.5X50MM 2CPOLY
16418555P	RELIN MAS SCREW+LS, 8.5X55MM 2CPOLY
16418560P	RELIN MAS SCREW+LS, 8.5X60MM 2CP.I.
16418565P	RELIN MAS SCREW+LS, 8.5X65MM 2CP.I.



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
16418570P	RELIN MAS SCREW+LS, 8.5X70MM 2CP.I.
16418580P	RELIN MAS SCREW+LS, 8.5X80MM 2CP.I.
16418590P	RELIN MAS SCREW+LS, 8.5X90MM 2CP.I.
16419510P	RELIN MAS SCREW+LS, 9.5X100MM 2CP.I.
16419535P	RELIN MAS SCREW+LS, 9.5X35MM 2CPOLY
16419540P	RELIN MAS SCREW+LS, 9.5X40MM 2CPOLY
16419545P	RELIN MAS SCREW+LS, 9.5X45MM 2CPOLY
16419550P	RELIN MAS SCREW+LS, 9.5X50MM 2CPOLY
16419555P	RELIN MAS SCREW+LS, 9.5X55MM 2CPOLY
16419560P	RELIN MAS SCREW+LS, 9.5X60MM 2CP.I.
16419570P	RELIN MAS SCREW+LS, 9.5X70MM 2CP.I.
16419580P	RELIN MAS SCREW+LS, 9.5X80MM 2CP.I.
16419590P	RELIN MAS SCREW+LS, 9.5X90MM 2CP.I.
16514525P	RELIN MAS SCREW+LS, 4.5X25MM TRACPOLY
16514530P	RELIN MAS SCREW+LS, 4.5X30MM TRACPOLY
16514535P	RELIN MAS SCREW+LS, 4.5X35MM TRACPOLY
16514540P	RELIN MAS SCREW+LS, 4.5X40MM TRACPOLY
16514545P	RELIN MAS SCREW+LS, 4.5X45MM TRACPOLY
16515530P	RELIN MAS SCREW+LS, 5.5X30MM TRACPOLY
16515535P	RELIN MAS SCREW+LS, 5.5X35MM TRACPOLY
16515540P	RELIN MAS SCREW+LS, 5.5X40MM TRACPOLY
16515545P	RELIN MAS SCREW+LS, 5.5X45MM TRACPOLY
16515550P	RELIN MAS SCREW+LS, 5.5X50MM TRACPOLY
16516530P	RELIN MAS SCREW+LS, 6.5X30MM TRACPOLY
16516535P	RELIN MAS SCREW+LS, 6.5X35MM TRACPOLY
16516540P	RELIN MAS SCREW+LS, 6.5X40MM TRACPOLY
16516545P	RELIN MAS SCREW+LS, 6.5X45MM TRACPOLY
16516550P	RELIN MAS SCREW+LS, 6.5X50MM TRACPOLY
16516555P	RELIN MAS SCREW+LS, 6.5X55MM TRACPOLY
16517535P	RELIN MAS SCREW+LS, 7.5X35MM TRACPOLY
16517540P	RELIN MAS SCREW+LS, 7.5X40MM TRACPOLY
16517545P	RELIN MAS SCREW+LS, 7.5X45MM TRACPOLY
16517550P	RELIN MAS SCREW+LS, 7.5X50MM TRACPOLY
16517555P	RELIN MAS SCREW+LS, 7.5X55MM TRACPOLY
16518540P	RELIN MAS SCREW+LS, 8.5X40MM TRACPOLY
16518545P	RELIN MAS SCREW+LS, 8.5X45MM TRACPOLY
16518550P	RELIN MAS SCREW+LS, 8.5X50MM TRACPOLY
16519540P	RELIN MAS SCREW+LS, 9.5X40MM TRACPOLY
16519545P	RELIN MAS SCREW+LS, 9.5X45MM TRACPOLY



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8051100	ARM15T LOCK SCREW, ILIAC
8052025	ARM15T ROD. 25MM CONVENTIONAL PRE-BENT
8052030	ARM15T ROD. 30MM CONVENTIONAL PRE-BENT
8052035	ARM15T ROD. 35MM CONVENTIONAL PRE-BENT
8052040	ARM15T ROD. 40MM CONVENTIONAL PRE-BENT
8052045	ARM15T ROD, 45MM CONVENTIONAL PRE-BENT
8052050	ARM15T ROD. 50MM CONVENTIONAL PRE-BENT
8052055	ARM15T ROD. 55MM CONVENTIONAL PRE-BENT
8052060	ARM15T ROD. 60MM CONVENTIONAL PRE-BENT
8052065	ARM15T ROD. 65MM CONVENTIONAL PRE-BENT
8052070	ARM15T ROD, 70MM CONVENTIONAL PRE-BENT
8052075	ARM15T ROD. 75MM CONVENTIONAL PRE-BENT
8052080	ARM15T ROD, 80MM CONVENTIONAL PRE-BENT
8052085	ARM15T ROD, 85MM CONVENTIONAL PRE-BENT
8052090	ARM15T ROD, 90MM CONVENTIONAL PRE-BENT
8052095	ARM15T ROD, 95MM CONVENTIONAL PRE-BENT
8052100	ARM15T ROD, 100MM CONVENTIONAL PRE-BENT
8052110	ARM15T ROD. 110MM CONVENTIONAL PRE-BENT
8052120	ARM15T ROD. 120MM CONVENTIONAL PRE-BENT
8052140	ARM15T ROD. 140MM CONVENTIONAL PRE-BENT
8052160	ARM15T ROD, 160MM CONVENTIONAL PRE-BENT
8052180	ARM15T ROD. 180MM CONVENTIONAL PRE-BENT
8052200	ARM15T ROD. 200MM CONVENTIONAL PRE-BENT
8052300	ARM15T ROD, 300MM CONVENTIONAL PRE-BENT
8053500	ARM15T ROD. 500MM PB LORDOTIC KYPHOTIC
8056510	ARM15T SCREW, 6.5X100MM ILIAC
8056560	ARM15T SCREW, 6.5X60MM ILIAC
8057510	ARM15T SCREW, 7.5X100MM ILIAC
8058510	ARM15T SCREW, 8.5X100MM ILIAC
8058570	ARM15T SCREW, 8.5X70MM ILIAC
8058580	ARM15T SCREW, 8.5X80MM ILIAC
8058590	ARM15T SCREW, 8.5X90MM ILIAC
8151100	ARM15S LOCK SCREW, CLOSED HEAD
8252004	ASF15C LATERAL OFFSET ROD
8252040	ARM15T COCR ROD. 40MM STRAIGHT
8252045	ARM15T COCR ROD, 45MM STRAIGHT
8252050	ARM15T COCR ROD. 50MM STRAIGHT
8252055	ARM15T COCR ROD. 55MM STRAIGHT
8252060	ARM15T COCR ROD. 60MM STRAIGHT



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8252065	ARM15T COCR ROD. 65MM STRAIGHT
8252070	ARM15T COCR ROD. 70MM STRAIGHT
8252075	ARM15T COCR ROD, 75MM STRAIGHT
8252080	ARM15T COCR ROD. 80MM STRAIGHT
8252085	ARM15T COCR ROD. 85MM STRAIGHT
8252100	ARM15T COCR ROD, 100MM STRAIGHT
8252200	ARM15T COCR ROD, 200MM STRAIGHT
8252250	ARM15T COCR ROD, 250MM STRAIGHT
8252300	ARM15T COCR ROD. 300MM STRAIGHT
8252500	ARM15T COCR ROD. 500MM STRAIGHT
8252600	ARM15T COCR ROD, 600MM STRAIGHT
8253025	ARM15T COCR ROD, 25MM CONVENTIONAL PB
8253030	ARM15T COCR ROD, 30MM CONVENTIONAL PB
8253035	ARM15T COCR ROD. 35MM CONVENTIONAL PB
8253040	ARM15T COCR ROD. 40MM CONVENTIONAL PB
8253045	ARM15T COCR ROD. 45MM CONVENTIONAL PB
8253050	ARM15T COCR ROD, 50MM CONVENTIONAL PB
8253055	ARM15T COCR ROD. 55MM CONVENTIONAL PB
8253060	ARM15T COCR ROD, 60MM CONVENTIONAL PB
8253065	ARM15T COCR ROD. 65MM CONVENTIONAL PB
8253070	ARM15T COCR ROD. 70MM CONVENTIONAL PB
8253075	ARM15T COCR ROD. 75MM CONVENTIONAL PB
8253080	ARM15T COCR ROD. 80MM CONVENTIONAL PB
8253085	ARM15T COCR ROD, 85MM CONVENTIONAL P6
8253090	ARM15T COCR ROD. 90MM CONVENTIONAL PB
8253095	ARM15T COCR ROD. 95MM CONVENTIONAL PB
8253100	ARM15T COCR ROD. 100MM CONVENTIONAL PB
8253110	ARM15T COCR ROD. 110MM CONVENTIONAL PB
8253120	ARM15T COCR ROD, 120MM CONVENTIONAL PB
8254025	ARM15T SCREW, 4X25MM PLS
8254030	ARM15T SCREW, 4X30MM PLS
8254035	ARM15T SCREW, 4X35MM PLS
8254040	ARM15T SCREW, 4X40MM PLS
8254045	ARM15T SCREW, 4X45MM PLS
8254525	ARM15T SCREW, 4.5X25MM PLS
8254530	ARM15T SCREW, 4.5X30MM PLS
8254535	ARM15T SCREW, 4.5X35MM PLS
8254540	ARM15T SCREW, 4.5X40MM PLS
8254545	ARM15T SCREW, 4.5X45MM PLS



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8255025	ARM15T SCREW, 5X25MM PLS
8255030	ARM15T SCREW, 5X30MM PLS
8255035	ARM15T SCREW, 5X35MM PLS
8255040	ARM15T SCREW, 5X40MM PLS
8255045	ARM15T SCREW, 5X45MM PLS
8255050	ARM15T SCREW, 5X50MM PLS
8255525	ARM15T SCREW, 5.5X25MM PLS
8255530	ARM15T SCREW, 5.5X30MM PLS
8255535	ARM15T SCREW, 5.5X35MM PLS
8255540	ARM15T SCREW, 5.5X40MM PLS
8255545	ARM15T SCREW, 5.5X45MM PLS
8255550	ARM15T SCREW, 5.5X50MM PLS
8255555	ARM15T SCREW, 5.5X55MM PLS
8256525	ARM15T SCREW, 6.5X25MM PLS
8256530	ARM15T SCREW, 6.5X30MM PLS
8256535	ARM15T SCREW, 6.5X35MM PLS
8256540	ARM15T SCREW, 6.5X40MM PLS
8256545	ARM15T SCREW, 6.5X45MM PLS
8256550	ARM15T SCREW, 6.5X50MM PLS
8256555	ARM15T SCREW, 6.5X55MM PLS
8256560	ARM15T SCREW, 6.5X60MM PLS
8257530	ARM15T SCREW, 7.5X30MM PLS
8257535	ARM15T SCREW, 7.5X35MM PLS
8257540	ARM15T SCREW, 7.5X40MM PLS
8257545	ARM15T SCREW, 7.5X45MM PLS
8257550	ARM15T SCREW, 7.5X50MM PLS
8257555	ARM15T SCREW, 7.5X55MM PLS
8257560	ARM15T SCREW, 7.5X60MM PLS
8274030	ARM15T SCREW, 4X30MM REDUCTION
8274035	ARM15T SCREW, 4X35MM REDUCTION
8274040	ARM15T SCREW, 4X40MM REDUCTION
8274045	ARM15T SCREW, 4X45MM REDUCTION
8274525	ARM15T SCREW, 4.5X25MM REDUCTION
8274530	ARM15T SCREW, 4.5X30MM REDUCTION
8274535	ARM15T SCREW, 4.5X35MM REDUCTION
8274540	ARM15T SCREW, 4.5X40MM REDUCTION
8274545	ARM15T SCREW, 4.5X45MM REDUCTION
8274550	ARM15T SCREW, 4.5X50MM REDUCTION
8275025	ARM15T SCREW, 5X25MM REDUCTION



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8275030	ARM15T SCREW, 5X30MM REDUCTION
8275035	ARM15T SCREW, 5X35MM REDUCTION
8275040	ARM15T SCREW, 5X40MM REDUCTION
8275045	ARM15T SCREW, 5X45MM REDUCTION
8275050	ARM15T SCREW, 5X50MM REDUCTION
8275525	ARM15T SCREW, 5.5X25MM REDUCTION
8275530	ARM15T SCREW, 5.5X30MM REDUCTION
8275535	ARM15T SCREW, 5.5X35MM REDUCTION
8275540	ARM15T SCREW, 5.5X40MM REDUCTION
8275545	ARM15T SCREW, 5.5X45MM REDUCTION
8275550	ARM15T SCREW, 5.5X50MM REDUCTION
8275555	ARM15T SCREW, 5.5X55MM REDUCTION
8275560	ARM15T SCREW, 5.5X60MM REDUCTION
8276530	ARM15T SCREW, 6.5X30MM REDUCTION
8276535	ARM15T SCREW, 6.5X35MM REDUCTION
8276540	ARM15T SCREW, 6.5X40MM REDUCTION
8276545	ARM15T SCREW, 6.5X45MM REDUCTION
8276550	ARM15T SCREW, 6.5X50MM REDUCTION
8276555	ARM15T SCREW, 6.5X55MM REDUCTION
8276560	ARM15T SCREW, 6.5X60MM REDUCTION
8276565	ARM15T SCREW, 6.5X65MM REDUCTION
8277530	ARM15T SCREW, 7.5X30MM REDUCTION
8277535	ARM15T SCREW, 7.5X35MM REDUCTION
8277540	ARM15T SCREW, 7.5X40MM REDUCTION
8277545	ARM15T SCREW, 7.5X45MM REDUCTION
8277550	ARM15T SCREW, 7.5X50MM REDUCTION
8277555	ARM15T SCREW, 7.5X55MM REDUCTION
8277560	ARM15T SCREW, 7.5X60MM REDUCTION
8277565	ARM15T SCREW, 7.5X65MM REDUCTION
8352200	ARM15S ROD. 200MM ULTRA STRENGTH
8352300	ARM15S ROD. 300MM ULTRA STRENGTH
8352500	ARM15S ROD, 500MM ULTRA STRENGTH
8353300	ARM15C ROD. 300MM PB LORDOTIC KYPHOTIC
8353500	ARM15C ROD. 500MM PB LORDOTIC KYPHOTIC
8354025	ARM15S SCREW, 4X25MM PLS
8354030	ARM15S SCREW, 4X30MM PLS
8354035	ARM15S SCREW, 4X35MM PLS
8354040	ARM15S SCREW, 4X40MM PLS
8354045	ARM15S SCREW, 4X45MM PLS



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8354525	ARM15S SCREW, 4.5X25MM PLS
8354530	ARM15S SCREW, 4.5X30MM PLS
8354535	ARM15S SCREW, 4.5X35MM PLS
8354540	ARM15S SCREW, 4.5X40MM PLS
8354545	ARM15S SCREW, 4.5X45MM PLS
8355025	ARM15S SCREW, 5X25MM PLS
8355030	ARM15S SCREW, 5X30MM PLS
8355035	ARM15S SCREW, 5X35MM PLS
8355040	ARM15S SCREW, 5X40MM PLS
8355045	ARM15S SCREW, 5X45MM PLS
8355050	ARM15S SCREW, 5X50MM PLS
8355525	ARM15S SCREW, 5.5X25MM PLS
8355530	ARM15S SCREW, 5.5X30MM PLS
8355535	ARM15S SCREW, 5.5X35MM PLS
8355540	ARM15S SCREW, 5.5X40MM PLS
8355545	ARM15S SCREW, 5.5X45MM PLS
8355550	ARM15S SCREW, 5.5X50MM PLS
8355555	ARM15S SCREW, 5.5X55MM PLS
8356530	ARM15S SCREW, 6.5X30MM PLS
8356535	ARM15S SCREW, 6.5X35MM PLS
8356540	ARM15S SCREW, 6.5X40MM PLS
8356545	ARM15S SCREW, 6.5X45MM PLS
8356550	ARM15S SCREW, 6.5X50MM PLS
8356555	ARM15S SCREW, 6.5X55MM PLS
8364525	ARM16S SCREW, 4.5X25MM PLS
8364530	ARM16S SCREW, 4.5X30MM PLS
8364535	ARM16S SCREW, 4.5X35MM PLS
8364540	ARM16S SCREW, 4.5X40MM PLS
8364545	ARM16S SCREW, 4.5X45MM PLS
8365530	ARM16S SCREW, 5.5X30MM PLS
8365535	ARM16S SCREW, 5.5X35MM PLS
8365540	ARM16S SCREW, 5.5X40MM PLS
8365545	ARM16S SCREW, 5.5X45MM PLS
8366530	ARM16S SCREW, 6.5X30MM PLS
8366535	ARM16S SCREW, 6.5X35MM PLS
8366540	ARM16S SCREW, 6.5X40MM PLS
8366545	ARM16S SCREW, 6.5X45MM PLS
8374525	ARM15S SCREW, 4.5X25MM REDUCTION
8374530	ARM15S SCREW, 4.5X30MM REDUCTION



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8374535	ARM15S SCREW, 4.5X35MM REDUCTION
8374540	ARM15S SCREW, 4.5X40MM REDUCTION
8374545	ARM15S SCREW, 4.5X45MM REDUCTION
8375530	ARM15S SCREW, 5.5X30MM REDUCTION
8375535	ARM15S SCREW, 5.5X35MM REDUCTION
8375540	ARM15S SCREW, 5.5X40MM REDUCTION
8375545	ARM15S SCREW, 5.5X45MM REDUCTION
8375550	ARM15S SCREW, 5.5X50MM REDUCTION
8375555	ARM15S SCREW, 5.5X55MM REDUCTION
8376530	ARM15S SCREW, 6.5X30MM REDUCTION
8376535	ARM15S SCREW, 6.5X35MM REDUCTION
8376540	ARM15S SCREW, 6.5X40MM REDUCTION
8376545	ARM15S SCREW, 6.5X45MM REDUCTION
8376550	ARM15S SCREW, 6.5X50MM REDUCTION
8376555	ARM15S SCREW, 6.5X55MM REDUCTION
8377530	ARM15S SCREW, 7.5X30MM REDUCTION
8377535	ARM15S SCREW, 7.5X35MM REDUCTION
8377540	ARM15S SCREW, 7.5X40MM REDUCTION
8377545	ARM15S SCREW, 7.5X45MM REDUCTION
8377550	ARM15S SCREW, 7.5X50MM REDUCTION
8377555	ARM15S SCREW, 7.5X55MM REDUCTION
8384530	ARM16S SCREW, 4.5X30MM REDUCTION
8384535	ARM16S SCREW, 4.5X35MM REDUCTION
8384545	ARM16S SCREW, 4.5X45MM REDUCTION
8385040	ARM16S SCREW, 5X40MM REDUCTION
8385045	ARM16S SCREW, 5X45MM REDUCTION
8385050	ARM16S SCREW, 5X50MM REDUCTION
8385530	ARM16S SCREW, 6.35X30MM REDUCTION
8385540	ARM16S SCREW, 6.35X40MM REDUCTION
8385545	ARM16S SCREW, 6.35X45MM REDUCTION
8386540	ARM16S SCREW, 6.5X40MM REDUCTION
8386545	ARM16S SCREW, 6.5X45MM REDUCTION
8386550	ARM16S SCREW, 6.5X50MM REDUCTION
8387540	ARM16S SCREW, 7.5X40MM REDUCTION
8387545	ARM16S SCREW, 7.5X45MM REDUCTION
8387550	ARM16S SCREW, 7.5X50MM REDUCTION
8431510	ARM15T SCREW, 10.5X100MM ILILAC SL-CLOSED
8431511	ARM15T SCREW, 10.5X105MM ILIAC SL-CLOSED
8431512	ARM15T SCREW, 10.5X110MM ILIAC SL-CLOSED



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8431514	ARM15T SCREW, 10.5X120MM ILIAC SL-CLOSED
8431560	ARM15T SCREW, 10.5X60MM ILIAC SL-CLOSED
8431570	ARM15T SCREW, 10.5X70MM ILIAC SL-CLOSED
8431580	ARM15T SCREW, 10.5X80MM ILIAC SL-CLOSED
8431590	ARM15T SCREW, 10.5X90MM ILIAC SL-CLOSED
8436510	ARM15T SCREW, 6.5X100MM ILIAC SL-CLOSED
8436560	ARM15T SCREW, 6.5X60MM ILIAC SL-CLOSED
8436570	ARM15T SCREW, 6.5X70MM ILIAC SL-CLOSED
8436580	ARM15T SCREW, 6.5X80MM ILIAC SL-CLOSED
8436590	ARM15T SCREW, 6.5X90MM ILIAC SL-CLOSED
8437510	ARM15T SCREW, 7.5X100MM ILIAC SL-CLOSED
8437512	ARM15T SCREW, 7.5X110MM ILIAC SL-CLOSED
8437560	ARM15T SCREW, 7.5X60MM ILIAC SL-CLOSED
8437570	ARM15T SCREW, 7.5X70MM ILIAC SL-CLOSED
8437580	ARM15T SCREW, 7.5X80MM ILIAC SL-CLOSED
8437590	ARM15T SCREW, 7.5X90MM ILIAC SL-CLOSED
8438510	ARM15T SCREW, 8.5X100MM ILIAC SL-CLOSED
8438511	ARM15T SCREW, 8.5X105MM ILIAC SL-CLOSED
8438512	ARM15T SCREW, 8.5X110MM ILIAC SL-CLOSED
8438513	ARM15T SCREW, 8.5X115MM ILIAC SL-CLOSED
8438514	ARM15T SCREW, 8.5X120MM ILIAC SL-CLOSED
8438560	ARM15T SCREW, 8.5X60MM ILIAC SL-CLOSED
8438570	ARM15T SCREW, 8.5X70MM ILIAC SL-CLOSED
8438580	ARM15T SCREW, 8.5X80MM ILIAC SL-CLOSED
8438590	ARM15T SCREW, 8.5X90MM ILIAC SL-CLOSED
8439510	ARM15T SCREW, 9.5X100MM ILIAC SL-CLOSED
8439511	ARM15T SCREW, 9.5X105MM ILIAC SL-CLOSED
8439512	ARM15T SCREW, 9.5X110MM ILIAC SL-CLOSED
8439513	ARM15T SCREW, 9.5X115MM ILIAC SL-CLOSED
8439514	ARM15T SCREW, 9.5X120MM ILIAC SL-CLOSED
8439560	ARM15T SCREW, 9.5X60MM ILIAC SL-CLOSED
8439570	ARM15T SCREW, 9.5X70MM ILIAC SL-CLOSED
8439580	ARM15T SCREW, 9.5X80MM ILIAC SL-CLOSED
8439590	ARM15T SCREW, 9.5X90MM ILIAC SL-CLOSED
8441110	ARM15T SCREW, 11.5X100MM ILIAC SL-OPEN
8441111	ARM15T SCREW, 11.5X105MM ILIAC SL-OPEN
8441170	ARM15T SCREW, 11.5X70MM ILIAC SL-OPEN
8441180	ARM15T SCREW, 11.5X80MM ILIAC SL-OPEN
8441190	ARM15T SCREW, 11.5X90MM LILAC SL-OPEN



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8441210	ARM15T SCREW, 12.5X100MM ILIAC SL-OPEN
8441270	ARM15T SCREW, 12.5X70MM ILIAC SL-OPEN
8441280	ARM15T SCREW, 12.5X80MM ILIAC SL-OPEN
8441290	ARM15T SCREW, 12.5X90MM ILIAC SL-OPEN
8441510	ARM15T SCREW, 10.5X100MM ILIAC SL-OPEN
8441512	ARM15T SCREW, 10.5X110MM ILIAC SL-OPEN
8445514	ARM15T SCREW, 10.5X120MM ILIAC SL-OPEN
8441560	ARM15T SCREW, 10.5X60MM ILIAC SL-OPEN
8441570	ARM15T SCREW, 10.5X70MM ILIAC SL-OPEN
8441580	ARM15T SCREW, 10.5X80MM ILIAC SL-OPEN
8441590	ARM15T SCREW, 10.5X90MM LLIAC SL-OPEN
8446510	ARM15T SCREW, 6.5X100MM ILIAC SL-OPEN
8446560	ARM15T SCREW, 6.5X60MM LLIAC SL-OPEN
8446565	ARM15T SCREW, 6.5X65MM LLIAC SL-OPEN
8446570	ARM15T SCREW, 6.5X70MM LLIAC SL-OPEN
8446575	ARM15T SCREW, 6.5X75MM LLIAC SL-OPEN
8446580	ARM15T SCREW, 6.5X80MM LLIAC SL-OPEN
8446585	ARM15T SCREW, 6.5X85MM LLIAC SL-OPEN
8446590	ARM15T SCREW, 6.5X90MM LLIAC SL-OPEN
8446595	ARM15T SCREW, 6.5X95MM LLIAC SL-OPEN
8447510	ARM15T SCREW, 7.5X100MM LLIAC SL-OPEN
8447511	ARM15T SCREW, 7.5X105MM LLIAC SL-OPEN
8447512	ARM15T SCREW, 7.5X110MM LLIAC SL-OPEN
8447560	ARM15T SCREW, 7.5X60MM LLIAC SL-OPEN
8447565	ARM15T SCREW, 7.5X65MM ILIAC SL-OPEN
8447570	ARM15T SCREW, 7.5X70MM LLIAC SL-OPEN
8447575	ARM15T SCREW, 7.5X75MM LLIAC SL-OPEN
8447580	ARM15T SCREW, 7.5X80MM LLIAC SL-OPEN
8447585	ARM15T SCREW, 7.5X85MM LLIAC SL-OPEN
8447590	ARM15T SCREW, 7.5X90MM LLIAC SL-OPEN
8447595	ARM15T SCREW, 7.5X95MM LLIAC SL-OPEN
8448510	ARM15T SCREW, 8.5X100MM LLIAC SL-OPEN
8448511	ARM15T SCREW, 8.5X105MM LLIAC SL-OPEN
8448512	ARM15T SCREW, 8.5X110MM LLIAC SL-OPEN
8448514	ARM15T SCREW, 8.5X120MM LLIAC SL-OPEN
8448560	ARM15T SCREW, 8.5X60MM LLIAC SL-OPEN
8448565	ARM15T SCREW, 8.5X65MM LLIAC SL-OPEN
8448570	ARM15T SCREW, 8.5X70MM LLIAC SL-OPEN
8448575	ARM15T SCREW, 8.5X75MM ILIAC SL-OPEN



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8448580	ARM15T SCREW, 8.5X80MM LLIAC SL-OPEN
8448585	ARM15T SCREW, 8.5X85MM LLIAC SL-OPEN
8448590	ARM15T SCREW, 8.5X90MM LLIAC SL-OPEN
8448595	ARM15T SCREW, 8.5X95MM LLIAC SL-OPEN
8449510	ARM1 5T SCREW, 9.5X100MM LLIAC SL-OPEN
8449511	ARM15T SCREW, 9.5X105MM LLIAC SL-OPEN
8449512	ARM15T SCREW, 9.5X110MM ILIAC SL-OPEN
8449514	ARM15T SCREW, 9.5X120MM LLIAC SL-OPEN
8449560	ARM15T SCREW, 9.5X60MM LLIAC SL-OPEN
8449565	ARM15T SCREW, 9.5X65MM LLIAC SL-OPEN
8449570	ARM15T SCREW, 9.5X70MM LLIAC SL-OPEN
8449575	ARM15T SCREW, 9.5X75MM LLIAC SL-OPEN
8449580	ARM15T SCREW, 9.5X80MM LLIAC SL-OPEN
8449585	ARM15T SCREW, 9.5X85MM I LIAC SL-OPEN
8449590	ARM15T SCREW, 9.5X90MM LLIAC SL-OPEN
8449595	ARM15T SCREW, 9.5X95MM LLIAC SL-OPEN
8450120	ARMADA TI CONNECTOR. 20MM CROSS FIXED
8450122	ARMADA TI CONNECTOR. 22.5MM CROSS FIXED
8450125	ARMADA TI CONNECTOR. 25MM CROSS FIXED
8450127	ARMADA TI CONNECTOR. 27.5MM CROSS FIXED
8450130	ARMADA TI CONNECTOR. 30MM CROSS FIXED
8450132	ARMADA TI CONNECTOR, 32.5MM CROSS FIXED
8450135	ARMADA TI CONNECTOR, 35MM CROSS FIXED
8450137	ARMADA TI CONNECTOR, 37.5MM CROSS FIXED
8450140	ARMADA TI CONNECTOR, 40MM CROSS FIXED
8450145	ARMADA TI CONNECTOR, 45MM CROSS FIXED
8450150	ARMADA TI CONNECTOR. 50MM CROSS FIXED
8450155	ARMADA TI CONNECTOR, 55MM CROSS FIXED
8450160	ARMADA TI CONNECTOR, 60MM CROSS FIXED
8450165	ARMADA TI CONNECTOR, 65MM CROSS FIXED
8450170	ARMADA TI CROSS CONNECTOR, 70MM FIXED
8450175	ARMADA TI CONNECTOR, 75MM CROSS FIXED
8450180	ARMADA TI CROSS CONNECTOR, 80MM FIXED
8450245	ARM15T CROSS CONNECTOR, 45-50MM
8450250	ARM15T CROSS CONNECTOR, 50-60MM
8450260	ARM15T CROSS CONNECTOR. 60-70MM
8450270	ARM15T CROSS CONNECTOR. 70-85MM
8450302	ARM15T CONNECTOR. 5.5MM IN-LINE ROD
8450325	ARM15T LP ADJ. CROSS CONNECTOR 25MM



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8450330	ARM15T LP ADJ. CROSS CONNECTOR 30MM
8450335	ARM15T LP ADJ. CROSS CONNECTOR 35MM
8450340	ARM15T LP ADJ. CRASS CONNECTOR 40MM
8450350	ARM15T LP ADJ. CROSS CONNECTOR 50MM
8450365	ARM15T LP ADJ. CROSS CONNECTOR 55MM
8450402	ARM15T CONNECTOR. 5.5MM PARALLEL ROD
8451100	ARM15T LOCK SCREW
8451302	ASF15T CONNECTOR, IN-LANE ROD
8451402	ARM15T CONNECTOR. PARALLEL ROD
8451530	ARM15T SCREW, 10.5X30MM POLYAXIAL
8451535	ARM15T SCREW, 10.5X35MM POLYAXIAL
8451540	ARM15T SCREW, 10.5X40MM POLYAXIAL
8451545	ARM15T SCREW, 10.5X45MM POLYAXIAL
8451550	ARM15T SCREW, 10.5X50MM POLYAXIAL
8451555	ARM15T SCREW, 10.5X55MM POLYAXIAL
8451602	ARM15T CONNECTOR, 5.5MM 4-HOLE ROD
8451701	ASF ROD-ROD. 5.5MM OPEN 6.35MM OPEN
8452004	ASF15T LATERAL OFFSET ROD
8452025	ARM15T ROD, 25MM
8452030	ARM15T ROD, 30MM
8452035	ARM15T ROD. 35MM
8452040	ARM15T ROD 40MM
8452045	ARM15T ROD, 45MM
8452050	ARM15T ROD, 50MM
8452055	ARM15T ROD. 55MM
8452060	ARM15T ROD. 60MM
8452065	ARM15T ROD, 65MM
8452070	ARM15T ROD. 70MM
8452075	ARM15T ROD. 75MM
8452080	ARM15T ROD. 80MM
8452085	ARM15T ROD. 85MM
8452090	ARM15T ROD 90MM
8452095	ARM15T ROD. 95MM
8452100	ARM15T ROD, 100MM
8452120	ARM15T ROD. 120MM
8452140	ARM15T ROD, 140MM
8452152	ASF15T ROD-ROD, 5.5MM O 5.5MM O
8452156	ASF15T ROD-ROD, 5.5MM O 6.35MM O
8452166	ASF15T ROD-ROD, 6.35MM O 6.35MM O



República de Colombia  
Ministerio de Salud y Protección Social  
Instituto Nacional de Vigilancia de Medicamentos y Alimentos – INVIMA

**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**  
**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
8452200	ARM15T ROD, 200MM
8452250	ARM15T ROD, 250MM
8452255	ASF15T ROD-ROD, 5.5MM O 5.5MM O WIDE
8452256	ASF15T ROD-ROD. 5.5MM O 6.35MM O WIDE
8452260	ARM15T ROD, 260MM
8452266	ASF15T ROD-ROD, 6.35MM O 6.35MM O WIDE
8452300	ARM15T ROD. 300MM
8452355	ASF ROD-ROD. 5.5MMOPEN 5.5MMCLOSED 90°
8452356	ASF ROD-ROD. 5.5MMOPEN 6.35MMCLOSED 90°
8452365	ASF ROD-ROD. 6.35MMOPEN 5.5MMCLOSED 90°
8452366	ASF ROD-ROD. 6.35MMOPEN 6.35MMCLOSED 90°
8452412	ASF ROD CONNECTOR. 5.5MM 15X200MM LH
8452422	ASF ROD CONNECTOR, 5.5MM 25X200MM LH
8452425	ASF25T CONN 2-HOIE INLINE 5.5-5.5
8452426	ASF26T CONN 2 HOLE LNLIN 6.35-5.5
8452432	ASF ROD CONNECTOR, 5.5MM 35X200MM LH
8452442	ASF ROD CONNECTOR, 5.5MM 40X200MM LH
8452445	ASF25T CONN 4-HOLE INLINE 5.5-5.5
8452446	ASF26T CONN 4-HOLE INLINE 6.35-5.5
8452500	ARM15T ROD 500MM
8452512	ASF ROD CONNECTOR 5.5MM 15X200MM RH
8452522	ASF ROD CONNECTOR 5.5MM 25X200MM RH
8452532	ASF ROD CONNECTOR 5.5MM 35X200MM RH
8452542	ASF ROD CONNECTOR 5.5MM 40X200MM RH
8452555	ASF CONN ROD-ROD 5.5-O TO 5.5MM-O
8452556	ASF CONN ROD-ROD 6.35-O TO 5.5MM-O
8452600	ARM15T ROD. 600MM
8452755	ASF CONN ROD-ROD 5.5-O TO 5.5MM-C
8452766	ASF CONN ROD-ROD 6.35-O TO 6.35MM-C
8452855	ASF25T CONN ROD-ROD 5.5H-5.5O
8452856	ASF26T CONN ROD-ROD 6.35H-5.5O
8564525	ARM16S SCREW 4.5X25MM POLYAXIAL
8564530	ARM16S SCREW 4.5X30MM POLYAXIAL
8564535	ARM16S SCREW 4.5X35MM POLYAXIAL
8564540	ARM16S SCREW 4.5X40MM POLYAXIAL
8564545	ARM16S SCREW 4.5X45MM POLYAXIAL
8565530	ARM16S SCREW 5.5X30MM POLYAXIAL
8565535	ARM16S SCREW 5.5X35MM POLYAXIAL
8565540	ARM16S SCREW 5.5X40MM POLYAXIAL



República de Colombia  
Ministerio de Salud y Protección Social  
Instituto Nacional de Vigilancia de Medicamentos y Alimentos – INVIMA

**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**  
**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8565545	ARM16S SCREW 5.5X45MM POLYAXIAL
8565550	ARM16S SCREW 5.5X50MM POLYAXIAL
8565555	ARM16S SCREW 5.5X55MM POLYAXIAL
8566535	ARM16S SCREW 6.5X35MM POLYAXIAL
8566540	ARM16S SCREW 6.5X40MM POLYAXIAL
8566545	ARM16S SCREW 6.5X45MM POLYAXIAL
8566550	ARM16S SCREW 6.5X50MM POLYAXIAL
8566555	ARM16S SCREW 6.5X55MM POLYAXIAL
8567535	ARM16S SCREW 7.5X35MM POLYAXIAL
8567540	ARM16S SCREW 7.5X40MM POLYAXIAL
8567545	ARM16S SCREW 7.5X45MM POLYAXIAL
8567550	ARM16S SCREW 7.5X50MM POLYAXIAL
8567555	ARM16S SCREW 7.5X55MM POLYAXIAL
8567570	ARM16S SCREW 7.5X70MM POLYAXIAL
8567580	ARM16S SCREW 7.5X80MM POLYAXIAL
8567590	ARM16S SCREW 7.5X90MM POLYAXIAL
8568535	ARM16S SCREW 8.5X35MM POLYAXIAL
8568540	ARM16S SCREW 8.5X40MM POLYAXIAL
8568545	ARM16S SCREW 8.5X45MM POLYAXIAL
8568550	ARM16S SCREW 8.5X50MM POLYAXIAL
8568580	ARM16S SCREW 8.5X80MM POLYAXIAL
8568590	ARM16S SCREW 8.5X90MM POLYAXIAL
8576510	ARM16S SCREW, 6.5X100MM LLIAC SL-CLOSED
8576570	ARM16S SCREW, 6.5X70MM LLIAC SL-CLOSED
8576580	ARM16S SCREW, 6.5X80MM LLIAC SL-CLOSED
8576590	ARM16S SCREW, 6.5X90MM LLIAC SL-CLOSED
8577510	ARM16S SCREW, 7.5X100MM LLIAC SL-CLOSED
8577570	ARM16S SCREW, 7.5X70MM ILIAC SL-CLOSED
8577580	ARM16S SCREW, 7.5X80MM ILIAC SL-CLOSED
8577590	ARM16S SCREW, 7.5X90MM LLIAC SL-CLOSED
8578510	ARM16S SCREW, 8.5X100MM LLIAC SL-CLOSED
8578570	ARM16S SCREW, 8.5X70MM ILIAC SL-CLOSED
8578580	ARM16S SCREW, 8.5X80MM LLIAC SL-CLOSED
8578590	ARM16S SCREW, 8.5X90MM LLIAC SL-CLOSED
8634035	ARM15T SCREW, 4.0X35MM RED UNIPLANAR
8634040	ARM15T SCREW, 4.0X40MM RED UNIPLANAR
8634045	ARM15T SCREW, 4.0X45MM RED UNIPLANAR
8634525	ARM15T SCREW, 4.5X25MM RED UNIPLANAR
8634530	ARM15T SCREW, 4.5X30MM RED UNIPLANAR



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8634535	ARM15T SCREW, 4.5X35MM RED UNIPLANAR
8634540	ARM15T SCREW, 4.5X40MM RED UNIPLANAR
8634545	ARM15T SCREW, 4.5X45MM RED UNIPLANAR
8635025	ARM15T SCREW, 5.0X25MM RED UNIPLANAR
8635030	ARM15T SCREW, 5.0X30MM RED UNIPLANAR
8635035	ARM15T SCREW, 5.0X35MM RED UNIPLANAR
8635040	ARM15T SCREW, 5.0X40MM RED UNIPLANAR
8635045	ARM15T SCREW, 5.0X45MM RED UNIPLANAR
8635050	ARM15T SCREW, 5.0X50MM RED UNIPLANAR
8635530	ARM15T SCREW, 5.5X30MM RED UNIPLANAR
8635535	ARM15T SCREW, 5.5X35MM RED UNIPLANAR
8635540	ARM15T SCREW, 5.5X40MM RED UNIPLANAR
8635545	ARM15T SCREW, 5.5X45MM RED UNIPLANAR
8635550	ARM15T SCREW, 5.5X50MM RED UNIPLANAR
8635555	ARM15T SCREW, 5.5X55MM RED UNIPLANAR
8636030	ARM15T SCREW, 6.0X30MM RED UNIPLANAR
8636035	ARM15T SCREW, 6.0X35MM RED UNIPLANAR
8636040	ARM15T SCREW, 6.0X40MM RED UNIPLANAR
8636045	ARM15T SCREW, 6.0X45MM RED UNIPLANAR
8636050	ARM15T SCREW, 6.0X50MM RED UNIPLANAR
8636055	ARM15T SCREW, 6.0X55MM RED UNIPLANAR
8644025	ARM15T SCREW, 4X25MM UNIPLANAR
8644030	ARM15T SCREW, 4X30MM UNIPLANAR
8644035	ARM15T SCREW, 4X35MM UNIPLANAR
8644040	ARM15T SCREW, 4X40MM UNIPLANAR
8644525	ARM15T SCREW, 4.5X25MM UNIPLANAR
8644530	ARM15T SCREW, 4.5X30MM UNIPLANAR
8644535	ARM15T SCREW, 4.5X35MM UNIPLANAR
8644540	ARM15T SCREW, 4.5X40MM UNIPLANAR
8644545	ARM15T SCREW, 4.5X45MM UNIPLANAR
8645025	ARM15T SCREW, 5X25MM UNIPLANAR
8645030	ARM15T SCREW, 5X30MM UNIPLANAR
8645035	ARM15T SCREW, 5X35MM UNIPLANAR
8645040	ARM15T SCREW, 5X40MM UNIPLANAR
8645045	ARM15T SCREW, 5X45MM UNIPLANAR
8645050	ARM15T SCREW, 5X50MM UNIPLANAR
8645530	ARM15T SCREW, 5.5X30MM UNIPLANAR
8646030	ARM15T SCREW, 6X30MM UNIPLANAR
8646035	ARM15T SCREW, 6X35MM UNIPLANAR



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8646040	ARM15T SCREW, 6X40MM UNIPLANAR
8646045	ARM15T SCREW, 6X45MM UNIPLANAR
8646050	ARM15T SCREW, 6X50MM UNIPLANAR
8646055	ARM15T SCREW, 6X55MM UNIPLANAR
8650520	ARM15T CONNECTOR. 20MM OPEN OFFSET
8650530	ARM15T CONNECTOR. 30MM OPEN OFFSET
8650560	ARM15T CONNECTOR, 60MM OPEN OFFSET
8651520	ARM15T CONNECTOR. 20MM CLOSED OFFSET
8651530	ARM15T CONNECTOR, 30MM CLOSED OFFSET
8651540	ARM15T CONNECTOR. 40MM CLOSED OFFSET
8651560	ARM15T CONNECTOR. 60MM CLOSED OFFSET
8653530	ARM15T CONN. 30MM LEFT OPEN ILIAC OFFSET
8653560	ARM15T CONN, 60MM LEFT OPEN ILIAC OFFSET
8653630	ARM15T CONN, 30MM RT OPEN ILIAC OFFSET
8653660	ARM15T CONN, 60MM RT OPEN ILIAC OFFSET
8654530	ARM15T SCREW, 4.5X30MM UNIAXIAL
8656540	ARM15T SCREW, 6.5X40MM UNIAXIAL
8745025	ARM15S SCREW, 5.0X25MM RED UNIPLANAR
8745030	ARM15S SCREW, 5.0X30MM RED UNIPLANAR
8745035	ARM15S SCREW, 5.0X35MM RED UNIPLANAR
8746030	ARM15S SCREW, 6.0X30MM RED UNIPLANAR
8746035	ARM15S SCREW, 6.0X35MM RED UNIPLANAR
8746040	ARM15S SCREW, 6.0X40MM RED UNIPLANAR
8746045	ARM15S SCREW, 6.0X45MM RED UNIPLANAR
8750520	ARM15S CONNECTOR, 20MM OPEN OFFSET
8750530	ARM15S CONNECTOR. 30MM OPEN OFFSET
8750560	ARM15S CONNECTOR, 60MM OPEN OFFSET
8751520	ARM15S CONNECTOR. 20MM CLOSED OFFSET
8751530	ARM15S CONNECTOR, 30MM CLOSED OFFSET
8751560	ARM15S CONNECTOR. 60MM CLOSED OFFSEL
8754525	ARM15S SCREW, 4.5X25MM UNIPLANAR
8754530	ARM15S SCREW, 4.5X30MM UNIPLANAR
8754535	ARM15S SCREW, 4.5X35MM UNIPLANAR
8754540	ARM15S SCREW, 4.5X40MM UNIPLANAR
8754545	ARM15S SCREW, 4.5X45MM UNIPLANAR
8755025	ARM15S SCREW, 5X25MM UNIAXIAL
8755030	ARM15S SCREW, 5X30MM UNIAXIAL
8755035	ARM15S SCREW, 5X35MM UNIAXIAL
8755040	ARM15S SCREW, 5X40MM UNIAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8755045	ARM15S SCREW, 5X45MM UNIAXIAL
8755050	ARM15S SCREW, 5X50MM UNIAXIAL
8756030	ARM15S SCREW, 6X30MM UNIAXIAL
8756035	ARM15S SCREW, 6X35MM UNIAXIAL
8756040	ARM15S SCREW, 6X40MM UNIAXIAL
8756045	ARM15S SCREW, 6X45MM UNIAXIAL
8756050	ARM15S SCREW, 6X50MM UNIAXIAL
8756055	ARM15S SCREW, 6X55MM UNIAXIAL
8760520	ARM16S CONNECTOR 20MM OPEN ONSET
8760530	ARM16S CONNECTOR 30MM OPEN OFFSET
8760560	ARM16S CONNECTOR 60MM OPEN OFFSET
8854025	ARM15T SCREW, 4X25MM FIXED
8854030	ARM15T SCREW, 4X30MM FIXED
8854035	ARM15T SCREW, 4X35MM FIXED
8854040	ARM15T SCREW, 4X40MM FIXED
8854045	ARM15T SCREW, 4X45MM FIXED
8854525	ARM15T SCREW, 4.5X25MM FIXED
8854530	ARM15T SCREW, 4.5X30MM FIXED
8854535	ARM15T SCREW, 4.5X35MM FIXED
8854540	ARM15T SCREW, 4.5X40MM FIXED
8854545	ARM15T SCREW, 4.5X45MM FIXED
8854550	ARM15T SCREW, 4.5X50MM FIXED
8855025	ARM15T SCREW, 5X25MM FIXED
8855030	ARM15T SCREW, 5X30MM FIXED
8855035	ARM15T SCREW, 5X35MM FIXED
8855040	ARM15T SCREW, 5X40MM FIXED
8855045	ARM15T SCREW, 5X45MM FIXED
8855050	ARM15T SCREW, 5X50MM FIXED
8855525	ARM15T SCREW, 5.5X25MM FIXED
8855530	ARM15T SCREW, 5.5X30MM FIXED
8855535	ARM15T SCREW, 5.5X35MM FIXED
8855540	ARM15T SCREW, 5.5X40MM FIXED
8855545	ARM15T SCREW, 5.5X45MM FIXED
8855550	ARM15T SCREW, 5.5X50MM FIXED
8855555	ARM15T SCREW, 5.5X55MM FIXED
8855560	ARM15T SCREW, 5.5X60MM FIXED
8855565	ARM15T SCREW, 5.5X65MM FIXED
8856525	ARM15T SCREW, 6.5X25MM FIXED
8856530	ARM15T SCREW, 6.5X30MM FIXED



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8856535	ARM15T SCREW, 6.5X35MM FIXED
8856540	ARM15T SCREW, 6.5X40MM FIXED
8856545	ARM15T SCREW, 6.5X45MM FIXED
8856550	ARM15T SCREW, 6.5X50MM FIXED
8856555	ARM15T SCREW, 6.5X55MM FIXED
8856560	ARM15T SCREW, 6.5X60MM FIXED
8856565	ARM15T SCREW, 6.5X65MM FIXED
8857525	ARM15T SCREW, 7.5X25MM FIXED
8857530	ARM15T SCREW, 7.5X30MM FIXED
8857535	ARM15T SCREW, 7.5X35MM FIXED
8857540	ARM15T SCREW, 7.5X40MM FIXED
8857545	ARM15T SCREW, 7.5X45MM FIXED
8857550	ARM15T SCREW, 7.5X50MM FIXED
8857555	ARM15T SCREW, 7.5X55MM FIXED
8857560	ARM15T SCREW, 7.5X60MM FIXED
8857565	ARM15T SCREW, 7.5X65MM FIXED
8858535	ARM15T SCREW, 8.5X35MM FIXED
8954025	ARM15S SCREW, 4X25MM FIXED
8954030	ARM15S SCREW, 4X30MM FIXED
8954035	ARM15S SCREW, 4X35MM FIXED
8954040	ARM15S SCREW, 4X40MM FIXED
8954045	ARM15S SCREW, 4X45MM FIXED
8954525	ARM15S SCREW, 4.5X25MM FIXED
8954530	ARM15S SCREW, 4.5X30MM FIXED
8954535	ARM15S SCREW, 4.5X35MM FIXED
8954540	ARM15S SCREW, 4.5X40MM FIXED
8954545	ARM15S SCREW, 4.5X45MM FIXED
8955025	ARM15S SCREW, 5X25MM FIXED
8955030	ARM15S SCREW, 5X30MM FIXED
8955035	ARM15S SCREW, 5X35MM FIXED
8955040	ARM15S SCREW, 5X40MM FIXED
8955045	ARM15S SCREW, 5X45MM FIXED
8955050	ARM15S SCREW, 5X50MM FIXED
8955525	ARM15S SCREW, 5.5X25MM FIXED
8955530	ARM15S SCREW, 5.5X30MM FIXED
8955535	ARM15S SCREW, 5.5X35MM FIXED
8955540	ARM15S SCREW, 5.5X40MM FIXED
8955545	ARM15S SCREW, 5.5X45MM FIXED
8955550	ARM15S SCREW, 5.5X50MM FIXED



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8955555	ARM15S SCREW, 5.5X55MM FIXED
8956530	ARM15S SCREW, 6.5X30MM FIXED
8956535	ARM15S SCREW, 6.5X35MM FIXED
8956540	ARM15S SCREW, 6.5X40MM FIXED
8956545	ARM15S SCREW, 6.5X45MM FIXED
8956550	ARM15S SCREW, 6.5X50MM FIXED
8956555	ARM15S SCREW, 6.5X55MM FIXED
8957530	ARM15S SCREW, 7.5X30MM FIXED
8957535	ARM15S SCREW, 7.5X35MM FIXED
8957540	ARM15S SCREW, 7.5X40MM FIXED
8957545	ARM15S SCREW, 7.5X45MM FIXED
8957550	ARM15S SCREW, 7.5X50MM FIXED
8957555	ARM15S SCREW, 7.5X55MM FIXED
8964525	ARM16S SCREW, 4.5X25MM FIXED
8964530	ARM16S SCREW, 4.5X30MM FIXED
8964535	ARM16S SCREW, 4.5X35MM FIXED
8964540	ARM16S SCREW, 4.5X40MM FIXED
8964545	ARM16S SCREW, 4.5X45MM FIXED
8965530	ARM16S SCREW, 5.5X30MM FIXED
8965535	ARM16S SCREW, 5.5X35MM FIXED
8965540	ARM16S SCREW, 5.5X40MM FIXED
8965545	ARM16S SCREW, 5.5X45MM FIXED
8966535	ARM16S SCREW, 6.5X35MM FIXED
8966540	ARM16S SCREW, 6.5X40MM FIXED
8966545	ARM16S SCREW, 6.5X45MM FIXED
8966550	ARM16S SCREW, 6.5X50MM FIXED
8967535	ARM16S SCREW, 7.5X35MM FIXED
8967540	ARM16S SCREW, 7.5X40MM FIXED
8967545	ARM16S SCREW, 7.5X45MM FIXED
8967550	ARM16S SCREW, 7.5X50MM FIXED
8052025006	ARMADA TI ROD, 5.5X25MM CONVENTIONAL PB
8052030006	ARMADA TI ROD, 5.5X30MM CONVENTIONAL PB
8052035006	ARMADA TI ROD, 5.5X35MM CONVENTIONAL PB
8052040006	ARMADA TI ROD, 5.5X40MM CONVENTIONAL PB
8052045006	ARMADA TI ROD, 5.5X45MM CONVENTIONAL PB
8052050006	ARMADA TI ROD, 5.5X50MM CONVENTIONAL PB
8052055006	ARMADA TI ROD, 5.5X55MM CONVENTIONAL PB
8052060006	ARMADA TI ROD. 5.5X60MM CONVENTIONAL PB
8052065006	ARMADA TI ROD. 5.5X65MM CONVENTIONAL PB



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8052070006	ARMADA TI ROD, 5.5X70MM CONVENTIONAL PB
8052075006	ARMADA TI ROD, 5.5X75MM CONVENTIONAL PB
8052080006	ARMADA TI ROD, 5.5X80MM CONVENTIONAL PB
8052085006	ARMADA TI ROD. 5.5X85MM CONVENTIONAL PB
8052090006	ARMADA TI ROD, 5.5X90MM CONVENTIONAL PB
8052095006	ARMADA TI ROD, 5.5X95MM CONVENTIONAL PB
8052100006	ARMADA TI ROD, 5.5X100MM CONVENTIONAL PB
8252500006	ARM15T COCR ROD, 500MM STRAIGHT
8254525006	ARMADA TI SCREW, 4.5X25MM PLS
8254530006	ARMADA TI SCREW, 4.5X30MM PLS
8254535006	ARMADA TI SCREW, 4.5X35MM PLS
8254540006	ARMADA TI SCREW, 4.5X40MM PLS
8254545006	ARMADA TI SCREW, 4.5X45MM PLS
8255525006	ARMADA TI SCREW, 5.5X25MM PLS
8255530006	ARMADA TI SCREW, 5.5X30MM PLS
8255535006	ARMADA TI SCREW, 5.5X35MM PLS
8255540006	ARMADA TI SCREW, 5.5X40MM PLS
8255545006	ARMADA TI SCREW, 5.5X45MM PLS
8256530006	ARMADA TI SCREW, 6.5X30MM PLS
8256535006	ARMADA TI SCREW, 6.5X35MM PLS
8256540006	ARMADA TI SCREW, 6.5X40MM PLS
8256545006	ARMADA TI SCREW, 6.5X45MM PLS
8256550006	ARMADA TI SCREW, 6.5X50MM PLS
8256555006	ARMADA TI SCREW, 6.5X55MM PLS
8274525006	ARM15T SCREW, 4.5X25MM REDUCTION
8274530006	ARM15T SCREW, 4.5X30MM REDUCTION
8275525006	ARM15T SCREW, 5.5X25MM REDUCTION
8275530006	ARM15T SCREW, 5.5X30MM REDUCTION
8275535006	ARM15T SCREW, 5.5X35MM REDUCTION
8275540006	ARM15T SCREW, 5.5X40MM REDUCTION
8275545006	ARM15T SCREW, 5.5X45MM REDUCTION
8275550006	ARM15T SCREW, 5.5X50MM REDUCTION
8275555006	ARM15T SCREW, 5.5X55MM REDUCTION
8276530006	ARM15T SCREW, 6.5X30MM REDUCTION
8276535006	ARN15T SCREW, 6.5X35MM REDUCTION
8276540006	ARM15T SCREW, 6.5X40MM REDUCTION
8276545006	ARM15T SCREW, 6.5X45MM REDUCTION
8276550006	ARM15T SCREW, 6.5X50MM REDUCTION
8276555006	ARM15T SCREW, 6.5X55MM REDUCTION



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8277530006	ARM15T SCREW. 7.5X30MM REDUCTION
8277535006	ARM15T SCREW. 7.5X35MM REDUCTION
8277540006	ARM15T SCREW, 7.5X40MM REDUCTION
8277545006	ARM15T SCREW, 7.5X45MM REDUCTION
8277550006	ARM15T SCREW, 7.5X50MM REDUCTION
8277555006	ARM15T SCREW. 7.5X55MM REDUCTION
8405535001	ARMADA SCREW. 5.5X35MM POLY FENESTRATED
8405540001	ARMADA SCREW, 5.5X40MM POLY FENESTRATED
8405545001	ARMADA SCREW, 5.5X45MM POLY FENESTRATED
8405550001	ARMADA SCREW, 5.5X50MM POLY FENESTRATED
8406535001	ARMADA SCREW. 6.5X35MM POLY FENESTRATED
8406540001	ARMADA SCREW. 6.5X40MM POLY FENESTRATED
8406545001	ARMADA SCREW. 6.5X45MM POLY FENESTRATED
8406550001	ARMADA SCREW, 6.5X50MM POLY FENESTRATED
8406555001	ARMADA SCREW, 6.5X55MM POLY FENESTRATED
8407540001	ARMADA SCREW, 7.5X40MM POLY FENESTRATED
8407545001	ARMADA SCREW, 7.5X45MM POLY FENESTRATED
8407550001	ARMADA SCREW. 7.5X50MM POLY FENESTRATED
8407555001	ARMADA SCREW, 7.5X55MM POLY FENESTRATED
8408540001	ARMADA SCREW, 8.5X40MM POLY FENESTRATED
8408545001	ARMADA SCREW. 8.5X45MM POLY FENESTRATED
8408550001	ARMADA SCREW, 8.5X50MM POLY FENESTRATED
8408555001	ARMADA SCREW, 8.5X55MM POLY FENESTRATED
8408560001	ARMADA SCREW, 8.5X60MM POLY FENESTRATED
8450120006	ARMADA TI CROSS CONNECTOR, 20MM FIXED
8450122006	ARMADA TI CROSS CONNECTOR. 22.5MM FIXED
8450125006	ARMADA TI CROSS CONNECTOR, 25MM FIXED
8450127006	ARMADA TI CROSS CONNECTOR, 27.5MM FIXED
8450130006	ARMADA TI CROSS CONNECTOR, 30MM FIXED
8450132006	ARMADA TI CROSS CONNECTOR, 32.5MM FIXED
8450135006	ARMADA TI CROSS CONNECTOR, 35MM FIXED
8450137006	ARMADA TI CROSS CONNECTOR, 37.5MM FIXED
8450140006	ARMADA TI CROSS CONNECTOR, 40MM FIXED
8450145006	ARMADA TI CROSS CONNECTOR, 45MM FIXED
8450150006	ARMADA TI CROSS CONNECTOR, 50MM FIXED
8450155006	ARMADA TI CROSS CONNECTOR. 55MM FIXED
8450160006	ARMADA TI CROSS CONNECTOR, 60MM FIXED
8450245006	ARMADA TI CONNECTOR, 45-50MM ADJ TRANS
8450250006	ARMADA TI CONNECTOR, 50-60MM ADJ TRANS



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
8450260006	ARMADA TI CONNECTOR. 60-70MM ADJ TRANS
8451302006	ASF CONNECTOR. IN-LINE ROD
8451402006	ASF CONNECTOR, PARALLEL ROD
8451602006	ASF CONNECTOR, 4-HOLE ROD
8452025006	ARMADA TI ROD. 5.5X25MM
8452030006	ARMADA TI ROD, 5.5X30MM
8452035006	ARMADA TI ROD, 5.5X35MM
8452040006	ARMADA TI ROD, 5.5X40MM
8452045006	ARMADA TI ROD, 5.5X45MM
8452050006	ARMADA TI ROD. 5.5X50MM
8452055006	ARMADA TI ROD, 5.5X55MM
8452060006	ARMADA TI ROD, 5.5X60MM
8452065006	ARMADA TI ROD. 5.5X65MM
8452070006	ARMADA TI ROD, 5.5X70MM
8452080006	ARMADA TI ROD, 5.5X80MM
8452090006	ARMADA TI ROD. 5.5X90MM
8452100006	ARMADA TI ROD. 5.5X100MM
8452140066	ARMADA TI ROD, 5.5X140MM
8452160006	ARMADA TI ROD. 5.5X160MM
8452200066	ARMADA TI ROD, 5.5X200MM
8452300006	ARMADA TI ROD. 5.5X300MM
8452500006	ARMADA TI ROD. 5.5X500MM
8454025006	ARM15T SCREW, 4.0X25MM POLYAXIAL
8454030006	ARM15T SCREW, 4.0X30MM POLYAXIAL
8454035006	ARM15T SCREW, 4.0X35MM POLYAXIAL
8454040006	ARM15T SCREW. 4.0X40MM POLYAXIAL
8454525006	ARM15T SCREW. 4.5X25MM POLYAXIAL
8454530005	ARM15T SCREW, 4.5X30MM POLYAXIAL
8454535006	ARM15T SCREW, 4.5X35MM POLYAXIAL
8454540006	ARM15T SCREW. 4.5X40MM POLYAXIAL
8454545006	ARM15T SCREW. 4.5X45MM POLYAXIAL
8455025006	ARM15T SCREW, 5.0X25MM POLYAXIAL
8455030006	ARM15T SCREW, 5.0X30MM POLYAXIAL
8455035006	ARM15T SCREW. 5.0X35MM POLYAXIAL
8455040006	ARM15T SCREW, 5.0X40MM POLYAXIAL
8455045006	ARM15T SCREW. 5.0X45MM POLYAXIAL
8455050006	ARM15T SCREW. 5.0X50MM POLYAXIAL
8455525006	ARM15T SCREW, 5.5X25MM POLYAXIAL
8455530006	ARM15T SCREW, 5.5X30MM POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8455535006	ARM15T SCREW, 5.5X35MM POLYAXIAL
8455540006	ARM15T SCREW, 5.5X40MM POLYAXIAL
8455545006	ARM15T SCREW. 5.5X45MM POLYAXIAL
8455550006	ARM15T SCREW, 5.5X50MM POLYAXIAL
8455555006	ARM15T SCREW, 5.5X55MM POLYAXIAL
8456510006	ARM15T SCREW. 6.5X100MM POLYAXIAL
8456530006	ARM15T SCREW, 6.5X30MM POLYAXIAL
8456535006	ARM15T SCREW, 6.5X35MM POLYAXIAL
8456540006	ARM15T SCREW, 6.5X40MM POLYAXIAL
8456545006	ARM15T SCREW. 6.5X45MM POLYAXIAL
8456550006	ARM15T SCREW, 6.5X50MM POLYAXIAL
8456555006	ARM15T SCREW. 6.5X55MM POLYAXIAL
8456560006	ARM15T SCREW, 6.5X60MM POLYAXIAL
8456570006	ARM15T SCREW, 6.5X70MM POLYAXIAL
8456580006	ARM15T SCREW. 6.5X80MM POLYAXIAL
8456590006	ARM15T SCREW. 6.5X90MM POLYAXIAL
8457510006	ARM15T SCREW, 7.5X100MM POLYAXIAL
8457530006	ARM15T SCREW. 7.5X30MM POLYAXIAL
8457535006	ARM15T SCREW, 7.5X35MM POLYAXIAL
8457540006	ARM15T SCREW, 7.5X40MM POLYAXIAL
8457545006	ARM15T SCREW, 7.5X45MM POLYAXIAL
8457550006	ARM15T SCREW, 7.5X50MM POLYAXIAL
8457555006	ARM15T SCREW, 7.5X55MM POLYAXIAL
8457560006	ARM15T SCREW. 7.5X60MM POLYAXIAL
8457570006	ARM15T SCREW. 7.5X70MM POLYAXIAL
8457580006	ARM15T SCREW, 7.5X80MM POLYAXIAL
8457590006	ARM15T SCREW. 7.5X90MM POLYAXIAL
8458510006	ARM15T SCREW, 8.5X100MM POLYAXIAL
8458530006	ARM15T SCREW. 8.5X30MM POLYAXIAL
8458535006	ARM15T SCREW. 8.5X35MM POLYAXIAL
8458540006	ARM15T SCREW. 8.5X40MM POLYAXIAL
8458545006	ARM15T SCREW, 8.5X45MM POLYAXIAL
8458550006	ARM15T SCREW, 8.5X50MM POLYAXIAL
8458555006	ARM15T SCREW. 8.5X55MM POLYAXIAL
8458560006	ARM15T SCREW, 8.5X60MM POLYAXIAL
8458570006	ARM15T SCREW, 8.5X70MM POLYAXIAL
8458580006	ARM15T SCREW, 8.5X80MM POLYAXIAL
8458590006	ARM15T SCREW, 8.5X90MM POLYAXIAL
8459510006	ARM15T SCREW. 9.5X100MM POLYAXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8459560006	ARM15T SCREW, 9.5X60MM POLYAXIAL
8459570006	ARM15T SCREW, 9.5X70MM POLYAXIAL
8459580006	ARM15T SCREW, 9.5X80MM POLYAXIAL
8459590006	ARM15T SCREW, 9.5X90MM POLYAXIAL
8461100006	ARMADA TI LOCK SCREW, OPEN
8464402006	ASF CONNECTOR, PARALLEL ROD
8644525006	ARMADA TI SCREW, 4.5X25MM UNIPLANAR
8644530006	ARMADA TI SCREW, 4.5X30MM UNIPLANAR
8644535006	ARMADA TI SCREW, 4.5X35MM UNIPLANAR
8644540006	ARMADA TI SCREW, 4.5X40MM UNIPLANAR
8644545006	ARMADA TI SCREW, 4.5X45MM UNIPLANAR
8645525006	ARMADA TI SCREW, 5.5X25MM UNIPLANAR
8645530006	ARMADA TI SCREW, 5.5X30MM UNIPLANAR
8645535006	ARMADA TI SCREW, 5.5X35MM UNIPLANAR
8645540006	ARMADA TI SCREW, 5.5X40MM UNIPLANAR
8645545006	ARMADA TI SCREW, 5.5X45MM UNIPLANAR
8646030006	ARMADA TI SCREW, 6.0X30MM UNIPLANAR
8646035006	ARMADA TI SCREW, 6.0X35MM UNIPLANAR
8646040006	ARMADA TI SCREW, 6.0X40MM UNIPLANAR
8646045006	ARMADA TI SCREW, 6.0X45MM UNIPLANAR
8646050006	ARMADA TI SCREW, 6.0X50MM UNIPLANAR
8646055006	ARMADA TI SCREW, 6.0X55MM UNIPLANAR
8650520006	ARM15T CONNECTOR. 20MM OFFSET OPEN
8650530006	ARM15T CONNECTOR, 30MM OFFSET OPEN
8650560006	ARM15T CONNECTOR, 60MM OFFSET OPEN
8651520006	ARMADA TI CONNECTOR, 20MM OFFSET CLOSED
8651530006	ARMADA TI CONNECTOR, 30MM OFFSET CLOSED
8651560006	ARMADA TI CONNECTOR, 60MM OFFSET CLOSED
8854025006	ARM15T SCREW, 4.0X25MM FIXED
8854030006	ARM15T SCREW, 4.0X30MM FIXED
8854035006	ARM15T SCREW, 4.0X35MM FIXED
8854040006	ARM15T SCREW. 4.0X40MM FIXED
8854045006	ARM15T SCREW, 4.0X45MM FIXED
8854525006	ARM15T SCREW, 4.5X25MM FIXED
8854530006	ARM15T SCREW. 4.5X30MM FIXED
8854535006	ARM15T SCREW. 4.5X35MM FIXED
8854540006	ARM15T SCREW, 4.5X40MM FIXED
8854545006	ARM15T SCREW, 4.5X45MM FIXED
8854550006	ARM15T SCREW, 4.5X50MM FIXED



República de Colombia  
Ministerio de Salud y Protección Social  
Instituto Nacional de Vigilancia de Medicamentos y Alimentos – INVIMA

**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**  
**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8855025006	ARM15T SCREW. 5.0X25MM FIXED
8855030006	ARM15T SCREW, 5.0X30MM FIXED
8855035006	ARM15T SCREW, 5.0X35MM FIXED
8855040006	ARM15T SCREW. 5.0X40MM FIXED
8855045006	ARM15T SCREW. 5.0X45MM FIXED
8855050006	ARM15T SCREW. 5.0X50MM FIXED
8855525066	ARM15T SCREW. 5.5X25MM FIXED
8855530006	ARM15T SCREW, 5.5X30MM FIXED
8855535006	ARM15T SCREW. 5.5X35MM FIXED
8855540006	ARM15T SCREW, 5.5X40MM FIXED
8855545006	ARM15T SCREW. 5.5X45MM FIXED
8855550006	ARM15T SCREW, 5.5X50MM FIXED
8855555006	ARM15T SCREW, 5.5X55MM FIXED
8855560006	ARM15T SCREW. 5.5X60MM FIXED
8856530006	ARM15T SCREW, 6.5X30MM FIXED
8856535006	ARM15T SCREW, 6.5X35MM FIXED
8856540006	ARM15T SCREW, 6.5X40MM FIXED
8856545006	ARM15T SCREW, 6.5X45MM FIXED
8856550006	ARM15T SCREW. 6.5X50MM FIXED
8856555006	ARM15T SCREW, 6.5X55MM FIXED
8856560006	ARM15T SCREW. 6.5X60MM FIXED
8857530006	ARM15T SCREW, 7.5X30MM FIXED
8857535006	ARM15T SCREW, 7.5X35MM FIXED
8857540006	ARM15T SCREW. 7.5X40MM FIXED
8857545008	ARM15T SCREW, 7.5X45MM FIXED
8857550006	ARM15T SCREW. 7.5X50MM FIXED
8857555006	ARM15T SCREW, 7.5X55MM FIXED
8052025P	ARM15T ROD, 25MM CONVENTIONAL PRE-BENT
8052030P	ARM15T ROD, 30MM CONVENTIONAL PRE-BENT
8052035P	ARM15T ROD, 35MM CONVENTIONAL PRE-BENT
8052040P	ARM15T ROD. 40MM CONVENTIONAL PRE-BENT
8052045P	ARM15T ROD, 45MM CONVENTIONAL PRE-BENT
8052050P	ARM15T ROD, 50MM CONVENTIONAL PRE-BENT
8052055P	ARM15T ROD, 55MM CONVENTIONAL PRE-BENT
8052060P	ARM15T ROD. 60MM CONVENTIONAL PRE-BENT
8052065P	ARM15T ROD. 65MM CONVENTIONAL PRE-BENT
8052070P	ARM15T ROD, 70MM CONVENTIONAL PRE-BENT
8052075P	ARM15T ROD, 75MM CONVENTIONAL PRE-BENT
8052080P	ARM15T ROD, 80MM CONVENTIONAL PRE-BENT



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8052085P	ARM15T ROD, 85MM CONVENTIONAL PRE-BENT
8052090P	ARM15T ROD. 90MM CONVENTIONAL PRE-BENT
8052095P	ARM15T RED, 95MM CONVENTIONAL PRE-BENT
8052100P	ARM15T ROD, 100MM CONVENTIONAL PRE-BENT
8052110P	ARM15T ROD, 110MM CONVENTIONAL PRE-BENT
8052120P	ARM15T ROD. 120MM CONVENTIONAL PRE-BENT
8252300P	ARM15T COCR ROD, 300MM STRAIGHT
8252500P	ARM15T COCR ROD, 500MM STRAIGHT
8255030P	ARM15T SCREW, 5.0X30MM PLS
8255035P	ARM15T SCREW, 5.0X35MM PLS
8255040P	ARM15T SCREW, 5.0X40MM PLS
8255045P	ARM15T SCREW. 5.0X45MM PLS
8255525P	ARM15T SCREW, 5.5X25MM PLS
8255530P	ARM15T SCREW. 5.5X30MM PLS
8255535P	ARM15T SCREW, 5.5X35MM PLS
8255540P	ARM15T SCREW, 5.5X40MM PLS
8255545P	ARM15T SCREW, 5.5X45MM PLS
8256535P	ARM15T SCREW. 6.5X35MM PLS
8256540P	ARM15T SCREW. 6.5X40MM PLS
8256545P	ARM15T SCREW, 6.5X45MM PLS
8256550P	ARM15T SCREW, 6.5X50MM PLS
8256555P	ARM15T SCREW, 6.5X55MM PLS
8275530P	ARM15T SCREW. 5.5X30MM REDUCTION
8275535P	ARM15T SCREW, 5.5X35MM REDUCTION
8275540P	ARM15T SCREW, 5.5X40MM REDUCTION
8275545P	ARM14T SCREW, 5.5X45MM REDUCTION
8276530P	ARM15T SCREW, 6.5X30MM REDUCTION
8276535P	ARM15T SCREW, 6.5X35MM REDUCTION
8276540P	ARM15T SCREW, 6.5X40MM REDUCTION
8276545P	ARM15T SCREW, 6.5X45MM REDUCTION
8276550P	ARM15T SCREW, 6.5X50MM REDUCTION
8276555P	ARM15T SCREW, 6.5X55MM REDUCTION
8277540P	ARM15T SCREW, 7.5X40MM REDUCTION
8277545P	ARM15T SCREW, 7.5X45MM REDUCTION
8405535001P	ARMADA SCREW, 5.5X35MM POLY FENESTRATED
8405540001P	ARMADA SCREW, 5.5X40MM POLY FENESTRATED
8405545001P	ARMADA SCREW. 5.5X45MM POLY FENESTRATED
8406540001P	ARMADA SCREW, 6.5X40MM POLY FENESTRATED
8406545001P	ARMADA SCREW, 6.5X45MM POLY FENESTRATED



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
8406550001P	ARMADA SCREW, 6.5X50MM POLY FENESTRATED
8406555001P	ARMADA SCREW, 6.5X55MM POLY FENESTRATED
8407545001P	ARMADA SCREW, 7.5X45MM POLY FENESTRATED
8407550001P	ARMADA SCREW, 7.5X50MM POLY FENESTRATED
8407555001P	ARMADA SCREW, 7.5X55MM POLY FENESTRATED
8450122P	ARMADA TI CONNECTOR, 22.5MM CROSS FIXED
8450125P	ARMADA TI CONNECTOR, 25MM CROSS FIXED
8450127P	ARMADA TI CONNECTOR, 27.5MM CROSS FIXED
8450130P	ARMADA TI CONNECTOR. 30MM CROSS FIXED
8450132P	ARMADA TI CONNECTOR, 32.5MM CROSS FIXED
8450135P	ARMADA TI CONNECTOR. 35MM CROSS FIXED
8450137P	ARMADA TI CONNECTOR. 37.5MM CROSS FIXED
8450140P	ARMADA TI CONNECTOR, 40MM CROSS FIXED
8450145P	ARMADA TI CONNECTOR. 45MM CROSS FIXED
8450150P	ARMADA TI CONNECTOR. 50MM CROSS FIXED
8450155P	ARMADA TI CONNECTOR. 55MM CROSS FIXED
8450160P	ARMADA TI CONNECTOR, 60MM CROSS FIXED
8450245P	ARM15T CROSS CONNECTOR. 45-50MM
8450250P	ARM15T CROSS CONNECTOR. 50-60MM
8450260P	ARM15T CROSS CONNECTOR. 60-70MM
8451302P	ASF15T CONNECTOR, IN-LINE ROD
8451402P	ARM15T CONNECTOR. PARALLEL ROD
8452025P	ARM15T RED, 25MM
8452030P	ARM15T ROD, 30MM
8452035P	ARM15T ROD, 35MM
8452040P	ARM15T ROD, 40MM
8452045P	ARM15T ROD, 45MM
8452050P	ARM15T ROD, 50MM
8452055P	ARM15T ROD, 55MM
8452060P	ARM15T ROD. 60MM
8452065P	ARM15T ROD, 65MM
8452070P	ARM15T ROD, 70MM
8452080P	ARM15T ROD, 80MM
8452090P	ARM15T ROD, 90MM
8452100P	ARM15T ROD, 100MM
8452120P	ARM15T ROD. 120MM
8452200P	ARM15T ROD, 200MM
8452300P	ARM15T ROD, 300MM
8452500P	ARM15T ROD, 500MM



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8455530P	ARM15T SCREW, 5.5X30MM POLYAXIAL
8455535P	ARM15T SCREW, 5.5X35MM POLYAXIAL
8455540P	ARM15T SCREW, 5.5X40MM POLYAXIAL
8455545P	ARM15T SCREW, 5.5X45MM POLYAXIAL
8455550P	ARM15T SCREW, 5.5X50MM POLYAXIAL
8455555P	ARM15T SCREW. 5.5X55MM POLYAXIAL
8456530P	ARM15T SCREW, 6.5X30MM POLYAXIAL
8456535P	ARM15T SCREW, 6.5X35MM POLYAXIAL
8456540P	ARM15T SCREW, 6.5X40MM POLYAXIAL
8456545P	ARM15T SCREW. 6.5X45MM POLYAXIAL
8456550P	ARM15T SCREW. 6.5X50MM POLYAXIAL
8456555P	ARM15T SCREW, 6.5X55MM POLYAXIAL
8457530P	ARM15T SCREW, 7.5X30MM POLYAXIAL
8457535P	ARM15T SCREW, 7.5X35MM POLYAXIAL
8457540P	ARM15T SCREW, 7.5X40MM POLYAXIAL
8457545P	ARM15T SCREW, 7.5X45MM POLYAXIAL
8457550P	ARM15T SCREW, 7.5X50MM POLYAXIAL
8457555P	ARM15T SCREW, 7.5X55MM PALYAXIAL
8458530P	ARM15T SCREW. 8.5X30MM POLYAXIAL
8458535P	ARM15T SCREW. 8.5X35MM POLYAXIAL
8458540P	ARM15T SCREW. 8.5X40MM POLYAXIAL
8458545P	ARM15T SCREW. 8.5X45MM POLYAXIAL
8458550P	ARM15T SCREW. 8.5X50MM POLYAXIAL
8461100P	ARMADA TI LOCK SCREW, OPEN TULIP
8461302P	ASF16T ROD CONNECTOR, 5.5MM IN-LINE
8461402P	ASF CONNECTOR, 5.5-6.25MM PARALLEL
8635030P	ARM15T SCREW, 5.0X30MM RED UNIPLANAR
8635035P	ARM15T SCREW. 5.0X35MM RED UNIPLANAR
8635040P	ARM15T SCREW. 5.0X40MM RED UNIPLANAR
8636035P	ARM15T SCREW, 6.0X35MM RED UNIPLANAR
8636040P	ARM15T SCREW, 6.0X40MM RED UNIPLANAR
8636045P	ARM15T SCREW. 6.0X45MM RED UNIPLANAR
8645030P	ARM15T SCREW. 5.0X30MM UNIPLANAR
8645035P	ARM15T SCREW. 5.0X35MM UNIPLANAR
8645040P	ARM15T SCREW. 5.0X40MM UNIPLANAR
8645045P	ARM15T SCREW. 5.0X45MM UNIPLANAR
8646030P	ARM15T SCREW. 6.0X30MM UNIPLANAR
8646035P	ARM15T SCREW, 6.0X35MM UNIPLANAR
8646040P	ARM15T SCREW. 6.0X40MM UNIPLANAR



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8646045P	ARM15T SCREW, 6.0X45MM UNIPLANAR
8650520P	ARM15T CONNECTOR. 20MM OPEN OFFSET
8650530P	ARM15T CONNECTOR, 30MM OPEN OFFSET
8650560P	ARM15T CONNECTOR, 60MM OPEN OFFSET
8855530P	ARM15T SCREW. 5.5X30MM FIXED
8855535P	ARM15T SCREW, 5.5X35MM FIXED
8855540P	ARM15T SCREW, 5.5X40MM FIXED
8855545P	ARM15T SCREW, 5.5X45MM FIXED
8856530P	ARM15T SCREW, 6.5X30MM FIXED
8856535P	ARM15T SCREW. 6.5X35MM FIXED
8856540P	ARM15T SCREW, 6.5X40MM FIXED
8856545P	ARM15T SCREW, 6.5X45MM FIXED
8856550P	ARM15T SCREW, 6.5X50MM FIXED
8856555P	ARM15T SCREW. 6.5X55MM FIXED
8857535P	ARM15T SCREW. 7.5X35MM FIXED
8857540P	ARM15T SCREW, 7.5X40MM FIXED
8857545P	ARM15T SCREW, 7.5X45MM FIXED
8857550P	ARM15T SCREW, 7.5X50MM FIXED
8857555P	ARM15T SCREW, 7.5X55MM FIXED
1518897	VP SET SCREW. OFFSET & OCCIPITAL CONN
1520967	VP SET SCREW. CROSS CONNECTOR
1525190	VP SET SCREW. PARALLEL ROD CONNECTOR
1525218	VP SET SCREW. INLINE ROD CONNECTOR
1632920	VP IL SET SCREW. CROSS CONNECTOR
7901310	VUEPOINT SCREW. 3.5X10MM POLYAXIAL TULIP
7901312	VUEPOINT SCREW. 3.5X12MM POLYAXIAL TULIP
7901314	VUEPOINT SCREW. 3.5X14MM POLYAXIAL TULIP
7901316	VUEPOINT SCREW. 3.5X16MM POLYAXIAL TULIP
7901318	VUEPOINT SCREW. 3.5X18MM POLYAXIAL TULIP
7901320	VUEPOINT SCREW. 3.5X20MM POLYAXIAL TULIP
7901322	VUEPOINT SCREW. 3.5X22MM POLYAXIAL TULIP
7901324	VUEPOINT SCREW. 3.5X24MM POLYAXIAL TULIP
7901326	VUEPOINT SCREW. 3.5X26MM POLYAXIAL TULIP
7901328	VUEPOINT SCREW. 3.5X28MM POLYAXIAL TULIP
7901330	VUEPOINT SCREW. 3.5X30MM POLYAXIAL TULIP
7901332	VUEPOINT SCREW. 3.5X32MM POLYAXIAL TULIP
7901334	VUEPOINT SCREW. 3.5X34MM POLYAXIAL TULIP
7901526	VUEPOINT SCREW. 4.5X26MM PALYAXIAL TULIP
7901528	VUEPOINT SCREW. 4.5X28MM POLYAXIAL TULIP



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
7901530	VUEPOINT SCREW. 4.5X30MM POLYAXIAL TULIP
7901532	VUEPOINT SCREW. 4.5X32MM POLYAXIAL TULIP
7901534	VUEPOINT SCREW. 4.5X34MM POLYAXIAL TULIP
7901536	VUEPOINT SCREW. 4.5X36MM POLYAXIAL TULIP
7901538	VUEPOINT SCREW. 4.5X38MM POLYAXIAL TULIP
7901540	VUEPOINT SCREW. 4.5X40MM POLYAXIAL TULIP
7902310	VUEPOINT SCREW. 3.5X10MM FAV ANGLE TULIP
7902312	VUEPOINT SCREW. 3.5X12MM FAV ANGLE TULIP
7902314	VUEPOINT SCREW. 3.5X14MM FAV ANGLE TULIP
7902316	VUEPOINT SCREW. 3.5X16MM FAV ANGLE TULIP
7902318	VUEPOINT SCREW. 3.5X18MM FAV ANGLE TULIP
7902320	VUEPOINT SCREW. 3.5X20MM FAY ANGLE TULIP
7902322	VUEPOINT SCREW. 3.5X22MM FAV ANGLE TULIP
7902324	VUEPOINT SCREW. 3.5X24MM FAY ANGLE TULIP
7902326	VUEPOINT SCREW. 3.5X26MM FAV ANGLE TULIP
7902410	VUEPOINT SCREW. 4X10MM FAV ANGLE TULIP
7902412	VUEPOINT SCREW. 4X12MM FAV ANGLE TULIP
7902414	VUEPOINT SCREW. 4X14MM FAV ANGLE TULIP
7902416	VUEPOINT SCREW. 4X16MM FAV ANGLE TULIP
7902418	VUEPOINT SCREW. 4X18MM FAV ANGLE TULIP
7902420	VUEPOINT SCREW. 4X20MM FAV ANGLE TULIP
7902422	VUEPOINT SCREW. 4X22MM FAV ANGLE TULIP
7902424	VUEPOINT SCREW. 4X24MM FAV ANGLE TULIP
7902426	VUEPOINT SCREW. 4X26MM FAV ANGLE TULIP
7902428	VUEPOINT SCREW. 4X28MM FAV ANGLE TULIP
7902430	VUEPOINT SCREW. 4X30MM FAV ANGLE TULIP
7902432	VUEPOINT SCREW. 4X32MM FAV ANGLE TULIP
7902434	VUEPOINT SCREW. 4X34MM FAV ANGLE TULIP
7903001	VUEPOINT CONNECTOR. OFFSET OPEN
7903002	VUEPOINT CONNECTOR. EYELET OPTICAL
7903003	VUEPOINT CONNECTOR. OFFSET LONG
7903006	VUEPOINT SET SCREW. TULIP
7903126	VUEPOINT CROSS CONN. 3.5X26MM FIXED ROD
7903128	VUEPOINT CROSS CONN. 3.5X28MM FIXED ROD
7903130	VUEPOINT CROSS CONN. 3.5X30MM FIXED ROD
7903132	VUEPOINT CROSS CONN. 3.5X32MM FIXED ROD
7903134	VUEPOINT CROSS CONN. 3.5X34MM FIXED ROD
7903136	VUEPOINT CROSS CONN. 3.5X36MM FIXED ROD
7903138	VUEPOINT CROSS CONN. 3.5X38MM FIXED ROD



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
7903140	VUEPOINT CROSS CONN. 3.5X40MM FIXED ROD
7903142	VUEPOINT CROSS CONN. 3.5X42MM FIXED ROD
7903144	VUEPOINT CROSS CONN. 3.5X44MM FIXED ROD
7903146	VUEPOINT CROSS CONN. 3.5X46MM FIXED ROD
7903148	VUEPOINT CROSS CONN. 3.5X48MM FIXED ROD
7903150	VUEPOINT CROSS CONN. 3.5X50MM FIXED ROD
7903324	VUEPOINT SCREW. 3.5X24MM LAG FAV ANGLE
7903326	VUEPOINT SCREW. 3.5X26MM LAG FAV ANGLE
7903328	VUEPOINT SCREW. 3.5X28MM LAG FAV ANGLE
7903330	VUEPOINT SCREW. 3.5X30MM LAG FAV ANGLE
7903332	VUEPOINT SCREW. 3.5X32MM LAG FAV ANGLE
7903334	VUEPOINT SCREW. 3.5X34MM LAG FAV ANGLE
7904035	VUEPOINT ROD CONN. 3.5MM TO 3.5MM INLINE
7904045	VUEPOINT ROD CONN. 3.5MM TO 4.5MM INLINE
7904055	VUEPOINT ROD CONN. 3.5MM TO 5.5MM INLINE
7904060	VUEPOINT ROD. 3.5X60MM
7904100	VUEPOINT PLATE. 35MM KEEL
7904120	VUEPOINT ROD. 3.5X120MM
7904135	VUEPOINT ROD CONN. 3.5MM-3.5MM OFFSET
7904145	VUEPOINT ROD CONN. 3.5TM-4.5MM OFFSET
7904200	VUEPOINT PLATE. 40MM KEEL
7904240	VUEPOINT ROD. 3.5X240MM
7904300	VUEPOINT PLATE. 45MM KEEL
7904406	VUEPOINT SCREW. 4.5X6MM CORT OCCIPITAL
7904408	VUEPOINT SCREW. 4.5X8MM CORT OCCIPITAL
7904410	VUEPOINT SCREW. 4.5X10MM CORT OCCIPITA
7904412	VUEPOINT SCREW. 4.5X12MM CORT OCCIPITAL
7904414	VUEPOINT SCREW. 4.5X14MM CORT OCCIPITAL
7904506	VUEPOINT SCREW. 5X6MM CORT OCCIPITAL
7904508	VUEPOINT SCREW. 5X8MM CORT OCCIPITAL
7904510	VUEPOINT SCREW. 5X10MM CORT OCCIPITAL
7904512	VUEPOINT SCREW. 5X12MM CORT OCCIPITAL
7904514	VUEPOINT SCREW. 5X14MM CORT OCCIPITAL
7905308	VUEPOINT SCREW. 3.5X8MM MULTI AXIAL
7905310	VUEPOINT SCREW. 3.5X10MM MULTI AXIAL
7905312	VUEPOINT SCREW. 3.5X12MM MULTI AXIAL
7905314	VUEPOINT SCREW. 3.5X14MM MULTI AXIAL
7905316	VUEPOINT SCREW. 3.5X16MM MULTI AXIAL
7905318	VUEPOINT SCREW. 3.5X18MM MULTI AXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
7905320	VUEPOINT SCREW. 3.5X20MM MULTI AXIAL
7905322	VUEPOINT SCREW. 3.5X22MM MULTI AXIAL
7905324	VUEPOINT SCREW. 3.5X24MM MULTI AXIAL
7905326	VUEPOINT SCREW. 3.5X26MM MULTI AXIAL
7905328	VUEPOINT SCREW. 3.5X28MM MULTI AXIAL
7905330	VUEPOINT SCREW. 3.5X30MM MULTI AXIAL
7905332	VUEPOINT SCREW. 3.5X32MM MULTI AXIAL
7905334	VUEPOINT SCREW. 3.5X34MM MULTI AXIAL
7905336	VUEPOINT SCREW. 3.5X36MM MULTI AXIAL
7905338	VUEPOINT SCREW. 3.5X38MM MULTI AXIAL
7905340	VUEPOINT SCREW. 3.5X40MM MULTI AXIAL
7905342	VUEPOINT SCREW. 3.5X42MM MULTI AXIAL
7905344	VUEPOINT SCREW. 3.5X44MM MULTI AXIAL
7905346	VUEPOINT SCREW. 3.5X46MM MULTI AXIAL
7905348	VUEPOINT SCREW. 3.5X48MM MULTI AXIAL
7905350	VUEPOINT SCREW. 3.5X50MM MULTI AXIAL
7905410	VUEPOINT SCREW. 4X10MM MULTI AXIAL
7905412	VUEPOINT SCREW. 4X12MM MULTI AXIAL
7905414	VUEPOINT SCREW. 4X14MM MULTI AXIAL
7905416	VUEPOINT SCREW. 4X16MM MULTI AXIAL
7905418	VUEPOINT SCREW. 4X18MM MULTI AXIAL
7905420	VUEPOINT SCREW. 4X20MM MULTI AXIAL
7905422	VUEPOINT SCREW. 4X22MM MULTI AXIAL
7905424	VUEPOINT SCREW. 4X24MM MULTI AXIAL
7905426	VUEPOINT SCREW. 4X26MM MULTI AXIAL
7905428	VUEPOINT SCREW. 4X28MM MULTI AXIAL
7905430	VUEPOINT SCREW. 4X30MM MULTI AXIAL
7905432	VUEPOINT SCREW. 4X32MM MULTI AXIAL
7905434	VUEPOINT SCREW. 4X34MM MULTI AXIAL
7905436	VUEPOINT SCREW. 4X36MM MULTI AXIAL
7905438	VUEPOINT SCREW. 4X38MM MULTI AXIAL
7905440	VUEPOINT SCREW. 4X40MM MULTI AXIAL
7905442	VUEPOINT SCREW. 4X42MM MULTI AXIAL
7908444	VUEPOINT SCREW. 4X44MM MULTI AXIAL
7905446	VUEPOINT SCREW. 4X46MM MULTI AXIAL
7905448	VUEPOINT SCREW. 4X48MM MULTI AXIAL
7905450	VUEPOINT SCREW. 4X50MM MULTI AXIAL
7905508	VUEPOINT SCREW. 3.5X8MM FAV ANGLE MA
7905510	VUEPOINT SCREW. 3.5X10MM FAV ANGLE MA



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
7905512	VUEPOINT SCREW. 3.5X12MM FAV ANGLE MA
7905514	VUEPOINT SCREW. 3.5X14MM FAV ANGLE MA
7905516	VUEPOINT SCREW. 3.5X16MM FAV ANGLE MA
7905518	VUEPOINT SCREW. 3.5X18MM FAV ANGLE MA
7905520	VUEPOINT SCREW. 3.5X20MM FAV ANGLE MA
7905522	VUEPOINT SCREW. 3.5X22MM FAV ANGLE MA
7905524	VUEPOINT SCREW. 3.5X24MM FAV ANGLE MA
7905526	VUEPOINT SCREW. 3.5X26MM FAV ANGLE MA
7905528	VUEPOINT SCREW. 3.5X28MM FAV ANGLE MA
7905530	VUEPOINT SCREW. 3.5X30MM FAV ANGLE MA
7905532	VUEPOINT SCREW. 3.5X32MM FAY ANGLE MA
7905534	VUEPOINT SCREW. 3.5X34MM FAV ANGLE MA
7905536	VUEPOINT SCREW. 3.5X36MM FAV ANGLE MA
7905538	VUEPOINT SCREW. 3.5X38MM FAV ANGLE MA
7905540	VUEPOINT SCREW. 3.5X40MM FAV ANGLE MA
7905542	VUEPOINT SCREW. 3.5X42MM FAV ANGLE MA
7905544	VUEPOINT SCREW. 3.5X44MM FAV ANGLE MA
7905546	VUEPOINT SCREW. 3.5X46MM FAV ANGLE MA
7905548	VUEPOINT SCREW. 3.5X48MM FAV ANGLE MA
7905550	VUEPOINT SCREW. 3.5X50MM FAV ANGLE MA
7905610	VUEPOINT SCREW. 4X10MM FAV ANGLE MA
7905612	VUEPOINT SCREW. 4X12MM FAV ANGLE MA
7905614	VUEPOINT SCREW. 4X14MM FAV ANGLE MA
7905616	VUEPOINT SCREW. 4X16MM FAV ANGLE MA
7905618	VUEPOINT SCREW. 4X18MM FAV ANGLE MA
7905620	VUEPOINT SCREW. 4X20MM FAV ANGLE MA
7905622	VUEPOINT SCREW. 4X22MM FAV ANGLE MA
7905624	VUEPOINT SCREW. 4X24MM FAV ANGLE MA
7905626	VUEPOINT SCREW. 4X26MM FAV ANGLE MA
7905628	VUEPOINT SCREW. 4X28MM FAV ANGLE MA
7905630	VUEPOINT SCREW. 4X30MM FAV ANGLE MA
7905632	VUEPOINT SCREW. 4X32MM FAV ANGLE MA
7905634	VUEPOINT SCREW. 4X34MM FAV ANGLE MA
7905636	VUEPOINT SCREW. 4X36MM FAV ANGLE MA
7905638	VUEPOINT SCREW. 4X38MM FAV ANGLE MA
7905640	VUEPOINT SCREW. 4X40MM FAV ANGLE MA
7905642	VUEPOINT SCREW. 4X42MM FAV ANGLE MA
7905644	VUEPOINT SCREW. 4X44MM FAV ANGLE MA
7905646	VUEPOINT SCREW. 4X46MM FAV ANGLE MA



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
7905648	VUEPOINT SCREW. 4X48MM FAV ANGLE MA
7905650	VUEPOINT SCREW. 4X50MM FAV ANGLE MA
7905710	VUEPOINT SCREW. 4.5X10MM MULTI AXIAL
7905712	VUEPOINT SCREW. 4.5X12MM MULTI AXIAL
7905714	VUEPOINT SCREW. 4.5X14MM MULTI AXIAL
7905716	VUEPOINT SCREW. 4.5X16MM MULTI AXIAL
7905718	VUEPOINT SCREW. 4.5X18MM MULTI AXIAL
7905720	VUEPOINT SCREW. 4.5X20MM MULTI AXIAL
7905722	VUEPOINT SCREW. 4.5X22MM MULTI AXIAL
7905724	VUEPOINT SCREW. 4.5X24MM MULTI AXIAL
7905726	VUEPOINT SCREW. 4.5X26MM MULTI AXIAL
7905728	VUEPOINT SCREW. 4.5X28MM MULTI AXIAL
7905730	VUEPOINT SCREW. 4.5X30MM MULTI AXIAL
7905732	VUEPOINT SCREW. 4.5X32MM MULTI AXIAL
7905734	VUEPOINT SCREW. 4.5X34MM MULTI AXIAL
7905736	VUEPOINT SCREW. 4.5X36MM MULTI AXIAL
7905738	VUEPOINT SCREW. 4.5X38MM MULTI AXIAL
7905740	VUEPOINT SCREW. 4.5X40MM MULTI AXIAL
7905742	VUEPOINT SCREW. 4.5X42MM MULTI AXIAL
7905744	VUEPOINT SCREW. 4.5X44MM MULTI AXIAL
7905746	VUEPOINT SCREW. 4.5X46MM MULTI AXIAL
7905748	VUEPOINT SCREW. 4.5X48MM MULTI AXIAL
7905750	VUEPOINT SCREW. 4.5X50MRN MULTI AXIAL
7905824	VP SCREW. 3.5.X24MM FAV ANG PART THREAD
7905826	VP SCREW. 3.5.X26MM FAV ANG PART THREAD
7905828	VP SCREW. 3.5.X28MM FAV ANG PART THREAD
7905830	VP SCREW. 3.5.X30MM FAV ANG PART THREAD
7905832	VP SCREW. 3.5.X32MM FAV ANG PART THREAD
7905834	VP SCREW. 3.5.X34MM FAV ANG PART THREAD
7905836	VP SCREW. 3.5.X36MM FAV ANG PART THREAD
7905924	VUEPOINT OCT SCREW. 3.5X24/10MM P.T.
7905925	VUEPOINT OCT SCREW. 3.5X24/12MM P.T
7905926	VUEPOINT OCT SCREW. 3.5X26/10MM P.T
7905927	VUEPOINT OCT SCREW. 3.5X26/12MM P.T.
7905928	VUEPOINT OCT SCREW. 3.5X28/10MM P.T.
7905929	VUEPOINT OCT SCREW. 3.5X28/12MM P.T
7905930	VUEPOINT OCT SCREW. 3.5X30/10MM PT.
7905931	VUEPOINT OCT SCREW. 3.5X30/12MM P.T.
7905932	VUEPOINT OCT SCREW. 3.5X32/16MM P.T.



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
7905933	VUEPOINT OCT SCREW. 3.5X32/12MM P.T.
7905934	VUEPOINT OCT SCREW. 3.5X34/10MM P.T.
7905935	VUEPOINT OCT SCREW. 3.5X34/12MM PT.
7905936	VUEPOINT OCT SCREW. 3.5X36/10MM P.T.
7905937	VUEPOINT OCT SCREW. 3.5X36/12MM P.T.
7906008	VUEPOINT OFFSET CONN. 11MM OPEN
7906022	VUEPOINT OFFSET CONN. 25MM OPEN
7906025	VUEPOINT OFFSET CONN. ACUTE ANGLE
7906035	VUEPOINT OCT CONN. 3.5-3.5MM INLINE ROD
7906045	VUEPOINT OCT CONN. 3.5-4.5MM INLINE ROD
7906055	VUEPOINT CCT CONN. 3.5-5.5MM INLINE ROD
7906065	VUEPOINT OCT CONN. 3.5-6.25MM INLINE ROD
7906120	VUEPOINT OCT ROD. 3.5X120MM
7906126	VUEPOINT OCT CONNECTOR. 26MM FIXED CROSS
7906128	VUEPOINT OCT CONNECTOR. 28MM FIXED CROSS
7906130	VUEPOINT OCT CONNECTOR. 30MM FIXED CROSS
7906132	VUEPOINT OCT CONNECTOR. 32MM FIXED CROSS
7906134	VUEPOINT OCT CONNECTER. 34MM FIXED CROSS
7906136	VUEPOINT OCT CONNECTOR. 36MM FIXED CROSS
7906138	VUEPOINT OCT CONNECTER. 38MM FIXED CROSS
7906140	VUEPOINT OCT CONNECTOR. 40MM FIXED CROSS
7906142	VUEPOINT OCT CONNECTOR. 42MM FIXED CROSS
7906144	VUEPOINT OCT CONNECTOR. 44MM FIXED CROSS
7906146	VUEPOINT OCT CONNECTOR. 46MM FIXED CROSS
7906148	VUEPOINT OCT CONNECTOR. 48MM FIXED CROSS
7906150	VUEPOINT OCT CONNECTOR. 50MM FIXED CROSS.
7906160	VUEPOINT OCT ROD. 3.5X60MM
7906240	VUEPOINT OCT ROD. 3.5X240MM
7906300	VUEPOINT OCT ROD. 240MM 45° PRECURVED
7906302	VUEPOINT OCT CONN. OCCIPITAL EYELET
7906306	VP SCREW. 4.5X6MM CORTICAL OCCIPITAL
7906308	VP SCREW. 4.5X8MM CORTICAL OCCIPITAL
7906310	VP SCREW. 4.5X10MM CORTICAL OCCIPITAL
7906312	VP SCREW. 4.5X12MM CORTICAL OCCIPITAL
7906314	VP SCREW. 4.5X14MM CORTICAL OCCIPITAL
7906335	VUEPOINT OCT PLATE. 35MM KEEL
7906340	VUEPOINT OCT PLATE. 40MM KEEL
7906345	VUEPOINT OCT PLATE. 45MM KEEL
7906355	VUEPOINT OCT ROD. 3.5-5.5MM TRANSITION



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
7906356	VP SCREW. 5X6MM CORTICAL OCCIPITAL
7906358	VP SCREW. 5X8MM CORTICAL OCCIPITAL
7906360	VP SCREW. 5X10MM CORTICAL OCCIPITAL
7906362	VP TI SCREW. 5X12MM CORTICAL OCCIPITAL
7906364	VP TI SCREW. 5X14MM CORTICAL OCCIPITAL
7906435	VP CONNECTOR. 3.5 TO 3.5MM ROD OFFSET
7906445	VP CONNECTOR. 3.5 TO 4.5MM ROD OFFSET
7906455	VP CONNECTOR. 3.5 TO 5.5MM ROD OFFSET
7906465	VP CONNECTOR. 3.5 TO 6.25MM ROD OFFSET
7907760	VP SET SCREW. OFFSET & OCCIPITAL CONN
7907765	VP TI SET SCREW. PARALLEL ROD CONN
7907770	VP TI SET SCREW. IN-LINE ROD CONN
7907775	VP TI SET SCREW. CROSS CONN
7908026	VUEPOINT CROSS CONN. 26MM GRIT BLASTED
7908028	VUEPOINT CROSS CONN. 28MM GRIT BLASTED
7908030	VUEPOINT CROSS CONN. 30MM GRIT BLASTED
7908032	VUEPOINT CROSS CONN. 32MM GRIT BLASTED
7908034	VUEPOINT CROSS CONN. 34MM GRIT BLASTED
7908036	VUEPOINT CROSS CONN. 36MM GRIT BLASTED
7908038	VUEPOINT CROSS CONN. 38MM GRIT BLASTED
7908040	VUEPOINT CROSS CONN. 40MM GRIT BLASTED
7908042	VUEPOINT CROSS CONN. 42MM GRIT BLASTED
7908044	VUEPOINT CROSS CONN. 44MM GRIT BLASTED
7908046	VUEPOINT CROSS CONN. 46MM GRIT BLASTED
7908048	VUEPOINT CROSS CONN. 48MM GRIT BLASTED
7908050	VUEPOINT CROSS CONN. 50MM GRIT BLASTED
7908052	VP TI SET SCREW. CROSS CONNECTOR GR
7908240	VUEPOINT COCR ROD. 3.5X240MM
7908300	VP COCR ROD. 240MM PRE-CURVED 45°
7908345	VUEPOINT COCR ROD. 3.5 TO 4.5 TAPERED
7908355	VUEPOINT COCR ROD. 3.5 TO 5.5 TAPERED
7908360	VUEPOINT COCR ROD. Ø3.5-6.0X420MM TRANS
7908362	VP COCR ROD. Ø3.5-6.25X420MM TRANS
8790000	VUEPOINT SET SCREW. TULP AND KEEL PLATE
8790406	VUEPOINT SCREW. 4.5X6MM OCCIPITAL
8790408	VUEPOINT SCREW. 4.5X8MM OCCIPITAL
8790410	VUEPOINT SCREW. 4.5X10MM OCCIPITAL
8790412	VUEPOINT SCREW. 4.5X12MM OCCIPITAL
8791312	VUEPOINT SCREW. 3.5X12MM MULTI AXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8791314	VUEPOINT SCREW. 3.5X14MM MULTI AXIAL
8791414	VUEPOINT SCREW. 4.0X14MM MULTI AXIAL
8791420	VUEPOINT SCREW. 4.0X20MM MULTI AXIAL
8792360	VUEPOINT CROSS CONNECTOR. ADJUSTABLE
8793001	VUEPOINT ROD. HINGED
8793330	VUEPOINT SCREW. 3.5X30MM FAVORED ANGLE
8794524	VUEPOINT SCREW. 4.5X24MM MULTI AXIAL
8795001	VUEPOINT OCCIPITAL PLATE. ADJUSTABLE
8945308	MAS PCF SCREW. 3.5X8MM
8945310	MAS PCF SCREW. 3.5X10MM
8945312	MAS PCF SCREW. 3.5X12MM
8945314	MAS PCF SCREW. 3.5X14MM
8945314	MAS PCF SCREW. 3.5X14MM
8945316	MAS PCF SCREW. 3.5X16MM
8945316	MAS PCF SCREW. 3.5X16MM
8945318	MAS PCF SCREW. 3.5X18MM
8945320	MAS PCF SCREW. 3.5X20MM
8945322	MAS PCF SCREW. 3.5X22MM
8945324	MAS PCF SCREW. 3.5X24MM
8945326	MAS PCF SCREW. 3.5X26MM
8945328	MAS PCF SCREW. 3.8X28MM
8945330	MAS PCF SCREW. 3.5X30MM
8945332	MAS PCF SCREW. 3.5X32MM
8945334	MAS PCF SCREW. 3.5X34MM
8945336	MAS PCF SCREW. 3.5X36MM
8945338	MAS PCF SCREW. 3.5X38MM
8945410	MAS PCF SCREW. 4.0X10MM
8945412	MAS PCF SCREW. 4.0X12MM
8945414	MAS PCF SCREW. 4.0X14MM
8945416	MAS PCF SCREW. 4.0X16MM
8972360	VUEPOINT IL CROSS CONN. ADJ TULIP-TULIP
8973535	VUEPOINT IL TI ROD. Ø3.5X35MM CURVED
8973545	VUEPOINT IL TI ROD. Ø3.5X45MM CURVED
8973555	VUEPOINT IL TI ROD Ø3.5X55MM CURVED
8973565	VUEPOINT LL TI ROD Ø3.5X65MM CURVED
8973575	VUEPOINT IL TI ROD. Ø3.5X75MM CURVED
8975000	VUEPOINT IL KEEL PLATE. ADJUSTABLE
8975115	VUEPOINT IL SCREW. 5.0X15MM MA
8975120	VUEPOINT II SCREW. 5.0X20MM MA



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8975122	VUEPOINT LL SCREW 5.0X22MM MA
8975124	VUEPOINT IL SCREW. 5.0X24MM MA
8975125	VUEPOINT LL SCREW. 5.0X25MM MA
8975126	VUEPOINT IL SCREW. 5.0X26MM MA
8975128	VUEPOINT IL SCREW. 5.0X28MM MA
8975130	VUEPOINT IL SCREW. 5.0X30MM MA
8975132	VUEPOINT II SCREW. 5.0X32MM MA
8975134	VUEPOINT IL SCREW. 5.0X34MM MA
8975135	VUEPOINT IL SCREW. 5.0X35MM MA
8975136	VUEPOINT IL SCREW 5.0X36MM MA
8975138	VUEPOINT IL SCREW.5.0X38MM MA
8975140	VUEPOINT IL SCREW. 5.0X40MM MA
8975142	VUEPOINT LL SCREW. 5.0X42MM MA
8975144	VUEPOINT IL SCREW. 5.0X44MM MA
8975145	VUEPOINT IL SCREW. 5.0X45MM MA
8975146	VUEPOINT LL SCREW. 5.0X46MM MA
8975148	VUEPOINT IL SCREW. 5.0X48MM MA
8975150	VUEPOINT IL SCREW. 5.0X50MM MA
8975215	VUEPOINT IL SCREW. 5.5X15MM MA
8975220	VUEPOINT IL SCREW. 5.5X20MM MA
8975222	VUEPOINT IL SCREW. 5.5X22MM MA
8975224	VUEPOINT IL SCREW. 5.5X24MM MA
8975225	VUEPOINT IL SCREW. 5.5X25MM MA
8975226	VUEPOINT IL SCREW. 5.5X26MM MA
8975228	VUEPOINT IL SCREW. 5.5X28MM MA
8975230	VUEPOINT IL SCREW. 5.5X30MM MA
8975232	VUEPOINT IL SCREW. 5.5X32MM MA
8975234	VUEPOINT IL SCREW. 5.5X34MM MA
8975235	VUEPOINT IL SCREW. 5.5X35MM MA
8975236	VUEPOINT IL SCREW. 5.5X36MM MA
8975238	VUEPOINT IL SCREW. 5.5X38MM MA
8975240	VUEPOINT IL SCREW. 5.5X40MM MA
8975242	VUEPOINT II SCREW. 5.5X42MM MA
8975244	VUEPOINT IL SCREW. 5.5X44MM MA
8975245	VUEPOINT IL SCREW. 5.5X45MM MA
8975246	VUEPOINT II SCREW. 5.5X46MM MA
8975248	VUEPOINT IL SCREW. 5.5X48MM MA
8975250	VUEPOINT II SCREW. 5.5X50MM MA
8975308	VUEPOINT IL SCREW. 3.5X8MM MA



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8975310	VUEPOINT IL SCREW. 3.5X 10MM MA
8975312	VUEPOINT LI SCREW. 3.5X12MM MA
8975314	VUEPOINT IL SCREW. 3.6X14MM MA
8975316	VUEPOINT IL SCREW. 3.5X16MM MA
8975318	VUEPOINT IL SCREW. 3.5X18MM MA
8975320	VUEPOINT IL SCREW. 3.5X20MM MA
8975322	VUEPOINT IL SCREW. 3.5X22MM MA
8975324	VUEPOINT IL SCREW. 3.5X24MM MA
8975326	VUEPOINT IL SCREW. 3.5X26MM MA
8975328	VUEPOINT IL SCREW. 3.5X28MM MA
8975330	VUEPOINT IL SCREW. 3.5X30MM MA
8975332	VUEPOINT IL SCREW. 3.5X32MM MA
8975334	VUEPOINT II SCREW. 3.5X34MM MA
8975336	VUEPOINT IL SCREW. 3.5X36MM MA
8975338	VUEPOINT II SCREW. 3.5X38MM MA
8975340	VUEPOINT II SCREW. 3.5X40MM MA
8975410	VUEPOINT II SCREW. 4.0X10MM MA
8975412	VUEPOINT II SCREW. 4.0X12MM MA
8975414	VUEPOINT IL SCREW. 4.0X14MM MA
8975416	VUEPOINT LI SCREW. 4.0X 16MM MA
8975418	VUEPOINT IL SCREW. 4.0X18MM MA
8975420	VUEPOINT II SCREW. 4.0X20MM MA
8975422	VUEPOINT IL SCREW. 4.0X22MM MA
8975424	VUEPOINT IL SCREW. 4.0X24MM MA
8975426	VUEPOINT II SCREW. 4.0X26MM MA
8975428	VUEPOINT II SCREW. 4.0X28MM MA
8975430	VUEPOINT IL SCREW. 4.0X30MM MA
8975432	VUEPOINT IL SCREW. 4.0X32MM MA
8975434	VUEPOINT IL SCREW. 4.0X34MM MA
8975436	VUEPOINT II SCREW. 4.0X36MM MA
8975438	VUEPOINT IL SCREW. 4.0X38MM MA
8975440	VUEPOINT IL SCREW. 4.0X40MM MA
8975510	VUEPOINT II SCREW. 3.5X10MM FA
8975512	VUEPOINT IL SCREW. 3.5X12MM FA
8975514	VUEPOINT IL SCREW. 3.5X14MM FA
8975516	VUEPOINT IL SCREW. 3.5X16MM FA
8975518	VUEPOINT IL SCREW. 3.5X18MM FA
8975520	VUEPOINT IL SCREW. 3.5X20MM FA
8975522	VUEPOINT IL SCREW. 3.5X22MM FA



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8975524	VUEPOINT LL SCREW. 3..5X24MM FA
8975526	VUEPOINT LI SCREW. 3.5X26MM FA
8975528	VUEPOINT II SCREW. 3.5X28MM FA
8975530	VUEPOINT II SCREW. 3.5X30MM FA
8975532	VUEPOINT II SCREW. 3.5X32MM FA
8975534	VUEPOINT IL SCREW. 3.5X34MM FA
8975610	VUEPOINT IL SCREW. 4.0X10MM FA
8975612	VUEPOINT IL SCREW. 4.0X12MM FA
8975614	VUEPOINT LL SCREW. 4.0X14MM FA
8975616	VUEPOINT II SCREW. 4.0X16MM FA
8975618	VUEPOINT IL SCREW. 4.0X18MM FA
8975620	VUEPOINT IL SCREW. 4.0X20MM FA
8975622	VUEPOINT IL SCREW. 4.0X22MM FA
8975624	VUEPOINT LL SCREW. 4.0X24MM FA
8975626	VUEPOINT LL SCREW. 4.0X26MM FA
8975628	VUEPOINT LL SCREW. 4.0X28MM FA
8975630	VUEPOINT IL SCREW. 4.0X30MM FA
8975632	VUEPOINT II SCREW. 4.0X32MM FA
8975634	VUEPOINT IL SCREW. 4.0X34MM FA
8975636	VUEPOINT IL SCREW. 4.0X36MM FA
8975638	VUEPOINT IL SCREW. 4.0X38MM FA
8975640	VUEPOINT IL SCREW. 4.0X40MM FA
8975710	VUEPOINT IL SCREW. 4.5X10MM MA
8975712	VUEPOINT IL SCREW. 4.5X12MM MA
8975714	VUEPOINT IL SCREW. 4.5X14MM MA
8975716	VUEPOINT IL SCREW. 4.5X16MM MA
8975718	VUEPOINT II SCREW. 4.5X18MM MA
8975720	VUEPOINT II SCREW. 4.5X20MM MA
8975722	VUEPOINT IL SCREW. 4.5X22MM MA
8975724	VUEPOINT IL SCREW. 4.5X24MM MA
8975726	VUEPOINT IL SCREW. 4.5X26MM MA
8975728	VUEPOINT IL SCREW. 4.5X28MM MA
8975730	VUEPOINT IL SCREW. 4.5X30MM MA
8975732	VUEPOINT IL SCREW. 4.5X32MM MA
8975734	VUEPOINT IL SCREW. 4.5X34MM MA
8975736	VUEPOINT II SCREW. 4.5X36MM MA
8975738	VUEPOINT IL SCREW. 4.5X38MM MA
8975740	VUEPOINT IL SCREW. 4.5X40MM MA
8975820	VP II SCREW. 3.5X20MM FA PART THREAD



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
8975822	VP II SCREW. 3.5X22MM FA PART THREAD
8975824	VP II SCREW. 3.5X24MM FA PART THREAD
8975826	VP IL SCREW. 3.5X26MM FA PART THREAD
8975828	VP IL SCREW. 3.5X28MM FA PART THREAD
8975830	VP IL SCREW. 3.5X30MM FA PART THREAD
8975832	VP II SCREW. 3.5X32MM FA PART THREAD
8975834	VP II SCREW. 3.5X34MM FA PART THREAD
8975836	VP IL SCREW. 3.5X36MM FA PART THREAD
8975838	VP IL SCREW. 3.5X38MM FA PART THREAD
8975840	VP IL SCREW. 3.5X40MM FA PART THREAD
8976011	VUEPOINT IL CONNECTOR. 11MM OFFSET OPEN
8976025	VP IL CONN. 25MM OFFSET OPEN LONG
8976030	VP IL CONN. 45° ACUTE ANGLE OFFSET
8976035	VP IL CONN. 3.5-3.5MM INLINE ROD
8976040	VUEPOINT IL COCR ROD. 3.5X400MM STRAIGHT
8976055	VP II CONN. 3.5-5.5MM INLINE ROD
8976078	VUEPOINT II COCR ROD. Ø3.5X60MM STRAIGHT
8976079	VUEPOINT IL COCR ROD Ø3.5X65MM STRAIGHT
8976080	VUEPOINT IL COCR ROD. Ø3.5X70MM STRAIGHT
8976081	VUEPOINT IL COCR ROD. Ø3.5X75MM STRAIGHT
8976082	VUEPOINT IL COCR ROD. Ø3.5X80MM STRAIGHT
8976083	VUEPOINT IL COCR ROD. Ø3.5X90MM STRAIGHT
8976084	VUEPOINT IL COCR ROD. Ø3.5X100MM STR
8976085	VUEPOINT IL COCR ROD. Ø3.5X110MM STR
8976086	VUEPOINT LL COCR ROD. Ø3.5X30MM CURVED
8976087	VUEPOINT IL COCR ROD. Ø3.5X35MM CURVED
8976088	VUEPOINT II COCR ROD. Ø3.5X40MM CURVED
8976089	VUEPOINT II COCR ROD Ø3.5X45MM CURVED
8976090	VUEPOINT II COCR ROD. Ø3.5X50MM CURVED
8976091	VUEPOINT IL COCR ROD. Ø3.5X55MM CURVED
8976092	VUEPOINT IL COCR ROD. Ø3.5X60MM CURVED
8976093	VUEPOINT LL COCR ROD. Ø3.5X65MM CURVED
8976094	VUEPOINT II COCR ROD. Ø3.5X70MM CURVED
8976095	VUEPOINT IL COCR ROD. Ø3.5X75MM CURVED
8976096	VUEPOINT IL COCR ROD. Ø3.5X80MM CURVED
8976097	VUEPOINT LL COCR ROD. Ø3.5X90MM CURVED
8976098	VUEPOINT IL COCR ROD. Ø3.5X100MM CURVED
8976099	VUEPOINT II COCR ROD. Ø3.5X110MM CURVED
8976100	VUEPOINT IL CONN. 9MM OFFSET



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8976102	VUEPOINT IL CONN. 15MM OFFSET
8976103	VUEPOINT IL CONN. 25MM OFFSET
8976104	VUEPOINT II ROD. Ø3.5X40MM
8976105	VUEPOINT II ROD. Ø3.5X50MM
8976106	VUEPOINT IL ROD. 3.5X60MM
8976107	VUEPOINT IL ROD. Ø3.5X70MM
8976108	VUEPOINT II ROD. Ø3.5X80MM STRAIGHT
8976109	VUEPOINT IL ROD. Ø3.5X90MM
8976110	VUEPOINT IL ROD. Ø3.5X 100MM CONTOURED
8976111	VUEPOINT LL ROD. Ø3.5X110MM
8976112	VUEPOINT W ROD. 3.5X120MM
8976124	VUEPOINT IL ROD. 3.5X240MM
8976128	VUEPOINT II CONNECTOR. 28MM FIXED CROSS
8976130	VUEPOINT IL CONNECTOR. 30MM FIXED CROSS
8976132	VUEPOINT IL CONNECTOR. 32MM FIXED CROSS
8976134	VUEPOINT IL CONNECTOR. 34MM FIXED CROSS
8976136	VUEPOINT IL CONNECTOR. 36MM FIXED CROSS
8976138	VUEPOINT IL CONNECTOR. 38MM FIXED CROSS
8976140	VUEPOINT IIL TI ROD. 3.5X400MM STRAIGHT
8976141	VUEPOINT IL TI ROD. Ø3.5X110MM CURVED
8976142	VUEPOINT IL TI ROD. Ø3.5X120MM CURVED
8976143	VUEPOINT IL TI ROD. Ø3.5X30MM CURVED
8976144	VUEPOINT IL TI ROD. Ø3.5X40MM CURVED
8976145	VUEPOINT II TI ROD. Ø3.5X50MM CURVED
8976146	VUEPOINT II TI ROD. Ø3.5X60MM CURVED
8976147	VUEPOINT IL TI ROD. Ø3.5X70MM CURVED
8976149	VUEPOINT II TI ROD. Ø3.5X90MM CURVED
8976165	VUEPOINT IL CONN. EXTRA ROD
8976180	VUEPOINT II ROD. 80MM PRE-CONTOURED
8976240	VUEPOINT IL ROD. 240MM 45° PRECURVED
8976300	VUEPOINT IL ROD. HINGED
8976305	VUEPOINT IL EYELETS. LOW PROFILE
8976345	VUEPOINT II TI ROD. 3.5-4.5MM TRANSITION
8976347	VP II TI ROD. Ø3.5-4.5X650MM TRANS
8976348	VP IL TI ROD. Ø3.5-5.0X650MM TRANS
8976349	VP IL COCR ROD. Ø3.5-4.5X650MM TRANS
8976350	VP IL COCR ROD. Ø3.5-5.0X650MM TRANS
8976351	VP IL COCR ROD. Ø3.5-4.5X600MM TRANS
8976352	VP IL COCR ROD. Ø3.5-5.5X600MM TRANS



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
8976353	VP IL COCR ROD. Ø3.5-6.0X600MM TRANS
8976354	VP IL COCR ROD. Ø3.5-6.25X600MM TRANS.
8976355	VUEPOINT II ROD. 3.5MM-5.5MM TRANSITION
8976356	VP IL COCR ROD. Ø3.5-4.5X700MM TRANS
8976357	VP IL COCR ROD. Ø3.5-5.5X700MM TRANS
8976358	VP IL COCR ROD. Ø3.5-6.0X700MM TRANS
8976359	VP IL COCR ROD. Ø3.5-6.25X700MM TRANS
8976406	VP IL SCREW. 4.5X6MM OCCIPITAL
8976408	VP II SCREW. 4.5X8MM OCCIPITAL
8976410	VP II SCREW. 4.5X10MM OCCIPITAL
8976412	VP IL SCREW. 4.5X12MM OCCIPITAL
8976414	VP IL SCREW. 4.5X14MM OCCIPITAL
8976416	VP IL SCREW. 4.5X16MM OCCIPITAL
8976418	VUEPOINT IL SCREW. 4.5X18MM OCCIPITAL
8976420	VUEPOINT IL SCREW. 4.5X20MM OCCIPITAL
8976435	VUEPOINT IL CONN. 3.5MM-3.5MM OFFSETROD
8976455	VUEPOINT IL CONN. 3.5MM-5.5MM OFFSETROD
8976506	VP IL SCREW. 5.0X6MM OCCIPITAL
8976508	VP IL SCREW. 5.0X8MM OCCIPITAL
8976510	VP IL SCREW. 5.0X10MM OCCIPITAL
8976512	VP IL SCREW. 5.0X12MM OCCIPITAL
8976514	VP IL SCREW. 5.0X14MM OCCIPITAL
8976516	VP IL SCREW. 5.0X16MM OCCIPITAL
8976518	VUEPOINT IL SCREW. 5.0X18MM OCCIPITAL
8976520	VUEPOINT IL SCREW. 5.0X20MM OCCIPITAL
8977760	VP II SET SCREW. OFFSET OCCIPITAL CONN
8977765	VP IL SET SCREW. ROD-ROD CONN
8977770	VP IL SET SCREW. STATIC CROSS CONNECTOR
7905310006	VUEPOINT SCREW. 3.5X10MM MULTI-AXIAL
7905312006	VUEPOINT SCREW. 3.5X12MM MULTI-AXIAL
7905314006	VUEPOINT SCREW. 3.5X14MM MULTI-AXIAL
7905316006	VUEPOINT SCREW. 3.5X 16MM MULTI-AXIAL
7905318006	VUEPOINT SCREW. 3.5X18MM MULTI-AXIAL
7905320006	VUEPOINT SCREW. 3.5X20MM MULTI-AXIAL
7905322006	VUEPOINT SCREW. 3.5X22MM MULTI-AXIAL
7905324006	VUEPOINT SCREW. 3.5X24MM MULTI-AXIAL
7905326006	VUEPOINT SCREW. 3.5X26MM MULTI-AXIAL
7905328006	VUEPOINT SCREW. 3.5X28MM MULTI-AXIAL
7905330006	VUEPOINT SCREW. 3.5X30MM MULTI-AXIAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
7905332006	VUEPOINT SCREW. 3.5X32MM MULTI-AXIAL
7905334006	VUEPOINT SCREW. 3.5X34MM MULTI-AXIAL
7905336006	VUEPOINT SCREW. 3.5X36MM MULTI-AXIAL
7905338006	VUEPOINT SCREW. 3.5X38MM MULTI-AXIAL
7905340006	VUEPOINT SCREW. 3.5X40MM MULTI-AXIAL
7905342006	VUEPOINT SCREW. 3.5X42MM MULTI-AXIAL
7905344006	VUEPOINT SCREW. 3.5X44MM MULTI-AXIAL
7905346006	VUEPOINT SCREW. 3.5X46MM MULTI-AXIAL
7905348006	VUEPOINT SCREW. 3.5X48MM MULTI-AXIAL
7905350006	VUEPOINT SCREW. 3.5X50MM MULTI-AXIAL
7905410006	VUEPOINT SCREW. 4.0X10MM MULTI-AXIAL
7905412006	VUEPOINT SCREW. 4.0X12MM MULTI-AXIAL
7905414006	VUEPOINT SCREW. 4.0X14MM MULTI-AXIAL
7905416006	VUEPOINT SCREW. 4.0X16MM MULTI-AXIAL
7905418066	VUEPOINT SCREW. 4.0X18MM MULTI-AXIAL
7905420006	VUEPOINT SCREW. 4.0X20MM MULTI-AXIAL
7905422006	VUEPOINT SCREW. 4.0X22MM MULTI-AXIAL
7905424006	VUEPOINT SCREW. 4.0X24MM MULTI-AXIAL
7905426006	VUEPOINT SCREW. 4.0X26MM MULTI-AXIAL
7905428006	VUEPOINT SCREW. 4.0X28MM MULTI-AXIAL
7905430006	VUEPOINT SCREW. 4.0X30MM MULTI-AXIAL
7905432006	VUEPOINT SCREW. 4.0X32MM MULTI-AXIAL
7905434006	VUEPOINT SCREW. 4.0X34MM MULTI-AXIAL
7905436006	VUEPOINT SCREW. 4.0X36MM MULTI-AXIAL
7905438006	VUEPOINT SCREW. 4.0X38MM MULTI-AXIAL
7905440006	VUEPOINT SCREW. 4.0X40MM MULTI-AXIAL
7905442006	VUEPOINT SCREW. 4.0X42MM MULTI-AXIAL
7905444006	VUEPOINT SCREW. 4.0X44MM MULTI-AXIAL
7905446006	VUEPOINT SCREW. 4.0X46MM MULTI-AXIAL
7905448006	VUEPOINT SCREW. 4.0X48MM MULTI-AXIAL
7905450006	VUEPOINT SCREW. 4.0X50MM MULTI-AXIAL
7905510006	VUEPOINT SCREW. 3.5X10MM FAVOURED ANGLE
7905512006	VUEPOINT SCREW. 3.5X12MM FAVOURED ANGLE
7905514006	VUEPOINT SCREW. 3.5X14MM FAVOURED ANGLE
7905516006	VUEPOINT SCREW. 3.5X16MM FAVOURED ANGLE
7905518006	VUEPOINT SCREW. 3.5X18MM FAVOURED ANGLE
7905520006	VUEPOINT SCREW. 3.5X20MM FAVOURED ANGLE
7905522006	VUEPOINT SCREW. 3.5X22MM FAVOURED ANGLE
7905524006	VUEPOINT SCREW. 3.5X24MM FAVOURED ANGLE



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
7905526006	VUEPOINT SCREW. 3.5X26MM FAVOURED ANGLE
7905528006	VUEPOINT SCREW. 3.6X28MM FAVOURED ANGLE
7905530006	VUEPOINT SCREW. 3.5X30MM FAVOURED ANGLE
7905532006	VUEPOINT SCREW. 3.5X32MM FAVOURED ANGLE
7905534006	VUEPOINT SCREW. 3.5X34MM FAVOURED ANGLE
7905610006	VUEPOINT SCREW. 4.0X10MM FAVOURED ANGLE
7905612006	VUEPOINT SCREW. 4.0X12MM FAVOURED ANGLE
7905614006	VUEPOINT SCREW. 4.0X14MM FAVOURED ANGLE
7905616006	VUEPOINT SCREW. 4.0X16MM FAVOURED ANGLE
7905618006	VUEPOINT SCREW. 4.0X18MM FAVOURED ANGLE
7905620006	VUEPOINT SCREW. 4.0X20MM FAVOURED ANGLE
7905622006	VUEPOINT SCREW. 4.0X22MM FAVOURED ANGLE
7905624006	VUEPOINT SCREW. 4.0X24MM FAVOURED ANGLE
7905626006	VUEPOINT SCREW. 4.0X26MM FAVOURED ANGLE
7905628006	VUEPOINT SCREW. 4.0X28MM FAVOURED ANGLE
7905630006	VUEPOINT SCREW. 4.0X30MM FAVOURED ANGLE
7905632006	VUEPOINT SCREW. 4.0X32MM FAVOURED ANGLE
7905634006	VUEPOINT SCREW. 4.0X34MM FAVOURED ANGLE
7905636006	VUEPOINT SCREW. 4.0X36MM FAVOURED ANGLE
7905638006	VUEPOINT SCREW. 4.0X38MM FAVOURED ANGLE
7905640006	VUEPOINT SCREW. 4.0X40MM FAVOURED ANGLE
7905710006	VUEPOINT SCREW. 4.5X10MM MULTI-AXIAL
7905712006	VUEPOINT SCREW. 4.5X12MM MULTI-AXIAL
7905714006	VUEPOINT SCREW. 4.5X14MM MULTI-AXIAL
7905716006	VUEPOINT SCREW. 4.5X16MM MULTI-AXIAL
7905718006	VUEPOINT SCREW. 4.5X18MM MULTI-AXIAL
7905720006	VUEPOINT SCREW. 4.5X20MM MULTI-AXIAL
7905722006	VUEPOINT SCREW. 4.5X22MM MULTI-AXIAL
7905724006	VUEPOINT SCREW. 4.5X24MM MULTI-AXIAL
7905726006	VUEPOINT SCREW. 4.5X26MM MULTI-AXIAL
7905728006	VUEPOINT SCREW. 4.5X28MM MULTI-AXIAL
7905730006	VUEPOINT SCREW. 4.5X30MM MULTI-AXIAL
7905732006	VUEPOINT SCREW. 4.5X32MM MULTI-AXIAL
7905734006	VUEPOINT SCREW. 4.5X34MM MULTI-AXIAL
7905736006	VUEPOINT SCREW. 4.5X36MM MULTI-AXIAL
7905738006	VUEPOINT SCREW. 4.5X38MM MULTI-AXIAL
7905740006	VUEPOINT SCREW. 4.5X40MM MULTI-AXIAL
7905824006	VUEPOINT SCREW. 3.5X24MM FA PART THREAD
7905826006	VUEPOINT SCREW. 3.5X26MM FA PART THREAD



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
7905828006	VUEPOINT SCREW. 3.5X28MM FA PART THREAD
7905830006	VUEPOINT SCREW. 3.5X30MM FA PART THREAD
7905832006	VUEPOINT SCREW. 3.5X32MM FA PART THREAD
7905834006	VUEPOINT SCREW. 3.5X34MM FA PART THREAD
7905836006	VUEPOINT SCREW. 3.5X36MM FA PART THREAD
7905924006	VP SCREW. 3.5X24/10MM FA PART THREAD
7905926006	VP SCREW. 3.5X26/10MM FA PART THREAD
7905928006	VP SCREW. 3.5X28/10MM FA PART THREAD
7905930006	VP SCREW. 3.5X30/10MM FA PART THREAD
7905932006	VP SCREW. 3.5X32/10MM FA PART THREAD
7905934006	VP SCREW. 3.5X34/10MM FA PART THREAD
7905936006	VP SCREW. 3.5X36/10MM FA PART THREAD
7906008006	VUEPOINT CONN. 11MM OFFSET
7906022006	VUEPOINT CONN. 25MM OFFSET
7906025006	VUEPOINT CONN. ACUTE ANGLE OFFSET
7906035006	VUEPOINT CONN. 3.5-3.5MM INLINE ROD
7906055006	VUEPOINT CONN. 3.5-5.5MM INLINE ROD
7906120006	VUEPOINT TI ROD. 3.5X120MM
7906126006	VUEPOINT CONNECTOR. 26MM FIXED CROSS
7906128006	VUEPOINT CONNECTOR.28MM FIXED CROSS
7906130006	VUEPOINT CONNECTOR.30MM FIXED CROSS
7906132006	VUEPOINT CONNECTOR.32MM FIXED CROSS
7906134006	VUEPOINT CONNECTOR.34MM FIXED CROSS
7906136006	VUEPOINT CONNECTOR. 36MM FIXED CROSS
7906138006	VUEPOINT CONNECTOR.38MM FIXED CROSS
7906140006	VUEPOINT CONNECTOR.40MM FIXED CROSS
7906142006	VUEPOINT CONNECTOR.42MM FIXED CROSS
7906144006	VUEPOINT CONNECTOR.44MM FIXED CROSS
7906146006	VUEPOINT CONNECTOR.46MM FIXED CROSS
7906148006	VUEPOINT CONNECTOR. 48MM FIXED CROSS
7906150006	VUEPOINT CONNECTOR.50MM FIXED CROSS
7906160006	VUEPOINT ROD. 3.5X60MM
7906240006	VUEPOINT ROD. 3.5X240MM
7906300006	VUEPOINT ROD. 240MM 45° PRECURVED
7906302006	VUEPOINT CONNECTOR. OCCIPITAL EYELET
7906306006	VUEPOINT SCREW. 4.6X6MM CORTICAL OCCIP
7906308006	VUEPOINT SCREW. 4.5X8MM CORTICAL OCCIP
7906310006	VUEPOINT SCREW. 4.5X10MM CORTICAL OCCIP
7906312006	VUEPOINT SCREW. 4.5X12MM CORTICAL OCCIP



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
7906314006	VUEPOINT SCREW. 4.5X14MM CORTICAL OCCIP
7906316006	VUEPOINT SCREW. 4.5X16MM CORTICAL OCCIP
7906318006	VUEPOINT SCREW. 4.5X18MM CORTICAL OCCIP
7906320006	VUEPOINT SCREW. 4.5X20MM CORTICAL OCCIP
7906335006	VUEPOINT PLATE. 35MM KEEL
7906340006	VUEPOINT PLATE. 40MM KEEL
7906345006	VUEPOINT PLATE. 45MM KEEL
7906355006	VUEPOINT ROD. 3.5-5.5MM TRANSITION
7906356006	VUEPOINT SCREW. 5.0X6MM CORTICAL OCCIP
7906358006	VUEPOINT SCREW. 5.0X8MM CORTICAL OCCIP
7906360006	VUEPOINT SCREW. 5.0X10MM CORTICAL OCCIP
7906362006	VUEPOINT SCREW. 5.0X12MM CORTICAL OCCIP
7906364006	VUEPOINT SCREW. 5.0X 14MM CORTICAL OCCIP
7906366006	VUEPOINT SCREW. 5.0X16MM CORTICAL OCCIP
7906368006	VUEPOINT SCREW. 5.0X18MM CORTICAL OCCIP
7906370006	VUEPOINT SCREW. 5.0X20MM CORTICAL OCCIP
7906435006	VUEPOINT CONN. 3.5 TO 3.5MM ROD OFFSET
7906455006	VUEPOINT CONN. 3.5 TO 5.5MM ROD OFFSET
7907760006	VP SET SCREW. OFFSET & OCCIPITAL CONN
7907765006	VUEPOINT SET SCREW. PARALLEL ROD CONN
7907770006	VUEPOINT SET SCREW. IN-LINE ROD CONN
7907775006	VUEPOINT SET SCREW. CROSS CONNECTER
7908240P	VUEPOINT COCR ROD. 3.5X240MM
7908300P	VP COCR ROD. 240MM PRE-CURVED 45°
7908355P	VUEPOINT COCR ROD. 3.5 TO 5.5 TAPERED
7908360P	VUEPOINT COCR ROD. Ø3.5-6.0X420MM TRANS
8972360P3	VUEPOINT IL CROSS CONN. ADJ TULIP-TULIP
8973535P	VUEPOINT II TI ROD. Ø3.5X35MM CURVED
8973545P	VUEPOINT II TI ROD. Ø3.5X45MM CURVED
8973555P	VUEPOINT II TI ROD. Ø3.5X55MM CURVED
8973565P	VUEPOINT II TI ROD. Ø3.5X65MM CURED
8973575P	VUEPOINT II TI ROD. Ø3.5X75MM CURVED
8975000P	VUEPOINT IL KEEL PLATE. ADJUSTABLE
8976011P3	VUEPOINT II CONN. 11MMX0° LATERAL OFFSET
8976015P3	VUEPOINT IL CONN. 15MMX0° LATERAL OFFSET
8976020P3	VUEPOINT II CONN. 20MMX0° LATERAL OFFSET
8976025P3	VUEPOINT IL CONN. 25MMX0° LATERAL OFFSET
8976030P3	VUEPOINT IL CONN. 30MMX45° LATERALOFFSET
8976031P3	VUEPOINT IL CONN. 30MMX35° LATERALOFFSET



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8976032P3	VUEPOINT IL CONN. 30MMX55° LATERALOFFSET
8976033P3	VUEPOINT II CONN. 30MMX65° LATERALOFFSET
8976035P3	VUEPOINT IL CONN. 3.5-3.5MM INLINE ROD
8976045P3	VUEPOINT IL CONN. 3.5-4.5MM INLINE ROD
8976055P3	VUEPOINT IL CONN. 3.5-5.5MM INLINE ROD
8976065P3	VUEPOINT IL CONN. 3.5-6.35MM INLINE ROD
8976078P	VUEPOINT I COCR ROD. Ø3.5X60MM STRAIGHT
8976079P	VUEPOINT IL COCR ROD. Ø3.5X65MM STRAIGHT
8976080P	VUEPOINT II COCR ROD. Ø3.5X70MM STRAIGHT
8976081P	VUEPOINT II COCR ROD. Ø3.5X75MM STRAIGHT
8976082P	VUEPOINT II COCR ROD. Ø3.5X80MM STRAIGHT
8976083P	VUEPOINT II COCR ROD. Ø3.5X90MM STRAIGHT
8976084P	VUEPOINT IL COCR ROD. Ø3.5X100MM STR
8976085P	VUEPOINT IL COCR ROD. Ø3.5X110MM STR
8976086P	VUEPOINT IL COCR ROD. Ø3.5X30MM CURVED
8976087P	VUEPOINT IL COCR ROD. Ø3.5X35MMN CURVED
8976088P	VUEPOINT II COCR ROD. Ø3.5X40MM CURVED
8976089P	VUEPOINT IL COCR ROD. Ø3.5X45MM CURVED
8976090P	VUEPOINT IL COCR ROD. Ø3.5X50MM CURVED
8976091P	VUEPOINT II COCR ROD. Ø3.5X55MM CURVED
8976092P	VUEPOINT II COCR ROD. Ø3.5X60MM CURVED
8976093P	VUEPOINT II COCR ROD. Ø3.5X65MM CURVED
8976094P	VUEPOINT II COCR ROD. Ø3.5X70MM CURVED
8976095P	VUEPOINT IL COCR ROD. Ø3.5X75MM CURVED
8976096P	VUEPOINT II COCR ROD. Ø3.5X80MM CURVED
8976097P	VUEPOINT II COCR ROD. Ø3.5X90MM CURVED
8976098P	VUEPOINT II COCR ROD. Ø3.5X100MM CURVED
8976099P	VUEPOINT II COCR ROD. Ø3.5X110MM CURVED
8976104P	VUEPOINT IL ROD. Ø3.5X40MM
8976105P	VUEPOINT IL ROD. Ø3.5X50MM
8976106P	VUEPOINT IL ROD. 3.5X60MM
8976107P	VUEPOINT IL ROD. Ø3.5X70MM
8976108P	VUEPOINT IL TI ROD. Ø3.5X80MM STRAIGHT
8976109P	VUEPOINT IL ROD. Ø3.5X90MM
8976110P	VUEPOINT IL ROD. Ø3.5X100MM CONTOURED
8976111P	VUEPOINT IL ROD. Ø3.5X110MM
8976112P	VUEPOINT IL ROD. 3.5X120MM
8976124P	VUEPOINT IL ROD. 3.5X240MM
8976126P3	VUEPOINT IL CONNECTOR. 26MM FIXED CROSS



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8976128P3	VUEPOINT IL CONNECTOR. 28MM FIXED CROSS
8976130P3	VUEPOINT IL CONNECTOR. 30MM FIXED CROSS
8976132P3	VUEPOINT IL CONNECTOR. 32MM FIXED CROSS
8976134P3	VUEPOINT IL CONNECTOR. 34MM FIXED CROSS
8976136P3	VUEPOINT LI CONNECTOR. 36MM FIXED CROSS
8976138P3	VUEPOINT CONNECTOR. 38MM FIXED CROSS
8976141P	VUEPOINT II TI ROD. Ø3.5X110MM CURVED
8976142P	VUEPOINT II TI ROD. Ø3.5X120MM CURVED
8976143P	VUEPOINT II TI ROD. Ø3.5X30MM CURVED
8976144P	VUEPOINT IL TI ROD. Ø3.5X40MM CURVED
8976145P	VUEPOINT IL TI ROD. Ø3.5X50MM CURVED
8976146P	VUEPOINT IL TI ROD. Ø3.5X60MM CURVED
8976147P	VUEPOINT II TI ROD. Ø3.5X70MM CURVED
8976149P	VUEPOINT IL TI ROD. Ø3.5X90MM CURVED
8976180P	VUEPOINT IL ROD. 80MM PRE-CANTOURED
8976240P	VUEPOINT IL ROD 240MM 45° PRECURVED
8976300P	VUEPOINT II ROD. HINGED
8976305P3	VUEPOINT II EYELETS. LOW PROFILE
8976345P	VUEPOINT LL TI ROD. 3.5-4.5MM TRANSITION
8976355P	VUEPOINT II ROD. 3.5MM-5.5MM TRANSITION
8976406P3	VUEPOINT II SCREW. 4.5X6MM OCCIPITAL
8976408P3	VUEPOINT LI SCREW. 4.5X8MM OCCIPITAL
8976410P3	VUEPOINT LI SCREW. 4.5X10MM OCCIPITAL
8976412P3	VUEPOINT IL SCREW. 4.5X12MM OCCIPITAL
8976414P3	VUEPOINT IL SCREW. 4.5X14MM OCCIPITAL
8976416P3	VUEPOINT IL SCREW. 4.5X16MM OCCIPITAL
8976418P3	VUEPOINT IL SCREW. 4.5X18MM OCCIPITAL
8976420P3	VUEPOINT II SCREW. 4.5X20MM OCCIPITAL
8976435P3	VUEPOINT II CONN. 3.5MM-3.5MM OFFSETROD
8976455P2	VUEPOINT II CONN. 3.5MM-5 5MM OFFSETROD
8976455P3	VUEPOINT IL CONN. 3.5MM-5.5MM OFFSETROD
8976506P3	VUEPOINT IL SCREW. 5.0X6MM OCCIPITAL
8976508P3	VUEPOINT IL SCREW. 5.0X8MM OCCIPITAL
8976510P3	VUEPOINT IL SCREW. 5.0X10MM OCCIPITAL
8976512P3	VUEPOINT IL SCREW. 5.0X12MM OCCIPITAL
8976514P3	VUEPOINT IL SCREW. 5.0X14MM OCCIPITAL
8976516P3	VUEPOINT IL SCREW. 6.0X16MM OCCIPITAL
8976518P3	VUEPOINT IL SCREW. 5.0X18MM OCCIPITAL
8976520P3	VUEPOINT IL SCREW. 5.0X20MM OCCIPITAL



República de Colombia  
Ministerio de Salud y Protección Social  
Instituto Nacional de Vigilancia de Medicamentos y Alimentos – INVIMA

**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**  
**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

CÓDIGO, MODELO O REFERENCIA	DESCRIPCIÓN
8977760P3	VUEPOINT IL SET SCREW. OFFSET OCCIPITAL
8977765P3	VUEPOINT IL SET SCREW. ROD-ROD CONN
8977770P3	VUEPOINT IL SET SCREW. STATIC CROSS CONN
8972360P	VUEPOINT II CROSS CONN. ADJ TULIP-TULIP
8976011P	VUEPOINT IL CONN. 11MMX0° LATERAL OFFSET
8976015P	VUEPOINT II CONN. 15MMX0° LATERAL OFFSET
8976020P	VUEPOINT IL CONN. 20MMX0° LATERAL OFFSET
8976025P	VUEPOINT IL CONN. 25MMX0° LATERAL OFFSET
8976030P	VUEPOINT IL CONN. 30MMX45° LATERALOFFSET
8976031P	VUEPOINT IL CONN. 30MMX35° LATERALOFFSET
8976032P	VUEPOINT IL CONN. 30MMX55° LATERALOFFSET
8976033P	VUEPOINT IL CONN. 30MMX65° LATERALOFFSET
8976035P	VUEPOINT IL CONN. 3.5-3.5MM INLINE ROD
8976055P	VUEPOINT IL CONN. 3.5-5.5MM INLINE ROD
8976126P	VUEPOINT IL CONNECTOR. 26MM FIXED CROSS
8976128P	VUEPOINT II CONNECTOR. 28MM FIXED CROSS
8976130P	VUEPOINT IL CONNECTOR. 30MM FIXED CROSS
8976132P	VUEPOINT IL CONNECTOR. 32MM FIXED CROSS
8976134P	VUEPOINT IL CONNECTOR. 34MM FIXED CROSS
8976136P	VUEPOINT IL CONNECTOR. 36MM FIXED CROSS
8976138P	VUEPOINT IL CONNECTOR. 38MM FIXED CROSS
8976305P	VUEPOINT II EYELETS. LOW PROFILE
8976406P	VUEPOINT IL SCREW. 4.5X6MM OCCIPITAL
8976408P	VUEPOINT SCREW. 4.5X8MM OCCIPITAL
8976410P	VUEPOINT II SCREW. 4.5X10MM OCCIPITAL
8976412P	VUEPOINT II SCREW. 4.5X12MM OCCIPITAL
8976414P	VUEPOINT IL SCREW. 4.5X14MM OCCIPITAL
8976416P	VUEPOINT II SCREW. 4.5X16MM OCCIPITAL
8976418P	VUEPOINT II SCREW. 4.5X18MM OCCIPITAL
8976420P	VUEPOINT LI SCREW. 4.5X20MM OCCIPITAL
8976435P	VUEPOINT II CONN. 3.5MM-3.5MM OFFSETROD
8976455P	VUEPOINT IL CONN. 3.5MM-5.5MM OFFSETROD
8976506P	VUEPOINT IL SCREW. 5.0X6MM OCCIPITAL
8976508P	VUEPOINT II SCREW. 5.0X8MM OCCIPITAL
8976510P	VUEPOINT II SCREW. 5.0X10MM OCCIPITAL
8976512P	VUEPOINT II SCREW. 5.0X12MM OCCIPITAL
8976514P	VUEPOINT LL SCREW. 5.0X14MM OCCIPITAL
8976516P	VUEPOINT LI SCREW. 5.0X16MM OCCIPITAL
8976518P	VUEPOINT IL SCREW. 5.0X18MM OCCIPITAL



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8976520P	VUEPOINT IL SCREW. 5.0X20MM OCCIPITAL
8977760P	VUEPOINT IL SET SCREW. OFFSET OCCIPITAL
8977765P	VUEPOINT IL SET SCREW. ROD-ROD CONN
8977770P	VUEPOINT IL SET SCREW. STATIC CROSS CONN
8980000P	VUEPOINT IL KEEL PLATE + SET SCREW
8980008P3	VUEPOINT IL SCREW/SET SCREW 3.5X8MM MA
8980010P3	VUEPOINT II SCREW/SET SCREW. 3.5X 10MM MA
8980012P3	VUEPOINT IL SCREW/SET SCREW. 3.5X12MM MA
8980014P3	VUEPOINT IL SCREW/SET SCREW. 3.5X14MM MA
8980016P3	VUEPOINT LL SCREW/SET SCREW. 3.5X16MM MA
8980018P3	VUEPOINT IL SCREW/SET SCREW. 3.5X18MM MA
8980020P3	VUEPOINT IL SCREW/SET SCREW. 3.5X20MM MA
8980022P3	VUEPOINT LL SCREW/SET SCREW. 3.5X22MM MA
8980024P3	VUEPOINT IL SCREW/SET SCREW. 3.5X24MM MA
8980026P3	VUEPOINT IL SCREW/SET SCREW. 3.5X26MM MA
8980028P3	VUEPOINT LL SCREW/SET SCREW. 3.5X28MM MA
8980030P3	VUEPOINT IL SCREW/SET SCREW. 3.5X30MM MA
8980032P3	VUEPOINT LL SCREW/SET SCREW. 3.5X32MM MA
8980034P3	VUEPOINT IL SCREW/SET SCREW. 3.5X34MM MA
8980036P3	VUEPOINT IL SCREW/SET SCREW. 3.5X36MM MA
8980038P3	VUEPOINT IL SCREW/SET SCREW. 3.5X38MM MA
8980040P3	VUEPOINT IL SCREW/SET SCREW. 3.5X40MM MA
8980042P3	VUEPOINT IL SCREW/SET SCREW. 3.5X42MM MA
8980044P3	VUEPOINT IL SCREW/SET SCREW. 3.5X44MM MA
8980046P3	VUEPOINT II SCREW/SET SCREW. 3.5X46MM MA
8980048P3	VUEPOINT IL SCREW/SET SCREW. 3.5X48MM MA
8980050P3	VUEPOINT IL SCREW/SET SCREW. 3.6X50MM MA
8980108P3	VUEPOINT IL SCREW/SET SCREW. 4.6X8MM MA
8980110P3	VUEPOINT IL SCREW/SET SCREW. 4.0X10MM MA
8980112P3	VUEPOINT IL SCREW/SET SCREW. 4.0X12MM MA
8980114P3	VUEPOINT IL SCREW/SET SCREW. 4.0X14MM MA
8980116P3	VUEPOINT IL SCREW/SET SCREW. 4.0X16MM MA
8980118P3	VUEPOINT IL SCREW/SET SCREW. 4.0X18MM MA
8980120P3	VUEPOINT IL SCREW/SET SCREW. 4.0X20MM MA
8980122P3	VUEPOINT IL SCREW/SET SCREW. 4.0X22MM MA
8980124P3	VUEPOINT LL SCREW/SET SCREW. 4.0X24MM MA
8980126P3	VUEPOINT LI SCREW/SET SCREW. 4.0X26MM MA
8980128P3	VUEPOINT II SCREW/SET SCREW. 4.0X28MM MA
8980130P3	VUEPOINT LI SCREW/SET SCREW. 4.0X30MM MA



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8980132P3	VUEPOINT LI SCREW/SET SCREW. 4.0X32MM MA
8980134P3	VUEPOINT II SCREW/SET SCREW. 4.0X34MM MA
8980136P3	VUEPOINT II SCREW/SET SCREW. 4.0X36MM MA
8980138P3	VUEPOINT IL SCREW/SET SCREW. 4.0X38MM MA
8980140P3	VUEPOINT II SCREW/SET SCREW. 4.0X40MM MA
8980142P3	VUEPOINT IL SCREW/SET SCREW. 4.0X42MM MA
8980144P3	VUEPOINT II SCREW/SET SCREW. 4.0X44MM MA
8980146P3	VUEPOINT IL SCREW/SET SCREW. 4.0X46MM MA
8980148P3	VUEPOINT LL SCREW/SET SCREW. 4.0X48MM MA
8980150P3	VUEPOINT IL SCREW/SET SCREW. 4.0X50MM MA
8980208P3	VUEPOINT IL SCREW/SET SCREW. 3.5X8MM FA
8980210P3	VUEPOINT IL SCREW/SET SCREW. 3.5X10MM FA
8980212P3	VUEPOINT IL SCREW/SET SCREW. 3.5X12MM FA
8980214P3	VUEPOINT II SCREW/SET SCREW. 3.5X14MM FA
8980216P3	VUEPOINT II SCREW/SET SCREW. 3.5X 16MM FA
8980218P3	VUEPOINT IL SCREW/SET SCREW. 3.5X18MM FA
8980220P3	VUEPOINT IL SCREW/SET SCREW. 3.5X20MM FA
8980222P3	VUEPOINT IL SCREW/SET SCREW. 3.5X22MM FA
8980224P3	VUEPOINT II SCREW/SET SCREW. 3.5X24MM FA
8980226P3	VUEPOINT II SCREW/SET SCREW. 3.5X26MM FA
8980228P3	VUEPOINT II SCREW/SET SCREW. 3.5X28MM FA
8980230P3	VUEPOINT IL SCREW/SET SCREW. 3.5X30MM FA
8980232P3	VUEPOINT IL SCREW/SET SCREW. 3.5X32MM FA
8980234P3	VUEPOINT IL SCREW/SET SCREW. 3.5X34MM FA
8980236P3	VUEPOINT IL SCREW/SET SCREW. 3.5X36MM FA
8980308P3	VUEPOINT IL SCREW/SET SCREW. 4.0X8MM FA
8980310P3	VUEPOINT IL SCREW/SET SCREW. 4.0X10MM FA
8980312P3	VUEPOINT IL SCREW/SET SCREW. 4.0X12MM FA
8980314P3	VUEPOINT IL SCREW/SET SCREW. 4.0X14MM FA
8980316P3	VUEPOINT IL SCREW/SET SCREW. 4.0X16MM FA
8980318P3	VUEPOINT IL SCREW/SET SCREW. 4.0X18MM FA
8980320P3	VUEPOINT IL SCREW/SET SCREW. 4.0X20MM FA
8980322P3	VUEPOINT II SCREW/SET SCREW. 4.0X22MM FA
8980324P3	VUEPOINT IL SCREW/SET SCREW. 4.0X24MM FA
8980326P3	VUEPOINT IL SCREW/SET SCREW. 4.0X26MM FA
8980328P3	VUEPOINT IL SCREW/SET SCREW. 4.0X28MM FA
8980330P3	VUEPOINT IL SCREW/SET SCREW. 4.0X30MM FA
8980332P3	VUEPOINT IL SCREW/SET SCREW. 4.0X32MM FA
8980334P3	VUEPOINT LL SCREW/SET SCREW. 4.0X34MM FA



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8980336P3	VUEPOINT IL SCREW/SET SCREW. 4.0X36MM FA
8980338P3	VUEPOINT IL SCREW/SET SCREW. 4.0X38MM FA
8980340P3	VUEPOINT II SCREW/SET SCREW. 4.0X40MM FA
8980342P3	VUEPOINT LI SCREW/SET SCREW. 4.0X42MM FA
8980344P3	VUEPOINT IL SCREW/SET SCREW. 4.0X44MM FA
8980346P3	VUEPOINT II SCREW/SET SCREW. 4.0X46MM FA
8980408P3	VUEPOINT IL SCREW/SET SCREW. 4.5X8MM MA
8980410P3	VUEPOINT IL SCREW/SET SCREW. 4.5X10MM MA
8980412P3	VUEPOINT IL SCREW/SET SCREW. 4.5X12MM MA
8980414P3	VUEPOINT IL SCREW/SET SCREW. 4.5X14MM MA
8980416P3	VUEPOINT IL SCREW/SET SCREW. 4.5X16MM MA
8980418P3	VUEPOINT LL SCREW/SET SCREW. 4.5X 18MM MA
8980420P3	VUEPOINT IL SCREW/SET SCREW. 4.5X20MM MA
8980422P3	VUEPOINT LL SCREW/SET SCREW. 4.5X22MM MA
8980424P3	VUEPOINT IL SCREW/SET SCREW. 4.5X24MM MA
8980426P3	VUEPOINT II SCREW/SET SCREW. 4.5X26MM MA
8980428P3	VUEPOINT IL SCREW/SET SCREW. 4.5X28MM MA
8980430P3	VUEPOINT II SCREW/SET SCREW. 4.5X30MM MA
8980432P3	VUEPOINT IL SCREW/SET SCREW. 4.5X32MM MA
8980434P3	VUEPOINT IL SCREW/SET SCREW. 4.5X34MM MA
8980436P3	VUEPOINT IL SCREW/SET SCREW. 4.5X36MM MA
8980438P3	VUEPOINT IL SCREW/SET SCREW. 4.5X38MM MA
8980440P3	VUEPOINT IL SCREW/SET SCREW. 4.5X40MM MA
8980442P3	VUEPOINT IL SCREW/SET SCREW. 4.5X42MM MA
8980444P3	VUEPOINT IL SCREW/SET SCREW. 4.5X44MM MA
8980446P3	VUEPOINT IL SCREW/SET SCREW. 4.5XD6MM MA
8980448P3	VUEPOINT LL SCREW/SET SCREW. 4.5X48MM MA
8980450P3	VUEPOINT IL SCREW/SET SCREW. 4.5X50MM MA
8980518P3	VP II SCREW/SET SCREW. 3.5X18MM FA P.T.
8980520P3	VP IL SCREW/SET SCREW. 3.5X20MM FA P.T.
8980522P3	VP II SCREW/SET SCREW. 3.5X22MM FA PT.
8980524P3	VP II SCREW/SET SCREW. 3.5X24MM FA PT.
8980526P3	VP IL SCREW/SET SCREW. 3.5X26MM FA PT
8980528P3	VP IL SCREW/SET SCREW. 3.5X28MM FA P.T.
8980530P3	VP II SCREW/SET SCREW. 3.5X30MM FA P.T.
8980532P3	VP II SCREW/SET SCREW. 3.5X32MM FA P.T.
8980534P3	VP IL SCREW/SET SCREW. 3.5X34MM FA P.T.
8980536P3	VP II SCREW/SET SCREW. 3.5X36MM FA P.T.
8980538P3	VP LL SCREW/SET SCREW. 3.5X38MM FA P.T.



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8980540P3	VP IL SCREW/SET SCREW. 3.5X40MM FA P.T.
8980700P3	VUEPOINT IL CONN/SET SCREW. EXTRA ROD
8980008P	VUEPOINT II SCREW/SET SCREW. 3.5X8MM MA
8980010P	VUEPOINT II SCREW/SET SCREW. 3.5X10MM MA
8980012P	VUEPOINT IL SCREW/SET SCREW. 3.5X12MM MA
8980014P	VUEPOINT II SCREW/SET SCREW. 3.5X14MM MA
8980016P	VUEPOINT II SCREW/SET SCREW. 3.5X16MM MA
8980018P	VUEPOINT II SCREW/SET SCREW. 3.5X18MM MA
8980020P	VUEPOINT IL SCREW/SET SCREW. 3.5X20MM MA
8980022P	VUEPOINT IL SCREW/SET SCREW. 3.5X22MM MA
8980024P	VUEPOINT IL SCREW/SET SCREW. 3.5X24MM MA
8980026P	VUEPOINT TL SCREW/SET SCREW. 3.5X26MM MA
8980028P	VUEPOINT IL SCREW/SET SCREW. 3.5X28MM MA
8980030P	VUEPOINT IL SCREW/SET SCREW. 3.5X30MM MA
8980032P	VUEPOINT IL SCREW/SET SCREW. 3.5X32MM MA
8980034P	VUEPOINT IL SCREW/SET SCREW. 3.5X34MM MA
8980036P	VUEPOINT IL SCREW/SET SCREW. 3.5X36MM MA
8980038P	VUEPOINT IL SCREW/SET SCREW. 3.5X38MM MA
8980040P	VUEPOINT IL SCREW/SET SCREW. 3.5X40MM MA
8980042P	VUEPOINT.II SCREW/SET SCREW. 3.5X42MM MA
8980044P	VUEPOINT II SCREW/SET SCREW. 3.5X44MM MA
8980046P	VUEPOINT IL SCREW/SET SCREW. 3.5X46MM MA
8980048P	VUEPOINT IL SCREW/SET SCREW. 3.5X48MM MA
8980050P	VUEPOINT IL SCREW/SET SCREW. 3.5X50MM MA
8980108P	VUEPOINT II SCREW/SET SCREW. 4.0X8MM MA
8980110P	VUEPOINT IL SCREW/SET SCREW. 4.0X10MM MA
8980112P	VUEPOINT IL SCREW/SET SCREW. 4.0X12MM MA
8980114P	VUEPOINT IL SCREW/SET SCREW. 4.0X14MM MA
8980116P	VUEPOINT IL SCREW/SET SCREW. 4.0X16MM MA
8980118P	VUEPOINT IL SCREW/SET SCREW. 4.0X18MM MA
8980120P	VUEPOINT IL SCREW/SET SCREW. 4.0X20MM MA
8980122P	VUEPOINT IL SCREW/SET SCREW. 4.0X22MM MA
8980124P	VUEPOINT IL SCREW/SET SCREW. 4.0X24MM MA
8980126P	VUEPOINT IL SCREW/SET SCREW. 4.0X26MM MA
8980128P	VUEPOINT IL SCREW/SET SCREW. 4.0X28MM MA
8980130P	VUEPOINT IL SCREW/SET SCREW. 4.0X30MM MA
8980132P	VUEPOINT IL SCREW/SET SCREW. 4.0X32MM MA
8980134P	VUEPOINT IL SCREW/SET SCREW. 4.0X34MM MA
8980136P	VUEPOINT II SCREW/SET SCREW. 4.0X36MM MA



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

<b>CÓDIGO, MODELO O REFERENCIA</b>	<b>DESCRIPCIÓN</b>
8980138P	VUEPOINT II SCREW/SET SCREW. 4.0X38MM MA
8980140P	VUEPOINT II SCREW/SET SCREW. 4.0X40MM MA
8980142P	VUEPOINT II SCREW/SET SCREW. 4.0X42MM MA
8980144P	VUEPOINT II SCREW/SET SCREW. 4.0X44MM MA
8980146P	VUEPOINT IL SCREW/SET SCREW. 4.0X46MM MA
8980148P	VUEPOINT IL SCREW/SET SCREW. 4.0X48MM MA
8980150P	VUEPOINT IL SCREW/SET SCREW. 4.0X50MM MA
8980208P	VUEPOINT II SCREW/SET SCREW. 3.5X8MM FA
8980210P	VUEPOINT IL SCREW/SET SCREW. 3.5X10MM FA
8980212P	VUEPOINT IL SCREW/SET SCREW. 3.5X12MM FA
8980214P	VUEPOINT IL SCREW/SET SCREW. 3.5X14MM FA
8980216P	VUEPOINT IL SCREW/SET SCREW. 3.5X16MM FA
8980218P	VUEPOINT II SCREW/SET SCREW. 3.5X18MM FA
8980220P	VUEPOINT IL SCREW/SET SCREW. 3.5X20MM FA
8980222P	VUEPOINT IL SCREW/SET SCREW. 3.5X22MM FA
8980224P	VUEPOINT II SCREW/SET SCREW. 3.5X24MM FA
8980226P	VUEPOINT II SCREW/SET SCREW. 3.5X26MM FA
8980228P	VUEPOINT IL SCREW/SET SCREW. 3.5X28MM FA
8980230P	VUEPOINT II SCREW/SET SCREW. 3.5X30MM FA
8980232P	VUEPOINT IL SCREW/SET SCREW. 3.5X32MM FA
8980234P	VUEPOINT IL SCREW/SET SCREW. 3.5X34MM FA
8980236P	VUEPOINT IL SCREW/SET SCREW. 3.5X36MM FA
8980308P	VUEPOINT II SCREW/SET SCREW. 4.0X8MM FA
8980310P	VUEPOINT IL SCREW/SET SCREW. 4.0X10MM FA
8980312P	VUEPOINT IL SCREW/SET SCREW. 4.0X12MM FA
8980314P	VUEPOINT II SCREW/SET SCREW. 4.0X14MM FA
8980316P	VUEPOINT IL SCREW/SET SCREW. 4.0X16MM FA
8980318P	VUEPOINT IL SCREW/SET SCREW. 4.0X18MM FA
8980320P	VUEPOINT IL SCREW/SET SCREW. 4.0X20MM FA
8980322P	VUEPOINT II SCREW/SET SCREW. 4.0X22MM FA
8980324P	VUEPOINT II SCREW/SET SCREW. 4.0X24MM FA
8980326P	VUEPOINT II SCREW/SET SCREW. 4.0X26MM FA
8980328P	VUEPOINT II SCREW/SET SCREW. 4.0X28MM FA
8980330P	VUEPOINT II SCREW/SET SCREW. 4.0X30MM FA
8980332P	VUEPOINT IL SCREW/SET SCREW. 4.0X32MM FA
8980334P	VUEPOINT II SCREW/SET SCREW. 4.0X34MM FA
8980336P	VUEPOINT IL SCREW/SET SCREW. 4.0X36MM FA
8980338P	VUEPOINT II SCREW/SET SCREW. 4.0X38MM FA
8980340P	VUEPOINT II SCREW/SET SCREW. 4.0X40MM FA



**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**

**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

**CÓDIGO, MODELO O  
REFERENCIA**

**DESCRIPCIÓN**

8980342P	VUEPOINT IL SCREW/SET SCREW. 4.0X42MM FA
8980344P	VUEPOINT II SCREW/SET SCREW. 4.0X44MM FA
8980346P	VUEPOINT II SCREW/SET SCREW. 4.0X46MM FA
8980408P	VUEPOINT II SCREW/SET SCREW. 4.5X8MM MA
8980410P	VUEPOINT LL SCREW/SET SCREW. 4.5X10MM MA
8980412P	VUEPOINT II SCREW/SET SCREW. 4.5X12MM MA
8980414P	VUEPOINT II SCREW/SET SCREW. 4.5X14MM MA
8980416P	VUEPOINT IISCREW/SET SCREW. 4.5X16MM MA
8980418P	VUEPOINT II SCREW/SET SCREW. 4.5X18MM MA
8980420P	VUEPOINT II SCREW/SET SCREW. 4.5X20MM MA
8980422P	VUEPOINT IL SCREW/SET SCREW 4.5X22MM MA
8980424P	VUEPOINT IL SCREW/SET SCREW. 4.5X24MM MA
8980426P	VUEPOINT IL SCREW/SET SCREW. 4.5X26MM MA
8980428P	VUEPOINT IL SCREW/SET SCREW. 4.5X28MM MA
8980430P	VUEPOINT LL SCREW/SET SCREW. 4.5X30MM MA
8980432P	VUEPOINT IL SCREW/SET SCREW. 4.5X32MM MA
8980434P	VUEPOINT IL SCREW/SET SCREW 4.5X34MM MA
8980436P	VUEPOINT IL SCREW/SET SCREW. 4.5X36MM MA
8980438P	VUEPOINT IL SCREW/SET SCREW. 4.5X38MM MA
8980440P	VUEPOINT IL SCREW/SET SCREW. 4.5X40MM MA
8980442P	VUEPOINT IL SCREW/SET SCREW. 4.5X42MM MA
8980444P	VUEPOINT IL SCREW/SET SCREW 4.5X44MM MA
8980446P	VUEPOINT IL SCREW/SET SCREW. 4.5X46MM MA
8980448P	VUEPOINT IL SCREW/SET SCREW. 4.5X48MM MA
8980450P	VUEPOINT IL SCREW/SET SCREW. 4.5X50MM MA
8980518P	VP IL SCREW/SET SCREW. 3.5X18MM FA P.T.
8980520P	VP IL SCREW/SET SCREW. 3.5X20MM FA P.T.
8980522P	VP IL SCREW/SET SCREW. 3.5X22MM FA P.T.
8980524P	VP II SCREW/SET SCREW. 3.5X24MM FA P.T.
8980526P	VP IL SCREW/SET SCREW. 3.5X26MM FA P.T.
8980528P	VP II SCREW/SET SCREW. 3.5X28MM FA P.T.
8980530P	VP II SCREW/SET SCREW. 3.5X30MM FA P.T.
8980532P	VP II SCREW/SET SCREW. 3.5X32MM FA P.T.
8980534P	VP II SCREW/SET SCREW. 3.5X34MM FA P.T.
8980536P	VP II SCREW/SET SCREW. 3.5X36MM FA P.T.
8980538P	VP II SCREW/SET SCREW. 3.5X38MM FA PT
8980540P	VP IL SCREW/SET SCREW. 3.5X40MM FA PT.
8980700P	VUEPOINT II CONN/SET SCREW. EXTRA ROD

VIDA UTIL:

NO APLICA



La salud  
es de todos

Minsalud

República de Colombia  
Ministerio de Salud y Protección Social  
Instituto Nacional de Vigilancia de Medicamentos y Alimentos – INVIMA

**RESOLUCIÓN No. 2022005281 DE 25 de Marzo de 2022**  
**Por la cual se concede un Registro Sanitario**

La Directora Técnica de Dispositivos Médicos y otras Tecnologías del Instituto Nacional de Vigilancia de Medicamentos y Alimentos INVIMA, en ejercicio de las facultades Legales conferidas en el Decreto 2078 de 2012, Decreto Reglamentario 4725 de 2005 modificado por el Decreto 582 de 2017, y ley 1437 de 2011.

EXPEDIENTE No.: 20203135  
RADICACIÓN No.: 20211100250  
FECHA RADICADO INICIAL: 24/05/2021

**ARTICULO SEGUNDO.** - Se aprueba etiqueta del fabricante y sticker del importador aportadas mediante radicado 20211254942 del 30 de noviembre de 2021.

**ARTICULO TERCERO.-** Contra la presente resolución procede únicamente el Recurso de Reposición, que deberá interponerse ante la Dirección de Dispositivos Médicos y otras Tecnologías del INVIMA, dentro de los DIEZ (10) días siguientes a su notificación, en los términos señalados en el Código de Procedimiento Administrativo y de lo Contencioso Administrativo.

**ARTÍCULO CUARTO.-** La presente resolución rige a partir de la fecha de su expedición.

**COMUNIQUESE, NOTIFIQUESE Y CUMPLASE**

Dada en Bogotá D.C. a los 25 de Marzo de 2022

**LUCIA AYALA RODRIGUEZ**  
**DIRECTORA TÉCNICA DE DISPOSITIVOS MÉDICOS Y OTRAS TECNOLOGÍAS**  
Proyectó: Legal: anaranjoa, Técnico: agomezs Revisó: cordina\_varios